

QUICK LINKS

SURGICAL GOWNS

MICROCOOL*	page 3
AERO BLUE*	page 4
Surgical Gowns	page 5

SURGICAL DRAPES & PACKS

Cardiac Catheterisation / Radiology Drapes	page 6
Cardiovascular Drapes	page 8
ENT Drapes & Packs	page 10
Equipment Covers	page 12
General Drapes & Packs:	
Accessories	page 13
Basic Packs	page 14
Basic Drapes	page 15
Breast / Chest Drapes	page 17
Fenestrated Drapes	page 18
Laparoscopy / Laparotomy Drapes	page 19
Minor Procedure Drapes	page 21
Thyroid Drapes	page 22
Universal Packs	page 23
Universal Drapes	page 26
Neuro / Spine Drape	page 28
Obstetrics and Gynaecology Drapes & Packs:	
C-Section Drapes	page 29
Laparoscopy Packs	page 30
Laparoscopy Drapes	page 31
Lithotomy Packs	page 32
Lithotomy Drapes	page 33
Under Buttocks Drapes	page 34
Ophthalmic Drapes & Packs	page 35
Orthopaedic Drapes & Packs:	
Ortho Packs	page 36
Extremity Drapes	page 38
Hand Drapes	page 40
Hip Drapes	page 41
Knee Arthroscopy Drapes	page 42
Shoulder Drapes	page 43
U & Bar Drapes	page 45
Urology Drapes & Packs	page 49
Microbial Sealant	page 50
Index	page 51


SURGICAL GOWNS

MICROCOOL* BREATHABLE HIGH PERFORMANCE SURGICAL GOWNS WITH SECURE-FIT*

- Innovative patented glove-gown interface
- Patented Breathable Fabric
- Low Linting Fabric
- Raglan Sleeves
- Colour-Coded Neck Binding – Red
- Adjustable Hook & Loop Neckline
- Book Fold
- 1 x Blue Huck Towel is included for all sterile codes
- AAMI Level 4 / EN 13795 High Performance



Code	Description	Packaging
13359	MICROCOOL* Breathable High Performance Surgical Gown with SECURE-FIT*, Sterile, Small	Items/Case = 30
92320	MICROCOOL* Breathable High Performance Surgical Gown with SECURE-FIT*, Sterile, Large	Items/Case = 30
72454	MICROCOOL* Breathable High Performance Surgical Gown with SECURE-FIT*, Non-Sterile, Large	Items/Case = 40
92321	MICROCOOL* Breathable High Performance Surgical Gown with SECURE-FIT*, Sterile, X-Large	Items/Case = 28
30734	MICROCOOL* Breathable High Performance Surgical Gown with SECURE-FIT*, Non-Sterile, X-Large	Items/Case = 384
92323	MICROCOOL* Breathable High Performance Surgical Gown with SECURE-FIT*, Sterile, XX-Large	Items/Case = 26
72448	MICROCOOL* Breathable High Performance Surgical Gown with SECURE-FIT*, Non-Sterile, XX-Large	Items/Case = 32



MICROCOOL* BREATHABLE HIGH PERFORMANCE SPECIALTY SURGICAL GOWNS WITH SECURE-FIT*

- Innovative patented glove-gown interface
- Patented Breathable Fabric
- Low Linting Fabric
- Patented A-Line Design for Seated Procedures
- Raglan Sleeves
- Colour-Coded Neck Binding – Red
- Adjustable Hook & Loop Neckline
- Book Fold
- 1 x Blue Huck Towel is included for all sterile codes
- AAMI Level 4 / EN 13795 High Performance



Code	Description	Packaging
92330	MICROCOOL* Breathable High Performance Specialty Surgical Gown with SECURE-FIT*, Sterile, Large	Items/Case = 28
92331	MICROCOOL* Breathable High Performance Specialty Surgical Gown with SECURE-FIT*, Sterile, X-Large	Items/Case = 26
92333	MICROCOOL* Breathable High Performance Specialty Surgical Gown with SECURE-FIT*, Sterile, XX-Large	Items/Case = 26

-  Non-Reinforced Fabric
-  Fabric-Reinforced
-  High Performance

All codes in **bold** indicate a sterile product.

NOTE: Please see package inserts for directions for use for each of these products.

SURGICAL GOWNS



AERO BLUE* PERFORMANCE SURGICAL GOWNS

- Patented Breathable Fabric
- Low Linting Fabric
- Raglan Sleeves
- Colour-Coded Neck Binding – Green
- Adjustable Hook & Loop Neckline
- Book Fold
- 1 x Blue Huck Towel is included for all sterile codes
- AAMI Level 3 / EN 13795 High Performance

Code	Description	Packaging
41870	AERO BLUE* Performance Surgical Gown, Sterile, Small	Items/Case = 30
43941	AERO BLUE* Performance Surgical Gown, Non-Sterile, Small	Items/Case = 36
41871	AERO BLUE* Performance Surgical Gown, Sterile, Large	Items/Case = 30
41739	AERO BLUE* Performance Surgical Gown, Non-Sterile, Large	Items/Case = 40
42620	AERO BLUE* Performance Surgical Gown, Sterile, Large, X-Long	Items/Case = 28
41725	AERO BLUE* Performance Surgical Gown, Non-Sterile, Large, X-Long	Items/Case = 36
41872	AERO BLUE* Performance Surgical Gown, Sterile, X-Large	Items/Case = 28
41740	AERO BLUE* Performance Surgical Gown, Non-Sterile, X-Large	Items/Case = 32
42621	AERO BLUE* Performance Surgical Gown, Sterile, X-Large, X-Long	Items/Case = 30
41727	AERO BLUE* Performance Surgical Gown, Non-Sterile, X-Large, X-Long	Items/Case = 32
41873	AERO BLUE* Performance Surgical Gown, Sterile, XX-Large	Items/Case = 26
41741	AERO BLUE* Performance Surgical Gown, Non-Sterile, XX-Large	Items/Case = 32
43940	AERO BLUE* Performance Surgical Gown, Sterile, XX-Large, X-Long	Items/Case = 28
43623	AERO BLUE* Performance Surgical Gown, Non-Sterile, XX-Large, X-Long	Items/Case = 30
41874	AERO BLUE* Performance Surgical Gown, Sterile, XXX-Large	Items/Case = 24



Non-Reinforced Fabric



Fabric-Reinforced




High Performance

All codes in **bold** indicate a sterile product.

NOTE: Please see package inserts for directions for use for each of these products.

SURGICAL GOWNS



SURGICAL GOWNS

- EVOLUTION* 4 Fabric
- Low Linting Fabric
- Set-In-Sleeves
- Colour-Coded Neck Binding – Yellow
- Adjustable Hook & Loop Neckline
- Book Fold
- 1 x Blue Huck Towel is included for all sterile codes
- EN 13795 Standard Performance

Code	Description	Packaging
90003	Surgical Gown, Sterile, Small To Medium	Items/Case = 36
90018	Surgical Gown, Sterile, Large	Items/Case = 36
72978	Surgical Gown, Non-Sterile, Large	Items/Case = 640
90019	Surgical Gown, Sterile, Large, X-Long	Items/Case = 34
90048	Surgical Gown, Sterile, X-Large	Items/Case = 32
72988	Surgical Gown, Non-Sterile, X-Large	Items/Case = 640
90049	Surgical Gown, Sterile, X-Large, X-Long	Items/Case = 34
70392	Surgical Gown, Non-Sterile, X- Large, X-Long	Items/Case = 60
90050	Surgical Gown, Sterile, XX-Large	Items/Case = 28
90138	Surgical Gowns Double Pack	Items/Case = 18
Contents:		
72978	1-Surgical Gown, Large	
72988	1-Surgical Gown, X-Large	



Non-Reinforced Fabric



Fabric-Reinforced




High Performance

All codes in **bold** indicate a sterile product.

NOTE: Please see package inserts for directions for use for each of these products.


CARDIAC CATHETERISATION / RADIOLOGY DRAPE



PACEMAKER DRAPE, 221cm x 343cm

- Two Clear Side Windows 55cm x 342cm
- CONTROL* Plus Fabric Reinforcement
- Conformable Breather Bar 35cm
- Two Adhesive Fenestrations 15cm x 20cm
- Bilateral Pockets 15cm x 81cm


Code	Description	Packaging
89368	Pacemaker Drape, Sterile	Items/Case = 11
79368	Pacemaker Drape, Non-Sterile	Items/Case = 12



BRACHIAL ANGIOGRAPHY DRAPE, 111cm x 78cm

- CONTROL* Plus Fabric Reinforcement
- 7cm Adhesive Circular Fenestration

Code	Description	Packaging
89471	Brachial Angiography Drape, Sterile	Items/Case = 50



RADIAL FEMORAL ANGIOGRAPHY DRAPE, 240cm x 380cm

- Two 63cm x 350cm Windows
- Two 12cm x 7cm Oval Radial Fenestrations
- Two 12cm Round Femoral Fenestrations

Code	Description	Packaging
89583	HALYARD BASICS* Radial Femoral Angiography Drape, Sterile	Items/Case = 12

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


CARDIAC CATHETERISATION / RADIOLOGY DRAPE



EXTRA LONG FEMORAL ANGIOGRAPHY DRAPE WITH WINDOW, 208cm x 343cm

- Two 9cm Adhesive Circular Fenestrations
- One 68cm x 315cm Windows
- CONTROL* Plus Fabric Reinforcement


Code	Description	Packaging
89709	Extra Long Femoral Angiography Drape with Window, Sterile	Items/Case = 14
79709	Extra Long Femoral Angiography Drape with Window, Non-Sterile	Items/Case = 16



EXTRA LONG FEMORAL ANGIOGRAPHY DRAPE WITH WINDOWS, 221cm x 343cm

- Two 9cm Adhesive Circular Fenestrations
- Two 68cm x 343cm Windows
- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
89711	Extra Long Femoral Angiography Drape with Dual Windows, Sterile	Items/Case = 16
13296	Extra Long Femoral Angiography Drape with Dual Windows, Non-Sterile	Items/Case = 180



EXTRA-LONG FEMORAL ANGIOGRAPHY DRAPE, 220cm x 342cm

- Extra-Long 220cm x 342cm Film Based Drape
- 81cm x 342cm Absorbent Fabric Reinforcement
- 78cm x 342cm Patient Comfort Layer
- Two 8cm Adhesive Circular Fenestrations
- Two 68cm x 342cm Clear Plastic Windows

Code	Description	Packaging
79749	Extra-Long Femoral Angiography Drape, Non-Sterile	Items/Case = 144

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

CARDIOVASCULAR DRAPE




187cm x 241cm

PERINEAL DRAPE, 187cm x 241cm

- Split 44cm x 95cm
- CONTROL* Plus Fabric Reinforcement
- Positioning Tape


Code	Description	Packaging
76356	Perineal Drape, Non-Sterile	Items/Case = 17



CVARTS* HEAVY-DUTY UNIVERSAL DRAPE SYSTEM

- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
89454	CVARTS* Heavy-Duty Universal Drape System, Sterile	Items/Case = 5
Contents:		
1-CVARTS* Top Drape, 254cm x 152cm		
1-CVARTS* Bottom Drape, 193cm x 152cm		
2-CVARTS* Side Drapes, 96cm x 213cm, Longitudinal Flaps with Hook and Loop		



CVARTS* HEAVY-DUTY CV SPLIT DRAPE, 274cm x 203cm x 315cm

- Armboard Covers
- 15cm x 213cm Split
- Longitudinal Flaps with Hook and Loop
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
89452	CVARTS* Heavy-Duty CV Split Drape, Sterile	Items/Case = 6
79452	CVARTS* Heavy-Duty CV Split Drape, Non-Sterile	Items/Case = 7

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

CARDIOVASCULAR DRAPES




CVARTS* BILATERAL UNDERLEG / PERINEAL DRAPE, 213cm x 203cm

- Perineal Cover
- Umbilicus Cover
- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
89453	CVARTS* Bilateral Underleg / Perineal Drape, Sterile	Items/Case = 14
79453	CVARTS* Bilateral Underleg / Perineal Drape, Non-Sterile	Items/Case = 18


All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.

ENT DRAPES & PACKS




ENT PACK III

Code	Description	Packaging	
88223	ENT Pack III	Items/Case = 6	
Sterile Component Code	Non-Sterile Component Code	Contents	
89603	79603	1-Back Table Cover, Zone-Reinforced, 152cm x 191cm	
89431	79430	1-Mayo Stand Cover 80cm x 145cm	
89352	79351	1-Suture Bag	
		1-Head Drape 91cm x 127cm:	
		• Turban	
		• Positioning Tape	
		1-ENT Split Drape 193cm x 315cm:	
		• with 7 x 71cm Split	
		• Positioning Tape	
		• CONTROL* Fabric Reinforcement	
		• Tube Holders	



EENT PACK I

Code	Description	Packaging	
88711	EENT Pack I	Items/Case = 5	
Sterile Component Code	Non-Sterile Component Code	Contents	
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm	
95111	71311	1-ULTRA Surgical Gown, with Towel, Large	
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm	
		1-Absorbent Towel, 38cm x 56cm	
		1-Suture Bag	
89431	79430	1-Head Drape, 127cm x 91cm	
89352	79351	1-EENT Split Drape with CONTROL* Reinforcement, 193cm x 315cm with 7cm x 71cm Split	



EENT PACK II

Code	Description	Packaging	
88716	EENT Pack II	Items/Case = 5	
Sterile Component Code	Non-Sterile Component Code	Contents	
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm	
95111	71311	1-ULTRA Surgical Gown In Overwrap, with Towel, Large	
95111	71311	1-ULTRA Surgical Gown with Towel, Large	
		1-Microsurgery Mayo Stand Cover	
		1-Suture Bag	
89431	79430	1-Head Drape, 127cm x 91cm	
89352	79351	1-EENT Split Drape with CONTROL* Reinforcement, 193cm x 315cm with 7cm x 71cm Split	

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

EENT DRAPES & PACKS



EENT PACK IV


Code	Description	Packaging
88719	EENT Pack IV	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm 1-Surgical Gown with Reinforced Sleeves, In Overwrap, with Towel, Large 1-Surgical Gown, with Reinforced Sleeves, Large
89701	79720	1-Absorbent Towel, 38cm x 56cm
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
89438	79438	1-Head Drape, 104cm x 81cm 1-EENT Split Drape with Instrument Pad Reinforcement, 193cm x 315cm with 8cm x 64cm Split



EENT SPLIT DRAPE, 193cm x 315cm

- CONTROL* Fabric Reinforcement
- 7cm x 71cm Split
- Tube Holders
- Instrument Pad


Code	Description	Packaging
89352	EENT Split Drape, Sterile	Items/Case = 18
79351	EENT Split Drape, Non-Sterile	Items/Case = 25



HEAD DRAPE, 127cm x 91cm

- Turban, 83cm x 50cm


Code	Description	Packaging
89431	Head Drape, Sterile	Items/Case = 38
79430	Head Drape, Non-Sterile	Items/Case = 60



HEAD DRAPE, 104cm x 81cm

- Detachable Turban, 91cm x 61cm

Code	Description	Packaging
89438	Head Drape, Sterile	Items/Case = 38



DENTAL DRAPE, 120cm x 150cm

- Teardrop Fenestration 10cm x 15cm; 40cm From The Top
- Absorbent Reinforcement 30cm x 40cm

Code	Description	Packaging
89543	HALYARD BASICS* Dental Drape, Sterile	Items/Case = 15

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

EQUIPMENT COVERS



MAYO STAND COVERS

Code	Description	Packaging
88665	Mayo Stand Cover, Reinforced HALYARD BASICS* 58cm x 137cm, Sterile	Items/Case = 80
78656	Mayo Stand Cover, Reinforced HALYARD BASICS* 58cm x 137cm, Non-Sterile	Items/Case = 140
89601	Mayo Stand Cover, Reinforced, 58cm x 137cm, Sterile	Items/Case = 54
79607	Mayo Stand Cover, Reinforced, 58cm x 137cm, Non-Sterile	Items/Case = 1200
89602	Mayo Stand Cover, Reinforced, X-Large, 76cm x 145cm, Sterile	Items/Case = 44
79602	Mayo Stand Cover, Reinforced, X-Large, 76cm x 145cm, Non-Sterile	Items/Case = 120
14007	Mayo Stand Cover, Padded, X-Large, 76cm x 145cm, Sterile	Items/Case = 44
76358	Mayo Stand Cover, Padded, X-Large, 76cm x 145cm, Non-Sterile	Items/Case = 68
89603	Mayo Stand Cover, Reinforced, XX-Large, 80cm x 145cm, Sterile	Items/Case = 44
79603	Mayo Stand Cover, Reinforced, XX-Large, 80cm x 145cm, Non-Sterile	Items/Case = 120



BACK TABLE COVERS

Code	Description	Packaging
88666	Back Table Cover, Zone Reinforced HALYARDS BASICS* 112cm x 229cm, Sterile	Items/Case = 45
78664	Back Table Cover, Zone Reinforced HALYARDS BASICS* 112cm x 229cm, Non Sterile	Items/Case = 60
89632	Back Table Cover, HALYARD BASICS*, 100cm x 150cm, Sterile	Items/Case = 50
79482	Back Table Cover, 115cm x 152cm, Non-Sterile, Standard Fold	Items/Case = 90
42217	Back Table Cover, Standard, 112cm x 229cm, Sterile	Items/Case = 28
42218	Back Table Cover, Standard, 112cm x 229cm, Non-Sterile, Unfolded	Items/Case = 640
42231	Back Table Cover, Standard, 112cm x 229cm, Non-Sterile, Fan Fold	Items/Case = 400
42224	Back Table Cover, Standard, 152cm x 229cm, Sterile	Items/Case = 36
42225	Back Table Cover, Standard, 152cm x 229cm, Non-Sterile, Unfolded	Items/Case = 480
42233	Back Table Cover, Standard, 152cm x 229cm, Non-Sterile, Fan Fold	Items/Case = 400
42306	Back Table Cover, Heavy Duty, 152cm x 229cm, Sterile	Items/Case = 28
42307	Back Table Cover, Heavy Duty, 152cm x 229cm, Non-Sterile, Unfolded	Items/Case = 400
42366	Back Table Cover, Heavy Duty, 152cm x 229cm, Non-Sterile, Fan Fold	Items/Case = 200
42309	Back Table Cover, Heavy Duty, 177cm x 279cm, Sterile	Items/Case = 15
42310	Back Table Cover, Heavy Duty, 177cm x 279cm, Non-Sterile, Unfolded	Items/Case = 280
42367	Back Table Cover, Heavy Duty, 177cm x 279cm, Non-Sterile, Fan Fold	Items/Case = 160
29258	Back Table Cover, Padded Heavy Duty, 152cm x 229cm, Sterile	Items/Case = 20
29259	Back Table Cover, Padded Heavy Duty, 152cm x 229cm, Non-Sterile, Fan Fold	Items/Case = 22
29269	Back Table Cover, Padded Heavy Duty, 152cm x 229cm, Non-Sterile, Unfolded	Items/Case = 300
37498	Back Table Cover, Two Tier Heavy Duty, 152cm, Sterile	Items/Case = 10

OVERHEAD / CV TABLE COVERS

Code	Description	Packaging
89621	Overhead Table Cover, 216cm x 193cm, 112cm x 152cm Skirt, Sterile	Items/Case = 14

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS


ACCESSORIES

DISPOSABLE HUCK TOWELS, 51cm x 64cm

Code	Description	Packaging
89021	Disposable Huck Towel, Sterile, Pack of 2	Items/Case = 142
79725	Disposable Huck Towel, Non-Sterile	Items/Case = 300

TAPE STRIPS, 10cm x 55cm


Code	Description	Packaging
89053	Tape Strips, Sterile, Pack of 2	Items/Case = 300



DIATHERMY POUCH, 38cm x 41cm

- Adhesive Backing Strip

Code	Description	Packaging
89945	Diathermy Pouch, Sterile	Items/Case = 150
79945	Diathermy Pouch, Non-Sterile	Items/Case = 220



UTILITY DRAPE WITH TAPE, 38cm x 64cm

Code	Description	Packaging
89731	Utility Drape with Tape, Sterile, Pack of 2	Items/Case = 100
79731	Utility Drape with Tape, Non-Sterile	Items/Case = 550

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS

BASIC PACKS

BASIC PACK VI

Code	Description	Packaging
88161	Basic Pack VI	Items/Case = 7
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
95111	71311	1-ULTRA Surgical Gown, In Overwrap with Towel, Large
95121	71321	1-ULTRA Surgical Gown, with Towel, X-Large
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm
89731	79731	4-Utility Drapes, with Tape
		1-Suture Bag
89101	79101	1-Half Drape, 152cm x 112cm
89121	79121	1-Large Drape, 152cm x 193cm

BASIC PACK XV

Code	Description	Packaging
88178	Basic Pack XV	Items/Case = 15
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm
		1-Suture Bag
89121	79121	1-Large Drape, 152cm x 193cm

SET-UP PACK I

Code	Description	Packaging
88171	Set-Up Pack I	Items/Case = 14
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
95111	71311	1-Ultra Surgical Gown In Overwrap with Towel, Large
95121	71321	1-ULTRA Surgical Gown, X-Large
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm
		1-Suture Bag

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


GENERAL DRAPES & PACKS

BASIC DRAPES




UTILITY DRAPE, 90cm x 90cm

Code	Description	Packaging
89353	Utility Drape, Sterile	Items/Case = 120
76323	Utility Drape, Non-Sterile	Items/Case = 200




HALF DRAPE, 152cm x 112cm

Code	Description	Packaging
89101	Half Drape, Sterile	Items/Case = 54
79106	Half Drape, Non-Sterile	Items/Case = 1000




MEDIUM DRAPE, 193cm x 112cm

Code	Description	Packaging
89111	Medium Drape, Sterile	Items/Case = 50
79112	Medium Drape, Non-Sterile	Items/Case = 800



MEDIUM DRAPE WITH TAPE, 193cm x 105cm

Code	Description	Packaging
89023	Medium Drape with Tape, Sterile	Items/Case = 50



LARGE DRAPE, 152cm x 193cm

Code	Description	Packaging
89121	Large Drape, Sterile	Items/Case = 36
79120	Large Drape, Non-Sterile	Items/Case = 600

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


GENERAL DRAPES & PACKS



THREE-QUARTER DRAPE, 152cm x 193cm


- CONTROL* Plus Fabric Reinforcement
- No Adhesive Tape

Code	Description	Packaging
89141	Three-Quarter Drape, Sterile	Items/Case = 20



X-LARGE DRAPE, 193cm x 254cm


Code	Description	Packaging
89131	X-Large Drape, Sterile	Items/Case = 20
79132	X-Large Drape, Non-Sterile	Items/Case = 24



SMALL ABSORBENT DRAPE WITH TAPE, 57cm x 60cm

- CONTROL* Plus Absorbent Fabric

Code	Description	Packaging
89828	Small Absorbent Drape with Tape, Sterile	Items/Case = 72
79828	Small Absorbent Drape with Tape, Non-Sterile	Items/Case = 108



STANDARD ABSORBENT DRAPE WITH TAPE, 91cm x 91cm

- CONTROL* Plus Absorbent Fabric

Code	Description	Packaging
89343	Standard Absorbent Drape with Tape, Sterile	Items/Case = 50
79343	Standard Absorbent Drape with Tape, Non-Sterile	Items/Case = 70

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS


BREAST / CHEST DRAPES



CHEST U-DRAPE, 193cm x 305cm

- CONTROL* Plus Fabric Reinforcement
- 30cm x 76cm Split
- Tube Holders

Code	Description	Packaging
79335	Chest U-Drape, Non-Sterile	Items/Case = 22



BREAST DRAPE, 254cm x 193cm x 315cm

- 15cm x 15cm Fenestration
- CONTROL* Fabric Reinforcement
- Hook and Loop Tube and Line Control
- Armboard Covers

Code	Description	Packaging
89236	Breast Drape, Sterile	Items/Case = 10

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS


FENESTRATED DRAPES



CV FULL BODY DRAPE, 173cm x 279cm

- CONTROL* Plus Fabric Reinforcement
- 10cm Circular Adhesive


Code	Description	Packaging
31445	CV Full Body Drape, Non-Sterile	Items/Case = 27



EPIDURAL DRAPE WITH TAPE, 50cm x 75cm

- 8cm x 10cm Oval Fenestration, Centred
- Spunbond Film
- Positioning Tape


Code	Description	Packaging
89636	HALYARD BASICS* Epidural Drape with Tape, Sterile	Items/Case = 40



VARIABLE FENESTRATED DRAPE, 75cm x 100cm

- Variable Fenestration
- Half Circle On Each Drape with 12cm Diameter
- Spunbond Film
- Positioning Tape

Code	Description	Packaging
89680	HALYARD BASICS* Variable Fenestrated Drape, Sterile	Items/Case = 30



CV LINE DRAPE, 160cm x 240cm

- CONTROL* Fabric Reinforcement
- 12cm x 12cm Fenestration

Code	Description	Packaging
89835	CV Line Drape, Sterile	Items/Case = 10
79835	CV Line Drape, Non-Sterile	Items/Case = 14

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS


LAPAROSCOPY / LAPAROTOMY DRAPES



LAPAROSCOPIC / CHOLECYSTECTOMY DRAPE, 264cm x 193cm x 305cm

- 30cm x 34cm Adjustable Fenestration, Lateral Offset
- Hook and Loop Tube and Line Control
- CONTROL* Instrument Pad Reinforcement
- Armboard Covers


Code	Description	Packaging
89229	Laparoscopic / Cholecystectomy Drape, Sterile	Items/Case = 12



LAPAROSCOPIC ABDOMINAL DRAPE, 254cm x 193cm x 315cm

- 30cm x 33cm Abdominal Fenestration
- CONTROL* Fabric Reinforcement
- Clear Instrument and Line Control Pockets
- Hook and Loop Tube and Line Control
- Armboard Covers

Code	Description	Packaging
89233	Laparoscopic Abdominal Drape, Sterile	Items/Case = 15
79233	Laparoscopic Abdominal Drape, Non-Sterile	Items/Case = 18



LAPAROSCOPIC DRAPE, 264cm x 89cm x 315cm

- 31cm x 33cm Fenestration
- CONTROL* Plus Fabric Reinforcement
- Clear Film Pouch
- Tube Holders
- Arm Covers
- Under Leg Wings

Code	Description	Packaging
29240	Laparoscopic Drape, Sterile	Items/Case = 15
29241	Laparoscopic Drape, Non-Sterile	Items/Case = 19

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


GENERAL DRAPES & PACKS



PAEDIATRIC LAPAROTOMY DRAPE, 152cm x 274cm

- 7cm x 20cm Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders


Code	Description	Packaging
89241	Paediatric Laparotomy Drape, Sterile	Items/Case = 15
79240	Paediatric Laparotomy Drape, Non-Sterile	Items/Case = 25



LAPAROTOMY DRAPE WITH INCISE, 254cm x 183cm x 305cm

- 20cm x 33cm Adhesive Incise Area
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Armboard Covers

Code	Description	Packaging
89282	Laparotomy Drape with Incise, Sterile	Items/Case = 15
79282	Laparotomy Drape with Incise, Non-Sterile	Items/Case = 15



ROBOTIC ABDOMINAL / PELVIC DRAPE, 264cm x 406cm

- 30cm x 10cm x 48cm Trapezoidal Abdominal Fenestration
- Clear Film Poly Pouches
- Tube Holders and Line Control
- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
29283	Robotic Abdominal/Pelvic Drape, Sterile	Items/Case = 8
29284	Robotic Abdominal/Pelvic Drape, Non-Sterile	Items/Case = 10

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS


MINOR PROCEDURE DRAPES



MINOR PROCEDURE FENESTRATED DRAPE, 76cm x 76cm

- 4cm x 5cm Fenestration
- CONTROL* Fabric Reinforcement
- Positioning Tape


Code	Description	Packaging
89206	Minor Procedure Fenestrated Drape, Sterile	Items/Case = 100
79206	Minor Procedure Fenestrated Drape, Non-Sterile	Items/Case = 250



MINOR PROCEDURE FENESTRATED DRAPE, 122cm x 152cm

- 7cm Diameter Fenestration
- CONTROL* Fabric Reinforcement
- Positioning Tape


Code	Description	Packaging
89207	Minor Procedure Fenestrated Drape, Sterile	Items/Case = 40
79207	Minor Procedure Fenestrated Drape, Non-Sterile	Items/Case = 60



MINOR PROCEDURE FENESTRATED DRAPE, 122cm x 152cm

- 12cm x 15cm Fenestration
- CONTROL* Fabric Reinforcement
- Positioning Tape

Code	Description	Packaging
89208	Minor Procedure Fenestrated Drape, Sterile	Items/Case = 40
79208	Minor Procedure Fenestrated Drape, Non-Sterile	Items/Case = 60



MINOR PROCEDURE FENESTRATED DRAPE, 193cm x 315cm

- 15cm x 15cm Fenestration
- CONTROL* Fabric Reinforcement
- Positioning Tape

Code	Description	Packaging
89209	Minor Procedure Fenestrated Drape, Sterile	Items/Case = 19
79209	Minor Procedure Fenestrated Drape, Non-Sterile	Items/Case = 23

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS

THYROID DRAPES



THYROID DRAPE, 193cm x 315cm

- 12cm x 12cm Diamond Fenestration
- CONTROL* Plus Fabric Reinforcement
- Instrument Pad
- Tube Holders


Code	Description	Packaging
89255	Thyroid Drape, Sterile	Items/Case = 13
79255	Thyroid Drape, Non-Sterile	Items/Case = 18

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.


GENERAL DRAPES & PACKS

UNIVERSAL PACKS

UNIVERSAL PACK – MAJOR WITH LONG SIDE DRAPES



Code	Description		Packaging
88759	Universal Pack – Major		Items/Case = 5
	Sterile Component Code	Non-Sterile Component Code	Contents
	89601	79601	1-Back Table Cover, Fully Reinforced, 137cm x 248cm 1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Absorbent Towel, 15cm x 56cm 1-Suture Bag
	34902	34903	1-Bottom Drape with CONTROL* Plus Reinforcement, 193cm x 195cm
	34906	34901	1-Top Drape with CONTROL* Plus Reinforcement, 273cm x 127cm
	34904	34905	2-Side Drapes with CONTROL* Plus Reinforcement, 195cm x 105cm



UNIVERSAL PACK – MINOR WITH SHORT SIDE DRAPES

Code	Description		Packaging
88769	Universal Pack – Minor		Items/Case = 14
	Sterile Component Code	Non-Sterile Component Code	Contents
	89611	79613	1-Back Table Cover, Fully-Reinforced, 112cm x 224cm
	89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
		79097	2-Side Drapes with CONTROL* Reinforcement, 114cm x 84cm
		79096	1-Bottom Drape with CONTROL* Reinforcement, 165cm x 193cm
		79095	1-Top Drape with CONTROL* Reinforcement, 276cm x 150cm

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS



UNIVERSAL PACK I

Code	Description	Packaging
88878	Universal Pack – I	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
88666	78658	1-HALYARD BASICS* Back Table Cover, Zone Reinforced, 112cm x 229cm
88665	78656	1-HALYARD BASICS* Mayo Stand Cover, Reinforced, 58cm x 137cm
	79756	2-HALYARD BASICS* Side Drape, 193 x 112cm
	79754	1-HALYARD BASICS* Top Drape, 274 x 127cm
	79755	1-HALYARD BASICS* Bottom Drape, 193 x 173cm



UNIVERSAL PACK III WITH MEDIUM SIDE DRAPES

Code	Description	Packaging
88774	Universal Pack – III	Items/Case = 8
Sterile Component Code	Non-Sterile Component Code	Contents
89601	79601	1-Back Table Cover, Fully-Reinforced, 133cm x 248cm
		1-Mayo Stand Cover, 58cm x 137cm
		1-Suture Bag
		2-Side Drapes with CONTROL* Reinforcement, 145cm x 99cm
		1-Bottom Drape with CONTROL* Reinforcement, 193cm x 193cm
		1-Top Drape with CONTROL* Reinforcement, 276cm x 150cm

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS



U-BAR* PACK

Code	Description	Packaging	
88334	U-BAR* Pack	Items/Case = 16	
Sterile Component Code	Non-Sterile Component Code	Contents	
89601	79601	1-Back Table Cover, Zone Reinforced, 112cm x 224cm 1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag 1-U-Drape with CONTROL* Reinforcement, 180cm x 315cm, with 10cm x 102cm Split 1-Bar Drape with CONTROL* Reinforcement, 105cm x 228cm	



U-BAR* PACK – IV

Code	Description	Packaging	
88773	U-BAR* Pack – IV	Items/Case = 10	
Sterile Component Code	Non-Sterile Component Code	Contents	
89601	79601	1-Back Table Cover, Reinforced, 133cm x 248cm 1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-U-Drape with CONTROL* Reinforcement, 192cm x 315cm, with 11cm x 51cm Split 1-Bar Drape with CONTROL* Reinforcement, 248cm x 145cm	

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS


UNIVERSAL DRAPES



TOP DRAPE, 273cm x 127cm

- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape
- From Universal Pack 88759


Code	Description	Packaging
34906	Top Drape, Sterile	Items/Case = 15
34901	Top Drape, Non-Sterile	Items/Case = 275



BOTTOM DRAPE, 193cm x 195cm

- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape
- From Universal Pack 88759

Code	Description	Packaging
34902	Bottom Drape, Sterile	Items/Case = 12
34903	Bottom Drape, Non-Sterile	Items/Case = 242



SIDE DRAPE, 195cm x 105cm

- CONTROL* Plus Fabric Reinforcement
- Positioning Tape
- From Universal Pack 88759

Code	Description	Packaging
34904	Side Drape, Sterile	Items/Case = 20
34905	Side Drape, Non-Sterile	Items/Case = 440

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


GENERAL DRAPES & PACKS



TOP DRAPE, 276cm x 150cm

- CONTROL* Fabric Reinforcement
- Tube Holders
- Positioning Tape
- From Universal Pack 88769


Code	Description	Packaging
79095	Top Drape, Non-Sterile	Items/Case = 60



BOTTOM DRAPE, 165cm x 193cm

- CONTROL* Fabric Reinforcement
- Tube Holders
- Positioning Tape
- From Universal Pack 88769

Code	Description	Packaging
79096	Bottom Drape, Non-Sterile	Items/Case = 70




SIDE DRAPE, 114cm x 84cm

- CONTROL* Fabric Reinforcement
- Positioning Tape
- From Universal Pack 88769

Code	Description	Packaging
79097	Side Drape, Non-Sterile	Items/Case = 200

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.

NEURO / SPINE DRAPE



CRANIOTOMY DRAPE WITH POUCH, 310cm x 188cm x 340cm


- CONTROL* Plus Fabric Reinforcement
- 20cm x 30cm Oval Fenestration
- Clear Gusseted Anesthesia Panels
- Screen and Suction Port

Code	Description	Packaging
89062	Craniotomy Drape with Pouch, Sterile	Items/Case = 9
79075	Craniotomy Drape with Pouch, Non-Sterile	Items/Case = 15

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.

OBSTETRICS & GYNAECOLOGY

DRAPES & PACKS




C-SECTION DRAPES

C- SECTION DRAPE WITH CLEAR SCREEN, 264cm x 193cm x 305m

- Triangular Fenestration
- 360° Fluid Collection Pouch
- Clear Screen
- Tube Holders
- Armboard Covers


Code	Description	Packaging
44968	C-Section Drape with Clear Screen, Sterile	Items/Case = 6
44967	C-Section Drape with Clear Screen, Non-Sterile	Items/Case = 100



C-SECTION DRAPE WITH POUCH AND FENESTRATION, 264cm x 193cm x 305cm

- 36cm x 31cm Adhesive Area with Fenestration
- 360° Fluid Collection Pouch
- CONTROL* Fabric Reinforcement
- Tube Holders
- Armboard Covers


Code	Description	Packaging
44978	C-Section Drape with Pouch and Fenestration, Sterile	Items/Case = 6
44977	C-Section Drape with Pouch and Fenestration, Non-Sterile	Items/Case = 100



C-SECTION DRAPE WITH FENESTRATED INCISE, 254cm x 183cm x 305cm

- 38cm x 33cm Adhesive Incise Area with Triangular Fenestration
- 270° Fluid Collection Pouch, Malleable Band, with Suction Port
- Tube Holders
- Armboard Covers

Code	Description	Packaging
89528	C-Section Drape with Fenestrated Incise, Sterile	Items/Case = 8
79528	C-Section Drape with Fenestrated Incise, Non-Sterile	Items/Case = 150



C-SECTION DRAPE WITH FULL INCISE, 254cm x 183cm x 305cm

- 38cm x 33cm Adhesive Incise Area
- 270° Fluid Collection Pouch, Malleable Band, with Suction Port
- Tube Holders
- Armboard Covers

Code	Description	Packaging
89529	C-Section Drape with Full Incise, Sterile	Items/Case = 8
79529	C-Section Drape with Full Incise, Non-Sterile	Items/Case = 150

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


OBSTETRICS & GYNAECOLOGY DRAPES & PACKS

LAPAROSCOPY PACKS

LAPAROSCOPY PACK I



Code	Description	Packaging
88088	Laparoscopy Pack I	Items/Case = 10
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm 1-Absorbent Towel, 38cm x 56cm
89411		1-Under Buttocks Drape, Plastic, 75cm x 112cm
89408	79408	2-Leggings with 15cm Cuff, 78cm x 122cm
89119	79117	1-Laparoscopy "T" Drape with CONTROL* Plus Reinforcement, 249cm x 104cm x 251cm • 15cm x 20cm Rectangular Abdominal Fenestration • 7cm x 15cm Perineal Fenestration with Flap • Tube Holders



LAPAROSCOPY / PELVISCOPY PACK

Code	Description	Packaging
88222	Laparoscopy/Pelviscopy Pack	Items/Case = 6
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
95111	71311	1-ULTRA Surgical Gown, In Overwrap, with Towel, Large
95111	71311	1-ULTRA Surgical Gown, with Towel, Large
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
89413		1-Under Buttocks Drape, Plastic-Reinforced, Large 101cm x 116cm
89408	79408	2-Leggings, with 15cm Cuff, Sterile Pair, 78cm x 122cm 1-Laparoscopy / Pelviscopy Drape with CONTROL* Plus Reinforcement (without Leggings), 264cm x 89cm x 264cm • 18cm x 23cm Trapezoidal Abdominal Fenestration • 8cm x 15cm Perineal Fenestration with Flap

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

OBSTETRICS & GYNAECOLOGY DRAPES & PACKS


LAPAROSCOPY DRAPES



LAPAROSCOPY / PELVISCOPY DRAPE WITH FLUID COLLECTION POUCH, 284cm x 76cm x 314cm

- 17cm x 23cm Trapezoidal Abdominal Fenestration
- 7cm x 20cm Perineal Fenestration with Flap
- Fluid Collection Pouch with Screen and Suction Port
- CONTROL* Plus Fabric Reinforcement
- Attached Leggings
- Hook and Loop Line Controls
- Armboard Covers

Code	Description	Packaging
89219	Laparoscopy/Pelviscopy Drape Drape with Pouch, Sterile	Items/Case = 10



LAPAROSCOPY / PELVISCOPY DRAPE, 254cm x 76cm x 312cm

- 17cm x 23cm Trapezoidal Abdominal Fenestration
- 7cm x 20cm Perineal Fenestration with Flap
- CONTROL* Plus Fabric Reinforcement
- Attached Leggings
- Hook and Loop Tube and Line Controls

Code	Description	Packaging
89224	Laparoscopy/Pelviscopy Drape, Sterile	Items/Case = 10
79224	Laparoscopy/Pelviscopy Drape, Non-Sterile	Items/Case = 10

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

OBSTETRICS & GYNAECOLOGY


DRAPES & PACKS

LITHOTOMY PACKS



LITHOTOMY PACK II

Code	Description	Packaging
88521	Lithotomy Pack II	Items/Case = 10
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm 1-Absorbent Towel, 38cm x 56cm
89413	79413	1-Under Buttocks Drape, Plastic-Reinforced, 102cm x 117cm
89408	79408	2-Leggings with 15cm Cuff, 78cm x 122cm
89101	79101	1-Half Drape, 152cm x 111cm



LITHOTOMY PACK V

Code	Description	Packaging
88551	Lithotomy Pack V	Items/Case = 12
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
89701	79720	1-Absorbent Towel, 38cm x 56cm
89408	79408	2-Leggings, with 15cm Cuffs
	79551	1-Lithotomy "T" Drape with CONTROL* Plus Reinforcement, 152cm x 78cm x 190cm with 8cm x 15cm Perineal Fenestration

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

OBSTETRICS & GYNAECOLOGY DRAPES & PACKS


LITHOTOMY DRAPES



LITHOTOMY DRAPE, 254cm x 76cm x 221cm

- 7cm x 15cm Perineal Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Attached Leggings


Code	Description	Packaging
89203	Lithotomy Drape, Sterile	Items/Case = 16



GYN / UROLOGY DRAPE WITH POUCH, 160cm x 284cm x 332cm

- Pouch with Mesh Screen and Port
- Attached Leggings
- CONTROL* Plus Fabric Reinforcement
- 7cm x 15cm Perineal Fenestration
- Tube Holders
- Armboard Covers


Code	Description	Packaging
89204	Gyn/Urology Drape with Pouch, Sterile	Items/Case = 10
79204	Gyn/Urology Drape with Pouch, Non-Sterile	Items/Case = 10



LITHOTOMY "T" DRAPE, 152cm x 78cm x 190cm

- 7cm x 15cm Perineal Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
79551	Lithotomy "T" Drape, Non-Sterile	Items/Case = 55



LEGGINGS, 78cm x 122cm, WITH 15cm CUFF

Code	Description	Packaging
89408	Leggings, Sterile Pair	Items/Case = 26
79408	Leggings, Non-Sterile	Items/Case = 70

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

OBSTETRICS & GYNAECOLOGY

DRAPES & PACKS


UNDER BUTTOCKS DRAPES



UNDER BUTTOCKS DRAPE, ZONE-REINFORCED, MEDIUM, 96cm x 112cm


- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
79411	Under Buttocks Drape, Zone-Reinforced, Medium, Non-Sterile	Items/Case = 115



UNDER BUTTOCKS DRAPE, PLASTIC REINFORCED, LARGE, 101cm x 116cm


Code	Description	Packaging
89413	Under Buttocks Drape, Plastic Reinforced, Large, Sterile	Items/Case = 50



UNDER BUTTOCKS DRAPE WITH FLUID COLLECTION POUCH, 101cm x 112cm

- CONTROL* Plus Absorbent Pad
- 72cm x 66cm Triangular Pouch

Code	Description	Packaging
89414	Under Buttocks Drape, with Fluid Collection Pouch, Sterile	Items/Case = 30
79414	Under Buttocks Drape, with Fluid Collection Pouch, Non-Sterile	Items/Case = 80



UNDER BUTTOCKS DRAPE WITH FLUID COLLECTION POUCH, 101cm x 112cm

- Screen and Suction Port
- Adhesive Pouch Closure
- 69cm x 61cm x 7 cm Triangular Pouch

Code	Description	Packaging
89415	Under Buttocks Drape with Fluid Collection Pouch, Sterile	Items/Case = 40
79415	Under Buttocks Drape with Fluid Collection Pouch, Non-Sterile	Items/Case = 100

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


OPHTHALMIC DRAPE & PACKS



BASIC OPHTHALMIC DRAPE, 132cm x 142cm

- Lightweight Fabric
- 12cm x 9cm Adhesive Incise Area
- Fluid Collection Pouch, with Eye Wick

Code	Description	Packaging
89304	Basic Ophthalmic Drape, Sterile	Items/Case = 40
79304	Basic Ophthalmic Drape, Non-Sterile	Items/Case = 100



OCUARTS* FULL BODY OPHTHALMIC DRAPE, 183cm x 292cm

- 11cm x 11cm Adhesive Incise Area
- 5cm x 3cm Oval Fenestration
- 12cm x 15cm Optional Incise Drape
- 270° Fluid Collection Pouch, Foam Band
- Breather Bar
- Tube Holders


Code	Description	Packaging
89356	OCUARTS* Full Body Ophthalmic Drape, Sterile	Items/Case = 10
79356	OCUARTS* Full Body Ophthalmic Drape, Non-Sterile	Items/Case = 20

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.


ORTHOPAEDIC DRAPES & PACKS

ORTHOPAEDIC PACKS

ORTHO PACK I – MAJOR



Code	Description	Packaging
88411	Ortho Pack I – Major	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
92321	30734	2-MICROCOOL* Breathable High Performance Surgical Gown, with Towel, X-Large
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
89331	79330	1-Split Drape, 152cm x 178cm with 10cm x 53cm Split 1-Ortho U-Drape with CONTROL* Plus Reinforcement, 223cm x 274cm with 10cm x 102cm Split
89381	79360	1-Ortho Bar Drape with CONTROL* Plus Reinforcement, 254cm x 152cm



ORTHOARTS* HIP PACK WITH POCKETS

Code	Description	Packaging
88452	ORTHOARTS* Hip Pack with Pockets	Items/Case = 3
Sterile Component Code	Non-Sterile Component Code	Contents
89906	79006	1-Back Table Cover, Heavy-Duty, 152cm x 228cm
95111	71311	1-ULTRA Surgical Gown, in Overwrap, with Towel, Large
92321	30734	2-MICROCOOL* Breathable High Performance Surgical Gown, with Towel, X-Large
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
89821	79515	1-Stockinette, Large
89121	79121	1-Large Drape, 152cm x 193cm
89333	79333	1-Split Drape, Large, 152cm x 205cm with 10cm x 79cm Split
89345	79345	1-ORTHOARTS* Hip Drape with CONTROL* Plus Reinforcement, 289cm x 223cm x 320cm, with 21cm x 17cm Conformable Teardrop Fenestration 1-Yankauer Suction Tip 1-304cm Suction Tubing

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

ORTHOPAEDIC DRAPES & PACKS



ORTHOARTS* LOWER EXTREMITY

Code	Description	Packaging
88471	ORTHOARTS* Lower Extremity Pack	Items/Case = 3
Sterile Component Code	Non-Sterile Component Code	Contents
89906	79006	1-Back Table Cover, Heavy-Duty, 152cm x 228cm
95111	71311	1-ULTRA Surgical Gown, In Overwrap, with Towel, Large
		1-Extremity Prep Pad
95221	74921	1-ULTRA Fabric-Reinforced Surgical Gown, with Towel, X-Large
92331		1-MICROCOOL* Breathable High Performance Specialty Surgical Gown, with Towel, X-Large
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm
		1-Suture Bag
		1-304cm Suction Tubing
89821	79515	1-Stockinette, Large
89331	79330	1-Split Drape, 152cm x 178cm with 10cm x 53cm Split
89278	13233	1-ORTHOARTS* Lower Extremity Drape with CONTROL* Plus Reinforcement, with 6cm Conformable Circular Fenestration, 3-Drape Tape Strips

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

ORTHOPAEDIC DRAPES & PACKS


EXTREMITY DRAPES



LOWER EXTREMITY DRAPE, 274cm x 203cm x 315cm

- Armboard Covers
- 6cm Conformable Circular Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders


Code	Description	Packaging
89016	Lower Extremity Drape, Sterile	Items/Case = 12



LOWER EXTREMITY DRAPE, 223cm x 332cm

- CONTROL* Plus Fabric Reinforcement
- 6cm Conformable Circular Fenestration
- Tube Holders


Code	Description	Packaging
89276	Lower Extremity Drape, Sterile	Items/Case = 14



ORTHOARTS* LOWER EXTREMITY DRAPE, 289cm x 223cm x 320cm

- Armboard Covers
- 6cm Conformable Circular Fenestration
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89278	ORTHOARTS* Lower Extremity Drape, Sterile	Items/Case = 7
13233	ORTHOARTS* Lower Extremity Drape, Non-Sterile	Items/Case = 84



ORTHOARTS* EXTRA LONG LOWER EXTREMITY DRAPE, 289cm x 223cm x 381cm

- Armboard Covers
- 6cm Conformable Circular Fenestration
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
29219	ORTHOARTS* Extra Long Lower Extremity Drape, Sterile	Items/Case = 8
29220	ORTHOARTS* Extra Long Lower Extremity Drape, Non-Sterile	Items/Case = 8

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


ORTHOPAEDIC DRAPES & PACKS



UNIVERSAL EXTREMITY DRAPE, 229cm x 333cm

- CONTROL* Plus Fabric Reinforcement
- 6cm Conformable Circular Fenestration


Code	Description	Packaging
89279	Universal Extremity Drape, Sterile	Items/Case = 14
79279	Universal Extremity Drape, Non-Sterile	Items/Case = 150



BILATERAL LIMB DRAPE, 193cm x 304cm

- CONTROL* Plus Fabric Reinforcement
- Two 6cm Conformable Circular Fenestrations
- Tube Holders


Code	Description	Packaging
89291	Bilateral Limb Drape, Sterile	Items/Case = 14



BILATERAL LIMB DRAPE WITH ARMBOARD COVERS, 295cm x 223cm x 366cm

- CONTROL* Plus Fabric Reinforcement
- Two 6cm Conformable Circular Fenestrations
- Tube Holders
- Armboard Covers

Code	Description	Packaging
89380	Bilateral Limb Drape with Armboard Covers, Sterile	Items/Case = 12



UPPER EXTREMITY DRAPE, 193cm x 305cm

- 3cm Circular Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89271	Upper Extremity Drape, Sterile	Items/Case = 14
79273	Upper Extremity Drape, Non-Sterile	Items/Case = 20

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

ORTHOPAEDIC DRAPES & PACKS


HAND DRAPES



HAND DRAPE, 289cm x 360cm

- 157cm x 147cm Hand Table Cover
- 3cm Conformable Circular Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
77450	Hand Drape, Sterile	Items/Case = 11
13227	Hand Drape, Non-Sterile	Items/Case = 108



ORTHOARTS* HAND DRAPE, 193CM X 381CM


- 157cm x 147cm Gusseted Hand Table Cover
- 3cm Conformable Fenestration
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89272	ORTHOARTS* Hand Drape, Sterile	Items/Case = 6

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.

ORTHOPAEDIC DRAPES & PACKS


HIP DRAPES



ORTHOARTS* HIP DRAPE WITH LEG POCKETS, 289cm x 223cm x 320cm

- Leg Pockets On Both Sides
- Armboard Covers
- 21cm x 17cm Conformable Teardrop Fenestration
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders


Code	Description	Packaging
89345	ORTHOARTS* Hip Drape with Leg Pockets, Sterile	Items/Case = 6
79345	ORTHOARTS* Hip Drape with Leg Pockets, Non-Sterile	Items/Case = 6



ORTHOARTS* HIP DRAPE, 289cm x 223cm x 320cm

- Armboard Covers
- 21cm x 17cm Conformable Teardrop Fenestration
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89346	ORTHOARTS* Hip Drape, Sterile	Items/Case = 7
79346	ORTHOARTS* Hip Drape, Non-Sterile	Items/Case = 7



HIP DRAPE WITH CLEAR LEG POCKETS, 294cm x 223cm x 345cm

- Clear Leg Pockets
- Armboard Covers
- 21cm x 17cm Conformable Teardrop Fenestration
- CONTROL* Plus Fabric Reinforcement
- Hook and Loop Tube and Line Control

Code	Description	Packaging
89349	Hip Drape with Clear Leg Pockets, Sterile	Items/Case = 9
79349	Hip Drape with Clear Leg Pockets, Non-Sterile	Items/Case = 10

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

ORTHOPAEDIC DRAPES & PACKS

KNEE ARTHROSCOPY DRAPES


KNEE ARTHROSCOPY DRAPE, 228cm x 315cm



- Fluid Collection Pouch with Suction Port
- 76cm x 99cm x 8cm Triangular Pouch
- 6cm Conformable Circular Fenestration
- Tube Holders

Code	Description	Packaging
89265	Knee Arthroscopy Drape, Sterile	Items/Case = 10
79265	Knee Arthroscopy Drape, Non-Sterile	Items/Case = 150


KNEE ARTHROSCOPY DRAPE WITH ARMBORD COVERS, 290cm x 223cm x 315cm



- Armboard Covers
- Fluid Collection Pouch with Suction Port
- 76cm x 99cm x 9cm Triangular Pouch
- 6cm Conformable Circular Fenestration
- Tube Holders

Code	Description	Packaging
89371	Knee Arthroscopy Drape with Armboard Covers, Sterile	Items/Case = 10
79371	Knee Arthroscopy Drape with Armboard Covers, Non-Sterile	Items/Case = 13


ORTHOARTS* KNEE ARTHROSCOPY DRAPE, 289cm x 223cm x 294cm



- Armboard Covers
- Fluid Collection Pouch with Reflux Shields and Suction Port
- 76cm x 99cm x 9cm Triangular Pouch
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- 6cm Conformable Circular Fenestration
- Tube Holders

Code	Description	Packaging
89267	ORTHOARTS* Knee Arthroscopy Drape, Sterile	Items/Case = 7
79267	ORTHOARTS* Knee Arthroscopy Drape, Non-Sterile	Items/Case = 7

ORTHOARTS* EXTRA LONG KNEE ARTHROSCOPY DRAPE, 294cm x 233cm x 447cm



- Armboard Covers
- Fluid Collection Pouch with Reflux Shields and Suction Port
- 95cm x 112cm Pentagonal Pouch
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- 6cm Conformable Circular Fenestration
- Tube Holders

Code	Description	Packaging
36930	ORTHOARTS* Extra Long Knee Arthroscopy Drape, Sterile	Items/Case = 10
42429	ORTHOARTS* Extra Long Knee Arthroscopy Drape, Non-Sterile	Items/Case = 10

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

ORTHOPAEDIC DRAPES & PACKS


SHOULDER DRAPES



ORTHOARTS* BEACH CHAIR SHOULDER ARTHROSCOPY DRAPE, 411cm x 261cm

- Attached Fluid Collection Pouch
- 95cm x 97cm Pentagonal Pouch
- Hook and Loop Tube and Line Control
- 6cm Conformable Tear Drop Fenestration
- Heavy Duty Fabric


Code	Description	Packaging
36932	ORTHOARTS* Beach Chair Shoulder Arthroscopy Drape, Sterile	Items/Case = 10
42426	ORTHOARTS* Beach Chair Shoulder Arthroscopy Drape, Non-Sterile	Items/Case = 10



BEACH CHAIR SHOULDER DRAPE WITH ATTACHABLE FLUID COLLECTION POUCH, 437cm x 254cm

- Attachable Fluid Collection Pouch with Suction Port, 89cm x 109cm
- CONTROL* Plus Fabric Reinforcement
- Hook and Loop Tube and Line Control
- Pouch Attachment Zones
- 90cm x 110cm Trapezoidal Pouch

Code	Description	Packaging
89066	Beach Chair Shoulder Drape with Attachable Fluid Collection Pouch, Sterile	Items/Case = 6
79049	Beach Chair Shoulder Drape without Attachable Fluid Collection Pouch, Non-Sterile	Items/Case = 11
79061	Attachable Fluid Collection Pouch, Non-Sterile	Items/Case = 44



BEACH CHAIR SHOULDER ARTHROSCOPY DRAPE II, 406cm x 259cm

- Attached Fluid Collection Pouch
- 103m x 63cm x 2cm Pentagonal Pouch
- Hook and Loop Tube and Line Control
- 6cm Conformable Tear Drop Fenestration
- Instrument Pad Fabric Reinforcement

Code	Description	Packaging
89070	Beach Chair Shoulder Arthroscopy Drape, Sterile	Items/Case = 10
79070	Beach Chair Shoulder Arthroscopy Drape, Non-Sterile	Items/Case = 12

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


ORTHOPAEDIC DRAPES & PACKS



2-PIECE SHOULDER ARTHROSCOPY DRAPE PACK

- Body Split Drape, 223cm x 294cm
- 10cm x 68cm Split
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Shoulder Arthroscopy Drape, 254cm x 157cm
- Fluid Collection Pouch with Suction Ports
- 55cm x 112cm Trapezoidal Pouch
- Tube Holders


Code	Description	Packaging
89490	2-Piece Shoulder Arthroscopy Drape Pack, Sterile	Items/Case = 12
79490	2-Piece Shoulder Arthroscopy Drape Pack, Non-Sterile	Items/Case = 12



ORTHOARTS* SHOULDER ARTHROSCOPY DRAPE, 254cm x 152cm

- Fluid Collection Pouch
- 79cm x 149cm x 48cm Trapezoidal Pouch
- Heavy-Duty Fabric

Code	Description	Packaging
89492	ORTHOARTS* Shoulder Arthroscopy Drape, Sterile	Items/Case = 12
79492	ORTHOARTS* Shoulder Arthroscopy Drape, Non-Sterile	Items/Case = 15



ORTHOARTS* SHOULDER SPLIT DRAPE, 284cm x 152cm

- 15cm x 20cm Dual Splits
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89274	ORTHOARTS* Shoulder Split Drape, Sterile	Items/Case = 12

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

ORTHOPAEDIC DRAPES & PACKS


U & BAR DRAPES



U-DRAPE, 193cm x 305cm

- 10cm x 101cm Split
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape


Code	Description	Packaging
89301	U-Drape, Sterile	Items/Case = 17
13231	U-Drape, Non-Sterile	Items/Case = 190



LARGE U-DRAPE, 223cm x 317cm

- 15cm x 137cm Split
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape

Code	Description	Packaging
89311	Large Ortho U-Drape, Sterile	Items/Case = 11
79363	Large Ortho U-Drape, Non-Sterile	Items/Case = 15



ORTHOARTS* U-DRAPE, 223CM X 274CM

- 10cm x 101cm Split
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
89347	ORTHOARTS* U-Drape, Sterile,	Items/Case = 10
13294	ORTHOARTS* U-Drape, Non-Sterile	Items/Case = 132

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


ORTHOPAEDIC DRAPES & PACKS



ORTHOARTS* BODY SPLIT DRAPE, 223CM X 294CM

- 10cm x 68cm Split
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape


Code	Description	Packaging
89275	ORTHOARTS* Body Split Drape, Sterile	Items/Case = 10
79275	ORTHOARTS* Body Split Drape, Non-Sterile	Items/Case = 12



SPLIT DRAPE, 152cm x 178cm

- Polyethylene Material
- 10cm x 53cm Split

Code	Description	Packaging
89331	Split Drape, Sterile	Items/Case = 42
79330	Split Drape, Non-Sterile	Items/Case = 115



LARGE SPLIT DRAPE, 152cm x 205cm

- Polyethylene Material
- 11cm x 78cm Split

Code	Description	Packaging
89333	Large Split Drape, Sterile	Items/Case = 34
79333	Large Split Drape, Non-Sterile	Items/Case = 100

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


ORTHOPAEDIC DRAPES & PACKS



ORTHO SPLIT "T" DRAPE, 297cm x 210cm x 342cm

- Armboard Covers
- 10cm x 53cm Split
- CONTROL* Fabric Reinforcement
- Tube Holders
- Positioning Tape


Code	Description	Packaging
36931	Ortho Split "T" Drape, Sterile	Items/Case = 10



ORTHO BAR DRAPE, 254cm x 152cm

- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89381	Ortho Bar Drape, Sterile	Items/Case = 18
79360	Ortho Bar Drape, Non-Sterile	Items/Case = 28



VERTICAL ISOLATION DRAPE, 228cm x 320cm

- Clear Film
- Adhesive Incise Area 40cm x 76cm
- Attachable Fluid Collection Pouch with Suction Port
- 42cm x 84cm x 17cm Trapezoidal Pouch
- Tube Holders
- Central Release Tape Strip For Horizontal Rod Attachment

Code	Description	Packaging
77723	Vertical Isolation Drape, Sterile	Items/Case = 14

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

ORTHOPAEDIC DRAPES & PACKS



ISOLATION SCREEN, 170cm x 290cm

- Clear Plastic Drape 170cm x 290cm
- Adhesive Backing Strips

Code	Description	Packaging
89093	Isolation Screen, Sterile	Items/Case = 24
79093	Isolation Screen, Non-Sterile	Items/Case = 30



STOCKINETTES

Code	Description	Packaging
89811	Stockinette, Sterile, Medium, 23cm x 112cm	Items/Case = 36
79505	Stockinette, Non-Sterile, Medium, 23cm x 112cm	Items/Case = 75
89821	Stockinette, Sterile, Large, 30cm x 122cm	Items/Case = 25
79515	Stockinette, Non-Sterile, Large, 30cm x 122cm	Items/Case = 50

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


UROLOGY DRAPES & PACKS



238 cm x 175 cm x 170 cm


T.U.R. PACK

Code	Description	Packaging
77885	T.U.R. Pack	Items/Case = 12
	Non-Sterile Component	
Code	Contents	
	1-Back Table Cover, Zone-Reinforced, 152cm x 191cm	
	1-Absorbent Towel, 38cm x 56cm	
76885	1-T.U.R. Drape, 289cm x 175cm x 173cm:	
	• Attached Leggings	
	• Fluid Collection Pouch with Screen, Suction Port and Tie	
	• 54cm x 75cm Triangular Pouch	
	• Finger Cot	
	• 7cm Fenestration	
	• 25cm x 20cm Fenestration	



CYSTOSCOPY PACK I

Code	Description	Packaging
88611	Cystoscopy Pack I	Items/Case = 12
	Sterile Component	
	Non-Sterile Component	
Code	Code	Contents
		1-Back Table Cover, Zone-Reinforced, 112cm x 198cm
		1-Absorbent Towel, 38cm x 56cm
89408	79408	2-Leggings, with 15cm Cuff, 78cm x 122cm
		1-Cystoscopy "T" Drape with CONTROL* Plus Reinforcement, 152cm x 76cm x 213cm



PERCUTANEOUS NEPHROLITHOTOMY DRAPE (PCNL), 299cm x 182cm

- CONTROL* Fabric Reinforcement
- Tube Holders
- Fluid Collection Pouch
- 67cm x 61cm x 7cm Trapezoidal Pouch
- 18cm x 28cm Fenestration

Code	Description	Packaging
42427	Percutaneous Nephrolithotomy Drape, Sterile	Items/Case = 12

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

MICROBIAL SEALANT



HALYARD* INTEGUSEAL* MICROBIAL SEALANT

- Ergonomic Single-Use Applicator
- Evenly Covers Incision Area in a Single Stroke
- Widely-Used Cyanoacrylate Forms a Safe, Fast-Drying, Breathable Barrier
- Compatible with Most Common Surgical Techniques and Materials
- Need Not be Removed for Skin Closure
- Can be Used with a Variety of Preps and Incise Drapes
- Nonflammable, Latex-Free, Non-Irritating and Non-Cytotoxic

Code	Description	Coverage Area (approx)	Example of Procedures	Packaging
33736	IS 100	10" x 10" 25cm x 25cm	Cranial and Spinal Procedures, Laparotomies, Hernia Repair, Thyroid Surgery, Cardiac Catheterisation	Applicators/Box = 10 Boxes/Case = 2 Applicators/Case = 20
33737	IS 200	20" x 10" 50cm x 25cm	Cholecystectomy, Open Heart, Lung Surgery, Total Joint Replacement	Applicators/Box = 10 Boxes/Case = 2 Applicators/Case = 20

INDEX

Code	Page	Code	Page	Code	Page	Code	Page	Code	Page
13227	40	42427	49	79335	17	88719	11	89353	15
13231	45	42429	42	79343	16	88759	23	89356	35
13233	38	42621	4	79345	41	88769	23	89368	6
13294	45	43623	4	79346	41	88774	24	89371	42
13296	7	43940	4	79349	41	88773	25	89380	39
13359	3	43941	4	79351	11	88878	24	89381	47
14007	12	44967	29	79356	35	89016	38	89408	33
29219	38	44968	29	79360	47	89021	13	89413	34
29220	38	44977	29	79363	45	89023	15	89414	34
29240	19	44978	29	79368	6	89053	13	89415	34
29241	19	70392	5	79371	42	89062	28	89431	11
29283	20	72448	3	79408	33	89066	43	89438	11
29284	20	72454	3	79411	34	89070	43	89452	8
29258	12	72978	5	79414	34	89093	48	89453	9
29259	12	72988	5	79415	34	89101	15	89454	8
29269	12	76323	15	79430	11	89111	15	89471	6
30734	3	76356	8	79452	8	89121	15	89490	44
31445	18	76358	12	79453	9	89131	16	89492	44
33736	50	77450	40	79482	12	89141	16	89528	29
33737	50	77723	47	79490	44	89203	33	89529	29
34901	26	77885	49	79492	44	89204	33	89543	11
34902	26	78656	12	79505	48	89206	21	89583	6
34903	26	78664	12	79515	48	89207	21	89601	12
34904	26	79049	43	79528	29	89208	21	89602	12
34905	26	79061	43	79529	29	89209	21	89603	12
34906	26	79070	43	79551	33	89219	31	89621	12
36930	42	79075	28	79602	12	89224	31	89632	12
36931	47	79093	48	79603	12	89229	19	89636	18
36932	43	79095	27	79607	12	89233	19	89680	18
37498	12	79096	27	79709	7	89236	17	89709	7
41725	4	79097	27	79725	13	89241	20	89711	7
41727	4	79106	15	79731	13	89255	22	89731	13
41739	4	79112	15	79749	7	89265	42	89811	48
41740	4	79120	15	79828	16	89267	42	89821	48
41741	4	79132	16	79835	18	89271	39	89828	16
41870	4	79204	33	79945	13	89272	40	89835	18
41871	4	79206	21	88088	30	89274	44	89945	13
41872	4	79207	21	88161	14	89275	46	90003	5
41873	4	79208	21	88171	14	89276	38	90018	5
41874	4	79209	21	88178	14	89278	38	90019	5
42217	12	79224	31	88222	30	89279	39	90048	5
42218	12	79233	19	88223	10	89282	20	90049	5
42224	12	79240	20	88334	25	89291	39	90050	5
42225	12	79255	22	88411	36	89301	45	90138	5
42231	12	79265	42	88452	36	89304	35	92320	3
42233	12	79267	42	88471	37	89311	45	92321	3
42306	12	79273	39	88521	32	89333	46	92323	3
42307	12	79275	46	88551	32	89343	16	92330	3
42309	12	79279	39	88611	49	89345	41	92331	3
42310	12	79282	20	88665	12	89346	41	92333	3
42366	12	79304	35	88666	12	89347	45		
42367	12	79330	46	88711	10	89349	41		
42426	43	79333	46	88716	10	89352	11		



Infection Prevention Product Solutions
Halyard* Education Foundation
Tools & Best Practice
Commitment to Excellence and Sustainability
Experienced Customer Support
Knowledgeable Sales Force
Innovation & Technology



Halyard Australia Pty. Limited

ABN No 25 100 844 488

Level 7, 52 Alfred Street, Milsons Point NSW 2061

Phone: 1800 101 021

Email: CustomerService_ANZ@hyh.com

Fax: 1800 730 687

Halyard New Zealand

Office 2–5, Building 10

666 Great South Road, Penrose, Auckland 1051

Phone: 0800 482 211

www.halyardhealth.com.au



Halyard Health makes no representation that its products described in this catalogue will provide complete protection from exposure to contaminants, fluids or micro-organisms or will eliminate the risk of contracting infectious diseases.

*Registered Trademark or Trademark of Halyard Health, Inc. or its affiliates. © 2014 HYH. All rights reserved.
SS1612:WFHY723