



HALYARD

FORMERLY KIMBERLY-CLARK HEALTH CARE

HEALTHCARE PRODUCT CATALOGUE 2015

AUSTRALIA &
NEW ZEALAND



QUICK LINKS

[CLINICAL TOWELS](#)

[DIGESTIVE HEALTH ENTERAL FEEDING](#)

[DIGESTIVE HEALTH INTRODUCER KITS](#)

[FACIAL PROTECTION](#)

[FLUID ABSORPTION](#)

[MEDICAL GLOVES](#)

[PATIENT CARE](#)

[PROTECTIVE APPAREL](#)

[RESPIRATORY HEALTH CATHETER KIT](#)

[RESPIRATORY HEALTH CLOSED SUCTION SYSTEMS](#)

[RESPIRATORY HEALTH ENDOTRACHEAL TUBES](#)

[RESPIRATORY HEALTH ORAL CARE](#)

[STERILE CONTAINMENT](#)

[SURGICAL GOWNS](#)

[SURGICAL PACKS & DRAPES](#)



HALYARD

FORMERLY KIMBERLY-CLARK HEALTH CARE

CLINICAL TOWELS
CATALOGUE 2015

AUSTRALIA &
NEW ZEALAND


CLINICAL TOWELS



ISOWIPE* BACTERICIDAL WIPES

- Isopropyl Alcohol
- Bactericidal Action
- Measured Dispensing System
- Disinfection Wipe for Clean Surfaces


Code	Size	Description	Packaging
6835	42cm x 14cm	ISOWIPE* Bactericidal Wipes Canister	Eaches/Canister = 75 Canisters/Case = 12 Eaches/Case = 900
6836	42cm x 14cm	ISOWIPE* Refill Pack	Eaches/Refill = 75 Refills/Case = 12 Eaches/Case = 900
6837	21cm x 14cm	ISOWIPE* Minis Canister	Eaches/Canister = 150 Canisters/Case = 12 Eaches/Case = 1800
6838	21cm x 14cm	ISOWIPE* Minis Refill Pack	Eaches/Refill = 150 Refills/Case = 12 Eaches/Case = 1800
4989		Stainless Steel Wall Bracket	1



HALYARD* WASHCLOTH EXTRA

- Single Patient Use
- Soft on Skin
- Strong – Wet or Dry – and Ultra Absorbent
- Resealable Bags

Code	Size	Description	Packaging
2711	35cm x 43cm	HALYARD* Washcloth Extra	Eaches/Pack = 20 Packs/Case = 30 Eaches/Case = 600




TERI* WIPERS

- Individually Folded
- Strong Wet Strength with Scrim Reinforcement
- Soft and Absorbent Even After Sterilisation

Code	Size	Description	Packaging
4465	32.5cm x 33cm	TERI* Wipers Small	Eaches/Pack = 75 Packs/Case = 6 Eaches/Case = 450
4467	32.5cm x 33cm	TERI* Wipers Small Bulk Pack	Eaches/Pack = 450
4466	32.5cm x 60cm	TERI* Wipers Large	Eaches/Pack = 75 Packs/Case = 6 Eaches/Case = 450
4907		Clear Acrylic Dispenser Suits Code 4467	1


CLINICAL TOWELS



TERI*HYDROKNIT* WIPERS

- Super Absorbency
- Soft, Cloth-Like Feel
- High Strength
- Tear-Resistant


Code	Size	Description	Packaging
4476	32.5cm x 33cm	TERI*HYDROKNIT Wipers Small	Eaches/Pack = 75 Packs/Case = 6 Eaches/Case = 450
4478	32.5cm x 60cm	TERI*HYDROKNIT Wipers Large	Eaches/Pack = 75 Packs/Case = 6 Eaches/Case = 450



VERSA*TOWEL

- Soft and Bulky
- Crepe Wadding Provides Excellent Absorbency
- Strong Even When Wet
- 100 Sheets on Perforated Roll

Code	Size	Description	Packaging
4210	24.5cm x 41.5cm Sheet Size	VERSA*Towel Small	Rolls/Case = 16
4220	49cm x 41.5cm Sheet Size	VERSA*Towel Large	Rolls/Case = 8
4921		White Plastic Dispenser Suits Code 4210	1
4922		White Plastic Dispenser Suits Code 4220	1



HALYARD* ABSORBENT TRAY LINERS

- High Absorbency
- High Strength
- Cut To Suit Sterilisation Trays

Code	Size	Description	Packaging
4330	21.5cm x 24cm	HALYARD* Absorbent Tray Liners Small	Eaches/Pack = 125 Packs/Case = 4 Eaches/Case = 500
4331	24cm x 36.5cm	HALYARD* Absorbent Tray Liners Medium	Eaches/Pack = 75 Packs/Case = 6 Eaches/Case = 450
4332	30cm x 50.5cm	HALYARD* Absorbent Tray Liners Large	Eaches/Pack = 75 Packs/Case = 6 Eaches/Case = 450



HALYARD* NEUTRAL HAND WIPES

- Neutral Ingredients To Clean Visibly Soiled Hands
- No Alcohol
- Contains Vitamin E
- Resealable Soft Pack

Code	Size	Description	Packaging
8835	20cm x 20cm	HALYARD* Neutral Hand Wipes	Eaches/Pack = 15 Packs/Case = 40 Eaches/Case = 600

NOTE: Please see package inserts for directions for use for each of these products.



HALYARD

FORMERLY KIMBERLY-CLARK HEALTH CARE

DIGESTIVE HEALTH
ENTERAL FEEDING
CATALOGUE 2015

AUSTRALIA &
NEW ZEALAND


ENTERAL FEEDING

HALYARD* MIC* PERCUTANEOUS ENDOSCOPIC GASTROSTOMY (PEG) FEEDING TUBES/KIT

- Medical Grade Silicone Construction
- SECUR-LOK* External Retention Ring
- Universal and Bolus Feeding Port Connectors
- Medication Port
- Collapsible Internal Retention Bumper
- Radiopaque Stripe and Bumper
- Tubing Clamp
- Available in PUSH (OTW) and PULL Placement Methods
- Ethylene Oxide (ETO) Sterilised

Deluxe 0630-XX Series PUSH PEG Kit Contents:


- 1- Traction Removable PUSH PEG Tube
- 1- 5 Fr Retrieval Snare (1.8mm x 260cm)
- 1- 0136-XX Bolus Feeding Adapter
- 1- 0137-XX Threaded Feeding Adapter
- 1- SECUR-LOK* External Retention Ring
- 1- 22ga x 1 1/2", 3.81cm Filter Needle
- 1- 25ga x 1.5875cm Needle
- 1- 25ga x 1.5875cm Needle
- 1- Fenestrated Drape with Adhesive Tabs
- 1- 12ml Syringe
- 1- Disposable Scalpel with #11 Blade
- 1- 14ga x 2.75", 6.985cm Introducer Needle
- 1- Guidewire (.038" x 260cm)
- 1- Tubing Clamp
- 4- 4" x 4", 10.16cm x 10.16cm Gauze Pads
- 2- 2" x 2", 5.08cm x 5.08cm in Slit Gauze Pads
- 1- Suture Strand
- 1- Directions for Use



PUSH

Deluxe 0640-XX Series PULL PEG Kit Contents:

- 1- Traction Removable PULL PEG Tube
- 1- 5 Fr Retrieval Snare (1.8mm x 260cm)
- 1- 0136-XX Bolus Feeding Adapter
- 1- 0137-XX Threaded Feeding Adapter
- 1- SECUR-LOK* External Retention Ring
- 1- 22ga x 1.5", 3.81cm Filter Needle
- 1- 19ga x 1.5", 3.81cm Filter Needle
- 1- 25ga x 1.5875cm Needle
- 1- Fenestrated Drape with Adhesive Tabs
- 1- 12ml Syringe
- 1- Disposable Scalpel with #11 Blade
- 1- 14ga x 2.75", 6.985cm Introducer Needle
- 1- Looped Placement Wire
- 1- Tubing Clamp
- 4- 4" x 4", 10.16cm x 10.16cm Gauze Pads
- 2- 2" x 2", 5.08cm x 5.08cm in Slit Gauze Pads
- 1- Suture Strand
- 1- Directions for Use




PULL

Code	Size (Outer Diameter)	Description	Packaging
0630-14	14 Fr	PUSH (Over-the-Wire) Method	1 Each
0630-20	20 Fr	PUSH (Over-the-Wire) Method	1 Each
0630-24	24 Fr	PUSH (Over-the-Wire) Method	1 Each
0640-14	14 Fr	PULL Method	1 Each
0640-20	20 Fr	PULL Method	1 Each
0640-24	24 Fr	PULL Method	1 Each

Deluxe 0644-XX Series Basic PULL PEG Kit Contents:

- 1- Traction Removable PULL PEG Tube
- 1- 0136-XX Bolus Feeding Adapter
- 1- 0137-XX Threaded Feeding Adapter
- 1- SECUR-LOK* External Retention Ring
- 1- Disposable Scalpel with #11 Blade
- 1- 14ga x 6.985cm (2.75") Introducer Needle
- 1- Looped Placement Wire
- 1- Tubing Clamp
- 1- Directions for Use



PULL


Code	Size (Outer Diameter)	Description	Packaging
0644-14	14 Fr	Basic PULL Peg Kit	1 Each
0644-20	20 Fr	Basic PULL Peg Kit	1 Each
0644-24	24 Fr	Basic PULL Peg Kit	1 Eac

NOTE: Please see package inserts for directions for use for each of these products.

ENTERAL FEEDING

HALYARD* MIC* PERCUTANEOUS ENDOSCOPIC GASTROSTOMY FEEDING TUBE ACCESSORIES

For use with MIC* Percutaneous Endoscopic Gastrostomy Feeding Tubes




Code	Description	Packaging
0135-14	Universal Feeding Adapter – For use with 14 Fr MIC* PEG	1 Each
0135-20	Universal Feeding Adapter – For use with 20 Fr MIC* PEG	1 Each
0135-24	Universal Feeding Adapter – For use with 24 Fr MIC* PEG	1 Each
0136-14	Bolus Feeding Adapter – For use with 14 Fr MIC* PEG	1 Each
0136-20	Bolus Feeding Adapter – For use with 20 Fr MIC* PEG	1 Each
0136-24	Bolus Feeding Adapter – For use with 24 Fr MIC* PEG	1 Each

NOTE: Please see package inserts for directions for use for each of these products.

ENTERAL FEEDING

HALYARD* MIC* GASTROSTOMY FEEDING TUBES


- Medical Grade Silicone Construction
- Silicone Internal Retention Balloon
- SECUR-LOK* External Retention Ring
- Universal Feeding Port Connector
- Medication Port
- Tapered Distal Tip
- Dual Exit Ports
- Radiopaque Stripe
- Gamma Sterilised



Code	Size (Outer Diameter)	Balloon Volume	Packaging
0112-12LV	12 Fr	3-5ml Balloon	1 Each
0112-14LV	14 Fr	3-5ml Balloon	1 Each
0112-16LV	16 Fr	3-5ml Balloon	1 Each
0112-18	18 Fr	7-10ml Balloon	1 Each
0112-20	20 Fr	7-10ml Balloon	1 Each
0112-22	22 Fr	7-10ml Balloon	1 Each
0112-24	24 Fr	7-10ml Balloon	1 Each

HALYARD* MIC* BOLUS GASTROSTOMY FEEDING TUBES

- Medical Grade Silicone Construction
- Silicone Internal Retention Balloon
- SECUR-LOK* External Retention Ring
- Tapered Distal Tip
- Dual Exit Ports
- Radiopaque Stripe
- Gamma Sterilised




Code	Size (Outer Diameter)	Balloon Volume	Packaging
0100-14LV	14 Fr	3-5ml Balloon	1 Each
0100-16LV	16 Fr	3-5ml Balloon	1 Each
0100-18	18 Fr	7-10ml Balloon	1 Each
0100-20	20 Fr	7-10ml Balloon	1 Each
0100-22	22 Fr	7-10ml Balloon	1 Each
0100-24	24 Fr	7-10ml Balloon	1 Each
0100-26	26 Fr	7-10ml Balloon	1 Each
0100-28	28 Fr	7-10ml Balloon	1 Each
0100-30	30 Fr	7-10ml Balloon	1 Each

NOTE: Please see package inserts for directions for use for each of these products.

ENTERAL FEEDING

HALYARD* MIC* GASTROSTOMY FEEDING TUBES

- Medical Grade Silicone Construction for Visibility and Drapability
- Clearly Marked Feeding, Medication and Balloon Ports for Ease of Identification
- Ventilated External Retention Ring Designed to Allow Air to Circulate Around the Stoma Site and to Reduce Pressure 360° from the Stoma
- Radiopaque Stripe to Aid in Catheter Visualisation
- Graduated Centimetre Markings for Tube Position Verification
- LUER-LOK* Balloon Port to Ensure Secure Syringe Connection (also Compatible with LUER-SLIP* Syringes)
- Colour-Coded Balloon Port According to Catheter French Size for Ease of Identification



Code	Size (Outer Diameter)	Balloon Volume	Packaging
0110-14LV	14 Fr	3-5ml Balloon	1 Each
0110-16LV	16 Fr	3-5ml Balloon	1 Each
0110-18	18 Fr	7-10ml Balloon	1 Each
0110-20	20 Fr	7-10ml Balloon	1 Each
0110-22	22 Fr	7-10ml Balloon	1 Each
0110-24	24 Fr	7-10ml Balloon	1 Each

NOTE: Please see package inserts for directions for use for each of these products.


ENTERAL FEEDING

HALYARD* MIC-KEY* LOW-PROFILE GASTROSTOMY FEEDING TUBES/KITS

- Medical Grade Silicone Construction
- Low-Profile Design
- Tapered Distal Tip
- Silicone Internal Retention Balloon
- Distal Tip Recessed at 5ml
- Proximal Anti-Reflux Valve
- SECUR-LOK* Extension Set Connector Mechanism
- Radiopaque Stripe
- Wide Variety of Extension Sets Available
- Gamma Sterilised

Kit Includes:

- 1- MIC-KEY* Low-Profile Feeding Tube
- 1- MIC-KEY* Extension Set with SECUR-LOK*
Right Angle Connector and 2 Port "Y" and Clamp – 12", 30.5cm Length
- 1- MIC-KEY* Bolus Extension Set with Cath Tip, SECUR-LOK*
Straight Connector and Clamp – 12", 30.5cm Length
- 1- 6ml Syringe
- 1- 35ml Catheter Tip Syringe
- 4- Gauze Pads




Code	Size (Outer Diameter)	Stoma Length	Balloon Volume	Packaging
0120-12-0.8	12 Fr	0.8cm	3-5ml Balloon	1 Each
0120-12-1.0	12 Fr	1cm	3-5ml Balloon	1 Each
0120-12-1.2	12 Fr	1.2cm	3-5ml Balloon	1 Each
0120-12-1.5	12 Fr	1.5cm	3-5ml Balloon	1 Each
0120-12-1.7	12 Fr	1.7cm	3-5ml Balloon	1 Each
0120-12-2.0	12 Fr	2cm	3-5ml Balloon	1 Each
0120-12-2.3	12 Fr	2.3cm	3-5ml Balloon	1 Each
0120-12-2.5	12 Fr	2.5cm	3-5ml Balloon	1 Each
0120-12-2.7	12 Fr	2.7cm	3-5ml Balloon	1 Each
0120-12-3.0	12 Fr	3cm	3-5ml Balloon	1 Each
0120-12-3.5	12 Fr	3.5cm	3-5ml Balloon	1 Each
0120-12-4.0	12 Fr	4cm	3-5ml Balloon	1 Each
0120-12-4.5	12 Fr	4.5cm	3-5ml Balloon	1 Each
0120-14-0.8	14 Fr	0.8cm	5ml Balloon	1 Each
0120-14-1.0	14 Fr	1.0cm	5ml Balloon	1 Each
0120-14-1.2	14 Fr	1.2cm	5ml Balloon	1 Each
0120-14-1.5	14 Fr	1.5cm	5ml Balloon	1 Each
0120-14-1.7	14 Fr	1.7cm	5ml Balloon	1 Each
0120-14-2.0	14 Fr	2.0cm	5ml Balloon	1 Each
0120-14-2.3	14 Fr	2.3cm	5ml Balloon	1 Each
0120-14-2.5	14 Fr	2.5cm	5ml Balloon	1 Each
0120-14-2.7	14 Fr	2.7cm	5ml Balloon	1 Each
0120-14-3.0	14 Fr	3.0cm	5ml Balloon	1 Each
0120-14-3.5	14 Fr	3.5cm	5ml Balloon	1 Each
0120-14-4.0	14 Fr	4.0cm	5ml Balloon	1 Each
0120-14-4.5	14 Fr	4.5cm	5ml Balloon	1 Each
0120-14-5.0	14 Fr	5cm	5ml Balloon	1 Each

Continued over page

ENTERAL FEEDING

HALYARD* MIC-KEY* LOW-PROFILE GASTROSTOMY FEEDING TUBES/KITS (cont'd)




Code	Size (Outer Diameter)	Stoma Length	Balloon Volume	Packaging
0120-16-0.8	16 Fr	0.8cm	5ml Balloon	1 Each
0120-16-1.0	16 Fr	1.0cm	5ml Balloon	1 Each
0120-16-1.2	16 Fr	1.2cm	5ml Balloon	1 Each
0120-16-1.5	16 Fr	1.5cm	5ml Balloon	1 Each
0120-16-1.7	16 Fr	1.7cm	5ml Balloon	1 Each
0120-16-2.0	16 Fr	2.0cm	5ml Balloon	1 Each
0120-16-2.3	16 Fr	2.3cm	5ml Balloon	1 Each
0120-16-2.5	16 Fr	2.5cm	5ml Balloon	1 Each
0120-16-2.7	16 Fr	2.7cm	5ml Balloon	1 Each
0120-16-3.0	16 Fr	3.0cm	5ml Balloon	1 Each
0120-16-3.5	16 Fr	3.5cm	5ml Balloon	1 Each
0120-16-4.0	16 Fr	4.0cm	5ml Balloon	1 Each
0120-16-4.5	16 Fr	4.5cm	5ml Balloon	1 Each
0120-16-5.0	16 Fr	5cm	5ml Balloon	1 Each
0120-18-0.8	18 Fr	0.8cm	5ml Balloon	1 Each
0120-18-1.0	18 Fr	1.0cm	5ml Balloon	1 Each
0120-18-1.2	18 Fr	1.2cm	5ml Balloon	1 Each
0120-18-1.5	18 Fr	1.5cm	5ml Balloon	1 Each
0120-18-1.7	18 Fr	1.7cm	5ml Balloon	1 Each
0120-18-2.0	18 Fr	2.0cm	5ml Balloon	1 Each
0120-18-2.3	18 Fr	2.3cm	5ml Balloon	1 Each
0120-18-2.5	18 Fr	2.5cm	5ml Balloon	1 Each
0120-18-2.7	18 Fr	2.7cm	5ml Balloon	1 Each
0120-18-3.0	18 Fr	3.0cm	5ml Balloon	1 Each
0120-18-3.5	18 Fr	3.5cm	5ml Balloon	1 Each
0120-18-4.0	18 Fr	4.0cm	5ml Balloon	1 Each
0120-18-4.5	18 Fr	4.5cm	5ml Balloon	1 Each
0120-18-5.0	18 Fr	5.0cm	5ml Balloon	1 Each

Continued over page

ENTERAL FEEDING

HALYARD* MIC-KEY* LOW-PROFILE GASTROSTOMY FEEDING TUBES/KITS (cont'd)



Code	Size (Outer Diameter)	Stoma Length	Balloon Volume	Packaging
0120-20-0.8	20 Fr	0.8cm	5ml Balloon	1 Each
0120-20-1.0	20 Fr	1.0cm	5ml Balloon	1 Each
0120-20-1.2	20 Fr	1.2cm	5ml Balloon	1 Each
0120-20-1.5	20 Fr	1.5cm	5ml Balloon	1 Each
0120-20-1.7	20 Fr	1.7cm	5ml Balloon	1 Each
0120-20-2.0	20 Fr	2.0cm	5ml Balloon	1 Each
0120-20-2.3	20 Fr	2.3cm	5ml Balloon	1 Each
0120-20-2.5	20 Fr	2.5cm	5ml Balloon	1 Each
0120-20-2.7	20 Fr	2.7cm	5ml Balloon	1 Each
0120-20-3.0	20 Fr	3.0cm	5ml Balloon	1 Each
0120-20-3.5	20 Fr	3.5cm	5ml Balloon	1 Each
0120-20-4.0	20 Fr	4.0cm	5ml Balloon	1 Each
0120-20-4.5	20 Fr	4.5cm	5ml Balloon	1 Each
0120-20-5.0	20 Fr	5.0cm	5ml Balloon	1 Each
0120-24-1.5	24 Fr	1.5cm	5ml Balloon	1 Each
0120-24-1.7	24 Fr	1.7cm	5ml Balloon	1 Each
0120-24-2.0	24 Fr	2.0cm	5ml Balloon	1 Each
0120-24-2.3	24 Fr	2.3cm	5ml Balloon	1 Each
0120-24-2.5	24 Fr	2.5cm	5ml Balloon	1 Each
0120-24-2.7	24 Fr	2.7cm	5ml Balloon	1 Each
0120-24-3.0	24 Fr	3.0cm	5ml Balloon	1 Each
0120-24-3.5	24 Fr	3.5cm	5ml Balloon	1 Each
0120-24-4.0	24 Fr	4.0cm	5ml Balloon	1 Each
0120-24-4.5	24 Fr	4.5cm	5ml Balloon	1 Each
0120-24-5.0	24 Fr	5.0cm	5ml Balloon	1 Each


NOTE: Please see package inserts for directions for use for each of these products.

ENTERAL FEEDING

HALYARD* MIC-KEY* FEEDING TUBE ACCESSORIES

For use with MIC-KEY* Low-Profile Gastrostomy and Jejunal Feeding Tubes

- Proximal Tip Types: 2 Port "Y", Bolus, and Threaded
- SECUR-LOK* "Key" Connection Mechanism
- Clear Tubing
- Medication Port on Non-Bolus Extension Sets
- Individually Packaged Extension Sets
- Tubing Clamp
- Extension Sets Non-Sterile
- Stoma Measuring Device Gamma Sterilised



DEHP-Free Formulation Code	Standard Formulation Code	Length	Description	Packaging
0121-12	0221-12	30.5cm 12"	MIC-KEY* Extension Set with SECUR-LOK* Right Angle Connector and 2 Port "Y" and Clamp	5 Each
0121-24	0221-24	60.9cm 24"	MIC-KEY* Extension Set with SECUR-LOK* Right Angle Connector and 2 Port "Y" and Clamp	5 Each
0122-02		5.1cm 2"	MIC-KEY* Medication Set with SECUR-LOK* Right Angle Connector and 2 Port "Y"	5 Each
0123-12	0223-12	30.5cm 12"	MIC-KEY* Bolus Extension Set with Cath Tip, SECUR-LOK* Straight Connector and Clamp	5 Each
0123-24		60.9cm 24"	MIC-KEY* Bolus Extension Set with Cath Tip, SECUR-LOK* Straight Connector and Clamp	5 Each
0124-12		30.5cm 12"	MIC-KEY* Bolus Extension Set with Cath Tip, SECUR-LOK* Right Angle Connector and Clamp	5 Each
0124-24		60.9cm 24"	MIC-KEY* Bolus Extension Set with Cath Tip, SECUR-LOK* Right Angle Connector and Clamp	5 Each
0126-12		30.5cm 12"	MIC-KEY* Threaded Extension Set With SECUR-LOK* Right Angle Connector and Clamp	5 Each
0126-24		60.9cm 24"	MIC-KEY* Threaded Extension Set With SECUR-LOK* Right Angle Connector and Clamp	5 Each
98460			MIC-KEY* Stoma Measuring Device	10 Each

NOTE: Please see package inserts for directions for use for each of these products.


ENTERAL FEEDING

HALYARD* MIC* TRANSGASTRIC-JEJUNAL FEEDING TUBES – ENDOSCOPIC / RADIOLOGY PLACEMENT

- Medical Grade Silicone Construction
- Silicone Internal Retention Balloon
- SECUR-LOK* External Retention Ring
- Universal Jejunal Feeding Port Connector
- Gastric Decompression Port
- Multiple Gastric & Jejunal Exit Ports
- Tungsten Weighted Jejunal Portion
- Tapered Distal Tip
- Radiopaque Stripe
- Ethylene Oxide (ETO) Sterilised

Kit Includes:

- 1- Transgastric Jejunal Feeding Tube
- 1- Silicone Coated Guidewire
- 1- LUER-SLIP* Syringe
- 1- 35ml Catheter Tip Syringe
- 1- SECUR-LOK* External Retention Ring
- 1- Directions for Use
- 1- Patient Care Guidelines



Code	Size (Outer Diameter)	Jejunal Length†	Balloon Volume	Packaging
0650-16-15	16 Fr	15cm	7-10ml Balloon	1 Each
0650-16-22	16 Fr	22cm	7-10ml Balloon	1 Each
0650-16-30	16 Fr	30cm	7-10ml Balloon	1 Each
0650-16	16 Fr	45cm	7-10ml Balloon	1 Each
0650-18-22	18 Fr	22cm	7-10ml Balloon	1 Each
0650-18-30	18 Fr	30cm	7-10ml Balloon	1 Each
0650-18	18 Fr	45cm	7-10ml Balloon	1 Each
0650-22	22 Fr	45cm	7-10ml Balloon	1 Each

† Jejunal length measured from the distal aspect of the retention balloon to the catheter tip.


ENTERAL FEEDING

HALYARD* MIC-KEY* LOW-PROFILE TRANSGASTRIC-JEJUNAL TUBES – ENDOSCOPIC/RADIOLOGY PLACEMENT KITS

- Unobtrusive and easy to conceal
- Multiple feeding exit ports to help improve flow and minimise clogging
- Tungsten weighted jejunal portion to maintain jejunal position
- Radiopaque stripe to help aid in catheter visualisation
- Clearly marked separate gastric and jejunal ports for ease of identification
- High-clarity medical grade silicone designed for visibility and drapability
- Tapered distal tip for non-traumatic tube insertion

Kit Includes:

- 1- MIC-KEY* Low-Profile Transgastric-Jejunal Feeding Tube
- 1- Introducer Cannula
- 1- MIC-KEY* Extension Set with SECUR-LOK* Right Angle Connector and 2 Port "Y" and Clamp – 12", 30.5cm Length
- 1- MIC-KEY* Bolus Extension Set with Cath Tip, SECUR-LOK* Straight Connector and Clamp – 12", 30.5cm Length
- 1- 6ml LUER-SLIP* Syringe
- 4- Gauze Pads
- 1- Directions For Use
- 1- Patient Care Guidelines




Code	Size (Outer Diameter)	Stoma Length	Jejunal Length†	Balloon Volume	Packaging
0270-16-1.0-15	16 Fr	1.0cm	15cm	3-5ml Balloon	1 Each
0270-16-1.0-22	16 Fr	1.0cm	22cm	3-5ml Balloon	1 Each
0270-16-1.2-15	16 Fr	1.2cm	15cm	3-5ml Balloon	1 Each
0270-16-1.2-22	16 Fr	1.2cm	22cm	3-5ml Balloon	1 Each
0270-16-1.2-30	16 Fr	1.2cm	30cm	3-5ml Balloon	1 Each
0270-16-1.5-15	16 Fr	1.5cm	15cm	3-5ml Balloon	1 Each
0270-16-1.5-22	16 Fr	1.5cm	22cm	3-5ml Balloon	1 Each
0270-16-1.5-30	16 Fr	1.5cm	30cm	3-5ml Balloon	1 Each
0270-16-1.5-45	16 Fr	1.5cm	45cm	7-10ml Balloon	1 Each
0270-16-1.7-15	16 Fr	1.7cm	15cm	3-5ml Balloon	1 Each
0270-16-1.7-22	16 Fr	1.7cm	22cm	3-5ml Balloon	1 Each
0270-16-1.7-30	16 Fr	1.7cm	30cm	3-5ml Balloon	1 Each
0270-16-1.7-45	16 Fr	1.7cm	45cm	7-10ml Balloon	1 Each
0270-16-2.0-22	16 Fr	2.0cm	22cm	3-5ml Balloon	1 Each
0270-16-2.0-30	16 Fr	2.0cm	30cm	3-5ml Balloon	1 Each
0270-16-2.0-45	16 Fr	2.0cm	45cm	7-10ml Balloon	1 Each
0270-16-2.3-30	16 Fr	2.3cm	30cm	3-5ml Balloon	1 Each
0270-16-2.3-45	16 Fr	2.3cm	45cm	7-10ml Balloon	1 Each
0270-16-2.5-30	16 Fr	2.5cm	30cm	3-5ml Balloon	1 Each
0270-16-2.7-30	16 Fr	2.7cm	30cm	3-5ml Balloon	1 Each

† Jejunal length measured from the distal aspect of the retention balloon to the catheter tip.

Continued over page

ENTERAL FEEDING

HALYARD* MIC-KEY* LOW-PROFILE TRANSGASTRIC-JEJUNAL TUBES – ENDOSCOPIC/RADIOLOGY PLACEMENT KITS (cont'd)



Code	Size (Outer Diameter)	Stoma Length	Jejunal Length†	Balloon Volume	Packaging
0270-18-1.2-22	18 Fr	1.2cm	22cm	3-5ml Balloon	1 Each
0270-18-1.2-30	18 Fr	1.2cm	30cm	3-5ml Balloon	1 Each
0270-18-1.5-22	18 Fr	1.5cm	22cm	3-5ml Balloon	1 Each
0270-18-1.5-30	18 Fr	1.5cm	30cm	3-5ml Balloon	1 Each
0270-18-1.5-45	18 Fr	1.5cm	45cm	7-10ml Balloon	1 Each
0270-18-1.7-22	18 Fr	1.7cm	22cm	3-5ml Balloon	1 Each
0270-18-1.7-30	18 Fr	1.7cm	30cm	3-5ml Balloon	1 Each
0270-18-1.7-45	18 Fr	1.7cm	45cm	7-10ml Balloon	1 Each
0270-18-2.0-22	18 Fr	2.0cm	22cm	3-5ml Balloon	1 Each
0270-18-2.0-30	18 Fr	2.0cm	30cm	3-5ml Balloon	1 Each
0270-18-2.0-45	18 Fr	2.0cm	45cm	7-10ml Balloon	1 Each
0270-18-2.3-30	18 Fr	2.3cm	30cm	3-5ml Balloon	1 Each
0270-18-2.3-45	18 Fr	2.3cm	45cm	7-10ml Balloon	1 Each
0270-18-2.5-30	18 Fr	2.5cm	30cm	3-5ml Balloon	1 Each
0270-18-2.5-45	18 Fr	2.5cm	45cm	7-10ml Balloon	1 Each
0270-18-2.7-30	18 Fr	2.7cm	30cm	3-5ml Balloon	1 Each
0270-18-2.7-45	18 Fr	2.7cm	45cm	7-10ml Balloon	1 Each
0270-18-3.0-45	18 Fr	3.0cm	45cm	7-10ml Balloon	1 Each
0270-18-3.5-45	18 Fr	3.5cm	45cm	7-10ml Balloon	1 Each
0270-22-1.5-45	22 Fr	1.5cm	45cm	7-10ml Balloon	1 Each
0270-22-1.7-45	22 Fr	1.7cm	45cm	7-10ml Balloon	1 Each
0270-22-2.0-45	22 Fr	2.0cm	45cm	7-10ml Balloon	1 Each
0270-22-2.3-45	22 Fr	2.3cm	45cm	7-10ml Balloon	1 Each
0270-22-2.5-45	22 Fr	2.5cm	45cm	7-10ml Balloon	1 Each
0270-22-2.7-45	22 Fr	2.7cm	45cm	7-10ml Balloon	1 Each
0270-22-3.0-45	22 Fr	3.0cm	45cm	7-10ml Balloon	1 Each
0270-22-3.5-45	22 Fr	3.5cm	45cm	7-10ml Balloon	1 Each

† Jejunal length measured from the distal aspect of the retention balloon to the catheter tip.

ENTERAL FEEDING

HALYARD* MIC* JEJUNAL FEEDING TUBES

Code	Size (Outer Diameter)	Balloon Volume	Packaging
0200-12LV	12 Fr	3-5ml Balloon	1 Each
0200-14	14 Fr	7-10ml Balloon	1 Each
0200-16	16 Fr	7-10ml Balloon	1 Each
0200-18	18 Fr	7-10ml Balloon	1 Each
0200-20	20 Fr	7-10ml Balloon	1 Each
0200-22	22 Fr	7-10ml Balloon	1 Each
0200-24	24 Fr	7-10ml Balloon	1 Each

HALYARD* MIC-KEY* LOW-PROFILE JEJUNAL FEEDING TUBES

Code	Size (Outer Diameter)	Stoma Length	Balloon Volume	Packaging
0230-14-0.8	14 Fr	0.8cm	5ml Balloon	1 Each
0230-14-1.0	14 Fr	1.0cm	5ml Balloon	1 Each
0230-14-1.2	14 Fr	1.2cm	5ml Balloon	1 Each
0230-14-1.5	14 Fr	1.5cm	5ml Balloon	1 Each
0230-14-1.7	14 Fr	1.7cm	5ml Balloon	1 Each
0230-14-2.0	14 Fr	2.0cm	5ml Balloon	1 Each
0230-14-2.3	14 Fr	2.3cm	5ml Balloon	1 Each
0230-14-2.5	14 Fr	2.5cm	5ml Balloon	1 Each
0230-14-2.7	14 Fr	2.7cm	5ml Balloon	1 Each
0230-14-3.0	14 Fr	3.0cm	5ml Balloon	1 Each
0230-14-3.5	14 Fr	3.5cm	5ml Balloon	1 Each
0230-14-4.0	14 Fr	4.0cm	5ml Balloon	1 Each
0230-14-4.5	14 Fr	4.5cm	5ml Balloon	1 Each
0230-18-0.8	18 Fr	0.8cm	5ml Balloon	1 Each
0230-18-1.0	18 Fr	1.0cm	5ml Balloon	1 Each
0230-18-1.2	18 Fr	1.2cm	5ml Balloon	1 Each
0230-18-1.5	18 Fr	1.5cm	5ml Balloon	1 Each
0230-18-1.7	18 Fr	1.7cm	5ml Balloon	1 Each
0230-18-2.0	18 Fr	2.0cm	5ml Balloon	1 Each
0230-18-2.3	18 Fr	2.3cm	5ml Balloon	1 Each
0230-18-2.5	18 Fr	2.5cm	5ml Balloon	1 Each
0230-18-2.7	18 Fr	2.7cm	5ml Balloon	1 Each
0230-18-3.0	18 Fr	3.0cm	5ml Balloon	1 Each
0230-18-3.5	18 Fr	3.5cm	5ml Balloon	1 Each
0230-18-4.0	18 Fr	4.0cm	5ml Balloon	1 Each
0230-18-4.5	18 Fr	4.5cm	5ml Balloon	1 Each

NOTE: Please see package inserts for directions for use for each of these products.


ENTERAL FEEDING

HALYARD* MIC* JEJUNOSTOMY FEEDING TUBE

- Medical Grade Silicone Construction
- Suture Wings
- 9 Fr Internal Diameter
- Radiopaque Stripe
- Gamma Sterilised

Kit Includes:

- 1- Jejunostomy Tube
- 1- Subcutaneous Extraducer
- 1- Feeding Port Adapter
- 1- LUER-LOK* Cap



Code	Outer Diameter / Inner Diameter	Description	Packaging
0301-14	14 Fr / 9 Fr	Jejunostomy Feeding Tube	1 Each

HALYARD* MIC* JEJUNOSTOMY FEEDING TUBE ACCESSORIES

For use with the MIC* Jejunostomy Feeding Tube

Code	Description	Packaging
0301-50	Jejunostomy Repair Kit – Includes Silicone Med A Medical Adhesive, Stent, Syringe, Blunt Tip Needle, Suture, Silicone Overcollar	1 Each
0301-55	Jejunostomy LUER-LOK* Caps	5 Each
0301-60	Jejunostomy Feeding Port Adapter	5 Each

ENTERAL FEEDING

HALYARD* MIC* FEEDING TUBE ACCESSORIES


For use with MIC* Gastrostomy, Bolus Gastrostomy, PEG, Jejunal, Gastro-Enteric, and Transgastric-Jejunal Feeding Tubes. DEHP-Free Formulation [Di-(2-Ethylhexyl)phthalate]



DEHP-Free Formulation Code	Standard Formulation Code	Length	Description	Packaging
0105-06	0205-06	15cm 6"	Extension Tubing with Bolus and Stepped Connectors at Opposite Ends	1 Each
0105-12	0205-12	30cm 12"	Extension Tubing with Bolus and Stepped Connectors at Opposite Ends	1 Each

HALYARD* FOREIGN BODY HOOD PROTECTOR

- Non-Sterile
- Disposable



Code	Description	Packaging
69392†	Bell-Shaped, Foreign Body Hood Protector	10 Each

† Caution: This product contains natural rubber latex which may cause allergic reactions



HALYARD

FORMERLY KIMBERLY-CLARK HEALTH CARE

DIGESTIVE HEALTH
INTRODUCER KITS
CATALOGUE 2015

AUSTRALIA &
NEW ZEALAND

INTRODUCER KITS

HALYARD* MIC* INTRODUCER KIT COMPONENTS

Code	Description	Packaging
98701	HALYARD* SAF-T-PEXY* T-Fasteners plus Curved Haemostat	3/Pack
98702	HALYARD* 16Fr Serial Dilator, single <i>Recommend using with 12Fr feeding tube</i>	1
98703	HALYARD* 18Fr Serial Dilator, single <i>Recommend using with 14Fr feeding tube</i>	1
98704	HALYARD* 20Fr Serial Dilator, single <i>Recommend using with 16Fr feeding tube</i>	1
98705	HALYARD* 22Fr Serial Dilator, single <i>Recommend using with 18Fr feeding tube</i>	1
98706	HALYARD* 24Fr Serial Dilator, single <i>Recommend using with 20Fr feeding tube</i>	1

NOTE: Please see package inserts for directions for use for each of these products.

INTRODUCER KITS

HALYARD* MIC*, MIC-KEY* INTRODUCER KITS

Kit Includes:

- 1- Gastrointestinal Anchor Set with 4 SAF-T-PEXY* T-Fasteners
 - 1- Haemostat
 - 1- 18 gauge Safety Introducer Needle
 - 1- .035", .9mm J-tip Guidewire
 - 1- #11 Scalpel
 - 1- 12ml Syringe
 - 1- Serial Dilator with Peel-Away Sheath
- Feeding tube not included in kit.



Code	Description	Packaging
98420	MIC* G -12 Introducer Kit Includes 16 Fr Dilator for use with MIC* 12 Fr Gastrostomy Feeding Tube	1 Kit/Case
98421	MIC* G -14 Introducer Kit Includes 18 Fr Dilator for use with MIC* 14 Fr Gastrostomy Feeding Tube	1 Kit/Case
98422	MIC* G -16 Introducer Kit Includes 20 Fr Dilator for use with MIC* 16 Fr Gastrostomy Feeding Tube	1 Kit/Case
98423	MIC* G -18 Introducer Kit Includes 22 Fr Dilator for use with MIC* 18 Fr Gastrostomy Feeding Tube	1 Kit/Case
98424	MIC* G -20 Introducer Kit Includes 24 Fr Dilator for use with MIC* 20 Fr Gastrostomy Feeding Tube	1 Kit/Case

HALYARD* MIC* J / TJ INTRODUCER KITS

Kit Includes:

- 1- Gastrointestinal Anchor set with 4 SAF-T-PEXY* T-Fasteners
 - 1- Haemostat
 - 1- 18 gauge Safety Introducer Needle
 - 1- Seeking Catheter 0.38", .97mm Guidewire Compatible
 - 1- #11 Scalpel
 - 1- 12ml Syringe
 - 1- Serial Dilator With Peel-Away Sheath
- Feeding tube not included in kit.



Code	Description	Packaging
98435	MIC* J / TJ - 16 Introducer Kit Includes 20 Fr Dilator for use with MIC* 16 Fr Jejunal / Transgastric Jejunal Feeding Tube	1 Kit/Case
98436	MIC* J / TJ - 18 Introducer Kit Includes 22 Fr Dilator for use with MIC* 18 Fr Jejunal / Transgastric Jejunal Feeding Tube	1 Kit/Case



HALYARD

FORMERLY KIMBERLY-CLARK HEALTH CARE


FACIAL PROTECTION
CATALOGUE 2015

AUSTRALIA &
NEW ZEALAND

FACIAL PROTECTION




FLUID PROTECTION FACE MASKS

	Code	Description	Packaging
 47147	47107	FLUIDSHIELD* Fog-Free Procedure Mask, Pleat-Style with Earloops, Foam Band, Orange	Eaches/Box = 40 Boxes/Case = 10 Eaches/Case = 400
	47147	FLUIDSHIELD* Fog-Free Procedure Mask with Wraparound SPLASHGUARD* Visor, Pleat-Style with Earloops, Foam Band, Orange	Eaches/Box = 25 Boxes/Case = 4 Eaches/Case = 100
 47147	48207	FLUIDSHIELD* Fog-Free Surgical Mask, Pleat-Style with Ties, Orange	Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300
	48208	FLUIDSHIELD* Surgical Mask, Pleat-Style with Ties, Orange	Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300
 48247	48227	FLUIDSHIELD* DUCKBILL* Surgical Mask, Pleat-Style with Ties, Orange	Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300
	48247	FLUIDSHIELD* Fog-Free Surgical Mask with Wraparound SPLASHGUARD* Visor, Pleat-Style with Ties, Orange	Eaches/Box = 25 Boxes/Case = 4 Eaches/Case = 100
	48248	FLUIDSHIELD* Fog-Free Surgical Mask with Wraparound SPLASHGUARD* Visor, Pleat-Style with Ties and Foam Band, Orange	Eaches/Box = 25 Boxes/Case = 4 Eaches/Case = 100
	48297	FLUIDSHIELD* Care Bear Surgical Mask, Pleat-Style with Ties, Pink And Blue Teddy Bear Print	Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300
	62113	THE PROTECTOR* Fog-Free Splash-Resistant Surgical Mask, Pleat-Style with Ties, Blue with Orange Diamond Print	Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300
	62114	The Protector* Fog-Free Splash-Resistant Surgical Mask with Wraparound Splashguard* Visor, Pleat-Style with Ties, Blue with Orange Diamond Print	Eaches/Box = 25 Boxes/Case = 4 Eaches/Case = 100

NOTE: Please see package inserts for directions for use for each of these products.

FACIAL PROTECTION

N95 RESPIRATOR FACE MASKS



Code	Description	Packaging
46727	FLUIDSHIELD* PFR95* N95 Particulate Filter Respirator and Surgical Mask, Pouch-Style with Polyurethane Headbands, Orange	Eaches/Box = 35 Boxes/Case = 6 Eaches/Case = 210
46827	FLUIDSHIELD* PFR95* N95 Particulate Filter Respirator and Surgical Mask, Pouch-Style with Polyurethane Headbands, Small Size, Orange	Eaches/Box = 35 Boxes/Case = 6 Eaches/Case = 210



N95 RESPIRATOR ACCESSORIES

- Contain Instructions and Materials to Conduct Fit Testing for HALYARD PFR95* N95 Respirators


Code	Description	Packaging
47950	Qualitative Fit Test Kit	Eaches/Kit = 1
99003	Saccharin Solutions	Eaches/Box = 2 Boxes/Case = 1
99009	BITREX* Solutions	Eaches/Box = 2 Boxes/Case = 1
99018	Nebulisers	Eaches/Box = 2 Boxes/Case = 1
99175	Ported FLUIDSHIELD* PFR95* Particulate Filter Respirator and Surgical Mask with Polyurethane Headbands, Regular Size, Orange	Eaches/Box = 100 Boxes/Case = 1 Eaches/Case = 100
99181	Ported FLUIDSHIELD* PFR95* Particulate Filter Respirator and Surgical Mask with Polyurethane Headbands, Small Size, Orange	Eaches/Box = 100 Boxes/Case = 1 Eaches/Case = 100

NOTE: Please see package inserts for directions for use for each of these products.

FACIAL PROTECTION

STANDARD FACE MASKS


SURGICAL MASKS

	Code	Description	Packaging
	6000	HALYARD* Surgical Mask with Ties	Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300
	48100	THE LITE ONE* Surgical Mask, Pleat-Style with Ties, Blue	Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300
	48201	CLASSIC* Surgical Mask, Pleat-Style with Ties, Blue	Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300
	48210	Nurse's Surgical Mask, Pleat-Style with Ties, Blue	Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300
	48220	DUCKBILL* Surgical Mask, Pouch-Style with Ties, Blue	Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300
	48390	SO SOFT* Surgical Mask, Pleat-Style with Ties, White	Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300
	49214	Fog-Free Surgical Mask, Pleat-Style with Ties, Blue	Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300
	49215	Anti-Fog* Surgical Mask with DERMA-TOUCH* Tape, Film Strip, Green	Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300
	49216	Fog-Free DUCKBILL* Surgical Mask, Pouch-Style with Ties, Blue	Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300
	49235	Anti-Fog* Surgical Mask, Pleat-Style with Ties, Green	Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300
	62467	Clean Room Classic* Mask, Polythylene Film with Soft Ties, Sterile, Individually Wrapped	Eaches/Bag = 20 Bags/Case = 10 Eaches/Case = 200

NOTE: Please see package inserts for directions for use for each of these products.

FACIAL PROTECTION

PROCEDURE MASKS




Code	Description	Packaging
47080	Procedure Mask, Pleat-Style with Earloops, Blue	Eaches/Box = 50 Boxes/Case = 10 Eaches/Case = 500




47127	Child's Face Mask, Pleat-Syle with Earloops, Multi-Colour	Eaches/Box = 75 Boxes/Case = 10 Eaches/Case = 750
-------	---	---



47295	Teddy Bear Procedure Mask, Pleat-Style with Earloops, Pink And Blue Teddy Bear Print	Eaches/Box = 50 Boxes/Case = 10 Eaches/Case = 500
-------	--	---



6001	HALYARD* Procedure* Mask With Ear Loops, Blue	Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300
------	---	--



62356	THE LITE ONE* Procedure Mask, Pleat-Style with Earloops, Blue	Eaches/Box = 50 Boxes/Case = 10 Eaches/Case = 500
-------	---	---




62363	SO SOFT* Fog-Free Procedure Mask, Pleat-Style with Earloops, White	Eaches/Box = 50 Boxes/Case = 10 Eaches/Case = 500
-------	--	---



62470	CR Classic* Mask, Polythylene Film and Ear Loops Sterile, Individually Wrapped	Eaches/Bag = 20 Bags/Case = 10 Eaches/Case = 200
-------	--	--


MASK DISPENSER WALL BRACKET




Code	Description	Packaging
4966	Wall Bracket, Width 100mm	1
4967	Wall Bracket, Width 125mm	1
4968	Wall Bracket, Width 190mm	1

FACIAL PROTECTION

FACE MASK ACCESSORIES



Code	Description	Packaging
41204	GUARDALL SHIELD* Full-Length Visor Only (No Mask), Clear	Eaches/Set = 1 Sets/Case = 40 Eaches/Case = 40



Code	Description	Packaging
41205	GUARDALL SHIELD* Three-Quarter Length Visor Only (No Mask), Clear	Eaches/Set = 1 Sets/Case = 40 Eaches/Case = 40

EYEWEAR



Code	Description	Packaging
SV50A	SAFEVIEW* Assembled Glasses	5 Pkg. of 10 Each Case = 50 Assembled Glasses
SV100F	SAFEVIEW* Frames / Assorted Colours	10 Pkg. of 10 Each Case = 100 Frames
SV250L	SAFEVIEW* Replacement Lenses	10 Pkg. of 25 Each Case = 250 Lenses

NOTE: Please see package inserts for directions for use for each of these products.




HALYARD

FORMERLY KIMBERLY-CLARK HEALTH CARE

FLUID ABSORPTION
CATALOGUE 2015

AUSTRALIA &
NEW ZEALAND


FLUID ABSORPTION



HALYARD* BEDSHEET

- Soft Yet Thick and Durable
- Strong and Hygienic Cover
- Continuous Roll


Code	Size	Description	Packaging
4260	54cm x 80m	HALYARD* Bedsheet Roll	Rolls/Case = 6
4924		Dispenser White Steel Suits Code 4260	1



HALYARD* BENCHROLL

- Laminated 3-Ply Paper Tissue, 1-Ply Polyethylene
- Strong And Soft
- Absorbent
- Waterproof
- Continuous Roll

Code	Size	Description	Packaging
2911	41.5cm x 91m	HALYARD* Benchroll 4-Ply Roll	Rolls/Case = 2




HALYARD* PROTECTAPAD

- 3-Ply Paper Tissue, 1-Ply Soft Polyethylene Laminated Together
- Waterproof
- Absorbent
- Ideal as Dental Bib or Pillow Protector

Code	Size	Description	Packaging
2706	28.5cm x 21.5cm	HALYARD* Protectapad Small	Eaches/Pack = 100 Packs/Case = 8 Eaches/Case = 800
2705	28.5cm x 43cm	HALYARD* Protectapad Large	Eaches/Pack = 100 Packs/Case = 4 Eaches/Case = 40


FLUID ABSORPTION



HALYARD* UNDERPAD

- Protects Bedsheets, Mattresses, Chairs and Benches
- Highly Absorbent
- Layers of Tissue with Polyethylene Backing

Code	Size	Description	Packaging
2740	39.5cm x 55cm	HALYARD* Underpads 4-Ply Pad	Eaches/Case = 300
2765	39.5cm x 55cm	HALYARD* Underpads 5-Ply Pad	Eaches/Case = 300
2760	39.5cm x 55cm	HALYARD* Underpads 8-Ply Pad	Eaches/Case = 250



HALYARD* BED PROTECTOR

- Unique Design incorporating Super Absorbent Material
- Double-Layered Base with Strong Poly-Coated Cover that Prevents Strike-Through

Code	Size	Description	Packaging
4261	81cm x 78cm	HALYARD* Bed Protector	Eaches/Case = 25



HALYARD

FORMERLY KIMBERLY-CLARK HEALTH CARE

MEDICAL GLOVES
CATALOGUE 2015

AUSTRALIA &
NEW ZEALAND

MEDICAL GLOVES

NITRILE EXAM GLOVES

HALYARD* PURPLE NITRILE* EXAM GLOVES

- Powder-Free
- Latex-Free
- Non-Sterile
- Single Use Only
- Length: 9.5" (242mm)
- Beaded Cuff
- Ambidextrous
- CE Marked
- Purple
- CE Marked
- Textured Fingertips to Enhance Wet and Dry Grip



Code	Description	Packaging
52000M	PURPLE NITRILE* Exam Gloves, X-Small	Eaches/Box = 100 Boxes/Case = 10
52001M	PURPLE NITRILE* Exam Gloves, Small	Eaches/Box = 100 Boxes/Case = 10
52002M	PURPLE NITRILE* Exam Gloves, Medium	Eaches/Box = 100 Boxes/Case = 10
52003M	PURPLE NITRILE* Exam Gloves, Large	Eaches/Box = 100 Boxes/Case = 10
52004M	PURPLE NITRILE* Exam Gloves, X-Large	Eaches/Box = 90 Boxes/Case = 10

HALYARD* PURPLE NITRILE-XTRA* HIGH RISK EXAM GLOVES

- Longer Cuff for Extra Protection
- Powder-Free
- Latex-Free
- Non-Sterile
- Ambidextrous
- Beaded Cuff
- Single Use Only
- Length: 12" (305mm)
- CE Marked
- Purple
- Textured Fingertips to Enhance Wet and Dry Grip
- Tested For Use With Chemotherapy Drugs. (Data On File)



Code	Description	Packaging
50600M	PURPLE NITRILE-XTRA* Exam Gloves, X-Small	Eaches/Box = 50 Boxes/Case = 10
50601M	PURPLE NITRILE-XTRA* Exam Gloves, Small	Eaches/Box = 50 Boxes/Case = 10
50602M	PURPLE NITRILE-XTRA* Exam Gloves, Medium	Eaches/Box = 50 Boxes/Case = 10
50603M	PURPLE NITRILE-XTRA* Exam Gloves, Large	Eaches/Box = 50 Boxes/Case = 10
50604M	PURPLE NITRILE-XTRA* Exam Gloves, X-Large	Eaches/Box = 50 Boxes/Case = 10

NOTE: Please see package inserts for directions for use for each of these products.

MEDICAL GLOVES

HALYARD* STERLING NITRILE* EXAM GLOVES

- Powder-Free
- Latex-Free
- Non-Sterile
- Ambidextrous
- Beaded Cuff
- Single Use Only
- Length: 9.5" (242mm)
- Silver
- CE Marked
- Fit, Feel and Tactile Sensitivity of Latex
- Textured Fingertips to Enhance Wet and Dry Grip
- Higher Box Count to Save Space



Code	Description	Packaging
13938	STERLING NITRILE* Exam Gloves, X-Small	Eaches/Box = 200 Boxes/Case = 10
13940	STERLING NITRILE* Exam Gloves, Small	Eaches/Box = 200 Boxes/Case = 10
13941	STERLING NITRILE* Exam Gloves, Medium	Eaches/Box = 200 Boxes/Case = 10
13942	STERLING NITRILE* Exam Gloves, Large	Eaches/Box = 200 Boxes/Case = 10
13943	STERLING NITRILE* Exam Gloves, X-Large	Eaches/Box = 170 Boxes/Case = 10

NOTE: Please see package inserts for directions for use for each of these products.

MEDICAL GLOVES

VINYL EXAM GLOVES


HALYARD* SYNTHETIC POWDER-FREE VINYL EXAM GLOVES

- Powder-Free
- Latex-Free
- Non-Sterile
- Smooth Finish
- Ambidextrous
- Beaded Cuff
- Single Use Only
- Polyvinyl Chloride
- Length: 9.5" (242mm)
- Clear
- CE Marked

Code	Description	Packaging
55031M	SYNTHETIC Powder-Free Vinyl Exam Gloves, Small	Eaches/Box = 100 Boxes/Case = 10
55032M	SYNTHETIC Powder-Free Vinyl Exam Gloves, Medium	Eaches/Box = 100 Boxes/Case = 10
55033M	SYNTHETIC Powder-Free Vinyl Exam Gloves, Large	Eaches/Box = 100 Boxes/Case = 10
55034M	SYNTHETIC Powder-Free Vinyl Exam Gloves, X-Large	Eaches/Box = 90 Boxes/Case = 10

GLOVE DISPENSERS

HALYARD* GLOVE DISPENSERS



Code	Description	Packaging
50201	Exam Glove Dispenser, Clear Plastic, Single	1
50202	Exam Glove Dispenser, Clear Plastic, Double	1
50203	Exam Glove Dispenser, Clear Plastic, Triple	1



HALYARD

FORMERLY KIMBERLY-CLARK HEALTH CARE

PATIENT CARE
CATALOGUE 2015

AUSTRALIA &
NEW ZEALAND

PATIENT CARE

SINGLE-USE LINENS

MOISTURE-PROOF SINGLE-USE LINENS

- Coated with Polyethylene Film
- Complies with Standards of Flammability Testing – Consumer Product Safety Commission (16cfr Part 1610)
- Complies with Standards of Flammability Testing – NATIONAL Fire Protection Association (nfpa) Standard 702, Rated as Class I

Code	Description	Packaging
67843	Draw Sheet, Poly-Coated, Blue, 40" x 72", 101.5cm x 183cm, Bulk	Usages/Case = 80

SINGLE-USE LINENS

- Complies with Standards of Flammability Testing – Consumer Product Safety Commission (16cfr Part 1610)
- Complies with Standards of Flammability Testing – NATIONAL Fire Protection Association (nfpa) Standard 702, Rated as Class I

Code	Description	Packaging
67783	Stretcher Sheet, Blue, 40" x 84", 101.5cm x 213cm, Bulk	Usages/Case = 80
67803	Pillowcase, Blue, 20" x 29", 50.5cm x 73.5cm, Bulk	Usages/Case = 200






ICE PACKS

Code	Description	Packaging
33495	STAY-DRY* Ice Pack, Small, 4 Ties, 5" x 11.8", 13cm x 30cm	Eaches/Package = 25 Packages/Case = 2 Usages/Case = 50
33500	STAY-DRY* Ice Pack, Large, 4 Ties, 6.5" x 14", 16.5cm x 36cm	Eaches/Package = 25 Packages/Case = 2 Usages/Case = 50
33150	EyeCare (EENT) Ice Pack, 1 Strap 4.5" x 10", 11.4cm x 25.4cm	Eaches/Box = 15 Boxes/Case = 2 Usages/Case = 30
33101	Bilateral Facial Ice Pack, 5" x 11.8", 13cm x 30cm	Eaches/Box = 12 Boxes/Case = 2 Usages/Case = 24

NOTE: Please see package inserts for directions for use for each of these products.

PATIENT CARE

I.V. SUPPORT


	Code	Description	Packaging
	29900	HAND-AID* I.V. Wrist Support, Adult, 6.5" x 2.75", 16.5cm x 7cm, 9", 23cm Straps	Eaches/Package = 10 Packages/Case = 2 Usages/Case = 20
	29920	Antecubital I.V. Support, 12.5" x 4.5", 32cm x 11.4cm, 19", 48cm Straps	Eaches/Package = 10 Packages/Case = 1 Usages/Case = 10
	29965	I.V. Support, Neonatal, 3.5" x 1.5", 9cm x 4cm, 5", 13cm Straps	Eaches/Package = 10 Packages/Case = 3 Usages/Case = 30
	29967	I.V. Support, Infant/Child, 6.25" x 2", 16cm x 5cm, 9", 23cm Straps	Eaches/Package = 15 Packages/Case = 2 Usages/Case = 30
	29970	HAND-AID* I.V. Wrist Support, Pediatric, 5.5" x 2", 14cm x 5cm, 9" 23cm Straps	Eaches/Package = 10 Packages/Case = 2 Usages/Case = 20
	29990	HAND-AID* Arterial Wrist Support, Paediatric, 6" x 2.5", 16cm x 7cm, 14.5", 23cm Straps	Eaches/Package = 10 Packages/Case = 2 Usages/Case = 20
	29980	HAND-AID* Arterial Wrist Support, Adult, 9" x 3.5", 23cm x 9cm, 14.5", 37cm Straps	Eaches/Package = 10 Packages/Case = 2 Usages/Case = 20
	29985	HAND-AID* Arterial Replacement Straps/Pad/Starter Block	Eaches/Package = 25 Packages/Case = 2 Usages/Case = 50

NOTE: Please see package inserts for directions for use for each of these products.

PATIENT CARE

O.R. PRODUCTS – NURSING SERVICES

STIRRUPS / STRAPS



Code	Description	Packaging
52712	Personal Stirrup Strap with Slip Ring, 19" x 3.5", 48cm x 9cm	Items/Package = 10 Packages/Case = 4 Usages/Case = 40
55299	TS-30* Personal Utility Strap, 29.5" x 3", 76cm x 8cm	Items/Package = 12 Packages/Case = 2 Usages/Case = 24

TOURNIQUET COVERS

Code	Description	Packaging
53600	O.R. Tourniquet Cover, Small 21" x 14.5", 53.5cm x 37cm	Items/Package = 15 Packages/Case = 2 Usages/Case = 30
53605	O.R. Tourniquet Cover, Large 30.5" x 14.5", 78cm x 37cm	Items/Package = 15 Packages/Case = 2 Usages/Case = 30



HALYARD



FORMERLY KIMBERLY-CLARK HEALTH CARE

PROTECTIVE APPAREL
CATALOGUE 2015



AUSTRALIA &
NEW ZEALAND

PPE FOR OR

FLUID PROTECTION FACE MASKS

	Code	Description	Packaging
 48207	48207	FLUIDSHIELD* Fog-Free Surgical Mask, Pleat-Style with Ties, Orange	Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300
	48208	FLUIDSHIELD* Surgical Mask, Pleat-Style with Ties, Orange	Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300
	48227	FLUIDSHIELD* DUCKBILL* Surgical Mask, Pleat-Style with Ties, Orange	Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300
 48247	48247	FLUIDSHIELD* Fog-Free Surgical Mask with Wraparound SPLASHGUARD* Visor, Pleat-Style with Ties, Orange	Eaches/Box = 25 Boxes/Case = 4 Eaches/Case = 100
	48297	FLUIDSHIELD* Care Bear Surgical Mask, Pleat-Style with Ties, Pink And Blue Teddy Bear Print	Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300
	62113	THE PROTECTOR* Fog-Free Splash-Resistant Surgical Mask, Pleat-Style with Ties, Blue with Orange Diamond Print	Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300
	62114	The Protector* Fog-Free Splash-Resistant Surgical Mask with Wraparound Splashguard* Visor, Pleat-Style with Ties, Blue with Orange Diamond Print	Eaches/Box = 25 Boxes/Case = 4 Eaches/Case = 100

FACE MASK ACCESSORIES

	Code	Description	Packaging
 41204	41204	GUARDALL SHIELD* Full-Length Visor Only (No Mask), Clear	Eaches/Set = 1 Sets/Case = 40 Eaches/Case = 40
 41205	41205	GUARDALL SHIELD* Three-Quarter Length Visor Only (No Mask), Clear	Eaches/Set = 1 Sets/Case = 40 Eaches/Case = 40

NOTE: Please see package inserts for directions for use for each of these products.

PPE FOR OR

HEADWEAR



BOUFFANT CAPS

- 1-Layer Spunbond Fabric
- No Potential Fluid Contact


Code	Description	Packaging
3501	Bouffant Cap, Spunbond Fabric, Blue, Regular	Eaches/Dispenser = 100 Dispensers/Case = 5 Eaches/Case = 500
3502	Bouffant Cap, Spunbond Fabric, White, Regular	Eaches/Dispenser = 100 Dispensers/Case = 5 Eaches/Case = 500



SURGICAL CAP

- Blue
- Universal Size

Code	Description	Packaging
69520	Surgical Cap, KAYCEL* Fabric, Blue, Universal	Eaches/Dispenser = 100 Dispensers/Case = 5 Eaches/Case = 500



SURGICAL HOOD

- Protective 3-Layer SMS Fabric
- Tie Neck
- Universal Size


Code	Description	Packaging
69110	Protective Surgical Hood, 3-Layer Fabric, Neck Tie, Blue, Universal	Eaches/Dispenser = 100 Dispensers/Case = 3 Eaches/Case = 300

PPE FOR ICU

FOOTWEAR

KIMLON* SOF-SHOE OVERSHOES

Code	Description	Packaging
6925	KIMLON* Sof-Shoe Overshoes, SMS Fabric, Blue, Regular	Eaches/Dispenser = 100 Dispensers/Case = 2 Eaches/Case = 200



BASIC NON-SKID SHOE COVERS

- Blue
- Protective 3-Layer SMS Fabric

Code	Description	Packaging
3603	Basic Non-Skid Shoe Cover, Spunbond Fabric, Regular	Eaches/Dispenser = 100 Dispensers/Case = 3 Eaches/Case = 300
3604	Basic Non-Skid Shoe Cover, SMS Fabric, X-Large	Eaches/Dispenser = 100 Dispensers/Case = 2 Eaches/Case = 200




X-TRA TRACTION* SHOE COVERS

- Blue
- Protective 3-Layer SMS Fabric
- X-tra Traction

Code	Description	Packaging
69252	X-TRA TRACTION* Shoe Cover, 3-Layer Fabric, Traction Strips, Universal	Eaches/Dispenser = 100 Dispensers/Case = 3 Eaches/Case = 300
69254	X-TRA TRACTION* Shoe Cover, 3-Layer Fabric, Traction Strips, X-Large	Eaches/Dispenser = 80 Dispensers/Case = 3 Eaches/Case = 240

ANKLE TOPS




ANKLE GUARD* SHOE COVER

- Over the Ankle Coverage
- Protective 3-Layer SMS Fabric
- Adhesive Traction Strips
- Moderate Fluid Contamination

Code	Description	Packaging
69353	ANKLE GUARD* Shoe Cover, 3-Layer Fabric, Adhesive Traction Strips, X-Large	Eaches/Dispenser = 90 Dispensers/Case = 3 Eaches/Case = 270

NOTE: Please see package inserts for directions for use for each of these products.

PROTECTIVE APPAREL



HALYARD* ISOLATION GOWNS


- Medium-Weight SMS Fabric
- Tape-Tab Neck Closure
- Tie Waist
- For Use when Expected Risk of Exposure to Fluid is between Low and Moderate

Code	Description	Packaging
69979	Kimberly-Clark Isolation Gown, 3-Layer Fabric, Elastic Cuff, Yellow, Large	Eaches/Bundle = 10 Bundles/Case = 10 Eaches/Case = 100
69988	Kimberly-Clark Isolation Gown, 3-Layer Fabric, Elastic Cuff, Yellow, X-Large	Eaches/Bundle = 10 Bundles/Case = 10 Eaches/Case = 100
69981	Kimberly-Clark Isolation Gown, 3-Layer Fabric, Elastic Cuff, Blue, Large	Eaches/Bundle = 10 Bundles/Case = 10 Eaches/Case = 100
69987	Kimberly-Clark Isolation Gown, 3-Layer Fabric, Elastic Cuff, Blue, X-Large	Eaches/Bundle = 10 Bundles/Case = 10 Usages/Case = 100

NOTE: Please see package inserts for directions for use for each of these products.

PROTECTIVE APPAREL


IMPERVIOUS GOWNS



BASIC IMPERVIOUS OPEN BACK GOWNS

- 3-Layer Spunbond/Laminate
- Potential for Heavy Fluid Contact in Front
- Elastic Cuffs


Code	Description	Packaging
3507	Impervious Open Back Gown, Elastic Cuff, Ties, Yellow, Universal	Eaches/Bundle = 10 Bundles/Case = 5 Eaches/Case = 50
3508	Impervious Open Back Gown, Elastic Cuff, Ties, Blue, Universal	Eaches/Bundle = 10 Bundles/Case = 5 Eaches/Case = 50



BASIC IMPERVIOUS FULL BACK GOWNS

- 3-Layer Spunbond/Laminate
- Potential for Heavy Fluid Contact in Front
- Elastic Cuffs

Code	Description	Packaging
3509	Impervious Full Back Gown, Elastic Cuff, Ties, Blue, Universal	Eaches/Bundle = 10 Bundles/Case = 5 Eaches/Case = 50
3510	Impervious Full Back Gown, Elastic Cuff, Ties, White, Universal	Eaches/Bundle = 10 Bundles/Case = 5 Eaches/Case = 50
3511	Impervious Full Back Gown, Elastic Cuff, Ties, Yellow, Universal	Eaches/Bundle = 10 Bundles/Case = 5 Eaches/Case = 50




IMPERVIOUS OPEN BACK GOWN

- 3-Layer Fabric with Polyethylene Coating
- Potential for Heavy Fluid Contact in Front
- Large
- Knit Cuffs
- Bulk Packaging

Code	Description	Packaging
69455	Impervious Open Back Gown, Knit Cuff, Blue, Large	Eaches/Bundle = 30 Bundles/Case = 2 Eaches/Case = 60

NOTE: Please see package inserts for directions for use for each of these products.

PROTECTIVE APPAREL



CHEMOTHERAPY GOWN

- Tested with the 7 Chemotherapy Drugs Required by ASTM Test Method F739
- Tested with 3 Additional Standard Widely Used Chemotherapy Drugs
- Polypropylene with Polyethylene Coating
- Potential for Soaking and Heavy Fluid Contact in Front

Code	Description	Packaging
69606	Impervious Comfort Gown, Knit Cuffs, Perforated Open Back, Blue, Universal	Eaches/Case = 100



THUMBS-UP* IMPERVIOUS FILM GOWNS

- Impervious Polyethylene Film
- Potential for Soaking and Heavy Fluid Contact in Front
- Thumbhooks

Code	Description	Packaging
7000	Thumbs-up* Impervious Film Gown, Thumbhooks, Open Back, Blue, Regular	Eaches/Dispenser = 15 Dispensers/Case = 5 Eaches/Case = 75
7001	Thumbs-up* Impervious Film Gown, Thumbhooks, Open Back, Blue, X-Large	Eaches/Dispenser = 15 Dispensers/Case = 5 Eaches/Case = 75
7002	Thumbs-up* Impervious Film Gown, Thumbhooks, Open Back, Yellow, Regular	Eaches/Dispenser = 15 Dispensers/Case = 5 Eaches/Case = 75
7003	Thumbs-up* Impervious Film Gown, Thumbhooks, Open Back, Yellow, X-Large	Eaches/Dispenser = 15 Dispensers/Case = 5 Eaches/Case = 75



IMPERVIOUS FULL BACK COMFORT GOWN


- 3-Layer Spunbond/Laminate
- Spunbond Back
- Potential for Heavy Fluid Contact in Front

Code	Description	Packaging
4511	Impervious Full Back Comfort Gown, Knit Cuffs, elcro Neck Secure, Yellow, Universal	Eaches/Case = 50

NOTE: Please see package inserts for directions for use for each of these products.

PROTECTIVE APPAREL


SCRUB SUITS



SCRUB SHIRTS

- Blue
- Protective 3-Layer SMS Fabric
- Potential for Light Fluid Contact
- Scoop Neck
- Reinforced Pocket

Code	Description	Packaging
69701	HALYARD* Scrub Shirt, 3-Layer Fabric, Blue, Medium	Eaches/Case = 48
69702	HALYARD* Scrub Shirt, 3-Layer Fabric, Blue, Large	Eaches/Case = 48
69703	HALYARD* Scrub Shirt, 3-Layer Fabric, Blue, X-Large	Eaches/Case = 48




SCRUB PANTS

- Blue
- Protective 3-Layer SMS Fabric
- Potential for Light Fluid Contact
- Drawstring Waist
- Side and Rear Pockets
- Reinforced In-Seam

Code	Description	Packaging
69701	HALYARD* Scrub Shirt, 3-Layer Fabric, Blue, Medium	Eaches/Case = 48

PROTECTIVE APPAREL

PATIENT APPAREL



PATIENT EXAM GOWN

Code	Description	Packaging
6816	Patient Exam Gown, Short Sleeve, Blue, Universal	Eaches/Package = 10 Packages/Case = 10 Usages/Case = 100


PROTECTAGOWN

Code	Description	Packaging
6812	ProtectaGown, Long Sleeve, SMS Fabric, White, Universal	Eaches/Package = 10 Packages/Case = 5 Usages/Case = 50

NOTE: Please see package inserts for directions for use for each of these products.

PROTECTIVE APPAREL

HEADWEAR



BOUFFANT CAPS

- 1-Layer Spunbond Fabric
- No Potential Fluid Contact

Code	Description	Packaging
3501	Bouffant Cap, Spunbond Fabric, Blue, Regular	Eaches/Dispenser = 100 Dispensers/Case = 5 Eaches/Case = 500
3502	Bouffant Cap, Spunbond Fabric, White, Regular	Eaches/Dispenser = 100 Dispensers/Case = 5 Eaches/Case = 500



SURGICAL CAP

- Blue
- Universal Size

Code	Description	Packaging
69520	Surgical Cap, KAYCEL* Fabric, Blue, Universal	Eaches/Dispenser = 100 Dispensers/Case = 5 Eaches/Case = 500



SURGICAL HOOD

- Protective 3-Layer SMS Fabric
- Tie Neck
- Universal Size

Code	Description	Packaging
69110	Protective Surgical Hood, 3-Layer Fabric, Neck Tie, Blue, Universal	Eaches/Dispenser = 100 Dispensers/Case = 3 Eaches/Case = 300


NOTE: Please see package inserts for directions for use for each of these products.

PROTECTIVE APPAREL

FOOTWEAR

KIMLON* SOF-SHOE OVERSHOES

Code	Description	Packaging
6925	KIMLON* Sof-Shoe Overshoes, SMS Fabric, Blue, Regular	Eaches/Dispenser = 100 Dispensers/Case = 2 Eaches/Case = 200



BASIC NON-SKID SHOE COVERS

- Blue
- Protective 3-Layer SMS Fabric

Code	Description	Packaging
3603	Basic Non-Skid Shoe Cover, Spunbond Fabric, Regular	Eaches/Dispenser = 100 Dispensers/Case = 3 Eaches/Case = 300
3604	Basic Non-Skid Shoe Cover, SMS Fabric, X-Large	Eaches/Dispenser = 100 Dispensers/Case = 2 Eaches/Case = 200




X-TRA TRACTION* SHOE COVERS

- Blue
- Protective 3-Layer SMS Fabric
- X-tra Traction

Code	Description	Packaging
69252	X-TRA TRACTION* Shoe Cover, 3-Layer Fabric, Traction Strips, Universal	Eaches/Dispenser = 100 Dispensers/Case = 3 Eaches/Case = 300
69254	X-TRA TRACTION* Shoe Cover, 3-Layer Fabric, Traction Strips, X-Large	Eaches/Dispenser = 80 Dispensers/Case = 3 Eaches/Case = 240

PROTECTIVE APPAREL

ANKLE TOPS



ANKLE GUARD* SHOE COVER

- Over the Ankle Coverage
- Protective 3-Layer SMS Fabric
- Adhesive Traction Strips
- Moderate Fluid Contamination

Code	Description	Packaging
69353	ANKLE GUARD* Shoe Cover, 3-Layer Fabric, Adhesive Traction Strips, X-Large	Eaches/Dispenser = 90 Dispensers/Case = 3 Eaches/Case = 270

KNEE HIGHS

HI GUARD* ULTRA FULL-COVERAGE BOOT

- Knee High Boot
- Plastic Film Coating over 3-Layer Fabric
- Two 1', 2.5cm Foam Traction Strips
- Potential For Heavy Fluid Contamination

Code	Description	Packaging
69353	ANKLE GUARD* Shoe Cover, 3-Layer Fabric, Adhesive Traction Strips, X-Large	Eaches/Dispenser = 90 Dispensers/Case = 3

NOTE: Please see package inserts for directions for use for each of these products.



HALYARD

FORMERLY KIMBERLY-CLARK HEALTH CARE

RESPIRATORY HEALTH
CATHETER KIT
CATALOGUE 2015

AUSTRALIA &
NEW ZEALAND

RESPIRATORY HEALTH – CATHETER KIT

HALYARD* NON-BRONCHOSCOPIC BRONCHOALVEOLAR LAVAGE CATHETER KIT



Code	Outer Diameter	Inner Diameter	Description	Packaging
141	16 Fr / 5.2mm	12 Fr / 4mm	Non-Bronchoscopic Bronchoalveolar Lavage Catheter	5/Case

NOTE: Please see package inserts for directions for use for each of these products.



HALYARD

FORMERLY KIMBERLY-CLARK HEALTH CARE

RESPIRATORY HEALTH
CLOSED SUCTION SYSTEMS
CATALOGUE 2015

AUSTRALIA &
NEW ZEALAND

RESPIRATORY HEALTH – CLOSED SUCTION SYSTEMS

CLOSED SUCTION SYSTEMS / ACCESSORIES

PAEDIATRIC ENDOTRACHEAL TUBES

HALYARD* MULTI-ACCESS PORT CLOSED SUCTION SYSTEM

- Rotating Manifold
- Single Use Cassette
- Turbulent Cleaning Chamber
- Sleeve Tether
- MDI Drug Delivery
- Replacement Catheters

Code	Description	Packaging
240128-5	HALYARD* Multi-Access Port Closed Suction System for Adults, Endotracheal, MDI, 12 F	20/Case
240148-5	HALYARD* Multi-Access Port Closed Suction System for Adults, Endotracheal, MDI, 14 F	20/Case
2401238-5	HALYARD* Multi-Access Port Closed Suction System for Adults, Tracheostomy, DSY, MDI, 12 F	20/Case
2401438-5	HALYARD* Multi-Access Port Closed Suction System for Adults, Tracheostomy, DSY, MDI, 14 F	20/Case
240125-5	HALYARD* Multi-Access Port Replacement Catheter for Adults, Endotracheal, 12 F	20/Case
240145-5	HALYARD* Multi-Access Port Replacement Catheter for Adults, Endotracheal, 14 F	20/Case
2401253-5	HALYARD* Multi-Access Port Replacement Catheter for Adults, Tracheostomy, 12 F	20/Case
2401453-5	HALYARD* Multi-Access Port Replacement Catheter for Adults, Tracheostomy, 14 F	20/Case
2401286-5	HALYARD* Wet Pak* Multi-Access Port Closed Suction System for Adults, Endotracheal, DSY, MDI, 12 F	20/Case
2401486-5	HALYARD* Wet Pak* Multi-Access Port Closed Suction System for Adults, Endotracheal, DSY, MDI, 14 F	20/Case
24013-5	HALYARD* Seal Cassette	20/Case


NOTE: Please see package inserts for directions for use for each of these products.

RESPIRATORY HEALTH – CLOSED SUCTION SYSTEMS

HALYARD* CLOSED SUCTION SYSTEMS

HALYARD* ADULT TRACH CARE* 24 HOUR DOUBLE SWIVEL ELBOWS

- Gamma Sterilised
- One-Way Lavage / Irrigation Port
- Dual Swivels Provides Flexibility without Connectors
- Tactile Sleeve for Better Feel
- Numbered Catheter Allowing the Clinician to Limit the Distance of Catheter Advance
- Tracheostomy Sizes




Code	Outer Diameter	Length	Description	Packaging
2160-5	12 Fr / 4mm	21.3" / 54cm	Double Swivel Elbow	20/Case
21603-5	12 Fr / 4mm	12" / 30.5cm	Double Swivel Elbow Tracheostomy	20/Case
21609-5	12 Fr / 4mm	22.2" / 56cm	Double Swivel Elbow Directional Tip	20/Case
2210-5	14 Fr / 4.6mm	21.3" / 54cm	Double Swivel Elbow	20/Case
22103-5	14 Fr / 4.6mm	21.3" / 30.5cm	Double Swivel Elbow Tracheostomy	20/Case
22107-5	14 Fr / 4.6mm	21.3" / 54cm	Double Swivel Elbow Dual Lumen	20/Case
221037-5	14 Fr / 4.6mm	21.3" / 30.5cm	Double Swivel Elbow Tracheostomy Dual Lumen	20/Case
221079-5	14 Fr / 4.6mm	22.2" / 56cm	Double Swivel Elbow Dual Lumen Directional Tip	20/Case
22108-5	14 Fr / 4.6mm	21.3" / 54cm	Double Swivel Elbow MDI Built-In	20/Case
221038-5	14 Fr / 4.6mm	12" / 30.5cm	Double Swivel Elbow Tracheostomy MDI Built-In	20/Case
22109-5	14 Fr / 4.6mm	22.2" / 56cm	Double Swivel Elbow Directional Tip	20/Case
2260-5	16 Fr / 5.3mm	21.3" / 54cm	Double Swivel Elbow	20/Case

NOTE: Please see package inserts for directions for use for each of these products.

RESPIRATORY HEALTH – CLOSED SUCTION SYSTEMS

HALYARD* ADULT TRACH CARE* 24 HOUR DOUBLE SWIVEL ELBOW WET PAK*

- Gamma Sterilised
- One-Way Lavage / Irrigation Port
- Dual Swivels Provides Flexibility without Connectors
- Tactile Sleeve for Better Feel
- Numbered Catheter Allowing the Clinician to Limit the Distance of Catheter Advance
- Tracheostomy Sizes




Code	Outer Diameter	Length	Description	Packaging
21606-5	12 Fr / 4mm	21.3" / 54cm	Double Swivel Elbow WET PAK*, includes 12 x 15ml Vials Saline	20/Case
216069-5	12 Fr / 4mm	22.2" / 56cm	Double Swivel Elbow WET PAK* Directional Tip, includes 12 x 15ml Vials Saline	20/Case
216036-5	12 Fr / 4mm	12" / 30.5cm	Double Swivel Elbow WET PAK*, includes 12 x 15ml Vials Saline Tracheostomy	20/Case
22106-5	14 Fr / 4.6mm	21.3" / 54cm	Double Swivel Elbow WET PAK*	20/Case
221069-5	14 Fr / 4.6mm	21.3" / 54cm	Double Swivel Elbow WET PAK* Directional Tip, includes 12 x 15ml Vials Saline	20/Case
221036-5	14 Fr / 4.6mm	12" / 30.5cm	Double Swivel Elbow WET PAK* Tracheostomy	20/Case
22606-5	16 Fr / 5.3mm	21.3" / 54cm	Double Swivel Elbow WET PAK*	20/Case

NOTE: Please see package inserts for directions for use for each of these products.

RESPIRATORY HEALTH – CLOSED SUCTION SYSTEMS

HALYARD* ADULT TRACH CARE* 24 HOUR T-PIECES

- Gamma Sterilised
- Lavage / Irrigation Port
- Locking Thumb Port
- T-Piece
- Multiple French Sizes



Code	Outer Diameter	Length	Description	Packaging
2105-5	10 Fr / 3.3mm	21.3" / 54cm	Clear T-Piece	20/Case
2155-5	12 Fr / 4mm	21.3" / 54cm	Clear T-Piece	20/Case
2205-5	14 Fr / 4.6mm	21.3" / 54cm	Clear T-Piece	20/Case
220135-5	14 Fr / 4.6mm	12" / 30.5cm	Clear T-Piece Tracheostomy Swivel Adapter	20/Case
2255-5	16 Fr / 5.3mm	21.3" / 54cm	Clear T-Piece	20/Case
25-2210-5	14 Fr / 4.6mm	21.3" / 54cm	24hr Double Swivel Elbow with Safety Lavage Port	20/Case

HALYARD* ADULT TRACH CARE* 24 HOUR T-PIECE WET PAK*

- Available in T-Piece, Elbow and Double Swivel Elbow
- Lavage / Irrigation Port
- Catheters Packaged Separate from Vials
- 12 Vials per Catheter
- Vials Designed for Respiratory Use



Code	Outer Diameter	Length	Description	Packaging
22056-5	14 Fr / 4.6mm	21.3" / 54cm	Clear T-Piece WET PAK	20/Case
220568-5	14 Fr / 4.6mm	21.3" / 54cm	Clear T-Piece WET PAK MDI Adapter	20/Case
22012356-5	14 Fr / 4.6mm	12" / 30.5cm	Clear T-Piece WET PAK Tracheostomy Swivel Adapter 15mm x 22mm Adapter	20/Case
220123568-5	14 Fr / 4.6mm	12" / 30.5cm	Clear T-Piece WET PAK Tracheostomy Swivel Adapter 15mm x 22mm Adapter MDI Adapter	20/Case
22556-5	16 Fr / 5.3mm	21.3" / 54cm	Clear T-Piece WET PAK	20/Case


NOTE: Please see package inserts for directions for use for each of these products.

RESPIRATORY HEALTH – CLOSED SUCTION SYSTEMS

TURBO ADULT CLOSED SUCTION SYSTEMS

HALYARD* ADULT TRACH CARE* 72 HOUR DOUBLE SWIVEL ELBOWS

- Gamma Sterilised
- One-Way Lavage / Irrigation Port
- Dual Swivels Provides Flexibility Without Connectors
- Tactile Sleeve for Better Feel
- Numbered Catheter Allowing the Clinician to Limit the Distance of Catheter Advance
- Tracheostomy Sizes



Code	Outer Diameter	Length	Description	Packaging
227-5	14 Fr / 4.6mm	21.3" / 54cm	Double Swivel Elbow Endotracheal	20/Case
22703-5	14 Fr / 4.6mm	12" / 30.5cm	Double Swivel Elbow Tracheostomy	20/Case
22716-5	12 Fr / 4mm	21.3" / 54cm	Double Swivel Elbow Endotracheal	20/Case
2271603-5	12 Fr / 4mm	12" / 30.5cm	Double Swivel Elbow Tracheostomy	20/Case


HALYARD* ADULT TRACH CARE* 72 HOUR DOUBLE SWIVEL ELBOWS WET PAK*

Code	Outer Diameter	Length	Description	Packaging
227166-5	12 Fr / 4mm	21.3" / 54cm	Double Swivel Elbow WET PAK*, includes 12 x 15ml Vials Saline	20/Case
22716036-5	12 Fr / 4mm	12" / 30.5cm	Double Swivel Elbow WET PAK*, includes 12 x 15ml Vials Saline Tracheostomy	20/Case
2276-5	14 Fr / 4.6mm	21.3" / 54cm	Double Swivel Elbow WET PAK*, includes 12 x 15ml Vials Saline	20/Case
227036-5	14 Fr / 4.6mm	12" / 30.5cm	Double Swivel Elbow WET PAK*, includes 12 x 15ml Vials Saline Tracheostomy	20/Case

NOTE: Please see package inserts for directions for use for each of these products.

RESPIRATORY HEALTH – CLOSED SUCTION SYSTEMS


HALYARD* ADULT TRACH CARE* 72 HOUR DOUBLE SWIVEL ELBOWS



Code	Outer Diameter	Length	Description	Packaging
22705-5	14 Fr / 4.6mm	21.3" / 54cm	T-Piece Endotracheal	20/Case
2270135-5	14 Fr / 4.6mm	12" / 30.5cm	T-Piece Tracheostomy	20/Case
22715-5	12 Fr / 4mm	21.3" / 54cm	T-Piece Endotracheal	20/Case
227135-5	12 Fr / 4mm	12" / 30.5cm	T-Piece Tracheostomy	20/Case

HALYARD* ADULT TRACH CARE* 72 HOUR T-PIECE WET PAK*

- Available in T-Piece, Elbow and Double Swivel Elbow
- Lavage / Irrigation Port
- Catheters Packaged Separate from Vials
- 12 Vials Per Catheter
- Vials Designed for Respiratory Use



Code	Outer Diameter	Length	Description	Packaging
227156-5	12 Fr / 4.0mm	21.3" / 54cm	Adult WET PAK* T-Piece Standard	20/Case
2271356-5	12 Fr / 4.0mm	12" / 30.5cm	Adult WET PAK* T-Piece Tracheostomy	20/Case
227056-5	14 Fr / 4.6mm	21.3" / 54cm	Adult WET PAK* T-Piece Standard	20/Case
22701356-5	14 Fr / 4.6mm	12" / 30.5cm	Adult WET PAK* T-Piece Tracheostomy	20/Case

NOTE: Please see package inserts for directions for use for each of these products.


RESPIRATORY HEALTH – CLOSED SUCTION SYSTEMS

SPECIALTY CLOSED SUCTION SYSTEMS

HALYARD* ADULT TRACH CARE* 24 HOUR DUAL LUMENS

- Gamma Sterilised
- Dual Lavage / Irrigation Ports
- Available in T-Piece and Double Swivel Elbow

- Available in Directional Tip
- Allows Lavage to Take Place at the Patient End of the Catheter



Code	Outer Diameter	Length	Description	Packaging
22107-5	14 Fr / 4.6mm	21.3" / 54cm	Double Swivel Elbow Dual Lumen	20/Case
221037-5	14 Fr / 4.6mm	12" / 30.5cm	Double Swivel Elbow Dual Lumen Tracheostomy	20/Case
221079-5	14 Fr / 4.6mm	22.2" / 56cm	Double Swivel Elbow Dual Lumen Directional Tip	20/Case

HALYARD* ADULT TRACH CARE* 24 HOUR DIRECTIONAL TIPS

- Gamma Sterilised
- Lavage / Irrigation Port
- Available in T-Piece, Elbow and Double Swivel Elbow
- Allows for Suctioning of the Right or Left Main Stem Bronchus

- Radiopaque Line Provides for Observation of Catheter Path Under Fluoroscopy
- Colored Center Line Allows Maintaining Manipulation of the Catheter Tip to Access Proper Location of the Right or Left Main Stem Bronchus


Code	Outer Diameter	Length	Description	Packaging
21559-5	12 Fr / 4mm	22.2" / 56cm	Clear T-Piece, Directional Tip	20/Case
22059	14 Fr / 4.6mm	22.2" / 56cm	Clear T-Piece, Directional Tip	20/Case
220569-5	14 Fr / 4.6mm	22.2" / 56cm	Clear T-Piece, WET PAK*, Directional Tip	20/Case
22559-5	16 Fr / 5.3mm	22.2" / 56cm	Clear T-Piece Directional Tip	20/Case
22109-5	14 Fr / 4.6mm	22.2" / 56cm	Double Swivel Elbow, Directional Tip	20/Case
221079-5	14 Fr / 4.6mm	22.2" / 56cm	Double Swivel Elbow, Directional Tip Dual Lumen	20/Case
225579-5	16 Fr / 5.3mm	22.2" / 56cm	Closed T-Piece, Directional Tip	20/Case

NOTE: Please see package inserts for directions for use for each of these products.

RESPIRATORY HEALTH – CLOSED SUCTION SYSTEMS

PAEDIATRIC / NEONATAL CLOSED SUCTION SYSTEMS

HALYARD* NEONATAL AND PAEDIATRIC TRACH CARE* ELBOWS



Code	Outer Diameter	Length	Description	Packaging
202-5	5 Fr / 1.66mm	12" / 30.5cm	Elbow, Endotracheal, Tube Adapters (2mm, 2.5mm)	20/Case
206-5	6 Fr / 2mm	12" / 30.5cm	Elbow, Endotracheal, Tube Adapters (2.5mm, 3mm, 3.5mm)	20/Case
207-5	7 Fr / 2.3mm	12" / 30.5cm	Elbow, Endotracheal, Tube Adapters (2.5mm, 3mm, 3.5mm, 4mm)	20/Case
208-5	8 Fr / 2.6mm	14" / 35.5cm	Elbow, Endotracheal, Tube Adapters (3mm, 3.5mm, 4mm)	20/Case
210-5	10 Fr / 3.3mm	16" / 40.5cm	Elbow	20/Case

HALYARD* NEONATAL AND PAEDIATRIC TRACH CARE* "Y"




Code	Outer Diameter	Length	Description	Packaging
195-5	5 Fr / 1.66mm	12" / 30.5cm	"Y" Endotracheal, Tube Adapters (2.0mm, 2.5mm)	20/Case
196-5	6 Fr / 2mm	12" / 30.5cm	"Y" Endotracheal, Tube Adapters (2.5mm, 3.0mm, 3.5mm)	20/Case
197-5	7 Fr / 2.3mm	12" / 30.5cm	"Y" Endotracheal, Tube Adapters (3.0mm, 3.5mm, 4mm)	20/Case
198-5	8 Fr / 2.6mm	12" / 30.5cm	"Y" Endotracheal, Tube Adapters (3.5mm, 4.0mm, 4.5mm)	20/Case
1910-5	10 Fr / 3.3mm	16" / 40.5cm	"Y" Endotracheal, Tube Adapters (4.5mm, 5.0mm, 5.5mm)	20/Case
1912-5	12 Fr / 4mm	16" / 40.5cm	"Y" Endotracheal, Tube Adapters (5.5mm, 6.0mm, 6.5mm)	20/Case

NOTE: Please see package inserts for directions for use for each of these products.

RESPIRATORY HEALTH – CLOSED SUCTION SYSTEMS

HALYARD* NEONATAL AND PAEDIATRIC MAC CATHETERS



Code	Outer Diameter	"Y" Size	Description	Packaging
1900	5 Fr / 1.66mm	2.5mm / 7.5mm, 3.0mm / 7.5mm, 3.5mm / 7.5mm	MAC Catheter	5/Case
1940	5 Fr / 1.66mm	4.0mm / 7.5mm	MAC Catheter	5/Case

HALYARD* TRACH CARE* PAEDIATRIC / NEONATAL COMPONENT KITS

- Gamma Sterilised
- Neonatal Sizes
- Available in T-Pieces, Elbows and Double Swivel Elbows
- Multiple Configurations
- Customer-Oriented Presentation Designed to Meet Unique Needs

Code	Outer Diameter	Length	Description	Packaging
1968127	6 Fr / 2mm	12" / 30.5cm	Neonatal "Y" WET PAK* (6 Vials) READY CARE* (5 Ea, 8 Fr Oral Catheters) 9', 274cm Suction Tubing	20/Case
1988128	8 Fr / 2.6mm	12" / 30.5cm	Neonatal "Y" WET PAK* (6 Vials) READY CARE* (5 Ea, 10 Fr Oral Catheters) 9', 274cm Suction Tubing	20/Case

NOTE: Please see package inserts for directions for use for each of these products.



HALYARD

FORMERLY KIMBERLY-CLARK HEALTH CARE

RESPIRATORY HEALTH
ENDOTRACHEAL TUBES
CATALOGUE 2015


AUSTRALIA &
NEW ZEALAND

RESPIRATORY HEALTH – ENDOTRACHEAL TUBES

ADULT & PAEDIATRIC ENDOTRACHEAL TUBES

HALYARD* MICROCUFF* ENDOTRACHEAL TUBE, ORAL/NASAL MAGILL

- Ultra-low Cuff Pressure can Reduce the Risk of Tracheal Trauma
- Short, Cylindrical Cuff Placed Near the Tracheal Tube Tip
- Four Precision Bands to Facilitate and Confirm Optimal Tube Placement
- No Murphy Eye
- Oral/Nasal Magill and Pre-Curved Shapes




Code	Inner Diameter	Description	Packaging
35111	3.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral/Nasal Magill	10/Case
35112	3.5mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral/Nasal Magill	10/Case
35113	4.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral/Nasal Magill	10/Case
35114	4.5mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral/Nasal Magill	10/Case
35115	5.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral/Nasal Magill	10/Case
35116	5.5mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral/Nasal Magill	10/Case
35117	6.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral/Nasal Magill	10/Case
35118	6.5mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral/Nasal Magill	10/Case
35119	7.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral/Nasal Magill	10/Case

NOTE: Please see package inserts for directions for use for each of these products.

RESPIRATORY HEALTH – ENDOTRACHEAL TUBES

HALYARD* MICROCUFF* ENDOTRACHEAL TUBE, ORAL CURVED

- Advanced Micro-Thin Polyurethane Cuff Material Designed to Reduce Micro-Aspiration
- Provides a Superior Tracheal Seal and Reduces Leakage Past the Cuff
- Seals at Low Pressure
- Oral/Nasal Magill Curve With Murphy Eye



Code	Inner Diameter	Description	Packaging
35161	3.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral Curved	10/Case
35162	3.5mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral Curved	10/Case
35163	4.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral Curved	10/Case
35164	4.5mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral Curved	10/Case
35165	5.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral Curved	10/Case
35166	5.5mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral Curved	10/Case
35167	6.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral Curved	10/Case
35168	6.5mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral Curved	10/Case
35169	7.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral Curved	10/Case


NOTE: Please see package inserts for directions for use for each of these products.

RESPIRATORY HEALTH – ENDOTRACHEAL TUBES

ADULT ENDOTRACHEAL TUBES

HALYARD* MICROCUFF* ORAL/NASAL MAGILL ENDOTRACHEAL TUBE

- Advanced Micro-Thin Polyurethane Cuff Material Designed to Reduce Micro-Aspiration
- Provides a Superior Tracheal Seal and Reduces Leakage Past the Cuff
- Seals at Low Pressure
- Oral/Nasal Magill Curve With Murphy Eye



Code	Inner Diameter	Description	Packaging
35210	5.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Oral/Nasal Magill, Murphy Eye	10/Case
35211	5.5mm	HALYARD* MICROCUFF* Endotracheal Tube, Oral/Nasal Magill, Murphy Eye	10/Case
35212	6.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Oral/Nasal Magill, Murphy Eye	10/Case
35213	6.5mm	HALYARD* MICROCUFF* Endotracheal Tube, Oral/Nasal Magill, Murphy Eye	10/Case
35214	7.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Oral/Nasal Magill, Murphy Eye	10/Case
35215	7.5mm	HALYARD* MICROCUFF* Endotracheal Tube, Oral/Nasal Magill, Murphy Eye	10/Case
35216	8.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Oral/Nasal Magill, Murphy Eye	10/Case
35217	8.5mm	HALYARD* MICROCUFF* Endotracheal Tube, Oral/Nasal Magill, Murphy Eye	10/Case
35218	9.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Oral/Nasal Magill, Murphy Eye	10/Case
35220	10.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Oral/Nasal Magill, Murphy Eye	10/Case

NOTE: Please see package inserts for directions for use for each of these products.



HALYARD

FORMERLY KIMBERLY-CLARK HEALTH CARE


RESPIRATORY HEALTH
ORAL CARE
CATALOGUE 2015

AUSTRALIA &
NEW ZEALAND

RESPIRATORY HEALTH – ORAL CARE

HALYARD* TRACH CARE* ACCESSORIES

- Complement the TRACH CARE* Closed System Line of Products
- Can be Sold Individually or with TRACH CARE*
- Designed for Specific Patient Needs




Code	Description	Packaging
11000	Manual Percussor for Neonatal Pulmonary Use	12/Case
111	Swivel Adapter 15mm Inner Diameter x 15mm Outer Diameter	50/Case
112	Fixed Adapter 15mm Outer Diameter x 22mm Inner Diameter	50/Case
114	Parallel "Y"	50/Case
118	Unit Dose Saline .9%, Pink	144/Case
119	Metered Dose Adapter	20/Case
19000	Transport Kit "Y" Plug, SUC CVC	144/Case

NOTE: Please see package inserts for directions for use for each of these products.


RESPIRATORY HEALTH – ORAL CARE

HALYARD* READY CARE* CLOSED SUCTION PATIENT COMPONENTS



Code	Description	Packaging
	Components	
12601	1- Oral Suction Catheter	25/Case
12602	1- Suction Toothbrush	25/Case
12603	1- Angled-Tip Suction Swab	25/Case
12604	1- Oral Suction Tip	25/Case

HALYARD* READY CARE* DENTASWAB* COMPONENTS



Code	Description	Packaging
	DENTASWAB* Single Wrapped	
12255	250- Single Wrapped POLYPLUS DENTASWAB* with Dentifrice	4/Case
12254	250- Single Wrapped POLYPLUS DENTASWAB*	4/Case
12245	250- Single Wrapped DENTASWAB* with Dentifrice	4/Case
12246	250- Single Wrapped DENTASWAB* with Dentifrice	2/Case
12240	250- Single Wrapped DENTASWAB*	4/Case
	DENTASWAB* Bulk Packaged	
12252	500- POLYPLUS DENTASWAB*	2/Case
12249	20 - DENTASWAB* with Dentifrice	50/Case

NOTE: Please see package inserts for directions for use for each of these products.



HALYARD

FORMERLY KIMBERLY-CLARK HEALTH CARE

STERILE CONTAINMENT
CATALOGUE 2015

AUSTRALIA &
NEW ZEALAND

STERILE CONTAINMENT

HALYARD* ONE-STEP* STERILIZATION WRAP – SIMULTANEOUS

HALYARD* H100 STERILIZATION WRAP

Code	Size	Unite per Case
12712	30cm x 30cm	480
12715	38cm x 38cm	480
12718	46cm x 46cm	480
12720	51cm x 51cm	480
12724	61cm x 61cm	240
12730	76cm x 76cm	144
12736	91cm x 91cm	144
12740	102cm x 102cm	120
12745	114cm x 114cm	120
12748	122cm x 122cm	120
12754	137cm x 137cm	96
12772	137cm x 183cm	48

HALYARD* H200 STERILIZATION WRAP

Code	Size	Unite per Case
62012	30cm x 30cm	480
62015	38cm x 38cm	480
62018	46cm x 46cm	480
62024	61cm x 61cm	240
62030	76cm x 76cm	144
62036	91cm x 91cm	144
62040	102cm x 102cm	120
62048	122cm x 122cm	48
62054	137cm x 137cm	48
62072	137cm x 183cm	48

HALYARD* H300 STERILIZATION WRAP

Code	Size	Unite per Case
12824	61cm x 61cm	240
12830	76cm x 76cm	144
12836	91cm x 91cm	144
12840	102cm x 102cm	96
12845	114cm x 114cm	48
12848	122cm x 122cm	96
12854	137cm x 137cm	48
12872	183cm x 183cm	48

NOTE: Please see package inserts for directions for use for each of these products.

STERILE CONTAINMENT

HALYARD* H400 STERILIZATION WRAP

Code	Size	Unite per Case
62224	61cm x 61cm	120
62236	91cm x 91cm	72
62240	102cm x 102cm	72
62245	114cm x 114cm	48
62248	122cm x 122cm	24
62254	137cm x 137cm	24
62272	137cm x 183cm	24

HALYARD* H500 STERILIZATION WRAP

Code	Size	Unite per Case
62118	46cm x 46cm	288
62124	61cm x 61cm	120
62130	76cm x 76cm	120
62136	91cm x 91cm	72
62145	114cm x 114cm	48
62148	122cm x 122cm	24
62154	137cm x 137cm	24
62172	137cm x 183cm	24
62190	137cm x 229cm	24
62160	152cm x 152cm	24

HALYARD* H600 STERILIZATION WRAP

Code	Size	Unite per Case
62636	91cm x 91cm	72
62640	102cm x 102cm	72
62645	114cm x 114cm	48
62648	122cm x 122cm	24
62654	137cm x 137cm	24
62672	137cm x 183cm	24

STERILE CONTAINMENT

HALYARD* STERILIZATION WRAP – SEQUENTIAL

HALYARD* H100 STERILIZATION WRAP

Code	Size	Unite per Case
10712	30cm x 30cm	1,000
10715	38cm x 38cm	1,000
10718	46cm x 46cm	1,000
10720	51cm x 51cm	1,000
10724	61cm x 61cm	500
10730	76cm x 76cm	300
10736	91cm x 91cm	300
10740	102cm x 102cm	250
10745	114cm x 114cm	250
10748	122cm x 122cm	250
10754	137cm x 137cm	100
10772	137cm x 183cm	100

HALYARD* H200 STERILIZATION WRAP

Code	Size	Unite per Case
68012	30cm x 30cm	1,000
68018	46cm x 46cm	1,000
68024	61cm x 61cm	500
68030	76cm x 76cm	300
68036	91cm x 91cm	300
68040	102cm x 102cm	250
68048	122cm x 122cm	100
68054	137cm x 137cm	100
68072	137cm x 183cm	100

HALYARD* H300 STERILIZATION WRAP

Code	Size	Unite per Case
10824	61cm x 61cm	500
10830	76cm x 76cm	300
10836	91cm x 91cm	300
10840	102cm x 102cm	250
10845	114cm x 114cm	100
10848	122cm x 122cm	100
10854	137cm x 137cm	100
10872	137cm x 183cm	100

NOTE: Please see package inserts for directions for use for each of these products.

STERILE CONTAINMENT

HALYARD* H400 STERILIZATION WRAP

Code	Size	Unite per Case
68218	46cm x 46cm	600
68224	61cm x 61cm	250
68230	76cm x 76cm	200
68236	91cm x 91cm	150
68240	102cm x 102cm	150
68245	114cm x 114cm	100
68248	122cm x 122cm	50
68254	137cm x 137cm	50
68272	137cm x 183cm	50

HALYARD* H500 STERILIZATION WRAP

Code	Size	Unite per Case
68118	46cm x 46cm	600
68124	61cm x 61cm	250
68130	76cm x 76cm	250
68136	91cm x 91cm	150
68145	114cm x 114cm	100
68148	122cm x 122cm	50
68154	137cm x 137cm	50
68172	137cm x 183cm	50
68190	137cm x 229cm	50
68160	152cm x 152cm	50

HALYARD* H600 STERILIZATION WRAP

Code	Size	Unite per Case
68636	91cm x 91cm	72
68640	102cm x 102cm	72
68645	114cm x 114cm	48
68654	122cm x 122cm	24
68672	137cm x 183cm	24
68636	137cm x 183cm	24

STERILE CONTAINMENT

OTHER STERILISATION PRODUCTS

HALYARD* CONTAINER FILTERS

Code	Description	Packaging
68507	HALYARD* Container Filter, Aesculap, 19cm Round	Items/Case = 1,000
68509	HALYARD* Container Filter, Genesis, 23cm Square	Items/Case = 1,000

NOTE: Please see package inserts for directions for use for each of these products.



HALYARD

FORMERLY KIMBERLY-CLARK HEALTH CARE

SURGICAL GOWNS
CATALOGUE 2015

AUSTRALIA &
NEW ZEALAND


SURGICAL GOWNS

MICROCOOL* BREATHABLE HIGH PERFORMANCE SURGICAL GOWNS WITH SECURE-FIT*

- Innovative patented glove-gown interface
- Patented Breathable Fabric
- Low Linting Fabric
- Raglan Sleeves
- Colour-Coded Neck Binding – Red
- Adjustable Hook & Loop Neckline
- Book Fold
- AAMI Level 4 / EN 13795 High Performance



Code	Description	Packaging
13359	MICROCOOL* Breathable High Performance Surgical Gown with SECURE-FIT*, Sterile, Small	Items/Case = 30
92320	MICROCOOL* Breathable High Performance Surgical Gown with SECURE-FIT*, Sterile, Large	Items/Case = 30
72454	MICROCOOL* Breathable High Performance Surgical Gown with SECURE-FIT*, Non-Sterile, Large	Items/Case = 40
92321	MICROCOOL* Breathable High Performance Surgical Gown with SECURE-FIT*, Sterile, X-Large	Items/Case = 28
30734	MICROCOOL* Breathable High Performance Surgical Gown with SECURE-FIT*, Non-Sterile, X-Large	Items/Case = 384
92323	MICROCOOL* Breathable High Performance Surgical Gown with SECURE-FIT*, Sterile, XX-Large	Items/Case = 26
72448	MICROCOOL* Breathable High Performance Surgical Gown with SECURE-FIT*, Non-Sterile, XX-Large	Items/Case = 32



MICROCOOL* BREATHABLE HIGH PERFORMANCE SPECIALTY SURGICAL GOWNS WITH SECURE-FIT*

- Innovative patented glove-gown interface
- Patented Breathable Fabric
- Low Linting Fabric
- Patented A-Line Design for Seated Procedures
- Raglan Sleeves
- Colour-Coded Neck Binding – Red
- Adjustable Hook & Loop Neckline
- Book Fold
- AAMI Level 4 / EN 13795 High Performance



Code	Description	Packaging
92330	MICROCOOL* Breathable High Performance Specialty Surgical Gown with SECURE-FIT*, Sterile, Large	Items/Case = 28
92331	MICROCOOL* Breathable High Performance Specialty Surgical Gown with SECURE-FIT*, Sterile, X-Large	Items/Case = 26
92333	MICROCOOL* Breathable High Performance Specialty Surgical Gown with SECURE-FIT*, Sterile, XX-Large	Items/Case = 26

-  Non-Reinforced Fabric
-  Fabric-Reinforced
-  High Performance

All codes in **bold** indicate a sterile product.

NOTE: Please see package inserts for directions for use for each of these products.

SURGICAL GOWNS



AERO BLUE* PERFORMANCE SURGICAL GOWNS

- Patented Breathable Fabric
- Low Linting Fabric
- Raglan Sleeves
- Colour-Coded Neck Binding – Red
- Adjustable Hook & Loop Neckline
- Book Fold
- AAMI Level 3 / EN 13795 High Performance

Code	Description	Packaging
41870	AERO BLUE* Performance Surgical Gown, Sterile, Small	Items/Case = 30
43941	AERO BLUE* Performance Surgical Gown, Non-Sterile, Small	Items/Case = 36
41871	AERO BLUE* Performance Surgical Gown, Sterile, Large	Items/Case = 30
41739	AERO BLUE* Performance Surgical Gown, Non-Sterile, Large	Items/Case = 40
42620	AERO BLUE* Performance Surgical Gown, Sterile, Large, X-Long	Items/Case = 28
41725	AERO BLUE* Performance Surgical Gown, Non-Sterile, Large, X-Long	Items/Case = 36
41872	AERO BLUE* Performance Surgical Gown, Sterile, X-Large	Items/Case = 28
41740	AERO BLUE* Performance Surgical Gown, Non-Sterile, X-Large	Items/Case = 32
42621	AERO BLUE* Performance Surgical Gown, Sterile, X-Large, X-Long	Items/Case = 30
41727	AERO BLUE* Performance Surgical Gown, Non-Sterile, X-Large, X-Long	Items/Case = 32
41873	AERO BLUE* Performance Surgical Gown, Sterile, XX-Large	Items/Case = 26
41741	AERO BLUE* Performance Surgical Gown, Non-Sterile, XX-Large	Items/Case = 32
43940	AERO BLUE* Performance Surgical Gown, Sterile, XX-Large, X-Long	Items/Case = 28
43623	AERO BLUE* Performance Surgical Gown, Non-Sterile, XX-Large, X-Long	Items/Case = 30
41874	AERO BLUE* Performance Surgical Gown, Sterile, XXX-Large	Items/Case = 24



Non-Reinforced Fabric



Fabric-Reinforced




High Performance

All codes in **bold** indicate a sterile product.

NOTE: Please see package inserts for directions for use for each of these products.

SURGICAL GOWNS



SURGICAL GOWNS

- EVOLUTION* 4 Fabric
- Low Linting Fabric
- Set-In-Sleeves
- Colour-Coded Neck Binding – Yellow
- Adjustable Hook & Loop Neckline
- Book Fold
- EN 13795 Standard Performance

Code	Description	Packaging
90003	Surgical Gown, Sterile, Small To Medium	Items/Case = 36
90018	Surgical Gown, Sterile, Large	Items/Case = 36
72978	Surgical Gown, Non-Sterile, Large,	Items/Case = 640
90019	Surgical Gown, Sterile, Large, X-Long	Items/Case = 34
90048	Surgical Gown, Sterile, X-Large	Items/Case = 32
72988	Surgical Gown, Non-Sterile, X-Large,	Items/Case = 640
90049	Surgical Gown, Sterile, X-Large, X-Long	Items/Case = 34
70392	Surgical Gown, Non-Sterile, X- Large, X-Long	Items/Case = 60
90050	Surgical Gown, Sterile, XX-Large	Items/Case = 28
90138	Surgical Gowns Double Pack	Items/Case = 18
	Contents:	
72878	1-Surgical Gown, Large	
72988	1-Surgical Gown, X-Large	



Non-Reinforced Fabric



Fabric-Reinforced



High Performance

All codes in **bold** indicate a sterile product.

NOTE: Please see package inserts for directions for use for each of these products.



HALYARD

FORMERLY KIMBERLY-CLARK HEALTH CARE


SURGICAL PACKS
& DRAPES
CATALOGUE 2015

AUSTRALIA &
NEW ZEALAND

CONTENTS

Cardiac Catheterisation / Radiology Drapes	page 6
Cardiovascular Drapes & Packs	page 8
ENT Drapes & Packs	page 10
Equipment Covers	page 12
General Drapes & Packs:	
Accessories	page 13
Basic Packs	page 14
Basic Drapes	page 15
Breast / Chest Drapes	page 17
Fenestrated Drapes	page 18
Laparoscopy / Laparotomy Drapes	page 19
Minor Procedure Drapes	page 20
Thyroid Drapes	page 21
Universal Packs	page 22
Universal Drapes	page 25
Neuro / Spine Drapes	page 27
Obstetrics and Gynaecology Drapes & Packs:	
C-Section Drapes	page 28
Laparoscopy Packs	page 29
Laparoscopy Drapes	page 30
Lithotomy Packs	page 31
Lithotomy Drapes	page 32
Under Buttocks Drapes	page 33
Ophthalmic Drapes & Packs	page 34
Orthopaedic Drapes & Packs:	
Ortho Packs	page 35
Extremity Drapes	page 37
Hand Drapes	page 39
Hip Drapes	page 40
Knee Arthroscopy Drapes	page 41
Shoulder Drapes	page 42
U & Bar Drapes	page 44
Urology Drapes & Packs	page 47
Microbial Sealant	page 48


CARDIAC CATHETERISATION / RADIOLOGY DRAPE



PACEMAKER DRAPE, 221cm x 343cm

- Two Clear Side Windows 55cm x 342cm
- CONTROL* Plus Fabric Reinforcement
- Conformable Breather Bar 35cm
- Two Adhesive Fenestrations 15cm x 20cm
- Bilateral Pockets 15cm x 81cm


Code	Description	Packaging
89368	Pacemaker Drape, Sterile	Items/Case = 11
79368	Pacemaker Drape, Non-Sterile	Items/Case = 12



BRACHIAL ANGIOGRAPHY DRAPE, 111cm x 78cm

- CONTROL* Plus Fabric Reinforcement
- 7cm Adhesive Circular Fenestration

Code	Description	Packaging
89471	Brachial Angiography Drape, Sterile	Items/Case = 50



RADIAL FEMORAL ANGIOGRAPHY DRAPE, 240cm x 380cm

- Two 63cm x 350cm Windows
- Two 12cm x 7cm Oval Radial Fenestrations
- Two 12cm Round Femoral Fenestrations

Code	Description	Packaging
89583	HALYARD BASICS* Radial Femoral Angiography Drape, Sterile	Items/Case = 12

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


CARDIAC CATHETERISATION / RADIOLOGY DRAPE



EXTRA LONG FEMORAL ANGIOGRAPHY DRAPE WITH WINDOW, 208cm x 343cm

- Two 9cm Adhesive Circular Fenestrations
- One 68cm x 315cm Windows
- CONTROL* Plus Fabric Reinforcement


Code	Description	Packaging
89709	Extra Long Femoral Angiography Drape with Window, Sterile	Items/Case = 14
79709	Extra Long Femoral Angiography Drape with Window, Non-Sterile	Items/Case = 16



EXTRA LONG FEMORAL ANGIOGRAPHY DRAPE WITH WINDOWS, 221cm x 343cm

- Two 9cm Adhesive Circular Fenestrations
- Two 68cm x 343cm Windows
- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
89711	Extra Long Femoral Angiography Drape with Dual Windows, Sterile	Items/Case = 16
13296	Extra Long Femoral Angiography Drape with Dual Windows, Non-Sterile	Items/Case = 180



IMPERVIOUS EXTRA-LONG FEMORAL ANGIOGRAPHY DRAPE, 220cm x 342cm

- Extra-Long 220cm x 342cm Film Based Drape
- 81cm x 342cm Absorbent Fabric Reinforcement
- 78cm x 342cm Patient Comfort Layer
- Two 8cm Adhesive Circular Fenestrations
- Two 68cm x 342cm Clear Plastic Windows

Code	Description	Packaging
79749	Impervious Extra-Long Femoral Angiography Drape, Non-Sterile	Items/Case = 144

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

CARDIOVASCULAR DRAPES & PACKS



CARDIOVASCULAR PACK II

Code	Description	Packaging
77748	Cardiovascular Pack II	Items/Case = 3
Sterile Component Code	Non-Sterile Component Code	Contents
89603	79603	1-Back Table Cover, Reinforced, 191cm x 201cm 1-Mayo Stand Cover, 80cm x 145cm 2-Foot Covers 3-Absorbent Towels, 138cm x 56cm 2-Tape Strips, 10cm x 55cm 1-Suture Bag 1-Utility Drape with Tape, 76cm x 76cm 1-Perineal Drape, 187cm x 241cm: 1-CV Incise Sheet, 264cm x 193cm x 363cm
89053	76356	• Armboard Covers • Pouches • CONTROL* Plus Fabric Reinforcement • Tube Holders • 40cm x 50cm Adhesive Fenestration



PERINEAL DRAPE, 187cm x 241cm

- Split 44cm x 95cm
- CONTROL* Plus Fabric Reinforcement
- Positioning Tape

Code	Description	Packaging
76356	Perineal Drape, Non-Sterile	Items/Case = 17

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


CARDIOVASCULAR DRAPES & PACKS



CVARTS* HEAVY-DUTY UNIVERSAL DRAPE SYSTEM

- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement


Code	Description	Packaging
89454	CVARTS* Heavy-Duty Universal Drape System, Sterile	Items/Case = 5
Contents:		
1-CVARTS* Top Drape, 254cm x 152cm		
1-CVARTS* Bottom Drape, 193cm x 152cm		
2-CVARTS* Side Drapes, 96cm x 213cm, Longitudinal Flaps with Hook and Loop		



CVARTS* HEAVY-DUTY CV SPLIT DRAPE, 274cm x 203cm x 315cm

- Armboard Covers
- 15cm x 213cm Split
- Longitudinal Flaps with Hook and Loop
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
89452	CVARTS* Heavy-Duty CV Split Drape, Sterile	Items/Case = 6
79452	CVARTS* Heavy-Duty CV Split Drape, Non-Sterile	Items/Case = 7



CVARTS* BILATERAL UNDERLEG / PERINEAL DRAPE, 213cm x 203cm

- Impervious Perineal Cover
- Umbilicus Cover
- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
89453	CVARTS* Bilateral Underleg / Perineal Drape, Sterile	Items/Case = 14
79453	CVARTS* Bilateral Underleg / Perineal Drape, Non-Sterile	Items/Case = 18

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

ENT DRAPES & PACKS




ENT PACK III

Code	Description	Packaging	
88223	ENT Pack III	Items/Case = 3	
Sterile Component Code	Non-Sterile Component Code	Contents	
89603	79603	1-Back Table Cover, Zone-Reinforced, 152cm x 191cm	
89431	79430	1-Mayo Stand Cover 80cm x 145cm	
89352	79351	1-Suture Bag	
		1-Head Drape 91cm x 127cm:	
		• Turban	
		• Positioning Tape	
		1-ENT Split Drape 193cm x 315cm:	
		• with 7 x 71cm Split	
		• Positioning Tape	
		• CONTROL* Fabric Reinforcement	
		• Tube Holders	



EENT PACK I

Code	Description	Packaging	
88711	EENT Pack I	Items/Case = 5	
Sterile Component Code	Non-Sterile Component Code	Contents	
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm	
95111	71311	1-ULTRA Surgical Gown, with Towel, Large	
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm	
		1-Absorbent Towel, 38cm x 56cm	
		1-Suture Bag	
89431	79430	1-Head Drape, 127cm x 91cm	
89352	79351	1-EENT Split Drape with CONTROL* Reinforcement, 193cm x 315cm with 7cm x 71cm Split	



EENT PACK II

Code	Description	Packaging	
88716	EENT Pack II	Items/Case = 5	
Sterile Component Code	Non-Sterile Component Code	Contents	
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm	
95111	71311	1-ULTRA Surgical Gown In Overwrap, with Towel, Large	
95111	71311	1-ULTRA Surgical Gown with Towel, Large	
		1-Microsurgery Mayo Stand Cover	
		1-Suture Bag	
89431	79430	1-Head Drape, 127cm x 91cm	
89352	79351	1-EENT Split Drape with CONTROL* Reinforcement, 193cm x 315cm with 7cm x 71cm Split	

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

EENT DRAPES & PACKS



EENT PACK IV


Code	Description	Packaging
88719	EENT Pack IV Items/Case = 5	
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm 1-Surgical Gown with Reinforced Sleeves, In Overwrap, with Towel, Large 1-Surgical Gown, with Reinforced Sleeves, Large
89701	79720	1-Absorbent Towel, 38cm x 56cm
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
89438	79438	1-Head Drape, 104cm x 81cm 1-EENT Split Drape with Instrument Pad Reinforcement, 193cm x 315cm with 8cm x 64cm Split



EENT SPLIT DRAPE, 193cm x 315cm

- CONTROL* Fabric Reinforcement
- 7cm x 71cm Split
- Tube Holders
- Instrument Pad


Code	Description	Packaging
89352	EENT Split Drape, Sterile	Items/Case = 18
79351	EENT Split Drape, Non-Sterile	Items/Case = 25



HEAD DRAPE, 127cm x 91cm

- Turban, 83cm x 50cm


Code	Description	Packaging
89431	Head Drape, Sterile	Items/Case = 38
79430	Head Drape, Non-Sterile	Items/Case = 60



HEAD DRAPE, 104cm x 81cm

- Detachable Turban, 91cm x 61cm

Code	Description	Packaging
89438	Head Drape, Sterile	Items/Case = 38



DENTAL DRAPE, 120cm x 150cm

- Teardrop Fenestration 10cm x 15cm; 40cm From The Top
- Absorbent Reinforcement 30cm x 40cm

Code	Description	Packaging
89543	HALYARD BASICS* Dental Drape, Sterile	Items/Case = 15

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

EQUIPMENT COVERS



MAYO STAND COVERS

Code	Description	Packaging
89601	Mayo Stand Cover, Reinforced, 58cm x 137cm, Sterile	Items/Case = 54
79607	Mayo Stand Cover, Reinforced, 58cm x 137cm, Non-Sterile	Items/Case = 1200
89602	Mayo Stand Cover, Reinforced, X-Large, 76cm x 145cm, Sterile	Items/Case = 44
79602	Mayo Stand Cover, Reinforced, X-Large, 76cm x 145cm, Non-Sterile	Items/Case = 120
14007	Mayo Stand Cover, Padded, X-Large, 76cm x 145cm, Sterile	Items/Case = 44
76358	Mayo Stand Cover, Padded, X-Large, 76cm x 145cm, Non-Sterile	Items/Case = 68
89603	Mayo Stand Cover, Reinforced, XX-Large, 80cm x 145cm, Sterile	Items/Case = 44
79603	Mayo Stand Cover, Reinforced, XX-Large, 80cm x 145cm, Non-Sterile	Items/Case = 120



BACK TABLE COVERS

Code	Description	Packaging
89632	Back Table Cover, HALYARD BASICS*, 100cm x 150cm	Items/Case = 50
79482	Back Table Cover, 115cm x 152cm, Non-Sterile, Standard Fold	Items/Case = 90
42217	Back Table Cover, Standard, 112cm x 229cm, Sterile	Items/Case = 28
42218	Back Table Cover, Standard, 112cm x 229cm, Non-Sterile, Unfolded	Items/Case = 640
42231	Back Table Cover, Standard, 112cm x 229cm, Non-Sterile, Fan Fold	Items/Case = 400
42224	Back Table Cover, Standard, 152cm x 229cm, Sterile	Items/Case = 36
42225	Back Table Cover, Standard, 152cm x 229cm, Non-Sterile, Unfolded	Items/Case = 480
42233	Back Table Cover, Standard, 152cm x 229cm, Non-Sterile, Fan Fold	Items/Case = 400
42306	Back Table Cover, Heavy Duty, 152cm x 229cm, Sterile	Items/Case = 28
42307	Back Table Cover, Heavy Duty, 152cm x 229cm, Non-Sterile, Unfolded	Items/Case = 400
42366	Back Table Cover, Heavy Duty, 152cm x 229cm, Non-Sterile, Fan Fold	Items/Case = 200
42309	Back Table Cover, Heavy Duty, 177cm x 279cm, Sterile	Items/Case = 15
42310	Back Table Cover, Heavy Duty, 177cm x 279cm, Non-Sterile, Unfolded	Items/Case = 280
42367	Back Table Cover, Heavy Duty, 177cm x 279cm, Non-Sterile, Fan Fold	Items/Case = 160
29258	Back Table Cover, Padded Heavy Duty, 152cm x 229cm, Sterile	Items/Case = 20
29259	Back Table Cover, Padded Heavy Duty, 152cm x 229cm, Non-Sterile, Fan Fold	Items/Case = 22
29269	Back Table Cover, Padded Heavy Duty, 152cm x 229cm, Non-Sterile, Unfolded	Items/Case = 300

OVERHEAD / CV TABLE COVERS

Code	Description	Packaging
89621	Overhead Table Cover, 216cm x 193cm, 112cm x 152cm Skirt, Sterile	Items/Case = 14

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS


ACCESSORIES

DISPOSABLE HUCK TOWELS, 51cm x 64cm

Code	Description	Packaging
89021	Disposable Huck Towel, Sterile, Pack of 2	Items/Case = 142
79725	Disposable Huck Towel, Non-Sterile	Items/Case = 300

TAPE STRIPS, 10cm x 55cm


Code	Description	Packaging
89053	Tape Strips, Sterile, Pack of 2	Items/Case = 300



DIATHERMY POUCH, 38cm x 41cm

- Adhesive Backing Strip

Code	Description	Packaging
89945	Diathermy Pouch, Sterile	Items/Case = 150
79945	Diathermy Pouch, Non-Sterile	Items/Case = 220



UTILITY DRAPE WITH TAPE, 38cm x 64cm

Code	Description	Packaging
89731	Utility Drape with Tape, Sterile, Pack of 2	Items/Case = 100
79731	Utility Drape with Tape, Non-Sterile	Items/Case = 550

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS

BASIC PACKS

BASIC PACK VI

Code	Description	Packaging	
88161	Basic Pack VI	Items/Case = 7	
	Sterile Component Code	Non-Sterile Component Code	Contents
	89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
	95111	71311	1-ULTRA Surgical Gown, In Overwrap with Towel, Large
	95121	71321	1-ULTRA Surgical Gown, with Towel, X-Large
	89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm
	89731	79731	4-Utility Drapes, with Tape
			1-Suture Bag
	89101	79101	1-Half Drape, 152cm x 112cm
	89121	79121	1-Large Drape, 152cm x 193cm

BASIC PACK XV

Code	Description	Packaging	
88178	Basic Pack XV	Items/Case = 15	
	Sterile Component Code	Non-Sterile Component Code	Contents
	89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
	89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm
			1-Suture Bag
	89121	79121	1-Large Drape, 152cm x 193cm

SET-UP PACK I

Code	Description	Packaging	
88171	Set-Up Pack I	Items/Case = 14	
	Sterile Component Code	Non-Sterile Component Code	Contents
	89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
	95111	71311	1-Ultra Surgical Gown In Overwrap with Towel, Large
	95121	71321	1-ULTRA Surgical Gown, X-Large
	89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm
			1-Suture Bag

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


GENERAL DRAPES & PACKS

BASIC DRAPES




UTILITY DRAPE, 90cm x 90cm

Code	Description	Packaging
89353	Utility Drape, Sterile	Items/Case = 120
76323	Utility Drape, Non-Sterile	Items/Case = 200




HALF DRAPE, 152cm x 112cm

Code	Description	Packaging
89101	Half Drape, Sterile	Items/Case = 54
79106	Half Drape, Non-Sterile	Items/Case = 1000




MEDIUM DRAPE, 193cm x 112cm

Code	Description	Packaging
89111	Medium Drape, Sterile	Items/Case = 50
79112	Medium Drape, Non-Sterile	Items/Case = 800



MEDIUM DRAPE WITH TAPE, 193cm x 105cm

Code	Description	Packaging
89023	Medium Drape with Tape, Sterile	Items/Case = 50



LARGE DRAPE, 152cm x 193cm

Code	Description	Packaging
89121	Large Drape, Sterile	Items/Case = 36
79120	Large Drape, Non-Sterile	Items/Case = 600

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


GENERAL DRAPE & PACKS



THREE-QUARTER DRAPE, 152cm x 193cm


- CONTROL* Plus Fabric Reinforcement
- No Adhesive Tape

Code	Description	Packaging
89141	Three-Quarter Drape, Sterile	Items/Case = 20



X-LARGE DRAPE, 193cm x 254cm


Code	Description	Packaging
89131	X-Large Drape, Sterile	Items/Case = 20
79132	X-Large Drape, Non-Sterile	Items/Case = 24



SMALL ABSORBENT DRAPE WITH TAPE, 57cm x 60cm

- CONTROL* Plus Absorbent Fabric

Code	Description	Packaging
89828	Small Absorbent Drape with Tape, Sterile	Items/Case = 72
79828	Small Absorbent Drape with Tape, Non-Sterile	Items/Case = 108



STANDARD ABSORBENT DRAPE WITH TAPE, 91cm x 91cm

- CONTROL* Plus Absorbent Fabric

Code	Description	Packaging
89343	Standard Absorbent Drape with Tape, Sterile	Items/Case = 50
79343	Standard Absorbent Drape with Tape, Non-Sterile	Items/Case = 70

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS


BREAST / CHEST DRAPES



CHEST U-DRAPE, 193cm x 305cm

- CONTROL* Plus Fabric Reinforcement
- 30cm x 76cm Split
- Tube Holders

Code	Description	Packaging
79335	Chest U-Drape, Non-Sterile	Items/Case = 22



BREAST DRAPE, 254cm x 193cm x 315cm


- 15cm x 15cm Fenestration
- CONTROL* Fabric Reinforcement
- Hook and Loop Tube and Line Control
- Armboard Covers

Code	Description	Packaging
89236	Breast Drape, Sterile	Items/Case = 10

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS


FENESTRATED DRAPES



EPIDURAL DRAPE WITH TAPE, 50cm x 75cm

- 8cm x 10cm Oval Fenestration, Centred
- Spunbond Film
- Positioning Tape


Code	Description	Packaging
89636	HALYARD BASICS* Epidural Drape with Tape, Sterile	Items/Case = 40



VARIABLE FENESTRATED DRAPE, 75cm x 100cm

- Variable Fenestration
- Half Circle On Each Drape with 12cm Diameter
- Spunbond Film
- Positioning Tape

Code	Description	Packaging
89680	HALYARD BASICS* Variable Fenestrated Drape, Sterile	Items/Case = 30



CV LINE DRAPE, 160cm x 240cm

- CONTROL* Fabric Reinforcement
- 12cm x 12cm Fenestration

Code	Description	Packaging
89835	CV Line Drape, Sterile	Items/Case = 10
79835	CV Line Drape, Non-Sterile	Items/Case = 14

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS


LAPAROSCOPY / LAPAROTOMY DRAPES



LAPAROSCOPIC / CHOLECYSTECTOMY DRAPE, 264cm x 193cm x 305cm

- 30cm x 34cm Adjustable Fenestration, Lateral Offset
- Hook and Loop Tube and Line Control
- CONTROL* Instrument Pad Reinforcement
- Armboard Covers


Code	Description	Packaging
89229	Laparoscopic / Cholecystectomy Drape, Sterile	Items/Case = 12



LAPAROSCOPIC ABDOMINAL DRAPE, 254cm x 193cm x 315cm

- 30cm x 33cm Abdominal Fenestration
- CONTROL* Fabric Reinforcement
- Clear Instrument and Line Control Pockets
- Hook and Loop Tube and Line Control
- Armboard Covers


Code	Description	Packaging
89233	Laparoscopic Abdominal Drape, Sterile	Items/Case = 15
79233	Laparoscopic Abdominal Drape, Non-Sterile	Items/Case = 18



PAEDIATRIC LAPAROTOMY DRAPE, 152cm x 274cm

- 7cm x 20cm Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89241	Paediatric Laparotomy Drape, Sterile	Items/Case = 15
79240	Paediatric Laparotomy Drape, Non-Sterile	Items/Case = 25



LAPAROTOMY DRAPE WITH INCISE, 254cm x 183cm x 305cm

- 20cm x 33cm Adhesive Incise Area
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Armboard Covers

Code	Description	Packaging
89282	Laparotomy Drape with Incise, Sterile	Items/Case = 15
79282	Laparotomy Drape with Incise, Non-Sterile	Items/Case = 15

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS


MINOR PROCEDURE DRAPES



MINOR PROCEDURE FENESTRATED DRAPE, 76cm x 76cm

- 4cm x 5cm Fenestration
- CONTROL* Fabric Reinforcement
- Positioning Tape


Code	Description	Packaging
89206	Minor Procedure Fenestrated Drape, Sterile	Items/Case = 100
79206	Minor Procedure Fenestrated Drape, Non-Sterile	Items/Case = 250



MINOR PROCEDURE FENESTRATED DRAPE, 122cm x 152cm

- 7cm Diameter Fenestration
- CONTROL* Fabric Reinforcement
- Positioning Tape


Code	Description	Packaging
89207	Minor Procedure Fenestrated Drape, Sterile	Items/Case = 40
79207	Minor Procedure Fenestrated Drape, Non-Sterile	Items/Case = 60



MINOR PROCEDURE FENESTRATED DRAPE, 122cm x 152cm

- 12cm x 15cm Fenestration
- CONTROL* Fabric Reinforcement
- Positioning Tape

Code	Description	Packaging
89208	Minor Procedure Fenestrated Drape, Sterile	Items/Case = 40
79208	Minor Procedure Fenestrated Drape, Non-Sterile	Items/Case = 60



MINOR PROCEDURE FENESTRATED DRAPE, 193cm x 315cm

- 15cm x 15cm Fenestration
- CONTROL* Fabric Reinforcement
- Positioning Tape

Code	Description	Packaging
89209	Minor Procedure Fenestrated Drape, Sterile	Items/Case = 19
79209	Minor Procedure Fenestrated Drape, Non-Sterile	Items/Case = 23

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS

THYROID DRAPES



THYROID DRAPE, 193CM X 315CM

- 12cm x 12cm Diamond Fenestration
- CONTROL* Plus Fabric Reinforcement
- Instrument Pad
- Tube Holders


Code	Description	Packaging
89255	Thyroid Drape, Sterile	Items/Case = 13
79255	Thyroid Drape, Non-Sterile	Items/Case = 18

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.


GENERAL DRAPES & PACKS

UNIVERSAL PACKS

UNIVERSAL PACK – MAJOR WITH LONG SIDE DRAPES



Code	Description		Packaging
88759	Universal Pack – Major		Items/Case = 5
	Sterile Component Code	Non-Sterile Component Code	Contents
	89601	79601	1-Back Table Cover, Fully Reinforced, 137cm x 248cm 1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Absorbent Towel, 15cm x 56cm 1-Suture Bag
	34902	34903	1-Bottom Drape with CONTROL* Plus Reinforcement, 193cm x 195cm
	34906	34901	1-Top Drape with CONTROL* Plus Reinforcement, 273cm x 127cm
	34904	34905	2-Side Drapes with CONTROL* Plus Reinforcement, 195cm x 105cm



UNIVERSAL PACK – MINOR WITH SHORT SIDE DRAPES

Code	Description		Packaging
88769	Universal Pack – Minor		Items/Case = 14
	Sterile Component Code	Non-Sterile Component Code	Contents
	89611	79613	1-Back Table Cover, Fully-Reinforced, 112cm x 224cm
	89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
		79097	2-Side Drapes with CONTROL* Reinforcement, 114cm x 84cm
		79096	1-Bottom Drape with CONTROL* Reinforcement, 165cm x 193cm
		79095	1-Top Drape with CONTROL* Reinforcement, 276cm x 150cm

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS



UNIVERSAL PACK III WITH MEDIUM SIDE DRAPES

Code	Description	Packaging	
88774	Universal Pack – III	Items/Case = 8	
Sterile Component Code	Non-Sterile Component Code	Contents	
89601	79601	1-Back Table Cover, Fully-Reinforced, 133cm x 248cm 1-Mayo Stand Cover, 58cm x 137cm 1-Suture Bag 2-Side Drapes with CONTROL* Reinforcement, 145cm x 99cm 1-Bottom Drape with CONTROL* Reinforcement, 193cm x 193cm 1-Top Drape with CONTROL* Reinforcement, 276cm x 150cm	



UNIVERSAL PACK III WITH GOWNS

Code	Description	Packaging	
88772	Universal Pack III with Gowns	Items/Case = 6	
Sterile Component Code	Non-Sterile Component Code	Contents	
89601	79601	1-Back Table Cover, Fully-Reinforced, 133cm x 248cm 1-Mayo Stand Cover, 58cm x 137cm	
90018	72978	1-Surgical Gown, Large with 79728 Huck Towel x 1 In 90cm x 90cm Wrap	
90048	72988	1-Surgical Gown, X-Large with 79728 Huck Towel x 1 1-Suture Bag 2-Side Drapes with CONTROL* Reinforcement, 145cm x 99cm 1-Bottom Drape with CONTROL* Reinforcement, 193cm x 193cm 1-Top Drape with CONTROL* Reinforcement, 276cm x 150cm	

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS



U-BAR* PACK

Code	Description	Packaging	
88334	U-BAR* Pack	Items/Case = 16	
Sterile Component Code	Non-Sterile Component Code	Contents	
89601	79601	1-Back Table Cover, Zone Reinforced, 112cm x 224cm 1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag 1-U-Drape with CONTROL* Reinforcement, 180cm x 315cm, with 10cm x 102cm Split 1-Bar Drape with CONTROL* Reinforcement, 105cm x 228cm	



U-BAR* PACK – IV

Code	Description	Packaging	
88773	U-BAR* Pack – IV	Items/Case = 10	
Sterile Component Code	Non-Sterile Component Code	Contents	
89601	79601	1-Back Table Cover, Reinforced, 133cm x 248cm 1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-U-Drape with CONTROL* Reinforcement, 192cm x 315cm, with 11cm x 51cm Split 1-Bar Drape with CONTROL* Reinforcement, 248cm x 145cm	

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

GENERAL DRAPES & PACKS


UNIVERSAL DRAPES



TOP DRAPE, 273cm x 127cm

- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape
- From Universal Pack 88759


Code	Description	Packaging
34906	Top Drape, Sterile	Items/Case = 15
34901	Top Drape, Non-Sterile	Items/Case = 275



BOTTOM DRAPE, 193cm x 195cm

- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape
- From Universal Pack 88759

Code	Description	Packaging
34902	Bottom Drape, Sterile	Items/Case = 12
34903	Bottom Drape, Non-Sterile	Items/Case = 242



SIDE DRAPE, 195cm x 105cm

- CONTROL* Plus Fabric Reinforcement
- Positioning Tape
- From Universal Pack 88759

Code	Description	Packaging
34904	Side Drape, Sterile	Items/Case = 20
34905	Side Drape, Non-Sterile	Items/Case = 440

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


GENERAL DRAPES & PACKS



TOP DRAPE, 276cm x 150cm

- CONTROL* Fabric Reinforcement
- Tube Holders
- Positioning Tape
- From Universal Pack 88769


Code	Description	Packaging
79095	Top Drape, Non-Sterile	Items/Case = 60



BOTTOM DRAPE, 165cm x 193cm

- CONTROL* Fabric Reinforcement
- Tube Holders
- Positioning Tape
- From Universal Pack 88769

Code	Description	Packaging
79096	Bottom Drape, Non-Sterile	Items/Case = 70




SIDE DRAPE, 114cm x 84cm

- CONTROL* Fabric Reinforcement
- Positioning Tape
- From Universal Pack 88769

Code	Description	Packaging
79097	Side Drape, Non-Sterile	Items/Case = 200

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.


NEURO / SPINE DRAPE



UNIVERSAL SPINE DRAPE, 284cm x 251cm x 350cm

- 10cm x 30cm Fenestration
- 30cm x 17cm Clear Instrument Pouches
- Hook and Loop Tube and Line Control
- Armboard Covers
- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
29221	Universal Spine Drape, Sterile	Items/Case = 8
29224	Universal Spine Drape, Non-Sterile	Items/Case = 9



CRANIOTOMY DRAPE WITH POUCH, 310cm x 188cm x 340cm

- CONTROL* Plus Fabric Reinforcement
- 20cm x 30cm Oval Fenestration
- Clear Gusseted Anesthesia Panels
- Screen and Suction Port


Code	Description	Packaging
89062	Craniotomy Drape with Pouch, Sterile	Items/Case = 9
79075	Craniotomy Drape with Pouch, Non-Sterile	Items/Case = 15

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.

OBSTETRICS & GYNAECOLOGY

DRAPES & PACKS


C-SECTION DRAPES



C-SECTION DRAPE WITH FENESTRATED INCISE, 254cm x 183cm x 305cm

- 38cm x 33cm Adhesive Incise Area with Triangular Fenestration
- 270° Fluid Collection Pouch, Malleable Band, with Suction Port
- Tube Holders
- Armboard Covers

Code	Description	Packaging
89528	C-Section Drape with Fenestrated Incise, Sterile	Items/Case = 8
79528	C-Section Drape with Fenestrated Incise, Non-Sterile	Items/Case = 150



C-SECTION DRAPE WITH FULL INCISE, 254cm x 183cm x 305cm

- 38cm x 33cm Adhesive Incise Area
- 270° Fluid Collection Pouch, Malleable Band, with Suction Port
- Tube Holders
- Armboard Covers


Code	Description	Packaging
89529	C-Section Drape with Full Incise, Sterile	Items/Case = 8
79529	C-Section Drape with Full Incise, Non-Sterile	Items/Case = 150

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.


OBSTETRICS & GYNAECOLOGY DRAPES & PACKS

LAPAROSCOPY PACKS

LAPAROSCOPY PACK I



Code	Description	Packaging
88088	Laparoscopy Pack I	Items/Case = 10
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm 1-Absorbent Towel, 38cm x 56cm
89411		1-Under Buttocks Drape, Plastic, 75cm x 112cm
89408	79408	2-Leggings with 15cm Cuff, 78cm x 122cm
89119	79117	1-Laparoscopy "T" Drape with CONTROL* Plus Reinforcement, 249cm x 104cm x 251cm • 15cm x 20cm Rectangular Abdominal Fenestration • 7cm x 15cm Perineal Fenestration with Flap • Tube Holders



LAPAROSCOPY / PELVISCOPY PACK

Code	Description	Packaging
88222	Laparoscopy/Pelviscopy Pack	Items/Case = 6
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
95111	71311	1-ULTRA Surgical Gown, In Overwrap, with Towel, Large
95111	71311	1-ULTRA Surgical Gown, with Towel, Large
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
89413		1-Under Buttocks Drape, Plastic-Reinforced, Large 101cm x 116cm
89408	79408	2-Leggings, with 15cm Cuff, Sterile Pair, 78cm x 122cm 1-Laparoscopy / Pelviscopy Drape with CONTROL* Plus Reinforcement (without Leggings), 264cm x 89cm x 264cm • 18cm x 23cm Trapezoidal Abdominal Fenestration • 8cm x 15cm Perineal Fenestration with Flap

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

OBSTETRICS & GYNAECOLOGY

DRAPES & PACKS


LAPAROSCOPY DRAPES



LAPAROSCOPY / PELVISCOPY DRAPE WITH FLUID COLLECTION POUCH, 284cm x 76cm x 314cm

- 17cm x 23cm Trapezoidal Abdominal Fenestration
- 7cm x 20cm Perineal Fenestration with Flap
- Fluid Collection Pouch with Screen and Suction Port
- CONTROL* Plus Fabric Reinforcement
- Attached Leggings
- Hook and Loop Line Controls
- Armboard Covers

Code	Description	Packaging
89219	Laparoscopy/Pelviscopy Drape Drape with Pouch, Sterile	Items/Case = 10



LAPAROSCOPY / PELVISCOPY DRAPE, 254cm x 76cm x 312cm

- 17cm x 23cm Trapezoidal Abdominal Fenestration
- 7cm x 20cm Perineal Fenestration with Flap
- CONTROL* Plus Fabric Reinforcement
- Attached Leggings
- Hook and Loop Tube and Line Controls

Code	Description	Packaging
89224	Laparoscopy/Pelviscopy Drape, Sterile	Items/Case = 10
79224	Laparoscopy/Pelviscopy Drape, Non-Sterile	Items/Case = 10

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

OBSTETRICS & GYNAECOLOGY


DRAPES & PACKS

LITHOTOMY PACKS



LITHOTOMY PACK II

Code	Description	Packaging
88521	Lithotomy Pack II	Items/Case = 10
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm 1-Absorbent Towel, 38cm x 56cm
89413	79413	1-Under Buttocks Drape, Plastic-Reinforced, 102cm x 117cm
89408	79408	2-Leggings with 15cm Cuff, 78cm x 122cm
89101	79101	1-Half Drape, 152cm x 111cm



LITHOTOMY PACK V

Code	Description	Packaging
88551	Lithotomy Pack V	Items/Case = 12
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
89701	79720	1-Absorbent Towel, 38cm x 56cm
89408	79408	2-Leggings, with 15cm Cuffs
	79551	1-Lithotomy "T" Drape with CONTROL* Plus Reinforcement, 152cm x 78cm x 190cm with 8cm x 15cm Perineal Fenestration

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

OBSTETRICS & GYNAECOLOGY DRAPES & PACKS


LITHOTOMY DRAPES



LITHOTOMY DRAPE, 254cm x 76cm x 221cm

- 7cm x 15cm Perineal Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Attached Leggings


Code	Description	Packaging
89203	Lithotomy Drape, Sterile	Items/Case = 16



GYN / UROLOGY DRAPE WITH POUCH, 160cm x 284cm x 332cm

- Pouch with Mesh Screen and Port
- Attached Leggings
- CONTROL* Plus Fabric Reinforcement
- 7cm x 15cm Perineal Fenestration
- Tube Holders
- Armboard Covers


Code	Description	Packaging
89204	Gyn/Urology Drape with Pouch, Sterile	Items/Case = 10
79204	Gyn/Urology Drape with Pouch, Non-Sterile	Items/Case = 10



LITHOTOMY "T" DRAPE, 152cm x 78cm x 190cm

- 7cm x 15cm Perineal Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
79551	Lithotomy "T" Drape, Non-Sterile	Items/Case = 55



LEGGINGS, 78cm x 122cm, WITH 15cm CUFF

Code	Description	Packaging
89408	Leggings, Sterile Pair	Items/Case = 26
79408	Leggings, Non-Sterile	Items/Case = 70

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

OBSTETRICS & GYNAECOLOGY

DRAPES & PACKS


UNDER BUTTOCKS DRAPES



UNDER BUTTOCKS DRAPE, ZONE-REINFORCED, MEDIUM, 96cm x 112cm


- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
79411	Under Buttocks Drape, Zone-Reinforced, Medium, Non-Sterile	Items/Case = 115



UNDER BUTTOCKS DRAPE, PLASTIC REINFORCED, LARGE, 101cm x 116cm


Code	Description	Packaging
89413	Under Buttocks Drape, Plastic Reinforced, Large, Sterile	Items/Case = 50



UNDER BUTTOCKS DRAPE WITH FLUID COLLECTION POUCH, 101cm x 112cm

- CONTROL* Plus Absorbent Pad

Code	Description	Packaging
89414	Under Buttocks Drape, with Fluid Collection Pouch, Sterile	Items/Case = 30
79414	Under Buttocks Drape, with Fluid Collection Pouch, Non-Sterile	Items/Case = 80



UNDER BUTTOCKS DRAPE WITH FLUID COLLECTION POUCH, 101cm x 112cm

- Screen and Suction Port
- Adhesive Pouch Closure

Code	Description	Packaging
89415	Under Buttocks Drape with Fluid Collection Pouch, Sterile	Items/Case = 40
79415	Under Buttocks Drape with Fluid Collection Pouch, Non-Sterile	Items/Case = 100

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


OPHTHALMIC DRAPE & PACKS



BASIC OPHTHALMIC DRAPE, 132cm x 142cm

- Lightweight Fabric
- 12cm x 9cm Adhesive Incise Area
- Fluid Collection Pouch, with Eye Wick

Code	Description	Packaging
89304	Basic Ophthalmic Drape, Sterile	Items/Case = 40
79304	Basic Ophthalmic Drape, Non-Sterile	Items/Case = 100



OCUARTS* FULL BODY OPHTHALMIC DRAPE, 183cm x 292cm

- 11cm x 11cm Adhesive Incise Area
- 5cm x 3cm Oval Fenestration
- 12cm x 15cm Optional Incise Drape
- 270° Fluid Collection Pouch, Foam Band
- Breather Bar
- Tube Holders


Code	Description	Packaging
89356	OCUARTS* Full Body Ophthalmic Drape, Sterile	Items/Case = 10
79356	OCUARTS* Full Body Ophthalmic Drape, Non-Sterile	Items/Case = 20

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.


ORTHOPAEDIC DRAPES & PACKS

ORTHOPAEDIC PACKS

ORTHO PACK I – MAJOR



Code	Description	Packaging
88411	Ortho Pack I – Major	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
92321	30734	2-MICROCOOL* Breathable High Performance Surgical Gown, with Towel, X-Large
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
89331	79330	1-Impervious Split Drape, 152cm x 178cm with 10cm x 53cm Split 1-Ortho U-Drape with CONTROL* Plus Reinforcement, 223cm x 274cm with 10cm x 102cm Split
89381	79360	1-Ortho Bar Drape with CONTROL* Plus Reinforcement, 254cm x 152cm



ORTHOARTS* HIP PACK WITH POCKETS

Code	Description	Packaging
88452	ORTHOARTS* Hip Pack with Pockets	Items/Case = 3
Sterile Component Code	Non-Sterile Component Code	Contents
89906	79006	1-Back Table Cover, Heavy-Duty, 152cm x 228cm
95111	71311	1-ULTRA Surgical Gown, In Overwrap, with Towel, Large
92321	30734	2-MICROCOOL* Breathable High Performance Surgical Gown, with Towel, X-Large
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
89821	79515	1-Impervious Stockinette, Large
89121	79121	1-Large Drape, 152cm x 193cm
89333	79333	1-Impervious Split Drape, Large, 152cm x 205cm with 10cm x 79cm Split
89345	79345	1-ORTHOARTS* Hip Drape with CONTROL* Plus Reinforcement, 289cm x 223cm x 320cm 1-Yankauer Suction Tip 1-304cm Suction Tubing

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

ORTHOPAEDIC DRAPES & PACKS



ORTHOARTS* LOWER EXTREMITY

Code	Description	Packaging
88471	ORTHOARTS* Lower Extremity Pack	Items/Case = 3
Sterile Component Code	Non-Sterile Component Code	Contents
89906	79006	1-Back Table Cover, Heavy-Duty, 152cm x 228cm
95111	71311	1-ULTRA Surgical Gown, In Overwrap, with Towel, Large
		1-Extremity Prep Pad
95221	74921	1-ULTRA Fabric-Reinforced Surgical Gown, with Towel, X-Large
92331		1-MICROCOOL* Breathable High Performance Specialty Surgical Gown, with Towel, X-Large
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm
		1-Suture Bag
		1-304cm Suction Tubing
89821	79515	1-Impervious Stockinette, Large
89331	79330	1-Impervious Split Drape, 152cm x 178cm with 10cm x 53cm Split
89278	13233	1-ORTHOARTS* Lower Extremity Drape with CONTROL* Plus Reinforcement, 3-Drape Tape Strips

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

ORTHOPAEDIC DRAPES & PACKS


EXTREMITY DRAPES



LOWER EXTREMITY DRAPE, 274cm x 203cm x 315cm

- Armboard Covers
- 6cm Conformable Circular Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders


Code	Description	Packaging
89016	Lower Extremity Drape, Sterile	Items/Case = 12



LOWER EXTREMITY DRAPE, 223cm x 332cm

- CONTROL* Plus Fabric Reinforcement
- 6cm Conformable Circular Fenestration
- Tube Holders


Code	Description	Packaging
89276	Lower Extremity Drape, Sterile	Items/Case = 14



ORTHOARTS* LOWER EXTREMITY DRAPE, 289cm x 223cm x 320cm

- Armboard Covers
- 6cm Conformable Circular Fenestration
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89278	ORTHOARTS* Lower Extremity Drape, Sterile	Items/Case = 7
13233	ORTHOARTS* Lower Extremity Drape, Non-Sterile	Items/Case = 84



ORTHOARTS* EXTRA LONG LOWER EXTREMITY DRAPE, 289cm x 223cm x 381cm

- Armboard Covers
- 6cm Conformable Circular Fenestration
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
29219	ORTHOARTS* Extra Long Lower Extremity Drape, Sterile	Items/Case = 8
29220	ORTHOARTS* Extra Long Lower Extremity Drape, Non-Sterile	Items/Case = 8

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


ORTHOPAEDIC DRAPE & PACKS



UNIVERSAL EXTREMITY DRAPE, 229cm x 333cm

- CONTROL* Plus Fabric Reinforcement


Code	Description	Packaging
89279	Universal Extremity Drape, Sterile	Items/Case = 14
79279	Universal Extremity Drape, Non-Sterile	Items/Case = 150



BILATERAL LIMB DRAPE, 193cm x 304cm

- CONTROL* Plus Fabric Reinforcement
- Two 6cm Conformable Circular Fenestrations
- Tube Holders

Code	Description	Packaging
89291	Bilateral Limb Drape, Sterile	Items/Case = 14



BILATERAL LIMB DRAPE WITH ARMBORD COVERS, 295cm x 223cm x 366cm


- CONTROL* Plus Fabric Reinforcement
- Two 6cm Conformable Circular Fenestrations
- Tube Holders
- Armboard Covers

Code	Description	Packaging
89380	Bilateral Limb Drape with Armboard Covers, Sterile	Items/Case = 12

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.

ORTHOPAEDIC DRAPES & PACKS


HAND DRAPES



HAND DRAPE, 289cm x 360cm

- 157cm x 147cm Hand Table Cover
- 3cm Conformable Circular Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
77450	Hand Drape, Sterile	Items/Case = 11
13227	Hand Drape, Non-Sterile	Items/Case = 108



ORTHOARTS* HAND DRAPE, 193CM X 381CM


- 157cm x 147cm Gusseted Hand Table Cover
- 3cm Conformable Fenestration
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89272	ORTHOARTS* Hand Drape, Sterile	Items/Case = 6

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.

ORTHOPAEDIC DRAPES & PACKS


HIP DRAPES



ORTHOARTS* HIP DRAPE WITH LEG POCKETS, 289cm x 223cm x 320cm

- Leg Pockets On Both Sides
- Armboard Covers
- 21cm x 17cm Conformable Teardrop Fenestration
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders


Code	Description	Packaging
89345	ORTHOARTS* Hip Drape with Leg Pockets, Sterile	Items/Case = 6
79345	ORTHOARTS* Hip Drape with Leg Pockets, Non-Sterile	Items/Case = 6



ORTHOARTS* HIP DRAPE, 289cm x 223cm x 320cm

- Armboard Covers
- 21cm x 17cm Conformable Teardrop Fenestration
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89346	ORTHOARTS* Hip Drape, Sterile	Items/Case = 7
79346	ORTHOARTS* Hip Drape, Non-Sterile	Items/Case = 7



HIP DRAPE WITH CLEAR LEG POCKETS, 294cm x 223cm x 345cm

- Clear Leg Pockets
- Armboard Covers
- 21cm x 17cm Conformable Teardrop Fenestration
- CONTROL* Plus Fabric Reinforcement
- Hook and Loop Tube and Line Control

Code	Description	Packaging
89349	Hip Drape with Clear Leg Pockets, Sterile	Items/Case = 9
79349	Hip Drape with Clear Leg Pockets, Non-Sterile	Items/Case = 10

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

ORTHOPAEDIC DRAPES & PACKS


KNEE ARTHROSCOPY DRAPES



KNEE ARTHROSCOPY DRAPE, 228cm x 315cm

- Fluid Collection Pouch with Suction Port
- 6cm Conformable Circular Fenestration
- Tube Holders


Code	Description	Packaging
89265	Knee Arthroscopy Drape, Sterile	Items/Case = 10
79265	Knee Arthroscopy Drape, Non-Sterile	Items/Case = 150



KNEE ARTHROSCOPY DRAPE WITH ARMBORD COVERS, 290cm x 223cm x 315cm

- Armboard Covers
- Fluid Collection Pouch with Suction Port
- 6cm Conformable Circular Fenestration
- Tube Holders


Code	Description	Packaging
89371	Knee Arthroscopy Drape with Armboard Covers, Sterile	Items/Case = 10
79371	Knee Arthroscopy Drape with Armboard Covers, Non-Sterile	Items/Case = 13



ORTHOARTS* KNEE ARTHROSCOPY DRAPE, 289cm x 223cm x 294cm

- Armboard Covers
- Fluid Collection Pouch with Reflux Shields and Suction Port
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- 6cm Conformable Circular Fenestration
- Tube Holders

Code	Description	Packaging
89267	ORTHOARTS* Knee Arthroscopy Drape, Sterile	Items/Case = 7
79267	ORTHOARTS* Knee Arthroscopy Drape, Non-Sterile	Items/Case = 7



ORTHOARTS* EXTRA LONG KNEE ARTHROSCOPY DRAPE, 294cm x 233cm x 447cm

- Armboard Covers
- Fluid Collection Pouch with Reflux Shields and Suction Port
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- 6cm Conformable Circular Fenestration
- Tube Holders

Code	Description	Packaging
36930	ORTHOARTS* Extra Long Knee Arthroscopy Drape, Sterile	Items/Case = 10
42429	ORTHOARTS* Extra Long Knee Arthroscopy Drape, Non-Sterile	Items/Case = 10

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

ORTHOPAEDIC DRAPES & PACKS


SHOULDER DRAPES



ORTHOARTS* BEACH CHAIR SHOULDER ARTHROSCOPY DRAPE, 411cm x 261cm

- Attached Fluid Collection
- Hook and Loop Tube and Line Control
- 6cm Conformable Tear Drop Fenestration
- Heavy Duty Fabric


Code	Description	Packaging
36932	ORTHOARTS* Beach Chair Shoulder Arthroscopy Drape, Sterile	Items/Case = 10
42426	ORTHOARTS* Beach Chair Shoulder Arthroscopy Drape, Non-Sterile	Items/Case = 10



BEACH CHAIR SHOULDER DRAPE WITH ATTACHABLE FLUID COLLECTION POUCH, 437cm x 254cm

- Attachable Fluid Collection Pouch with Suction Port, 89cm x 109cm
- CONTROL* Plus Fabric Reinforcement
- Hook and Loop Tube and Line Control
- Pouch Attachment Zones

Code	Description	Packaging
89066	Beach Chair Shoulder Drape with Attachable Fluid Collection Pouch, Sterile	Items/Case = 6
79049	Beach Chair Shoulder Drape without Attachable Fluid Collection Pouch, Non-Sterile	Items/Case = 11
79061	Attachable Fluid Collection Pouch, Non-Sterile	Items/Case = 44



BEACH CHAIR SHOULDER ARTHROSCOPY DRAPE II, 406cm x 259cm

- Attached Fluid Collection
- Hook and Loop Tube and Line Control
- 6cm Conformable Tear Drop Fenestration
- Instrument Pad Fabric Reinforcement

Code	Description	Packaging
89070	Beach Chair Shoulder Arthroscopy Drape, Sterile	Items/Case = 10
79070	Beach Chair Shoulder Arthroscopy Drape, Non-Sterile	Items/Case = 12

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


ORTHOPAEDIC DRAPES & PACKS



2-PIECE SHOULDER ARTHROSCOPY DRAPE PACK

- Body Split Drape, 223cm x 294cm
- 10cm x 68cm Split
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Shoulder Arthroscopy Drape, 254cm x 157cm
- Fluid Collection Pouch with Suction Ports
- Tube Holders


Code	Description	Packaging
89490	2-Piece Shoulder Arthroscopy Drape Pack, Sterile	Items/Case = 12
79490	2-Piece Shoulder Arthroscopy Drape Pack, Non-Sterile	Items/Case = 12



ORTHOARTS* SHOULDER ARTHROSCOPY DRAPE, 254cm x 152cm

- Fluid Collection Pouch
- Heavy-Duty Fabric

Code	Description	Packaging
89492	ORTHOARTS* Shoulder Arthroscopy Drape, Sterile	Items/Case = 12
79492	ORTHOARTS* Shoulder Arthroscopy Drape, Non-Sterile	Items/Case = 15



ORTHOARTS* SHOULDER SPLIT DRAPE, 284cm x 152cm

- 15cm x 20cm Dual Splits
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89274	ORTHOARTS* Shoulder Split Drape, Sterile	Items/Case = 12

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

ORTHOPAEDIC DRAPES & PACKS

U & BAR DRAPES


U-DRAPE, 193cm x 305cm



- 10cm x 101cm Split
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape

Code	Description	Packaging
89301	U-Drape, Sterile	Items/Case = 17
13231	U-Drape, Non-Sterile	Items/Case = 190


LARGE U-DRAPE, 223cm x 317cm



- 15cm x 137cm Split
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape

Code	Description	Packaging
89311	Large Ortho U-Drape, Sterile	Items/Case = 11
79363	Large Ortho U-Drape, Non-Sterile	Items/Case = 15


ORTHOARTS* U-DRAPE, 223CM X 274CM



- 10cm x 101cm Split
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
89347	ORTHOARTS* U-Drape, Sterile,	Items/Case = 10
13294	ORTHOARTS* U-Drape, Non-Sterile	Items/Case = 132

ORTHOARTS* BODY SPLIT DRAPE, 223CM X 294CM



- 10cm x 68cm Split
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape

Code	Description	Packaging
89275	ORTHOARTS* Body Split Drape, Sterile	Items/Case = 10
79275	ORTHOARTS* Body Split Drape, Non-Sterile	Items/Case = 12

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


ORTHOPAEDIC DRAPES & PACKS



IMPERVIOUS SPLIT DRAPE, 152cm x 178cm

- Polyethylene Material
- 10cm x 53cm Split


Code	Description	Packaging
89331	Impervious Split Drape, Sterile	Items/Case = 42
79330	Impervious Split Drape, Non-Sterile	Items/Case = 115



LARGE IMPERVIOUS SPLIT DRAPE, 152cm x 205cm

- Polyethylene Material
- 11cm x 78cm Split


Code	Description	Packaging
89333	Large Impervious Split Drape, Sterile	Items/Case = 34
79333	Large Impervious Split Drape, Non-Sterile	Items/Case = 100



ORTHO SPLIT "T" DRAPE, 297cm x 210cm x 342cm

- Armboard Covers
- 10cm x 53cm Split
- CONTROL* Fabric Reinforcement
- Tube Holders
- Positioning Tape

Code	Description	Packaging
36931	Ortho Split "T" Drape, Sterile	Items/Case = 10



ORTHO BAR DRAPE, 254cm x 152cm

- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89381	Ortho Bar Drape, Sterile	Items/Case = 18
79360	Ortho Bar Drape, Non-Sterile	Items/Case = 28

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


ORTHOPAEDIC DRAPES & PACKS



VERTICAL ISOLATION DRAPE, 228cm x 320cm

- Clear Film
- Adhesive Incise Area 40cm x 76cm
- Attachable Fluid Collection Pouch with Suction Port
- Tube Holders
- Central Release Tape Strip For Horizontal Rod Attachment

Code	Description	Packaging
77723	Vertical Isolation Drape, Sterile	Items/Case = 14



ISOLATION SCREEN, 170cm x 290cm

- Clear Plastic Drape 170cm x 290cm
- Adhesive Backing Strips

Code	Description	Packaging
89093	Isolation Screen, Sterile	Items/Case = 24
79093	Isolation Screen, Non-Sterile	Items/Case = 30



STOCKINETTES

Code	Description	Packaging
89811	Impervious Stockinette, Sterile, Medium, 23cm x 112cm	Items/Case = 36
79505	Impervious Stockinette, Non-Sterile, Medium, 23cm x 112cm	Items/Case = 75
89821	Impervious Stockinette, Sterile, Large, 30cm x 122cm	Items/Case = 25
79515	Impervious Stockinette, Non-Sterile, Large, 30cm x 122cm	Items/Case = 50

All codes in **bold** indicate a sterile product.

† Reinforcements include:


CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

UROLOGY DRAPEs & PACKS




T.U.R. PACK

Code	Description	Packaging
77885	T.U.R. Pack	Items/Case = 12
	Non-Sterile Component Code	Contents
		1-Back Table Cover, Zone-Reinforced, 152cm x 191cm
		1-Absorbent Towel, 38cm x 56cm
76885	1-T.U.R. Drape, 289cm x 175cm x 173cm:	
		• Attached Leggings
		• Fluid Collection Pouch with Screen, Suction Port and Tie
		• Finger Cot



CYSTOSCOPY PACK I

Code	Description	Packaging
88611	Cystoscopy Pack I	Items/Case = 12
	Sterile Component Code	Non-Sterile Component Code
		Contents
		1-Back Table Cover, Zone-Reinforced, 112cm x 198cm
		1-Absorbent Towel, 38cm x 56cm
89408	79408	2-Leggings, with 15cm Cuff, 78cm x 122cm
		1-Cystoscopy "T" Drape with CONTROL* Plus Reinforcement, 152cm x 76cm x 213cm



PERCUTANEOUS NEPHROLITHOTOMY DRAPE (PCNL), 299cm x 182cm

- CONTROL* Fabric Reinforcement
- Tube Holders
- Fluid Collection Pouch

Code	Description	Packaging
42427	Percutaneous Nephrolithotomy Drape, Sterile	Items/Case = 12

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

MICROBIAL SEALANT



HALYARD* INTEGUSEAL* MICROBIAL SEALANT

- Ergonomic Single-Use Applicator
- Evenly Covers Incision Area in a Single Stroke
- Widely-Used Cyanoacrylate Forms a Safe, Fast-Drying, Breathable Barrier
- Compatible with Most Common Surgical Techniques and Materials
- Need Not be Removed for Skin Closure
- Can be Used with a Variety of Preps and Incise Drapes
- Nonflammable, Latex-Free, Non-Irritating and Non-Cytotoxic

Code	Description	Coverage Area (approx)	Example of Procedures	Packaging
33736	IS 100	10" x 10" 25cm x 25cm	Cranial and Spinal Procedures, Laparotomies, Hernia Repair, Thyroid Surgery, Cardiac Catheterisation	Applicators/Box = 10 Boxes/Case = 2 Applicators/Case = 20
33737	IS 200	20" x 10" 50cm x 25cm	Cholecystectomy, Open Heart, Lung Surgery, Total Joint Replacement	Applicators/Box = 10 Boxes/Case = 2 Applicators/Case = 20



Infection Prevention Product Solutions
Halyard* Education Foundation
Tools & Best Practice
Commitment to Excellence and Sustainability
Experienced Customer Support
Knowledgeable Sales Force
Innovation & Technology



Halyard Australia Pty. Limited

ABN No 25 100 844 488

Level 7, 52 Alfred Street, Milsons Point NSW 2061

Phone: 1800 101 021

Email: CustomerService_ANZ@hyh.com

Fax: 1800 730 687

Halyard New Zealand

Office 2–5, Building 10

666 Great South Road, Penrose, Auckland 1051

Phone: 0800 482 211

www.halyardhealth.com.au



Halyard Health makes no representation that its products described in this catalogue will provide complete protection from exposure to contaminants, fluids or micro-organisms or will eliminate the risk of contracting infectious diseases.

*Registered Trademark or Trademark of Halyard Health, Inc. or its affiliates. © 2014 HYH. All rights reserved.
HYHM:INT2015:WFHY549