

DENTAL PRODUCTS CATALOG

HERE'S THE BATTLE PLAN
YOU NEED TO BE
VICTORIOUS IN THE
DENTAL MARKETPLACE



ABOUT HALYARD...

Halyard Health is a medical technology company focused on preventing infection, eliminating pain and speeding recovery. Just as a halyard fuels forward movement, Halyard Health's clinically-superior products and remarkable service help to advance health and healthcare worldwide. Halyard sells its recognized brands and products in more than 100 countries, and holds leading market positions in multiple categories across the portfolio.


Dental professionals are increasingly challenged to provide infection control for their patients, while protecting themselves. As one of the world's leading suppliers of infection control products, you can count on Halyard* for products you use every day in your practice including exam gloves, face masks, eyewear, sterilization wrap, dental surgical drapes and hand hygiene solutions.

This catalog details the wide array of products available to dental professionals.

For more information, visit
www.halyarddental.com

CONTENTS

2	Facial Protection
6	Protective Eyewear
10	Exam Gloves
14	Sterilization Wrap
17	Other Products



FACIAL PROTECTION



FACIAL PROTECTION



STANDARD PROCEDURE MASKS

- Standard masks for routine patient exams – no fluid exposure expected
- Pleat-style with earloops

Item Code	Description	Quantity
49700	Value Procedure Mask, Earloops, Blue	50 masks/box 10 boxes/case
47080	Standard Procedure Mask, Earloops, Blue	50 masks/box 10 boxes/case
47085	Standard Procedure Mask, Earloops, Green	50 masks/box 10 boxes/case
47090	Standard Procedure Mask, Earloops, White	50 masks/box 10 boxes/case
47117	Standard Procedure Mask, Earloops, Yellow	50 masks/box 10 boxes/case
62362	Fog-Free Procedure Mask, Foam Band, Earloops, Blue	50 masks/box 10 boxes/case
62363	So Soft* Procedure Mask, Fog-Free, Earloops, White	50 masks/box 10 boxes/case



FLUIDSHIELD* LEVEL 1 PROCEDURE & SURGICAL MASKS, ASTM F2100-11

- Meets ASTM F2100-11 Level 1 standards for fluid resistance, filtration and breathability
- Assured fluid protection for routine patient care when low level of exposure to fluids expected
- Pleat-style with inner layer of SO SOFT* fabric for comfort

Item Code	Description	Quantity
25867	FLUIDSHIELD* Level 1 Procedure Mask with SO SOFT* Lining, Earloops, Yellow	50 masks/box 10 boxes/case
25868	FLUIDSHIELD* Level 1 Procedure Mask with SO SOFT* Lining, Earloops, Lavender	50 masks/box 10 boxes/case
25869	FLUIDSHIELD* Level 1 Procedure Mask with SO SOFT* Lining, Earloops, Blue	50 masks/box 10 boxes/case
28802	FLUIDSHIELD* Level 1 Surgical Mask with SO SOFT*Lining, Ties, Blue	50 masks/box 6 boxes/case
28806	FLUIDSHIELD* Level 1 Fog-Free Surgical Mask with SO SOFT*Lining, Green	50 masks/box 6 boxes/case



FACIAL PROTECTION



FLUIDSHIELD* LEVEL 2 PROCEDURE & SURGICAL MASKS, ASTM F2100-11

- For protection from moderate risk of fluid splash and spray
- Meets ASTM F2100-11 Level 2 standards for fluid resistance, filtration and breathability
- Pleat-style available with earloops or ties



Item Code	Description	Quantity
28810	FLUIDSHIELD* Level 2 Procedure Mask, Earloops, Blue	50 masks/box 10 boxes/case
28821	FLUIDSHIELD* Level 2 Fog-Free Procedure Mask with Earloops, Wraparound Visor, Blue	25 masks/box 4 boxes/case
28808	FLUIDSHIELD* Level 2 Fog-Free Surgical Mask with Ties, Blue with Orange Cross Hatch	50 masks/box 6 boxes/case
28809	FLUIDSHIELD* Level 2 Fog-Free Surgical Mask with Ties and Wraparound Visor, Blue with Orange Cross Hatch	25 masks/box 4 boxes/case



FLUIDSHIELD* LEVEL 3 PROCEDURE & SURGICAL MASKS, ASTM F2100-11

- Protects from fluid splattering and splashing
- Meets ASTM F2100-11 Level 3 standards for fluid resistance, filtration and breathability
- Pleat-style available with earloops or ties



Item Code	Description	Quantity
28797	FLUIDSHIELD* Level 3 Fog-Free Procedure Mask with Earloops, Foam Band, Orange	40 masks/box 10 boxes/case
28800	FLUIDSHIELD* Level 3 Fog-Free Procedure Mask with Earloops and Wraparound Visor, Foam Band, Orange	25 masks/box 4 boxes/case
28820	FLUIDSHIELD* Level 3 Fog-Free Surgical Mask with Ties, Foam Band, Orange	50 masks/box 6 boxes/case
28804	FLUIDSHIELD* Level 3 Fog-Free Surgical Mask with Ties and Wraparound Visor, Foam Band, Orange	25 masks/box 4 boxes/case
00146	FLUIDSHIELD* Level 3 Fog-Free Procedure Mask with Earloops, Wraparound Splashguard* Visor, Blue	25 masks/box 4 boxes/case
00148	FLUIDSHIELD* Level 3 Fog-Free Procedure Mask with Earloops, Foam Band, Blue	40 masks/box 10 boxes/case



FACIAL PROTECTION



FLUIDSHIELD* LEVEL 3 MASKS WITH SO SOFT* LINING, ASTM F2100-11

- Protects from fluid splattering and splashing
- Meets ASTM F2100-11 Level 3 standards for fluid resistance, filtration and breathability
- Pleat-style available with earloops or ties
- Features SO SOFT* bio component fabric for soft comfort


Item Code	Description	Quantity
47107	FLUIDSHIELD* Level 3 Fog-Free Procedure Mask with SO SOFT* Lining, Foam Band, Earloops, Orange	40 masks/box 10 boxes/case
47147	FLUIDSHIELD* Level 3 Fog-Free Procedure Mask with SO SOFT* Lining, Wraparound Visor, Foam Band, Earloops, Orange	25 masks/box 4 boxes/case
48207	FLUIDSHIELD* Level 3 Fog-Free Surgical Mask with SO SOFT* Lining, Foam Band, Ties, Orange	50 masks/box 6 boxes/case
48247	FLUIDSHIELD* Level 3 Fog-Free Surgical Mask with SO SOFT* Lining, Wraparound Visor, Foam Band, Ties, Orange	25 masks/box 4 boxes/case



PARTICULATE FILTER RESPIRATORS

- NIOSH rated N95 particulate filter respirators, pouch-style
- For situations with airborne exposure risk
- Fluid resistant at 160 mm Hg. • Polyurethane headband

Item Code	Description	Quantity
46727	FLUIDSHIELD* N95 Particulate Filter Respirator and Surgical Mask, Regular Size, Orange	35 masks/box 6 boxes/case
46827	FLUIDSHIELD* N95 Particulate Filter Respirator and Surgical Mask, Small Size, Orange	35 masks/box 6 boxes/case



FACE SHIELDS

- Protect from splash, splatter and spray • Foam band • Fog-resistant

Item Code	Description	Quantity
41204	GUARDALL SHIELD*, Full Length Visor	40 shields/box 1 box/case
41205	GUARDALL SHIELD*, 3/4 Length Visor	40 shields/box 1 box/case

PROTECTIVE EYEWEAR



PROTECTIVE EYEWEAR



V30 NEMESIS* SAFETY EYEWEAR



- Sporty, flexible, lightweight design
- Soft-touch temples for added comfort
- Single-lens wraparound protection
- Every pair includes a neck cord
- Available in reader style lenses
- 99.9% uv PROTECTION
- Meets ANSI Z87.1+ Impact Standard
- Select eyewear is available with CSA and CE certification



Item Code	Description	Quantity
25676	Clear Lens, Black Frame	12 pairs/box
25679	Clear Anti-Fog Lens, Black Frame	12 pairs/box
25688	Smoke Mirror Lens, Black Frame	12 pairs/box
22475	Smoke Anti-Fog Lens, Black Frame	12 pairs/box
25685	Indoor/Outdoor Lens, Black Frame	12 pairs/box
14481	Blue Mirror Lens, Black Frame	12 pairs/box
25659	Amber Lens, Black Frame	12 pairs/box
22476	Amber Anti-Fog Lens, Black Frame	12 pairs/box
19639	Light Blue Lens, Blue Frame	12 pairs/box
19642	Copper Blue Shield Lens, Black Frame	12 pairs/box
22610	Amber Lens, Camo Frame	12 pairs/box
22611	Smoke Lens, Inferno/Red Frame	12 pairs/box
19644	Bronze Lens, Camo Frame	12 pairs/box
22608	Clear Anti-Fog Lens, Camo Frame	12 pairs/box
22609	Smoke Anti-Fog Lens, Camo Frame	12 pairs/box
25692	IRUV Shade 3.0 Lens, Black Frame	12 pairs/box
25671	IRUV Shade 5.0 Lens, Black Frame	12 pairs/box

PROTECTIVE EYEWEAR



V30 NEMESIS* S SAFETY EYEWEAR



- Sporty eyewear that fits the smaller, slimmer head without adjustments
- Soft touch temples for added comfort and fun temple tip colors for individuality
- Contemporary styling, available in a variety of frame colors and lens tints

Item Code	Description	Quantity
38474	Clear Hard Coat Lens, Black Frame with Black Tips	12 pairs/box
38476	Smoke Hard Coat Lens, Black Frame with Black Tips	12 pairs/box
38478	Clear Anti-Fog Lens, Black Frame with Pink Tips	12 pairs/box
38480	I/O Hard Coat Lens, Black Frame with Green Tips	12 pairs/box



JACKSON SAFETY* V60 NEMESIS* RX SAFETY EYEWEAR



- "Cheaters" style with diopters for vision assistance
- Larger lens for wraparound impact protection
- Enhanced nosepiece design to channel away sweat
- Every pair includes a neck cord

Item Code	Description	Quantity
28618	NEMESIS*-Rx Reader: Clear Lens, Black Frame, +1.0	6 pairs/box
28621	NEMESIS*-Rx Reader: Clear Lens, Black Frame, +1.5	6 pairs/box
28624	NEMESIS*-Rx Reader: Clear Lens, Black Frame, +2.0	6 pairs/box
28627	NEMESIS*-Rx Reader: Clear Lens, Black Frame, +2.5	6 pairs/box
28630	NEMESIS*-Rx Reader: Clear Lens, Black Frame, +3.0	6 pairs/box



PROTECTIVE EYEWEAR



SAFEVIEW* GLASSES – DISPOSABLE EYEWEAR

- Protection from splash, splatter or spray
- Multiple colors

Item Code	Description
SV50A	50 Assembled Glasses
SV100F	100 Replacement Frames
SV250L	250 Replacement Lenses



EXAM GLOVES



EXAM GLOVES



LAVENDER* NITRILE EXAM GLOVES

- For low-fluid, low-abrasion-risk procedures
- Perfect for routine patient care
- Ambidextrous
- Beaded cuffs
- Non-sterile
- Textured fingertips
- Thinner and lighter than standard nitrile
- High tactile sensitivity
- Higher count per box saves space and waste
- Excellent value

Item Code	Description	Quantity
52816	LAVENDER* Nitrile Exam Gloves X-Small	250 gloves/box 10 boxes/case
52817	LAVENDER* Nitrile Exam Gloves Small	250 gloves/box 10 boxes/case
52818	LAVENDER* Nitrile Exam Gloves Medium	250 gloves/box 10 boxes/case
52819	LAVENDER* Nitrile Exam Gloves Large	250 gloves/box 10 boxes/case
52820	LAVENDER* Nitrile Exam Gloves X-Large	230 gloves/box 10 boxes/case



STERLING* NITRILE EXAM GLOVES

- For procedures with low-to-moderate risk of exposure to fluid and abrasion
- Cleared for use in chemotherapy
- Flexible with high tactile sensitivity
- Non-sterile
- Ambidextrous
- Textured fingertips
- Beaded cuffs
- Provides nitrile protection while avoiding latex

Item Code	Description	Size	Quantity
50705	STERLING* Nitrile Exam Gloves	X-Small	200 gloves/box 10 boxes/case
50706	STERLING* Nitrile Exam Gloves	Small	200 gloves/box 10 boxes/case
50707	STERLING* Nitrile Exam Gloves	Medium	200 gloves/box 10 boxes/case
50708	STERLING* Nitrile Exam Gloves	Large	200 gloves/box 10 boxes/case
50709	STERLING* Nitrile Exam Gloves	X-Large	170 gloves/box 10 boxes/case



STERLING* NITRILE-XTRA* EXAM GLOVES

- 12" length for more protection against fluid and splash
- Flexible with high tactile sensitivity
- Nitrile protection while avoiding latex
- Extended, beaded cuff for better coverage over the wrist
- Certified to meet NFPA 1999-2008
- Cleared for use in chemotherapy

Item Code	Description	Size	Quantity
53137	STERLING* Nitrile-XTRA* Exam Gloves	X-Small	100 gloves/box 10 boxes/case
53138	STERLING* Nitrile-XTRA* Exam Gloves	Small	100 gloves/box 10 boxes/case
53139	STERLING* Nitrile-XTRA* Exam Gloves	Medium	100 gloves/box 10 boxes/case
53140	STERLING* Nitrile-XTRA* Exam Gloves	Large	100 gloves/box 10 boxes/case
53141	STERLING* Nitrile-XTRA* Exam Gloves	X-Large	100 gloves/box 10 boxes/case

EXAM GLOVES



STERLING* NITRILE STERILE PAIRS

- Sterile 9.5" Exam Glove • Cuffed and Walled
- For procedures with low-to-moderate risk of exposure to fluid and abrasion
- Flexible with high tactile sensitivity
- Provides nitrile protection while avoiding latex

Item Code	Description	Size	Quantity
53130	STERLING* Nitrile Sterile Pairs Exam Gloves	Small	50 pairs/box 4 boxes/case
53131	STERLING* Nitrile Sterile Pairs Exam Gloves	Medium	50 pairs/box 4 boxes/case
53132	STERLING* Nitrile Sterile Pairs Exam Gloves	Large	50 pairs/box 4 boxes/case



STERLING* NITRILE-XTRA* STERILE PAIRS

- Sterile 12" length for more protection against fluid and splash • Cuffed and Walled
- For procedures with low-to-moderate risk of exposure to fluid and abrasion
- Flexible with high tactile sensitivity
- Provides nitrile protection while avoiding latex

Item Code	Description	Size	Quantity
33022	STERLING* Nitrile-XTRA* Sterile Pairs Exam Gloves	X-Small	50 pairs/box 4 boxes/case
33023	STERLING* Nitrile-XTRA* Sterile Pairs Exam Gloves	Small	50 pairs/box 4 boxes/case
33024	STERLING* Nitrile-XTRA* Sterile Pairs Exam Gloves	Medium	50 pairs/box 4 boxes/case
33025	STERLING* Nitrile-XTRA* Sterile Pairs Exam Gloves	Large	50 pairs/box 4 boxes/case
33026	STERLING* Nitrile-XTRA* Sterile Pairs Exam Gloves	X-Large	50 pairs/box 4 boxes/case



STERLING* NITRILE-XTRA* STERILE SINGLES

- Sterile 12" length for more protection against fluid and splash • Cuffed and Walled
- For procedures with low-to-moderate risk of exposure to fluid and abrasion
- Flexible with high tactile sensitivity
- Provides nitrile protection while avoiding latex

Item Code	Description	Size	Quantity
33027	STERLING* Nitrile-XTRA* Sterile Singles Exam Gloves	X-Small	100 gloves/box 4 boxes/case
33028	STERLING* Nitrile-XTRA* Sterile Singles Exam Gloves	Small	100 gloves/box 4 boxes/case
33029	STERLING* Nitrile-XTRA* Sterile Singles Exam Gloves	Medium	100 gloves/box 4 boxes/case
33030	STERLING* Nitrile-XTRA* Sterile Singles Exam Gloves	Large	100 gloves/box 4 boxes/case
33031	STERLING* Nitrile-XTRA* Sterile Singles Exam Gloves	X-Large	100 gloves/box 4 boxes/case

EXAM GLOVES



PURPLE NITRILE* DENTAL EXAM GLOVES

- Proven barrier protection from bacteria, bloodborne viruses and chemicals during procedures where fluid exposure risk is moderate to high
- Available in both sterile and non-sterile versions
- 9.5" exam glove

Item Code	Description	Size	Quantity
53430	PURPLE NITRILE* Dental Exam Gloves	X-Small	100 gloves/box 10 boxes/case
53431	PURPLE NITRILE* Dental Exam Gloves	Small	100 gloves/box 10 boxes/case
53432	PURPLE NITRILE* Dental Exam Gloves	Medium	100 gloves/box 10 boxes/case
53433	PURPLE NITRILE* Dental Exam Gloves	Large	100 gloves/box 10 boxes/case
53434	PURPLE NITRILE* Dental Exam Gloves	X-Large	90 gloves/box 10 boxes/case



PURPLE NITRILE-XTRA* EXAM GLOVES

- Proven barrier protection from bacteria, bloodborne viruses and chemicals during procedures where fluid exposure risk is moderate to high
- Available in both sterile and non-sterile versions
- NFPA 1999-2008 Certified
- 12" exam glove

Item Code	Description	Size	Quantity
55090	PURPLE NITRILE-XTRA* Exam Gloves	X-Small	50 gloves/box 10 boxes/case
50601	PURPLE NITRILE-XTRA* Exam Gloves	Small	50 gloves/box 10 boxes/case
50602	PURPLE NITRILE-XTRA* Exam Gloves	Medium	50 gloves/box 10 boxes/case
50603	PURPLE NITRILE-XTRA* Exam Gloves	Large	50 gloves/box 10 boxes/case
50604	PURPLE NITRILE-XTRA* Exam Gloves	X-Large	50 gloves/box 10 boxes/case

STERILIZATION WRAP




STERILIZATION WRAP



HALYARD* SEQUENTIAL STERILIZATION WRAP H100 - COMPACT PACKAGING

- Strength and tear resistance
- Superior bacterial filtration efficiency (BFE) for airborne microbes
- Low linting
- Flame resistance
- Compatibility with steam and Ethylene Oxide (EO) sterilant; not appropriate for dry heat sterilization



Item Code	Size	Quantity
37159	9"x 9" / 22cm x 22cm	100 sheets/bag 10 bags/case 1000 sheets/case
37046	12"x 12" / 30cm x 30cm	100 sheets/bag 10 bags/case 1000 sheets/case
37047	15"x 15" / 38cm x 38cm	100 sheets/bag 10 bags/case 1000 sheets/case
37049	18"x 18" / 45cm x 45cm	100 sheets/bag 10 bags/case 1000 sheets/case
37050	20"x 20" / 50cm x 50cm	100 sheets/bag 10 bags/case 1000 sheets/case
37051	24"x 24" / 60cm x 60cm	100 sheets/bag 5 bags/case 500 sheets case

*Larger Sizes available

HALYARD* SEQUENTIAL STERILIZATION WRAP

- Sequential wrap
- Low-lint SMS fabric
- Steam, EO and vaporized hydrogen peroxide compatible
- Each application is imprinted with its name, size and lot code
- Made in the USA



Item Code	Size	Quantity
10712	12"x 12" / 30cm x 30cm	1000 sheets per case
10715	15"x 15" / 38cm x 38cm	1000 sheets per case
10718	18"x 18" / 45cm x 45cm	1000 sheets per case
10720	20"x 20" / 50cm x 50cm	1000 sheets per case
10724	24"x 24" / 60cm x 60cm	500 sheets per case
10730	30"x 30" / 76cm x 76cm	300 sheets per case
10736	36"x 36" / 91cm x 91cm	300 sheets per case
10740	40"x 40" / 101cm x 101cm	250 sheets per case
10745	45"x 45" / 114cm x 114cm	250 sheets per case
10748	48"x 48" / 121cm x 121cm	250 sheets per case
10754	54"x 54" / 137cm x 137cm	100 sheets per case
10772	54"x 72" / 137cm x 182cm	100 sheets per case

STERILIZATION WRAP



HALYARD* ONE-STEP* STERILIZATION WRAP H100

- Bonded simultaneous (non-sequential) wrap
- One application is equivalent to 2 sheets of wrap
- Low-lint SMS fabric
- Steam, EO and vaporized hydrogen peroxide compatible
- Protection of double wrapping in a quick, single step
- Each application is imprinted with its name, size and lot code
- Made in the USA



Item Code	Size	Quantity
12712	12"x 12" / 30cm x 30cm	480 applications per case
12715	15"x 15" / 38cm x 38cm	480 applications per case
12718	18"x 18" / 45cm x 45cm	480 applications per case
12720	20"x 20" / 50cm x 50cm	480 applications per case
12724	24"x 24" / 60cm x 60cm	240 applications per case
12730	30"x 30" / 76cm x 76cm	144 applications per case
12736	36"x 36" / 91cm x 91cm	144 applications per case
12740	40"x 40" / 101cm x 101cm	120 applications per case
12745	45"x 45" / 114cm x 114cm	120 applications per case
12748	48"x 48" / 121cm x 121cm	120 applications per case
12754	54"x 54" / 137cm x 137cm	96 applications per case

OTHER PRODUCTS



APPAREL



BASIC LAB COATS

- Protective 3-layer SMS fabric
- Protective in situations with potential for light fluid contact
- Traditional collar and knit cuff, snaps and 2 pockets

Item Code	Description	Size	Quantity
10120	Basic Lab Coat, White	Small	25 coats/case
10121	Basic Lab Coat, White	Medium	25 coats/case
10122	Basic Lab Coat, White	Large	25 coats/case
10123	Basic Lab Coat, White	X-large	25 coats/case
10124	Basic Lab Coat, White	XX-large	25 coats/case



LAB JACKETS

- Protective 3-layer SMS fabric
- Protective in situations with potential for moderate fluid contact
- Knit collar and cuffs
- 5 snaps and 3 pockets

Item Code	Description	Size	Quantity
10069	Universal Precautions Lab Jacket, White	Small	25 jackets/case
10071	Universal Precautions Lab Jacket, White	Medium	25 jackets/case
10072	Universal Precautions Lab Jacket, White	Large	25 jackets/case
10073	Universal Precautions Lab Jacket, White	X-large	25 jackets/case
10067	Universal Precautions Lab Jacket, Blue	Medium	25 jackets/case
10078	Universal Precautions Lab Jacket, Blue	Large	25 jackets/case
10087	Universal Precautions Lab Jacket, Blue	X-large	25 jackets/case

ICE PACKS



STAY-DRY* ICE PACKS

- Consistent cold therapy
- 3-layer construction
- Minimizes condensation and leaking
- Soft outer covering
- Clip closure and 4 ties

Item Code	Description	Size	Quantity
33495	STAY-DRY* Ice Pack Small, 4 Ties	5" x 12"	25 units/box 2 boxes/case
33500	STAY-DRY* Ice Pack Large, 4 Ties	6.5" x 14"	25 units/box 2 boxes/case



SOFT 'N COLD ICE PACKS

- Consistent cold therapy
- Inner plastic ice pack
- Soft knit outer cover
- Clip or LOC-TITE* closures

Item Code	Description	Size	Quantity
33625	Soft 'N Cold Ice Pack with Clip Closure, 4 ties, Small	4.5" x 10.5"	20 units/box 2 boxes/case
33630	Soft 'N Cold Ice Pack with Clip Closure, 4 ties, Large	6.25" x 9.5"	20 units/box 2 boxes/case



INSTANT ICE PACKS

Item Code	Description	Size	Quantity
59688	Instant Cold Pack, Large	6.25" x 8.5"	24 units/box 1 box/case
59681	Instant Cold Pack, Small	4" x 10"	24 units/box 1 box/case



BILATERAL FACIAL ICE PACKS

- Consistent cold therapy
- 3-layer construction
- Minimizes condensation
- Soft outer cover
- Single ice chamber for easy filing
- Conforms to upper face and eyes

Item Code	Description	Size	Quantity
33101	Bilateral Facial Ice Pack, 1 strap	5" x 12"	12 units/box 2 boxes/case

STERILIZATION GOWNS, DRAPES AND PACKS



STERILE GOWNS

Item Code	Description	Size
95111	HALYARD* ULTRA* Surgical Gown	Large
95121	HALYARD* ULTRA* Surgical Gown	X-Large
99284	HALYARD* BASICS* Surgical Gown	Large
99285	HALYARD* BASICS* Surgical Gown	X-Large



PACKS

Item Code	Description
88171	Set-Up Pack
88111	Basic Pack

STERILE DRAPES

Item Code	Description	Size
89534	HALYARD* BASICS* Universal Half Drape	60" x 44"
89141	HALYARD* BASICS* Universal Three-Quarter Drape	60" x 76"
89535	HALYARD* BASICS* Universal Large Drape	60" x 76"
89543	HALYARD* BASICS* Dental Drape	47" x 59"

TISSUE, TOWELS AND WIPERS



TISSUE

Item Code	Description	Quantity
28075	KLEENEX® BOUTIQUE Anti-Viral Facial Tissue White, 3-ply, 8.2" x 8.2"	75 sts. 27 boxes/case 2025 sts./case
12456	KLEENEX® COTTONELLE® Ultra Soft Bathroom Tissue White, 1-ply, 4.2" x 4.1"	200 sts. 48 boxes/case 9600 sts./case



TOWELS

Item Code	Description	Quantity
01701	KLEENEX® Hand Towels in a POP-UP* Box White, 1-ply, 9" x 10.5"	120 sts. 18 boxes/case 2160 sts./case
11268	KLEENEX® Ultra Soft Hand Towels in a POP-UP* Box White, 1-ply, 8.9" x 10"	70 sts. 18 boxes/case 2160 sts./case
04442	KLEENEX® SLIMFOLD* Towels, White, 1-ply, 7.5" x 11.6"	90 sts. 24 boxes/case 2160 sts./case
12388	SCOTT® SLIMROLL* Towels, White, 1-ply, White, 1-ply, 8" x 580'	580' per 6 rolls/case 3,480 ft./case pack
01960	SCOTT® SCOTTFOLD* M Towels, White, 1-ply, 7.84" x 12.4"	175 sts. 25 clips/case 4,375 sts./case



SURFACE DISINFECTION WIPES

Item Code	Description	Quantity
06411	KIMTECH® WETTASK® Wiping System White, 12" x 12.5"	90 sts. 6 rolls/case 540 sts./case



HAND HYGIENE SOLUTIONS



KLEENEX® ULTRA MOISTURIZING FOAM HAND SANITIZER (70%)

Item Code	Description	Size	Quantity
34604	Ultra Moisturizing Foam Hand Sanitizer	Bottle - 1.5 oz.	24/case
34610	Ultra Moisturizing Foam Hand Sanitizer	Bottle - 8 oz.	12/case
34618	Ultra Moisturizing Foam Hand Sanitizer	Bottle - 16 oz.	6/case
34643	Ultra Moisturizing Foam Hand Sanitizer	Electronic Cassette - 1.2 L	2/case
34700	Ultra Moisturizing Foam Hand Sanitizer	Manual Cassette - 1.0 L	6/case



KLEENEX® MOISTURIZING INSTANT HAND SANITIZER (62%)

Item Code	Description	Size	Quantity
31294	Moisturizing Instant Hand Sanitizer - Gel	Bottle - 8.0 oz.	12/case
33958	Moisturizing Instant Hand Sanitizer - Gel	Tear Pack - 0.05 oz.	1000/case
33976	Moisturizing Instant Hand Sanitizer - Gel	Bottle - 1.0 oz.	144/case
34089	Moisturizing Instant Hand Sanitizer - Foam	Bottle - 8.0 oz.	12/case
34091	Moisturizing Instant Hand Sanitizer - Foam	Bottle - 1.5 oz.	24/case
34767	Moisturizing Instant Hand Sanitizer - Foam	Bottle - 16.0 oz.	6/case
91560	Moisturizing Instant Hand Sanitizer - Foam	Manual Cassette - 1.0 L	6/case
91562	Moisturizing Instant Hand Sanitizer - Gel	Manual Cassette - 1.0 L	6/case
91590	Moisturizing Instant Hand Sanitizer - Foam	Electronic Cassette - 1.2 L	2/case
91600	Moisturizing Instant Hand Sanitizer - Gel	Electronic Cassette - 1.2 L	2/case
93050	Moisturizing Instant Hand Sanitizer - Gel	Bottle - 2.0 oz.	24/case



KLEENEX® GENERAL PURPOSE SOAP

Item Code	Description	Size	Quantity
91552	Foam Skin Cleanser with Moisturizers	Manual Cassette - 1.0 L	6/case
91592	Foam Skin Cleanser with Moisturizers	Electronic Cassette - 1.2 L	2/case
91556	Gentle Lotion Skin Cleanser	Manual Cassette - 1.0 L	6/case
91565	Green Certified Foam Skin Cleanser (Green Seal)	Manual Cassette - 1.0 L	6/case
91591	Green Certified Foam Skin Cleanser (Green Seal)	Electronic Cassette - 1.2 L	2/case

HAND HYGIENE SOLUTIONS



KLEENEX® ANTIBACTERIAL SOAP

Item Code	Description	Size	Quantity
91554	Antibacterial Foam Skin Cleanser (0.5% Triclosan)	Manual Cassette - 1.0 L	6/case
91594	Antibacterial Foam Skin Cleanser (0.5% Triclosan)	Electronic Cassette - 1.2 L	2/case
91559	Antibacterial Liquid Skin Cleanser (0.25% Triclosan)	Manual Cassette - 1.0 L	6/case
91555	E-2 Foam Skin Cleanser (1.75% PCMX)	Manual Cassette - 1.0 L	6/case
91595	E-2 Foam Skin Cleanser (1.75% PCMX)	Electronic Cassette - 1.2 L	2/case



KLEENEX® LOTION

Item Code	Description	Size	Quantity
35362	Kleenex® Moisturizing Hand and Body Lotion	Manual Cassette - 1.0 L	6/case
35363	Kleenex® Moisturizing Hand and Body Lotion	Bottle - 8 oz.	12/case



KLEENEX® DISPENSERS

Item Code	Description	Quantity
32499	Kleenex® MOD Electronic Touchless Skin Care Cassette, White	1/case
92144	Kleenex® Cassette Skin Care Dispenser, White	1/case



THE HEALTHY WORKPLACE PROJECT* DESK CADDY

THE HEALTHY WORKPLACE PROJECT* DESK CADDY

Item Code	Description	Size	Quantity
32181	Standard Desk Caddy (includes Products)	11" x 5.25" x 8.625"	4/case
32179	Premium Desk Caddy (includes products)	10.7" x 7.37" x 8.375"	4/case
30893	KLEENEX® Anti-Viral† Facial Tissue		68 sheets 27/case
31294	KLEENEX® Moisturizing Instant Hand Sanitizer	8 oz.	12/case
30760	SCOTT® Disinfectant Wipes		35 wipes 12/case

†In the tissue within 15 minutes. Virucidal against: Rhinoviruses Type 1A and 2 (rhinoviruses are the leading cause of the common cold), Influenza A and Influenza B (causes of the flu), Respiratory Syncytial Virus (RSV, the leading cause of lower respiratory infection in children).





KNOWLEDGE NETWORK* Accredited Education
In-Service Training and Technical Support
Credentialed Sales representatives
Tools & Best Practices
Clinical Research
Commitment to Excellence

For more information, please visit:

www.halyarddental.com

Call 1-844-HALYARD (1-844-425-9273)
in the United States and Canada.



NEMESIS, KLEENEX, KIMTECH, WETTASK, SCOTT, KLEENGUARD,
COTTONELLE and JACKSON SAFETY are Registered Trademarks
of Kimberly-Clark Worldwide, Inc. (c) 2015 KCWW.

*Registered Trademark or Trademark of Halyard Health, Inc. or its affiliates.
© 2015 HYH. All rights reserved. IP52 C1496