

# **HALYARD\* CHILD'S FACE MASK**

Featuring Disney® Characters



## THE FIRST MASK THAT MEETS FDA GUIDELINES FOR PEDIATRIC USE!

The Centers for Disease Control and Prevention protocols for respiratory etiquette recommend that patients with undiagnosed respiratory symptoms wear face masks in waiting rooms. That means active pediatric practices and hospital admission or waiting areas should have a child-size mask available for kids. We were the first company to design a mask specifically to fit small faces. Understanding that adult masks don't provide the fit, comfort and protection children need, Halyard Health now has the FIRST pediatric mask that meets FDA Guidelines for use with children. It's your clear choice for helping protect patients from colds, flu and other respiratory threats.

### **CHILD-SIZE FACE MASK WITH DISNEY® CHARACTERS MAKES PROTECTION FUN!**

Kid-friendly Disney characters add a sense of fun that helps them feel more comfortable with the idea of wearing a face mask.

#### **The HALYARD\* Child's Face Mask Features:**

- Colorful Disney characters on the mask and packaging
- Soft, lightweight materials for comfort
- Comfortable knitted earloops
- Formable nose wire for a custom fit on small faces

#### **The HALYARD\* Child's Face Mask is recommended for use by children who:**

- Are immuno-compromised
- Have flu-like symptoms
- Are visiting hospital patients



# HALYARD\* CHILD'S FACE MASK

Featuring Disney® Characters

## HELP YOUR PATIENTS UNDERSTAND RESPIRATORY ETIQUETTE

We can help your patients and parents understand and comply with the CDC Guidelines for Respiratory Etiquette with a variety of information accessories and products.

- Child Mask wall mounted bracket
- Downloadable pediatric Respiratory Etiquette poster
- Downloadable general Respiratory Etiquette poster

[www.halyardhealth.com](http://www.halyardhealth.com)



## HALYARD\* CHILD'S FACE MASK

Stock #	Descriptions	Masks/Box	Box/Case
32856	Child Face Mask	75	10
	PFE (0.1 µm)	BFE (3 µm)	ΔP (mm H <sub>2</sub> O/cm <sup>2</sup> )
	≥ 94%	≥ 96%	< 2.0



KNOWLEDGE NETWORK\* Accredited Education  
In-Service Training and Technical Support  
Credentialed Sales Representatives  
Tools & Best Practices  
Clinical Research  
Commitment to Excellence

For more information, please visit:

[www.halyardhealth.com](http://www.halyardhealth.com)

Call 1-844-HALYARD (1-844-425-9273)  
in the United States and Canada.

[www.halyardhealth.com](http://www.halyardhealth.com)



\*Registered Trademark or Trademark of Halyard Health, Inc. or its affiliates. © 2014 HYH. All rights reserved.  
© Kleenex is a registered trademark of Kimberly-Clark.  
Disney Elements ©DISNEY

IP52 C1486