



SweepZone® 310R Recessed Ultrasonic Cleaning Systems

Hinged Cover Design

- Provides additional work surface with flush lid design
- Unique recessed cover collects condensation
- Optional push button drain with built-in safety circuits shut off after two hours



Additional Features & Benefits

- SweepZone Technology insures optimum cleaning
- Digital Touch Pad Controls for accuracy and clean ability
- Internal Tank Dimensions: 15.5" L x 9" W x 6" D; 39.4cm x 22.9 cm x 15.2 cm
- Generator Dimensions: 16.25" L x 10.13" W x 6.25" H; 41.3cm x 25.7cm x 15.9cm
- Tank Assembly Surface Cut Out Dimensions: 16 1/4" W x 9 11/16" D Front to Back; 41.28 cm x 24.6 cm
- Tank Capacity: 3.25 Gallons, 12.3 Liters
- Control Panel Dimensions: 4 7/8"W x 9 3/4" H x 7"D including plug; 12.38 cm x 24.77 cm x 17.78 cm
- Shipping Weight: 26 lbs/11.80 kg

Model	Description
AG916	SweepZone 310R w/Remote Digital Timer, Drain, Mounting Rim
AG916LT	SweepZone 310R w/Hinged Lid, Remote Digital Timer, Drain, Mounting Rim
AG916LA	SweepZone 310R w/Hinged Lid, Remote Digital Timer, Auto Drain, Mounting Rim
26665	310/310R Cassette Rack – Holds 3 Large Cassettes
13518	310/310R Full Size Solid Stainless Steel Basket

Note: Models above shown with hinged lid upgrade

SweepZone® 310R Recessed Ultrasonic Cleaning Systems



310 Cassette Rack (PC 26665)



Instrument Cassette (PC 17752)



Accessory Kit (PC 535)



Solid Side Basket (PC 13518)

SweepZone 360R (PC AG915)

- Ideal for Orthodontics
- Timer, Drain & Tank Cover
- Accepts 5 large cassettes

Internal Tank Dimensions:

11.5" L x 9.25" W x 8.0" D

29.2 cm x 23.5 cm x 20.3 cm

Generator Dimensions:

16.25" L x 10.13" W x 6.25" H

41.3cm x 25.7cm x 15.9cm

Tank Capacity: 3.59 Gallons, 13.6 Liters

Shipping Weight: 22 lbs.

Product Code	Instrument Cassettes & Racks
17752	Large Cassette 310, 310R, 650, 650R
19368	Cassette Rack for 360R, Holds 3 Large Cassettes
17583	Accessory Box
17582	Bur Box (Holds 32 Burs)
26665	310, 310R Cassette Rack (Holds 3 Large Cassettes or 6 Small Cassettes)
15185	650, 650R Cassette Rack (Holds 3 Large Cassettes or 6 Small Cassettes)
	2 Racks Fit 1 Unit

Product Code	Stainless Steel Solid Side Draining Baskets
13518	310, 310R (Full Size)
18202	310, 310R (1/2 Size)
18676	360, 360R (Full Size)
17954	650, 650R (Full Size)

SweepZone 650R (PC AG917)

- Timer, Drain & Tank Cover
- Accepts 6 large cassettes

Internal Tank Dimensions:

19.75" L x 11.75" W x 6.5" D

50.2 cm x 29.8 cm x 16.5 cm

Generator Dimensions:

16.25" L x 10.13" W x 6.25" H

41.3cm x 25.7cm x 15.9cm

Tank Capacity: 6.50 Gallons, 24.6 Liters

Shipping Weight: 38 lbs.

Product Code	Accessory Kits
535	310, 310R - 1/2 Size Solid Side Basket, 1/2 Size Positioning Cover, 2-600 ML Beakers Complete, Bur Box, Bur Tray
1214	360, 360R - Solid Side Basket, Positioning Cover, 4-600 ML Beakers Complete, Bur Box, Bur Tray
1132	650, 650R - 1/2 Size Solid Side Basket, 1/2 Size Positioning Cover, 4-600 ML Beakers Complete, Bur Box, Bur Tray



L&R Manufacturing Company
577 Elm Street, P.O. Box 607
Kearny, NJ 07032-0607 USA
(P) 201.991.5331 (F) 201.991.5870
www.LRultrasonics.com info@LRultrasonics.com