

welcome to an

L&R CLEAN world

NEW

Barrier Milk

Ready to Use Spray

Instrument Protector/Lubricant



Reorder # 00139

Barrier Milk Ready to Use Spray helps prevent corrosion of instruments when used prior to sterilization. Protects and extends the life of instruments during autoclaving.

Simply spray the instrument directly before autoclaving with the unique water-based blend of lubricants and corrosion inhibitors. The fine mist works its way into hinges and box locks as well as covering all instruments with protection.

Specially formulated for use during Autoclaving with:

- Dental and Surgical instruments
- Hinges
- Box Locks
- General Equipment

Product Features:

- Convenient and Ready to use spray
- Corrosion inhibitor
- Lubricate
- Nitrate-Free
- Chromate-Free
- Phenolic-Free
- Non-Irritating
- Neutral PH
- Mild, Pleasant Odor



L&R MANUFACTURING COMPANY

577 Elm Street, PO Box 607, Kearny, NJ 07032-0607 USA Telephone: 201.991.5330 Fax: 201.991.5870
www.LRultrasonics.com info@LRultrasonics.com

