



Large IDC Cassette PC 17752

INSTRUMENT CASSETTES



IDC Large Cassette

Instrument Delivery Cassettes™

Our Instrument Delivery Cassette (IDC) System is a safe and efficient way to clean, store and sterilize dental/medical instruments. Protecting you and your staff from accidental punctures, the IDC system reduces the risk of cross-contamination while organizing and eliminating the need to handle contaminated instruments.

Features:

- Interlocking feet for compatibility with other major cassette brands, allowing the IDC to be added to your existing stock of cassettes.
- Stainless steel construction for enhanced cavitation and more effective cleaning.
- Detachable covers for added convenience.
- Instant snap-lock features for easier use.



IDC Accessory Box PC 17583



Bur Box PC 17582

Product Code	Description
17752	Large 18 Instrument Cassette for Quantrex & Sweepzone 310, 310R, 650, 650R
17583	Accessory Box holds small items to be ultrasonically cleaned and autoclaved—small enough to fit into most ultrasonic cleaning systems.
17582	Bur Box organizes and securely holds up to 32 burs for ultrasonic cleaning and sterilization. Fits inside the large cassette.
26665	310, 310R Cassette Rack accepts large, small & hygiene-type cassettes. Fits all Quantrex 310 and SweepZone 310 units. Utilizes entire internal tank.
19368	Cassette Rack 360R
26847	650, 650R Cassette Rack (Holds 3 Large Cassettes or 6 Small) 2 Racks Fit 1 Unit



310 Cassette Rack PC 26665