



because clean matters.



HOME

ABOUT US

APPLICATION BULLETINS

PRODUCT LITERATURE

PRODUCT LITERATURE

CONTACT US

DENTAL / MEDICAL ACCESSORIES

[Quantrex Systems](#)

[SweepZone Systems](#)

[Ultrasonic Cleaning Solutions](#)



[Accessory Kits](#)



**Stainless Steel
Auxiliary Pans**

Product Code	Stainless Steel Auxiliary Pans
11032	Quantrex® 90
10634	Quantrex 140
19366	SweepZone® 200
11093	Quantrex 210
12267	Quantrex 360 SweepZone 360R
13102	Quantrex 650 SweepZone 650R Full Size
17740	Quantrex 650 SweepZone 650R 1/2 Size

Instrument Delivery Cassettes™

Our Instrument Delivery Cassette (IDC) System is a safe and efficient way to clean, store and sterilize dental/medical instruments. Protecting you and your staff from accidental punctures, the IDC system reduces the risk of cross-contamination while organizing eliminating the need to handle contaminated instruments.

Features:

- Interlocking feet for compatibility with other major cassette brands, allowing the IDC to be added to your existing stock of cassettes
- Stainless steel construction for enhanced cavitation and more effective cleaning
- Detachable covers for added convenience
- Instant snap-lock features for easier use
- High-quality silicone rubber strips for holding instruments in place



Large IDC

18-instrument cassette works with the Quantrex and SweepZone 310, 310R and 650R.

Product Code	Model	Kit Contains
454	140	Solid Side Basket, Plastic Positioning Cover, 2-600 ML BeakersComplete, Bur Tray
535	310 310R	1/2 Size Solid Side Basket, 1/2 Size Positioning Cover, 2-600 ML BeakersComplete, Bur Box, Bur Tray
659	200	3/4 Size Solid Side Basket, Positioning Cover, 2-600 ML BeakersComplete, Bur Box, Bur Tray
735	PC-3	2-100 ML Beakers Complete 1-100 ML Positioning Cover
1132	650 650R	1/2 Size Solid Side Basket, 1/2 Size Positioning Cover, 4-600 ML BeakersComplete, Bur Box, Bur Tray
1214	360 360R	Solid Side Basket, Positioning Cover, 4-600 ML BeakersComplete, Bur Tray

Stainless Steel Draining

IDC Bur Box

Organizes and securely holds

Baskets

Product Code	Stainless Steel Solid Side Draining Baskets
10632	Quantrex 140
18670	SweepZone 200 (3/4 Size)
18666	SweepZone 200 (Full Size)
13518	Quantrex 310 SweepZone 310 & 310R (Full Size)
18202	Quantrex 310 SweepZone 310 & 310R (1/2 Size)
18676	Quantrex 360 SweepZone 360R (Full Size)
17954	Quantrex 650 SweepZone 650R (Full Size)
19620	SweepZone 1000 (Full Size)

Product Code	Stainless Steel Mesch Solid SideDraining Baskets
13329	PC-3
11086	Quantrex 90



up to 32 burs for ultrasonic cleaning and sterilization. Fits inside the large cassette.



IDC Accessory Box
Holds small items to be ultrasonically cleaned and autoclaved—small enough to fit into most ultrasonic cleaning systems.



310 Easy Access Stainless Steel Cassette Rack
Open sided design provides quick & easy access. Minimizes prep & cleaning time. Utilizes entire internal tank space. Fixed extended handles for a sure grip. Multiple cassette configurations possible. Accepts large, small & hygiene-type cassettes. Fits all Quantrex 310 & SweepZone 310 cleaners.

L&R Manufacturing Company
577 Elm Street, PO Box 607
Kearny, NJ 07032-0607 USA
(P) 201-991-5330 (F) 201-991-5870
www.LRultrasonics.com info@LRultrasonics.com