

# PROCART III™

## Mobile Treatment Console



**P**ROCART III MOBILE TREATMENT CONSOLE INTRODUCES THE LATEST ADVANCEMENTS IN SELF-CONTAINED DELIVERY SYSTEMS TECHNOLOGY INTO THE PROGRESSIVE DENTAL PRACTICE. PROCART III JOINS THE GROWING FAMILY OF OTHER SUCCESSFUL SELF-CONTAINED PRODUCTS FROM DNTLWORKS EQUIPMENT CORPORATION.

Where can you implement the new ProCart III? Increase your existing office's productivity by expanding into unused areas. Open a satellite office. For new dentists, open your first practice without costly plumbing and electrical wiring expenses.

ProCart III is perfectly suited for other dental treatment locations. Expand your vision to include hospitals, nursing homes, military installations, correctional facilities, community health centers, and other public health organizations.

### Advanced Features

- Built-in amalgam separator.
- Quick access to internal compartments.
- Rounded corners and edges to make cleaning easier.
- Locking rear casters.
- Superior sound reduction.



# PROCART III™

## Mobile Treatment Console

### Dimensions/Weight/Sound Level

Width:	24" (61 cm)
Depth:	26" (66 cm)
Height:	32" (81cm)
Weight:	112 lbs (50.8 Kg)
Sound Level:	44dB @ 3 feet (91.4 cm) Compressor and vacuum pump on



### Features

- Wet/Dry foot control disc.
- Moisture separator for the compressed air.
- Filtered exhaust air.
- Waste material is filtered twice.

### Electrical System

Voltage:	115 Volts / 60Hz (Model 2400) 220 Volts / 50Hz (Model 2410)
Wattage:	1200 - 1440 Watts

### Compressor System

Horsepower:	1/3 hp Oil-Free
Operating Pressure:	70 - 100 psi (4.8 Bar - 6.9 Bar)
Flow Rate:	2.5scfm @ 40 psi
Air Reservoir:	1 Gallon (3.8 Liters)

### Vacuum System

Horsepower:	1/3 hp oil-free
Operating Pressure:	12" Hg
Flow Rate:	4.6 scfm
Waste Receptacle:	2 qt (1.9 Liter)



### More Features

- Strong, light weight aluminum chassis.
- Two one-liter fresh water reservoirs.
- Three automatic handpiece positions.
- Two automatic evacuation positions.

### Still More Features

- Extra electrical outlets inside and outside.
- External air and water quick connects.
- Convenient side maneuvering handles.
- Easy-to-reach handpiece controls.

### Three separate oil-free pumps are built in.

- Compressor for handpieces and air/water syringe.
- Vacuum pump for HVE and saliva ejector.
- Waste purge pump for hands-free waste removal.



DNTLworks Equipment Corp.  
7300 South Tucson Way  
Centennial, CO 80112 USA

Toll-free: 1-800-847-0694 Phone: +303-693-1410 Fax: +303-693-6189 Email: [info@dntlworks.com](mailto:info@dntlworks.com)