PROCART ITM

SELF-CONTAINED PORTABLE OPERATORY



HE PROCART IMOFFERS POWERFUL PERFORMANCE IN A MODULAR, SELF-CONTAINED OPERATORY, MAKING IT A PREFERRED SYSTEM FOR TODAY'S PORTABLE DENTAL PROFESSIONALS.

Introducing the new generation ProCart I Model 2600 from DNTLworks Equipment corporation.

Suitable for use in offices, nursing homes, hospitals, or wherever on-site dentistry is needed, the ProCart I combines a powerful air compressor and separate vacuum pump for complete "four-handed" capability.

The ProCart I is a truly portable system. It offers exceptional performance and can easily disassemble into protective carrying cases, The air, water, and waste management are all contained inside the unit. Compatible with most common brands of handpeice, the ProCart I will also accommodate fiber optics.

It is also available in simple cart form, model 2500, without the compressor/vacuum module.

ProCart I is a registered trademark of DNTLworks Equipment Corporation © 2010





FEATURES INCLUDE:

- Oil-free 1/3 hp compressor
- Strong 1/3 hp vacuum motor
- Convenient two-piece design for portability
- Dual hose design for H.V.E and saliva ejector use
- Front mounted valve holders and controls for convenience
- Large waste container with automatic overflow shut-off
- 15 amp circuit protection
- 💗 Proudly made in the USA

Portable
Self-Contained
Lightweight
Easy to transport

DNTLworks Equipment Corp. 7300 South Tucson Way
Centennial, CO 80112 USA

Toll-free: 1-800-847-0694

Phone: +303-693-1410 Fax: +303-693-6189

E-mail: info@DNTLworks.com

PROCART ITM

SELF-CONTAINED PORTABLE OPERATORY

SPECIFICATIONS

DIMENSIONS/WEIGHT

Treatment Cart

Base: 23" x 23 (58.4 cm x 58.4 cm) Top: 15.5" x 14.5" (39.4 cm x 36.8 cm)

Height: 32" to 42" (adjustable)

(81.3 cm to 106.7 cm)

Weight: 38 lbs (17.3 kg)

Compressor/Vacuum/Water Module

Size: 11" D x 19" W x 10" H

(27.9 cm x 48.3 cm x 25.4 cm)

Built-in carrying handles

Quickly attaches to Treatment Cart

Weight: 48 lbs (21.8 kg)

ELECTRICAL

Power Supply: 115 vac, 60Hz

220 vac, 50Hz

Circuit Breaker: 15 amps

Hospital Grade Wiring: NEMA Type

5-15P Plug

COMPRESSOR/VACUUM/WATER MODULE

Compressor Motor: 1/3 hp, oil-free

Working Pressure: 56 psi

Flow Rate: 1.2 cfm @ 40 psi

Sound Level: 51 dB

Vacuum Motor: 1/3 hp, oil-free

Vacuum: 15" Hg Vacum Flow: 2.5 cfm

Waste Receptacle: 2 qt (1.9 liters)

Water Reservoir: 16.9 oz (500 ml)

Pressurized System: 40 psi Non-retracting Water System

Handpiece flush system

ProCart I is a registered trademark of DNTLworks Equipment Corporation © 2010





FEATURES INCLUDE

- Wet/Dry foot control
- Automatic controls for one, two and three handpieces
- Accommodates most brands of handpieces and fiber optics
- Non-retracting handpiece water coolant system
- 👸 Automatic handpiece flush system
- Air/water syringe
- F Adjustable high volume suction
- 🥉 Adjustable saliva ejector
- Smooth- design chassis and handpiece holders
- Side mounted controls

DNTLworks Equipment Corp. 7300 South Tucson Way Centennial, CO 80112 USA Toll-free: 1-800-847-0694 Phone: +303-693-1410 Fax: +303-693-6189

E-mail: info@DNTLworks.com

- Easily assembled, adjustable-height stand with lockable casters
- Fowder-coated, scratch resistant finish
- Air-filter moisture separator
- Self-contained water supply with 16.9 oz. (500 ml) reservoir
- * 2 quart (1.9 liters) waste container
- 💗 Repair/Tool kit
- Built-in carrying handles on Compressor Module
- Treatment Cart folds down into compact, easily transportable package
- Fragged carrying cases available

