

ULTRALITE™

Patient Chair with Scissor Base



DNTLworks™
Equipment Corporation

ECONOMY, PORTABILITY, AND BASIC DESIGN WERE THE GOALS DNTLworks EQUIPMENT CORPORATION ACHIEVED WHEN THE ULTRALITE™ PATIENT CHAIR WAS CREATED.

The UltraLite is manufactured with a light-weight frame and reinforced steel hinges. This chair features an adjustable knee break and a backrest that can be manually positioned from 70° upright to 180° supine. The seat height may also be changed for convenience prior to admitting the patient. Set-up and fold-down are achieved quickly and easily.

The UltraLite is ideal for the most rugged environments. They include missionary work, field hospitals, clinics, private home care, nursing homes, military, and many other field applications.

Options include armrest slings for elbow support, and a headrest cushion. A soft-sided carrying case is also available for convenient storage and transportation.

FEATURES INCLUDE:

- ♥ *Light-weight frame*
- ♥ *Comfortable rolled cushioning*
- ♥ *Easy-to-clean, smooth upholstery*
- ♥ *Knee-break*
- ♥ *Easy set-up/fold down*
- ♥ *Manual backrest and seat adjustment*
- ♥ *Optional armrest slings, headrest cushion, and soft-sided carrying case*



SPECIFICATIONS

- ♥ *Weight: 23 lbs (10.5 Kg)*
- ♥ *Folded Dimensions: 18.5" x 33.5" x 6.5"*
(47 cm x 85.1 cm x 16.5 cm)
- ♥ *Seat Adjustment: 17" to 22"*
(43.2 cm to 55.9 cm)
- ♥ *Backrest: 70° upright to 180° supine*

DNTLworks™ Equipment Corp.
15504 East Hinsdale Circle, Unit B
Englewood, CO 80112 U.S.A.
(303) 693-1410
FAX (303) 693-6189
(800) 847-0694

UltraLite is a trademark of
DNTLworks Equipment Corporation © 1998

