

PROBRITE™

Dental Operator & Headband Lights



DNTLworks™
Equipment Corporation

THE PROBRITE™ SERIES OF DENTAL OPERATORY AND HEADBAND LIGHTS COMBINES POWERFUL, LONG-LIFE HALOGEN BULBS IN A VARIETY OF STYLES AND MOUNTING CONFIGURATIONS.

ProBrite lights provide the healthcare professional with the perfect solution for difficult lighting requirements.

The ProBrite series includes direct-beam halogen bulb lights, fiber optic light transmission units and two types of headband lights.

FEATURES INCLUDE:

ProBrite Halogen Lights

- ♥ *Strong, flexible arm remains stationary where required.*
 - ♥ *Bright, even illumination.*
 - ♥ *Unobstructed view of field.*
- #### **Fiber Optic Lights**
- ♥ *Focusing, heat-free lamp.*
 - ♥ *Slimline, adjustable head.*
 - ♥ *Optional rugged, padded canvas carrying cases.*

Headband Lights

- ♥ *Wall plug-in or battery powered models.*
- ♥ *High quality, intense halogen illumination.*
- ♥ *Powerful, lightweight, comfortable, and easy-to-use.*



*Powerful
Lightweight
Packable
Transportable*

PROBRITE™

Dental Operator & Headband Lights



DNTLworks™
Equipment Corporation

SPECIFICATIONS

MOUNTING OPTIONS

-  *Wheeled floorstand*
-  *Chair-mount*
-  *Tripod stand*
-  *Surface mount*
-  *Post mount*

HALOGEN

ProBrite™ Halogen Lights come with a high-quality 40 inch (101.6 cm) flexible gooseneck and 78 inch (198.1 cm) power cord. Lightweight construction, weighs 6 lbs (2.7 Kg).

FIBER OPTIC

ProBrite™ Fiber Optic Lights come with a high-quality 34 inch (86.4 cm) flexible gooseneck and 78 inch (198.1 cm) power cord. The focusing sleeve adjusts the spot size from 4 inches (10.2 cm) to 10 inches (25.4 cm) at a distance of 16 inches (40.6 cm).

HEADBAND

ProBrite™ Headband Lights have an adjustable focusing lens with an articulating mount on a comfortable, soft plastic, adjustable headband.

Headband lights come in battery operated or wall plug-in configurations.

DNTLworks™ Equipment Corp.
15504 East Hinsdale Circle, Unit B
Englewood, CO 80112 U.S.A.
(303) 693-1410
FAX (303) 693-6189
(800) 847-0694

Note: All lights are available in either 110 volt or 220 volt models and include various plug options for international electrical requirements.

