



DNTLworks Equipment Corporation

Specifications

3040.00.00

NitrousVac Portable Vacuum System 120Vac



Amperage	6 Amps
Compressor Type	Oil-Free
Dimensions	22"x14-1/2"x9"
Hertz	60 Hz
Horse Power	1/3 Hp
Shipping Size	26"x16"x19"
Shipping Weight	43 lbs
Sound Level	44 db @ 3 feet
Vacuum Flow	5scfm
Vacuum inHg	7 inhg
Voltage	120Vac
Wattage	348
Weight	35 Lbs