

NEW

 Microbrush®

 **TRU**™
by Microbrush



TWO BENDS ARE BETTER THAN ONE.



NEW SIZE!



Ultrafine Tip Size

Available with a new 0.5-millimeter Ultrafine tip size, along with Regular (2.0mm), Fine (1.5mm) and Superfine (1.0mm) tip sizes for any clinical application.



Double Bend

The unique double-bending tip is engineered to give clinicians greater access to hard-to-reach areas of the mouth.

Longer Handle

The 12-centimeter handle is designed to give clinicians maximum tactile control.



NEW TO
TRU

REQUEST A FREE SAMPLE TODAY
AND READ TRU STORIES AT
MICROBRUSH.COM/TRU.

Advanced adhesive technology designed for greater consistency

MEET



A REVOLUTIONARY NEW MICRO-APPLICATOR

Microbrush® International has once again revolutionized clinical applicators with the introduction of the new TRU™ by Microbrush®.

TRU by Microbrush is the first in Microbrush International's extensive line of disposable applicators to feature a double-bending tip, a longer 12-centimeter handle and an Ultrafine 0.5-millimeter tip size option.

PRODUCT NUMBER	SIZE	COLOR
291159	Regular (2.0mm)	Blue
291160	Regular (2.0mm)	Green
291161	Fine (1.5mm)	Yellow
291162	Fine (1.5mm)	Pink
291163	Superfine (1.0mm)	Gold
291164	Superfine (1.0mm)	Silver
291165	Ultrafine (0.5mm)	Teal
291166	Ultrafine (0.5mm)	Orange

* Available in 80-count tubes with easy dispensing system



TACTILE CONTROL

12-centimeter handle with stay-put bending portion designed to provide maximum control



RELIABLE CONSISTENCY

Advanced adhesive technology designed for greater consistency



ULTIMATE ACCESS

Double-bending tip and new Ultrafine 0.5-millimeter head size designed for greater access

TRU
by Microbrush