


Trim® & Trim® II

Trim® and Trim® II are fast-setting, non-irritating PEMA temporary crown and bridge materials.

ORDERING INFORMATION

PROVISIONAL

0921900 Air Polishing Unit Kit (Red)



0921090 Air Polishing Unit Kit (Black)



New Hue	VITA	Trim	Trim II	Trim Plus	Ultra Trim
VW (Very White)	B0	*			
59	B1		*	*	
60	A1			*	*
61	B2	*	*	*	
62	A2		*	*	*
65	D3	*	*	*	
66	A3		*	*	
67	B3			*	
68			*	*	
69	C2		*	*	
77	B4		*	*	
81	A3.5		*	*	*
87	A4		*	*	
		clear*	incisal*		

*Clear and incisal colors are made available for characterization and for putting an incisal edge on temporaries.

RECOMMENDED FOR USE WITH PROVISIONAL MATERIALS

GLAZE®

Self Cure Resin Polish

GORILLA SEAL®

Ultra-clear light-cure Sealant

MASQUE®

Silicone Separator & Lubricant

PRESSURE POT

Hydraulic Water Press

TECHNICAL TIPS

Using Gorilla Seal® or Glaze® saves polishing time. Gorilla Seal will bond directly to the temporary with a halogen bulb.

Order Trim® VW (Very White) powder in addition to the standard kits to produce a temporary as white as bleached teeth



FEATURES AND BENEFITS

- 1 Trim II offers 6 VITA shades for greater selection
- 2 No distortion when withdrawn from undercuts
- 3 Powders can be mixed & matches together
- 4 Non-adhering to patient's natural dentition
- 5 PEMA formula suited for short-term temporaries

Trim[®] Plus

Trim[®]Plus is a PMMA formula chairside tooth shade resin for use in fabricating temporary crowns, bridges, partials, and splints.

ORDERING INFORMATION

PROVISIONAL MATERIALS

0921930 Standard Kit



FEATURES AND BENEFITS

- 1 Available in 6 VITA shades
- 2 Fast set time offers minimal chair time
- 3 Color stable creating a natural appearance
- 4 PMMA formula suited for long-term temporaries

Ultra Trim[™]

A cartridge dispensed temporary crown and bridge material. Ultra Trim[™] reduces in-office procedure time and provides excellent strength and greater flexibility without brittleness.


ORDERING INFORMATION

PROVISIONAL MATERIALS

0921803 Starter Kit A1 (Light)

0921804 Starter Kit A2 (Universal)

0921805 Starter Kit A3.5 (Dark)



0921809 E-Z Squeeze Universal Dispensing Gun

0921699 10:1 Mixing Tips (15)

0921806 50 ml A1 Cartridge

0921807 50 ml A2 Cartridge

0921808 50 ml A3.5 Cartridge



FEATURES AND BENEFITS

- 1 10:1 base to catalyst ratio
 - 2 Non-irritating formul, free of monomer odors
 - 3 Auto mixed
 - 4 Minimal shrinkage = excellent fit & marginal adaptation
- Chemically bonds to itself for quick repairs or if additional material is needed