

3M Science.
Applied to Life.™


3M™ Impregum™ Super Quick Polyether Impression Material

Capture every
detail in 2 minutes?
Yes, it can.



Quicker set.

With just a two minute setting time, the new Impregum Super Quick Material lives up to its name and allows for faster impressions in cases like single-unit crowns, implants, inlays/onlays and smaller bridges up to 3 units.



100 % Polyether. Trusted precision.

Impregum Super Quick Materials demonstrate excellent performance in moist situations and exhibit great flow behavior, leading to higher levels of accuracy and detailed, high-precision restorations.

The nature of polyether allows you to get the full 45 seconds of working time, avoiding inaccuracies associated with premature setting times common with VPS materials.




Make your patients feel comfortable.

An improved fresh, minty taste together with a short time in the mouth significantly increases your patients' comfort.

3M™ Impregum™ Super Quick Polyether Impression Material

- For monophase and 1-step impression techniques
- Perfect for impressions for single-unit crowns, implants, inlays/onlays and smaller bridges up to 3 units
- Impregum Super Quick Material is scannable with the current generation of scanners for a more efficient workflow.



The fastest way to polyether precision for smaller cases.

	Dispensing system	Impression technique	Single-unit crown	Small implant (≤ 2 units)	Inlay/onlay (≤ 2 units)	Small bridge (≤ 3 units)	Recommended max. working time (min:sec)	Intra-oral setting time at 37°C (min:sec)
3M™ Impregum™ Penta™ Super Quick (Medium Body)		Monophase or 1-step	✓	✓	✓	✓	0:45	2:00
3M™ Impregum™ Penta™ Super Quick (Heavy Body)		1-step	✓	✓	✓	✓	0:45	2:00
3M™ Impregum™ Super Quick (Light Body)		1-step	✓	✓	✓	✓	0:45	2:00

Ordering information



3M™ Impregum™ Penta™ Super Quick Medium Body Introductory Kit (69414)



3M™ Impregum™ Penta™ Super Quick Medium Body/Light Body Introductory Kit (69415)



3M™ Impregum™ Penta™ Super Quick Heavy Body/Light Body Introductory Kit (69407)

Item no. Product information

NEW

Introductory Kits for 3M™ Pentamix™ 3 Automatic Mixing Unit

69414

3M™ Impregum™ Penta™ Super Quick Medium Body

1 base paste MB – 300 ml, 1 catalyst MB – 60 ml, 13M™ Impregum™ Penta™ Super Quick MB cartridge for 3M™ Pentamix™ 3 Automatic Mixing Unit, 10 3M™ Penta™ Mixing Tips – Red, 13M™ Penta™ Elastomer syringe, 13M™ Polyether Tray Adhesive, 3 3M™ Astringent Retraction Paste capsules, 6 3M™ Dual Arch Trays (3 posterior, 3 anterior trays)

69415

3M™ Impregum™ Penta™ Super Quick Medium Body/Light Body

1 base paste MB – 300 ml, 1 catalyst MB – 60 ml, 2 3M™ Impregum™ Super Quick MB cartridges – 50 ml each, 13M™ Impregum™ Penta™ Super Quick MB cartridge for 3M™ Pentamix™ 3 Automatic Mixing Unit, 10 3M™ Penta™ Mixing Tips – Red, 5 3M™ Intra-oral Syringes – Purple, 5 3M™ Garant™ Mixing Tips – Purple, 5 3M™ Garant™ Intra-oral Tips – White, 13M™ Polyether Tray Adhesive, 3 3M™ Astringent Retraction Paste capsules, 6 3M™ Dual Arch Trays (3 posterior, 3 anterior trays)

69407

3M™ Impregum™ Penta™ Super Quick Heavy Body/Light Body

1 base paste HB – 300 ml, 1 catalyst HB – 60 ml, 2 3M™ Impregum™ Super Quick LB cartridges – 50 ml each, 13M™ Impregum™ Penta™ Super Quick HB cartridge for 3M™ Pentamix™ 3 Automatic Mixing Unit, 10 3M™ Penta™ Mixing Tips – Red, 5 3M™ Intra-oral Syringes – Purple, 5 3M™ Garant™ Mixing Tips – Purple, 5 3M™ Garant™ Intra-oral Tips – White, 13M™ Polyether Tray Adhesive, 3 3M™ Astringent Retraction Paste capsules, 6 3M™ Dual Arch Trays (3 posterior, 3 anterior trays)

NEW

Upgrade Kits for 3M™ Pentamix™ 3 Automatic Mixing Unit

69404

3M™ Impregum™ Penta™ Super Quick Medium Body

1 base paste MB – 300 ml, 1 catalyst MB – 60 ml, 10 3M™ Penta™ Mixing Tips – Red

69411

3M™ Impregum™ Penta™ Super Quick Heavy Body/Light Body

1 base paste HB – 300 ml, 1 catalyst HB – 60 ml, 13M™ Impregum™ Super Quick LB cartridge – 50 ml, 10 3M™ Penta™ Mixing Tips – Red, 5 3M™ Intra-oral Syringes – Purple, 5 3M™ Garant™ Mixing Tips – Purple, 5 3M™ Garant™ Intra-oral Tips – White

NEW

Refills

69385

3M™ Impregum™ Penta™ Super Quick Medium Body

2 base pastes – 300 ml each, 2 catalysts – 60 ml each

69379

3M™ Impregum™ Super Quick Light Body

4 cartridges – 50 ml each, 5 3M™ Intra-oral Syringes Purple, 5 3M™ Garant™ Mixing Tips purple

69406

3M™ Impregum™ Penta™ Super Quick Heavy Body

2 base pastes – 300 ml each, 2 catalysts – 60 ml each

NEW

Accessories

69387

3M™ Impregum™ Penta™ Super Quick Medium Body Cartridge

69401

3M™ Impregum™ Penta™ Super Quick Heavy Body Cartridge

71622

48 3M™ Dual Arch Impression Trays Posterior

71623

35 3M™ Dual Arch Impression Trays Anterior

69408

3M™ Polyether Contact Tray Adhesive

1 bottle – 17 ml



3M Oral Care
2510 Conway Avenue
St. Paul, MN 55144-1000
USA
1-800-634-2249

3M Canada
Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-888-363-3685

3M.com/Dental

3M, Garant, Impregum, Penta and Pentamix are trademarks of 3M Company or 3M Deutschland GmbH. Used under license in Canada. All other trademarks are owned by other companies.

© 3M 2018. All rights reserved.