



# Ethicon **Product Catalog**

# Products by Category

At Ethicon, we develop meaningful solutions that enhance surgical outcomes and improve the lives of patients. For more than 80 years, we have listened to the voice of our customers, working in partnership with clinicians, hospitals, payors, and other healthcare stakeholders to gain deep experience and perspective into ways to redefine surgery. Our global scope, along with our broad surgical portfolio, let us provide market-appropriate, value-driven solutions that help patients around the world live longer, healthier lives.

## Access

Ethicon offers you flexibility across surgical specialties, so you can choose the right trocar, with the right level of performance, when and where it is needed. The ENDOPATH® Trocar Portfolio brings value to every procedure.



3

## Surgical Stapling

Our surgical stapling portfolio addresses the challenges of tissue thickness, tissue movement during firing, staple line security, and tissue tension to provide you with greater stapling precision in surgery.



48

## Adjunctive Hemostats

The Ethicon portfolio of Adjunctive Hemostats is an area dedicated to aid in surgery by focusing on an array of technologies to minimize bleeding challenges. Solutions comprise innovative, biologically based products for surgical conditions that are often difficult or expensive to manage.



9

## Uterine & Pelvic Surgery

Ethicon has worked with gynecologists to create innovative new solutions for a greater choice of procedure. Enhanced technique combined with innovative technologies helps to improve surgical performance.



57

## Energy

With the broadest advanced energy portfolio, Ethicon provides you with unmatched choice for your laparoscopic and open surgical needs.



13

## Wound Closure

Ethicon has pioneered advances in skin closure that enable you to minimize tissue trauma, optimize healing, and improve patient outcomes. We offer a portfolio of innovative wound closure solutions designed for the specific needs of skin, fascia, and organ tissues.



61

## Hernia Repair & Fixation

Ethicon has been pioneering advances in hernia surgery to help restore patients back to their pre-hernia life. Our innovative value-led solutions are meant to provide strong and safe hernia repair and are designed for good long-term patient outcomes.



37

## Vascular Access & Infection Control

BIOPATCH® is the ONLY protective disk with a cleared indication and PROVEN to reduce the incidence of catheter-related bloodstream infections (CRBSIs), local infections and skin colonization in patients with central venous and arterial catheters.



63

## Ligation

The Ethicon portfolio of endoscopic and open ligation devices are designed to provide optimal patient surgical experiences. Our deep understanding of the interaction between surgical tools and living tissue has enabled us to design solutions for confident occlusion.



44

## Other Devices

Whatever your needs during the course of surgery, Ethicon is committed to providing you a solution. From laparoscopic hand instruments to gastric banding, drains, and reservoirs, you can count us to provide you with outstanding value, excellent service, and a broad range of product offerings.



70

## Products by Surgical Specialties

Bariatric Surgery

Colorectal Surgery

Gynecologic Surgery

Hernia Repair

Hepato-Biliary Surgery

Thoracic Surgery

# Access

Ethicon offers you flexibility across surgical specialties, so you can choose the right trocar, with the right level of performance, when and where it is needed. The ENDOPATH® Trocar Portfolio brings value to every procedure.

**Learn more about the portfolio on [Ethicon.com](https://www.ethicon.com)**



# Access

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## ENDOPATH XCEL® Bladeless Trocars with OPTIVIEW® Technology 5mm, 12mm



Product Code	Shaft Diameter	Shaft Length	Sleeve	EA/BX
<a href="#">2B5ST</a>	5mm	75mm	Stability	6
<a href="#">2B5LT</a>	5mm	100mm	Stability	6
<a href="#">2B5XT</a>	5mm	150mm	Stability	6
<a href="#">2B12LT</a>	12mm	100mm	Stability	6
<a href="#">2B12XT</a>	12mm	150mm	Stability	6

## ENDOPATH XCEL® Blunt Tip Trocar with OPTIVIEW® Technology 12mm




Product Code	Shaft Diameter	Shaft Length	Sleeve	EA/BX
<a href="#">2H12LP</a>	12mm	100mm	Stability	6

## ENDOPATH XCEL® Dilating Tip Trocar with OPTIVIEW® Technology 5mm, 12mm



Product Code	Shaft Diameter	Shaft Length	Sleeve	EA/BX
<a href="#">2D5ST</a>	5mm	75mm	Stability	6
<a href="#">2D5LT</a>	5mm	100mm	Stability	6
<a href="#">2D12LT</a>	12mm	100mm	Stability	6
<a href="#">2D12XT</a>	12mm	150mm	Stability	6

## ENDOPATH XCEL® Universal Sleeves with OPTIVIEW® Technology 5mm, 12mm



Product Code	Shaft Diameter	Shaft Length	Sleeve	EA/BX
<a href="#">2CB5ST</a>	5mm	75mm	Stability	6
<a href="#">2CB5LT</a>	5mm	100mm	Stability	6
<a href="#">2CB12LT</a>	12mm	100mm	Stability	6



# Access

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## ENDOPATH XCEL® Bladeless Trocars

5mm, 8mm, 11mm, 12mm, 15mm



Product Code	Shaft Diameter	Shaft Length	Sleeve	EA/BX
<a href="#">B5ST</a>	5mm	75mm	Stability	6
<a href="#">B5LT</a>	5mm	100mm	Stability	6
<a href="#">B5XT</a>	5mm	150mm	Stability	6
<a href="#">B8LT</a>	8mm	100mm	Stability	6
<a href="#">B11LT</a>	11mm	100mm	Stability	6
<a href="#">B11LP</a>	11mm	100mm	Smooth	6
<a href="#">B11LTH</a>	11mm	100mm, handled	Stability	6
<a href="#">B11LPH</a>	11mm	100mm, handled	Smooth	6
<a href="#">B12SRT</a>	12mm	75mm	Stability	6
<a href="#">B12LT</a>	12mm	100mm	Stability	6
<a href="#">B12LP</a>	12mm	100mm	Smooth	6
<a href="#">B12LTH</a>	12mm	100mm, handled	Stability	6
<a href="#">B12LPH</a>	12mm	100mm, handled	Smooth	6
<a href="#">B12XT</a>	12mm	150mm	Stability	6
<a href="#">B15LT</a>	15mm	100mm	Stability	6

## ENDOPATH XCEL® Dilating Tip Trocars

5mm, 11mm, 12mm



Product Code	Shaft Diameter	Shaft Length	Sleeve	EA/BX
<a href="#">D5ST</a>	5mm	75mm	Stability	6
<a href="#">D5LT</a>	5mm	100mm	Stability	6
<a href="#">D11LT</a>	11mm	100mm	Stability	6
<a href="#">D12LT</a>	12mm	100mm	Stability	6
<a href="#">D12XT</a>	12mm	150mm	Stability	6

# Access

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ENDOPATH XCEL® Blunt Tip Trocar  
12mm



Product Code	Shaft Diameter	Shaft Length	Sleeve	EA/BX
<a href="#">H12LP</a>	5mm	100mm	Smooth	6

ENDOPATH XCEL® Universal Sleeves  
5mm, 11mm, 12mm



Product Code	Shaft Diameter	Shaft Length	Sleeve	EA/BX
<a href="#">CB5ST</a>	5mm	75mm	Stability	6
<a href="#">CB5LT</a>	5mm	100mm	Stability	6
<a href="#">CB11LT</a>	11mm	100mm	Stability	6
<a href="#">CB12LT</a>	12mm	100mm	Stability	6

ENDOPATH BASX® Bladeless Trocars  
5mm, 11mm, 12mm



Product Code	Shaft Diameter	Shaft Length	Sleeve	EA/BX
<a href="#">TB5ST</a>	5mm	75mm	Stability	6
<a href="#">TB5LT</a>	5mm	100mm	Stability	6
<a href="#">TB11LT</a>	11mm	100mm	Stability	6
<a href="#">TB12ST</a>	12mm	75mm	Stability	6
<a href="#">TB12LT</a>	12mm	150mm	Stability	6

ENDOPATH BASX® Universal Sleeves  
5mm, 11mm, 12mm



Product Code	Shaft Diameter	Shaft Length	Sleeve	EA/BX
<a href="#">CTB5ST</a>	5mm	75mm	Stability	6
<a href="#">CTB5LT</a>	5mm	100mm	Stability	6
<a href="#">CTB11LT</a>	11mm	100mm	Stability	6
<a href="#">CTB12LT</a>	12mm	100mm	Stability	6

ENDOPATH® Mini/Micro  
Trocars Nonoptical  
2mm/3mm



Product Code	Shaft Diameter	Shaft Length	Sleeve	EA/BX
<a href="#">23NBS</a>	2mm/3mm	65mm	Smooth	6
<a href="#">23NBL</a>	2mm/3mm	100mm	Smooth	6

# Access

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## ENDOPATH® Non-Shielded Trocars 5mm



Product Code	Shaft Diameter	Shaft Length	Sleeve	EA/BX
<a href="#">355NS</a>	5mm	65mm	Non-shielded	6
<a href="#">35LNS</a>	5mm	100mm	Non-shielded	6

## ENDOPATH® Rigid Thoracic Trocar 10mm/12mm



Product Code	Shaft Diameter	Shaft Length	Sleeve	EA/BX
<a href="#">TT012</a>	10mm/12mm	75mm	Rigid	6

## FLEXIPATH® Flexible Thoracic Trocars



Product Code	Shaft Diameter	Shaft Length	Sleeve	EA/BX
<a href="#">FP007</a>	7mm	80mm	Flexible	6
<a href="#">FP015</a>	15mm	80mm	Flexible	6
<a href="#">FP020</a>	20mm	80mm	Flexible	6

## FLEXIPATH® Flexible Thoracic Trocar Pack



Product Code	Description	EA/BX
<a href="#">FPK02</a>	Contains: (1) 15mm obturator and (3) 15mm flexible sleeves	6

## Access – HALS/Insufflation Needles

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

### ENDOPATH DEXTRUS® Minimally Invasive Access System



Product Code	Description	Abdominal Wall Thickness	EA/BX
<a href="#">LFRO1</a>	Fixed-length access retractor-small	<4cm	6
<a href="#">FLRO2</a>	Fixed-length access retractor-medium	4-7cm	6
<a href="#">FLRO3</a>	Fixed-length access retractor-larger	>7cm	6
<a href="#">HAPO2</a>	Seal cap assembly with accessories	N/A	3

### ENDOPATH® Insufflation Needles Pneumoperitoneum



Product Code	Shaft Length	EA/BX
<a href="#">PN120</a>	120mm	12
<a href="#">PN150</a>	150mm	12

### ENDOPATH® Insufflation Needle Ultra Veress



Product Code	Diameter	EA/BX
<a href="#">UV120</a>	120mm	12

# Adjunctive Hemostat

The Ethicon portfolio of Adjunctive Hemostats is an area dedicated to aid in surgery by focusing on an array of technologies to minimize bleeding challenges. Solutions comprise innovative, biologically based products for surgical conditions that are often difficult or expensive to manage. Our Hemostasis Optimization Program is a systematic approach that enables improvement of adjunctive hemostat utilization. This program looks at solutions for surgical bleeding situations, and potentially reduces procedural costs and variability in product utilization.

**Learn more about the portfolio on [Ethicon.com](https://www.ethicon.com)**



# Adjunctive Hemostat

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## SURGICEL SNoW® Hemostat



Product Code	Description	EA/BX
<a href="#">2081</a>	1in x 2in	10
<a href="#">2082</a>	2in x 4in	10
<a href="#">2083</a>	4in x 4in	10

## SURGICEL® FIBRILLAR™ Hemostat



Product Code	Description	EA/BX
<a href="#">1961</a>	1in x 2in	10
<a href="#">1962</a>	2in x 4in	10
<a href="#">1963</a>	4in x 4in	10

## SURGICEL® NU-KNIT™ Hemostat



Product Code	Description	EA/BX
<a href="#">1943</a>	3in x 4in	24
<a href="#">1946</a>	6in x 9in	10

## SURGICEL® Original Hemostat



Product Code	Description	EA/BX
<a href="#">1951</a>	2in x 14in	24
<a href="#">1952</a>	4in x 8in	24
<a href="#">1953</a>	2in x 3in	24
<a href="#">1955</a>	0.5in x 2in	24

## SURGIFOAM® Absorbable Gelatin Powder



Product Code	Description	EA/BX
<a href="#">1978</a>	1.0g	6

# Adjunctive Hemostat

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

SURGIFOAM® Absorbable Gelatin  
Powder Kit



Product Code	Description	EA/BX
<a href="#">1979</a>	1.0g	6

SURGIFOAM® Absorbable Gelatin  
Sponge



Product Code	Description	EA/BX
<a href="#">4972</a>	2cm x 6cm x 7mm (12-7)	12
<a href="#">1973</a>	8cm x 6.25cm x 10mm (50)	4
<a href="#">1974</a>	8cm x 12.5cm x 10mm (100)	6
<a href="#">1975</a>	8cm x 12.5cm x 2mm (100c)	6

SURGIFOAM® Absorbable Gelatin  
Sponge Hemorrhoidectomy Sponge



Product Code	Description	EA/BX
<a href="#">1977</a>	8in x 3in	5

INSTAT® MCH Hemostat



Product Code	Description	EA/BX
<a href="#">1985</a>	1.0g	12

# Adjunctive Hemostat

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## EVARREST® Fibrin Sealant Patch



Product Code	Description	EA/BX
<a href="#">EVT5024</a>	EVARREST® Fibrin Sealant Patch 2"x4" (5.1cm x 10.2cm)	2

## SURGIFLO® Hemostatic Matrix



Product Code	Description	EA/BX
<a href="#">2993</a>	8mL SURGIFLO® Hemostatic Matrix Kit with thrombin	6
<a href="#">2994</a>	SURGIFLO® Hemostatic Matrix Kit with thrombin (fully sterile)	6
<a href="#">2991</a>	8mL SURGIFLO® Hemostatic Matrix	6
<a href="#">MS1995</a>	34cm SURGIFLO® Endoscopic Applicator	6

## EVICEL® Fibrin Sealant (Human)



Product Code	Description	EA/BX
<a href="#">3905</a>	EVICEL 5ML KIT	1
<a href="#">3902</a>	EVICEL 2ML KIT	1
<a href="#">3921S</a>	EVICEL AIRLESS PRAY ACC US EA	1
<a href="#">3920</a>	EVICEL 4CM CONTROL TIP (SET OF 10)	1
<a href="#">3909</a>	EVICEL 45CM FLEXIBLE TIP (SET OF 3)	1
<a href="#">3908</a>	EVICEL 35 CM APPLICATOR TIP	1
<a href="#">3905ASA</a>	10 EACH ES OF 3905 AND 1 BOX OF 39215	1
<a href="#">3902ASA</a>	10 EACH ES OF 3902 AND 1 BOX 39215	1
<a href="#">3901ASA</a>	10 EACH ES OF 3901 AND 1 BOX OF 39215	1
<a href="#">3901</a>	EVICEL 1ML KIT	1



# Energy

Ethicon provides you with unmatched choice for your Energy needs in laparoscopic and open surgery. HARMONIC® Ultrasonic devices offer precision, multi-functionality and strong vessel sealing. ENSEAL® Advanced Bipolar devices bring confidence and consistency to vessel sealing and deliver versatility in every device. MEGADYNE™ delivers an expansive suite of innovative Core Energy products.

**Learn more about the portfolio on [Ethicon.com](https://www.ethicon.com)**



# Advanced Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## HARMONIC® HD 1000i Shears



Product Code	Shaft Diameter	Shaft Length	Tip	EA/BX
<a href="#">HARHD20</a>	5mm	20cm	Curved	6
<a href="#">HARHD36</a>	5mm	36cm	Curved	6

## HARMONIC ACE®+7 Shears



Product Code	Shaft Diameter	Shaft Length	Tip	EA/BX
<a href="#">HARH23</a>	5mm	23cm	Curved	6
<a href="#">HARH36</a>	5mm	36cm	Curved	6
<a href="#">HARH45</a>	5mm	45cm	Curved	6

## HARMONIC ACE®+ Shears



Product Code	Shaft Diameter	Shaft Length	Tip	EA/BX
<a href="#">HAR23</a>	5mm	23cm	Curved	6
<a href="#">HAR36</a>	5mm	36cm	Curved	6

## HARMONIC FOCUS®+ Shears



Product Code	Length	Tip	EA/BX
<a href="#">HAR9F</a>	9cm	Curved	6

## HARMONIC FOCUS®+ Long Shears



Product Code	Length	Tip	EA/BX
<a href="#">HAR17F</a>	17cm	Curved	6

# Advanced Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## HARMONIC SYNERGY® Blades



Product Code	Shaft Length	Tip	EA/BX
<a href="#">SNGCB</a>	Adjustable 4-9cm	Curved	6
<a href="#">SNGHK</a>	Adjustable 4-9cm	Hook	6
<a href="#">SNGHK2</a>	Adjustable 4-9cm	Combination Hook	6

## HARMONIC® Blades



Product Code	Shaft Diameter	Shaft Length	Tip	EA/BX
<a href="#">HC325</a>	5mm	23cm	Curved	6
<a href="#">HDH05</a>	5mm	32cm	Hook	6
<a href="#">HBC05</a>	5mm	31cm	Ball	6

## HARMONIC® Hand Piece



Product Code	Compatible with	EA/BX
<a href="#">HP054</a>	HAR23, HAR36, HARH23, HARH36, HARH45, HBC05, HC325, HDH05	1

## HARMONIC® BLUE Hand Piece



Product Code	Compatible with	EA/BX
<a href="#">HPBLUE</a>	HAR9F, HAR17F, SNGCB, SNGHK, SNGHK2	1

# Advanced Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## HARMONIC® Accessories

Product Code	Description	EA/BX
<a href="#">HAOSSL</a>	HARMONIC ACE® 5mm open sleeve	6
<a href="#">HP54CTA</a>	HARMONIC® Hand Piece cleaning tool	1
<a href="#">HSA08</a>	Disposable hand switch adapter (only use with HPO54)	3
<a href="#">HSTOI</a>	Sterilization tray	1
<a href="#">HSTO2</a>	Test tip for HPO54	1
<a href="#">TLBOI</a>	Torque lock blade wrench	2
<a href="#">TTBLUE</a>	Test tip for HPBLUE	1

## ENSEAL® X1 Large Jaw



Product Code	Shaft Diameter	Shaft Length	Jaw	EA/BX
<a href="#">NSLX120L</a>	13mm	20cm	Curved	6

## ENSEAL® G2 Articulating Tissue Sealers



Product Code	Shaft Diameter	Shaft Length	Jaw	EA/BX
<a href="#">NSLG2C35A</a>	5mm	35cm	Curved	6
<a href="#">NSLG2C45A</a>	5mm	45cm	Curved	6
<a href="#">NSLG2S35A</a>	5mm	35cm	Straight	6
<a href="#">NSLG2S45A</a>	5mm	45cm	Straight	6

## ENSEAL® G2 Curved Tissue Sealers



Product Code	Shaft Diameter	Shaft Length	Jaw	EA/BX
<a href="#">NSLG2C14</a>	5mm	14cm	Curved	6
<a href="#">NSLG2C25</a>	5mm	25cm	Curved	6
<a href="#">NSLG2C35</a>	5mm	35cm	Curved	6
<a href="#">NSLG2C45</a>	5mm	45cm	Curved	6

# Advanced Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## ENSEAL® G2 Straight Tissue Sealers



Product Code	Shaft Diameter	Shaft Length	Jaw	EA/BX
<a href="#">NSLG2S14</a>	5mm	14cm	Straight	6
<a href="#">NSLG2S25</a>	5mm	25cm	Straight	6
<a href="#">NSLG2S35</a>	5mm	35cm	Straight	6
<a href="#">NSLG2S45</a>	5mm	45cm	Straight	6

## ENSEAL® Trio Tissue Sealers



Product Code	Shaft Diameter	Shaft Length	Jaw	EA/BX
<a href="#">ETRIO314H</a>	5mm	14cm	Curved	6
<a href="#">ETRIO325H</a>	5mm	25cm	Curved	6
<a href="#">ETRIO335H</a>	5mm	35cm	Curved	6
<a href="#">ETRIO345H</a>	5mm	45cm	Curved	6

## ENSEAL® Round Tip Tissue Sealers



Product Code	Shaft Diameter	Shaft Length	Jaw	EA/BX
<a href="#">NSEAL514RH</a>	5mm	14cm	Straight	6
<a href="#">NSEAL525RH</a>	5mm	25cm	Straight	6
<a href="#">NSEAL535RH</a>	5mm	35cm	Straight	6
<a href="#">NSEAL545RH</a>	5mm	45cm	Straight	6

## ETHICON GEN11 Generator



Product Code	Compatible with	EA/BX
<a href="#">GEN11</a>	All HARMONIC® and ENSEAL® devices	1

# Advanced Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## ETHICON GEN11 Generator Accessories



Product Code	Description	Compatible with	EA/BX
<a href="#">HGA11</a>	HARMONIC® connector	HPO54 and HPBLUE hand pieces	1
<a href="#">EGA11</a>	ENSEAL® connector	ENSEAL® Trio Tissue Sealer and ENSEAL® Round Tip Tissue Sealer	1
<a href="#">GEN11VK</a>	Output verification key	GEN11	1
<a href="#">FSW11</a>	Foot switch and cable	GEN11	1
<a href="#">CRT11</a>	Generator cart	GEN11	1
<a href="#">GEN11WAR1</a>	1-Year extended warranty		1
<a href="#">GEN11WAR2</a>	2-Year extended warranty		1

# Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## MEGADYNE™ MEGA POWER™ Generator



Product Code	Compatible with	EA/BX
<a href="#">1000</a>	Monopolar and Bipolar devices	1

## MEGADYNE™ MEGA POWER™ Generator accessories



Product Code	Product Description	EA/BX
<a href="#">1300SJ</a>	Mega Cart with Accessory Shelf	1
<a href="#">1450J</a>	Mega Power Bipolar Foot Switch	1
<a href="#">1400JJ</a>	Mega Power Monopolar Foot Switch	1
<a href="#">96007</a>	Monopolar Adapter	1

## MEGADYNE™ ACE BLADE™ 700 Soft Tissue Dissector



Product Code	Description	Tip Length	EA/BX
<a href="#">ME7251C</a>	Standard Blade C Connector	2.5"	6
<a href="#">ME7251E</a>	Standard Blade E Connector	2.5"	6
<a href="#">ME725M1C</a>	Modified Blade C Connector	2.5"	6
<a href="#">ME725M1E</a>	Modified Blade E Connector	2.5"	6

# Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## MEGADYNE™ E-Z CLEAN™ Coated Blade Electrodes



Product Code	Description	Tip Length	EA/BX
<a href="#">0012</a>	Standard Blade Electrode	2.5"	12
<a href="#">0012M</a>	Standard Blade Electrode Modified	2.5"	12
<a href="#">0012A</a>	X-Long Standard Blade Electrode	2.75"	12
<a href="#">0012AM</a>	X-Long Standard Blade Electrode Modified	2.75"	12
<a href="#">0014A</a>	Extended Blade Electrode	4"	12
<a href="#">0014AM</a>	Extended Blade Electrode Modified	4"	12
<a href="#">0014</a>	Extended Blade Electrode	6.5"	12
<a href="#">0014M</a>	Extended Blade Electrode Modified	6.5"	12

## MEGADYNE™ E-Z CLEAN™ Coated Needle Electrodes



Product Code	Description	Tip Length	EA/BX
<a href="#">0013</a>	Needle Electrode	2.75"	12
<a href="#">0013M</a>	Needle Electrode Modified	2.75"	12
<a href="#">0016A</a>	Extended Needle Electrode	4"	12
<a href="#">0016AM</a>	Needle Electrode Modified	4"	12
<a href="#">0016</a>	Extended Needle Electrode	6"	12
<a href="#">0016M</a>	Extended Needle Electrode Modified	6"	12

## MEGADYNE™ E-Z CLEAN™ Coated Blunt Needle Electrodes



Product Code	Description	Tip Length	EA/BX
<a href="#">0113</a>	Blunt Needle Electrode	2.75"	12
<a href="#">0113M</a>	Blunt Needle Electrode Modified	2.75"	12
<a href="#">0113A</a>	Blunt Needle Electrode	4"	12

## MEGADYNE™ E-Z CLEAN™ Coated Ball Electrodes




Product Code	Description	Tip Length	EA/BX
<a href="#">0015</a>	Ball Electrode	2"	12
<a href="#">0009</a>	Ball Electrode	5"	12



# Energy


Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## MEGADYNE™ E-Z CLEAN™ Coated Specialty Electrodes



Product Code	Description	Tip Length	EA/BX
<a href="#">0028</a>	Bayonet Needle	5.8"	12
<a href="#">0028M</a>	Bayonet Needle Modified	5.8"	12
<a href="#">0029M</a>	Bayonet Blade Modified	6.3"	12

## MEGADYNE™ MEGA FINE™ Needle Electrodes



Product Code	Description	Tip Length	EA/BX
<a href="#">0118</a>	Needle Electrode	2"	12
<a href="#">0118A</a>	Extended Needle Electrode	2.5"	12
<a href="#">0119</a>	Angled 45 Degree	2"	12
<a href="#">0119A</a>	Angled 45 Degree 3mm	2"	12
<a href="#">0120</a>	Angled 90 Degree	2"	12
<a href="#">0121</a>	Needle Electrode	6.5"	12

## MEGADYNE™ Precision Electrodes



Product Code	Description	Tip Length	EA/BX
<a href="#">0012AP</a>	Blade Electrode	2.75"	12
<a href="#">0012AMP</a>	Blade Electrode Modified	2.75"	12
<a href="#">0014AP</a>	Blade Electrode	4"	12
<a href="#">0014AMP</a>	Blade Electrode Modified	4"	12

# Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## MEGADYNE™ E-Z CLEAN™ Protective Nose Cone Electrodes



Product Code	Description	Tip Length	EA/BX
<a href="#">0012MD</a>	Nose Cone Blade	2.5"	24
<a href="#">0012AMD</a>	Nose Cone Blade	2.75"	24
<a href="#">0014AMD</a>	Nose Cone Blade	4"	24
<a href="#">0014MD</a>	Nose Cone Blade	6.5"	24
<a href="#">0013MD</a>	Nose Cone Needle	2.75"	24

## MEGADYNE™ E-Z CLEAN™ Laparoscopic Electrodes



Product Code	Tip Configuration	Compatible With	Length	EA/BX
<a href="#">0017</a>	Standard	Standard Pencils	33cm	6
<a href="#">0018</a>	Spatula	Standard Pencils	33cm	6
<a href="#">0018C</a>	Curved Spatula	Standard Pencils	33cm	6
<a href="#">0018CS</a>	Curved Spatula	Reusable Foot Control Cords	33cm	6
<a href="#">0019</a>	Curved	Standard Pencils	33cm	6
<a href="#">0019S</a>	Curved	Reusable Foot Control Cords	33cm	6
<a href="#">0019L</a>	Curved	Standard Pencils	45cm	6
<a href="#">0019LS</a>	Curved	Reusable Foot Control Cords	45cm	6
<a href="#">0020</a>	L-Hook	Standard Pencils	33cm	6
<a href="#">0020S</a>	L-Hook	Reusable Foot Control Cords	33cm	6
<a href="#">0020L</a>	L-Hook	Standard Pencils	45cm	6
<a href="#">0020LS</a>	L-Hook	Reusable Foot Control Cords	45cm	6
<a href="#">0021</a>	J-Hook	Standard Pencils	33cm	6
<a href="#">0021S</a>	J-Hook	Reusable Foot Control Cords	33cm	6
<a href="#">0021L</a>	J-Hook	Standard Pencils	45cm	6
<a href="#">0021LS</a>	J-Hook	Reusable Foot Control Cords	45cm	6
<a href="#">0022</a>	Needle	Standard Pencils	33cm	6
<a href="#">0022S</a>	Needle	Reusable Foot Control Cords	33cm	6
<a href="#">0024</a>	Ball	Standard Pencils	33cm	6
<a href="#">0026S</a>	O-Wire	Reusable Foot Control Cords	33cm	6

# Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## MEGADYNE™ E-Z CLEAN™ Laparoscopic Electrodes



Product Code	Tip Configuration	Compatible With	Length	EA/BX
<a href="#">0100</a>	L-Wire	Standard Pencils	33cm	6
<a href="#">0100S</a>	L-Wire	Reusable Foot Control Cords	33cm	6
<a href="#">0100L</a>	L-Wire	Standard Pencils	45cm	6
<a href="#">0100LS</a>	L-Wire	Reusable Foot Control Cords	45cm	6
<a href="#">0105S</a>	J-Wire	Reusable Foot Control Cords	35cm	6

## MEGADYNE™ Megatip Electrodes



Product Code	Tip Configuration	EA/BX
<a href="#">0600</a>	L-Wire	6
<a href="#">0600M</a>	L-Wire Modified	6
<a href="#">0605</a>	J-Wire	6
<a href="#">0618</a>	Spatula	6
<a href="#">0619</a>	Curved Blade	6
<a href="#">0620</a>	L-Hook	6
<a href="#">0620M</a>	L-Hook Modified	6
<a href="#">0621</a>	J-Hook	6
<a href="#">0621M</a>	J-Hook Modified	6

## MEGADYNE™ Reusable Indicator Shafts



Product Code	Description	Tip Length	EA/BX
<a href="#">0690</a>	Shaft	32cm	2
<a href="#">0690S</a>	Shaft with Split Stem	32cm	2
<a href="#">0695</a>	Shaft	38cm	2
<a href="#">0695S</a>	Shaft with Split Stem	38cm	2

## MEGADYNE™ Reusable Foot Control Cord



Product Code	Description	Length	EA/BX
<a href="#">0075</a>	Split stem	3m	1

# Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## MEGADYNE™ Rocker Switch Pencil





Product Code	Description	Cord Length	EA/BX
<a href="#">0030</a>	Pencil with E-Z CLEAN™ Electrode	10ft	50
<a href="#">0030L</a>	Pencil with E-Z CLEAN™ Electrode	15ft	25
<a href="#">0030H</a>	Pencil with E-Z CLEAN™ Electrode and holster	10ft	50
<a href="#">0030HL</a>	Pencil with E-Z CLEAN™ Electrode and holster	15ft	25
<a href="#">0036</a>	Pencil with E-Z CLEAN™ Electrode Modified	15ft	50
<a href="#">0036H</a>	Pencil with E-Z CLEAN™ Electrode Modified and holster	15ft	50

## MEGADYNE™ Button Switch Pencil



Product Code	Description	Cord Length	EA/BX
<a href="#">0035</a>	Pencil with E-Z CLEAN™ Electrode	10ft	50
<a href="#">0035L</a>	Pencil with E-Z CLEAN™ Electrode	15ft	25
<a href="#">0035H</a>	Pencil with E-Z CLEAN™ Electrode and holster	10ft	50
<a href="#">0035HL</a>	Pencil with E-Z CLEAN™ Electrode and holster	15ft	25
<a href="#">0037</a>	Pencil with E-Z CLEAN™ Electrode Modified	15ft	50
<a href="#">0037H</a>	Pencil with E-Z CLEAN™ Electrode Modified and holster	15ft	50

## MEGADYNE™ E-Z Pen Reusable Electrosurgical Pencil



Product Code	Description	Cord Length	EA/BX
<a href="#">3030H</a>	Rocker switch pencil with holster	10ft	10
<a href="#">3030HL</a>	Rocker switch pencil with holster	15ft	10
<a href="#">3035H</a>	Button switch pencil with holster	10ft	10
<a href="#">3035HL</a>	Button switch pencil with holster	15ft	10

# Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

MEGADYNE™ MEGA SOFT™ Universal  
Reusable Patient Return Electrode  
91cm x 51cm x 0.3cm



Product Code	Description	Electrodes per box	Cables per box
<a href="#">0845M2K01</a>	MEGA SOFT™ Universal with one M2K01 cable: Standard connector, 2.4m length	1	1
<a href="#">0845M2K02</a>	MEGA SOFT™ Universal with one M2K02 cable: Standard connector, 4.4m length	1	1
<a href="#">0845M2K03</a>	MEGA SOFT™ Universal with one M2K03 cable: Phone plug connector, 2.4m length	1	1
<a href="#">0845M2K04</a>	MEGA SOFT™ Universal with one M2K04 cable: Phone plug connector, 4.4m length	1	1
<a href="#">0845M2K05</a>	MEGA SOFT™ Universal with one M2K05 cable: Extended phone plug connector, 4.4m length	1	1
<a href="#">0845M2K06</a>	MEGA SOFT™ Universal with one M2K06 cable: Argon beam connector, 2.4m length	1	1
<a href="#">0845M2K08</a>	MEGA SOFT™ Universal with one M2K08 cable: Compatibility cable, 2.4m length	1	1
<a href="#">0845M2K09</a>	MEGA SOFT™ Universal with one M2K09 cable: Compatibility cable, 4.4m length	1	1

# Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

MEGADYNE™ MEGA SOFT™ Universal Dual Reusable Patient Return Electrode  
91cm x 50cm x 0.3175cm



Product Code	Description	Electrodes per box	Cables per box
<a href="#">0846M2K01</a>	MEGA SOFT™ Universal Dual with two M2K01 cables: Standard connector, 2.4m length	1	2
<a href="#">0846M2K02</a>	MEGA SOFT™ Universal Dual with two M2K02 cables: Standard connector, 4.4m length	1	2
<a href="#">0846M2K03</a>	MEGA SOFT™ Universal Dual with two M2K03 cables: Phone plug connector, 2.4m length	1	2
<a href="#">0846M2K04</a>	MEGA SOFT™ Universal Dual with two M2K04 cables: Phone plug connector, 4.4m length	1	2
<a href="#">0846M2K05</a>	MEGA SOFT™ Universal Dual with two M2K05 cables: Extended phone plug connector, 4.4m length	1	2
<a href="#">0846M2K06</a>	MEGA SOFT™ Universal Dual with two M2K06 cables: Argon beam connector, 2.4m length	1	2
<a href="#">0846M2K08</a>	MEGA SOFT™ Universal Dual with two M2K08 cables: Compatibility cable, 2.4m length	1	2
<a href="#">0846M2K09</a>	MEGA SOFT™ Universal Dual with two M2K09 cables: Compatibility cable, 4.4m length	1	2

# Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

MEGADYNE™ MEGA SOFT™ Reusable Patient Return Electrode  
117cm x 51cm x 1.25cm



Product Code	Description	Electrodes per box	Cables per box
<a href="#">0830M2K01</a>	MEGA SOFT™ with one M2K01 cable: Standard connector, 2.4m length	1	1
<a href="#">0830M2K02</a>	MEGA SOFT™ with one M2K02 cable: Standard connector, 4.4m length	1	1
<a href="#">0830M2K03</a>	MEGA SOFT™ with one M2K03 cable: Phone plug connector, 2.4m length	1	1
<a href="#">0830M2K04</a>	MEGA SOFT™ with one M2K04 cable: Phone plug connector, 4.4m length	1	1
<a href="#">0830M2K05</a>	MEGA SOFT™ with one M2K05 cable: Extended phone plug connector, 4.4m length	1	1
<a href="#">0830M2K06</a>	MEGA SOFT™ with one M2K06 cable: Argon beam connector, 2.4m length	1	1
<a href="#">0830M2K08</a>	MEGA SOFT™ with one M2K08 cable: Compatibility cable, 2.4m length	1	1
<a href="#">0830M2K09</a>	MEGA SOFT™ with one M2K09 cable: Compatibility cable, 4.4m length	1	1

MEGADYNE™ MEGA SOFT™ Dual Reusable Patient Return Electrode  
117cm x 51cm x 1.25cm



Product Code	Description	Electrodes per box	Cables per box
<a href="#">0835M2K01</a>	MEGA SOFT™ Dual with two M2K01 cables: Standard connector, 2.4m length	1	2
<a href="#">0835M2K02</a>	MEGA SOFT™ Dual with two M2K02 cables: Standard connector, 4.4m length	1	2
<a href="#">0835M2K03</a>	MEGA SOFT™ Dual with two M2K03 cables: Phone plug connector, 2.4m length	1	2
<a href="#">0835M2K04</a>	MEGA SOFT™ Dual with two M2K04 cables: Phone plug connector, 4.4m length	1	2
<a href="#">0835M2K05</a>	MEGA SOFT™ Dual with two M2K05 cables: Extended phone plug connector, 4.4m length	1	2
<a href="#">0835M2K06</a>	MEGA SOFT™ Dual with two M2K06 cables: Argon beam connector, 2.4m length	1	2
<a href="#">0835M2K08</a>	MEGA SOFT™ Dual with two M2K08 cables: Compatibility cable, 2.4m length	1	2
<a href="#">0835M2K09</a>	MEGA SOFT™ Dual with two M2K09 cables: Compatibility cable, 4.4m length	1	2

# Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

MEGADYNE™ MEGA SOFT™ Pediatric  
Reusable Patient Return Electrode  
66cm x 30.5cm x 1.3cm



Product Code	Description	Electrodes per box	Cables per box
<a href="#">0840M2K01</a>	MEGA SOFT™ Pediatric with one M2K01 cable: Standard connector, 2.4m length	1	1
<a href="#">0840M2K02</a>	MEGA SOFT™ Pediatric with one M2K02 cable: Standard connector, 4.4m length	1	1
<a href="#">0840M2K03</a>	MEGA SOFT™ Pediatric with one M2K03 cable: Phone plug connector, 2.4m length	1	1
<a href="#">0840M2K04</a>	MEGA SOFT™ Pediatric with one M2K04 cable: Phone plug connector, 4.4m length	1	1
<a href="#">0840M2K05</a>	MEGA SOFT™ Pediatric with one M2K05 cable: Extended phone plug connector, 4.4m length	1	1
<a href="#">0840M2K06</a>	MEGA SOFT™ Pediatric with one M2K06 cable: Argon beam connector, 2.4m length	1	1
<a href="#">0840M2K08</a>	MEGA SOFT™ Pediatric with one M2K08 cable: Compatibility cable, 2.4m length	1	1
<a href="#">0840M2K09</a>	MEGA SOFT™ Pediatric with one M2K09 cable: Compatibility cable, 4.4m length	1	1

MEGADYNE™ MEGA SOFT™ 2000  
Reusable Patient Return Electrode  
91cm x 51cm

Product Code	Description	Electrodes per box	Cables per box
<a href="#">0800M2K01</a>	MEGA SOFT™ 2000 with one M2K01 cable: Standard connector, 2.4m length	1	1
<a href="#">0800M2K02</a>	MEGA SOFT™ 2000 with one M2K02 cable: Standard connector, 4.4m length	1	1
<a href="#">0800M2K03</a>	MEGA SOFT™ 2000 with one M2K03 cable: Phone plug connector, 2.4m length	1	1
<a href="#">0800M2K04</a>	MEGA SOFT™ 2000 with one M2K04 cable: Phone plug connector, 4.4m length	1	1
<a href="#">0800M2K05</a>	MEGA SOFT™ 2000 with one M2K05 cable: Extended phone plug connector, 4.4m length	1	1
<a href="#">0800M2K06</a>	MEGA SOFT™ 2000 with one M2K06 cable: Argon beam connector, 2.4m length	1	1
<a href="#">0800M2K08</a>	MEGA SOFT™ 2000 with one M2K08 cable: Compatibility cable, 2.4m length	1	1
<a href="#">0800M2K09</a>	MEGA SOFT™ 2000 with one M2K09 cable: Compatibility cable, 4.4m length	1	1




# Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## MEGADYNE™ MEGA SOFT™ Reusable Patient Return Electrode

Product Code	Description	EA/BX
<a href="#">0825</a>	MEGA 2000™ Sheath	30

## MEGADYNE™ MEGA SOFT™ Cables and Accessories



Product Code	Description	Cord Length	EA/BX
<a href="#">M2KO1</a>	Standard Connector	8ft	1
<a href="#">M2KO2</a>	Standard Connector	14.4ft	1
<a href="#">M2KO3</a>	Phone Plug Connector, same as B-205 adapter	8ft	1
<a href="#">M2KO4</a>	Phone Plug Connector, same as B-205 adapter	14.4ft	1
<a href="#">M2KO5</a>	Extended Phone Plug Connector, same as B-210 adapter	14.4ft	1
<a href="#">M2KO6</a>	Argon Beam Connector, same as A-238 adapter	8ft	1
<a href="#">M2KO7</a>	Replacement Pigtail Cable	2ft	1
<a href="#">M2KO8</a>	Compatibility Cable	8ft	1
<a href="#">M2KO9</a>	Compatibility Cable	14.4ft	1
<a href="#">PKIT001</a>	MEGA SOFT™ Repair Kit		1

## MEGADYNE™ Mega Vac Smoke Evacuators and Accessories



Product Code	Description	EA/BX
<a href="#">2200J</a>	MEGADYNE™ Mega Vac PLUS Smoke Evacuator with Laparoscopic Mode	1
<a href="#">2100J</a>	MEGADYNE™ Mega Vac Smoke Evacuator	1
<a href="#">2211J</a>	ULPA Filter with Fluid Trap	10
<a href="#">2210J</a>	ULPA Filter with Fluid Trap	9
<a href="#">2220J</a>	Charcoal Filter	1
<a href="#">2250J</a>	MEGADYNE™ RF Sensor for Mega Vac	1
<a href="#">2160J</a>	MEGADYNE™ Mega Vac Smoke Evacuator Footswitch	1

# Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## MEGADYNE™ MINI VAC™ Smoke Evacuators and Accessories



Product Code	Description	EA/BX
<a href="#">ECVV120</a>	MINI VAC™ Smoke Evacuator 120v	1
<a href="#">ECVV220</a>	MINI VAC™ Smoke Evacuator 220v	1
<a href="#">MGEZLINK01</a>	MEGADYNE™ RF Sensor for MINI VAC™	
<a href="#">MGVS353</a>	MEGADYNE™ Micro Safe Filter 1 Pack	1
<a href="#">MGVS35302</a>	MEGADYNE™ Micro Safe Filter 2 Pack	2

<a href="#">MGVSFT10</a>	Fluid Trap for MINI VAC™	10
--------------------------	--------------------------	----

# Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## Smoke Evacuation Pencils and Accessories



Product Code	Description	EA/BX
<a href="#">252510</a>	MEGADYNE™ ZIP-PEN™ Smoke Evacuation Pencil	20
<a href="#">252510EC</a>	MEGADYNE™ ZIP-PEN™ Smoke Evacuation Pencil with 22mm connector	20
<a href="#">2540J</a>	MEGADYNE™ ZIP-PEN™ Extension Nozzle 2.7"	10
<a href="#">2560J</a>	MEGADYNE™ ZIP-PEN™ Extension Nozzle 5.2"	10
<a href="#">211010J</a>	MEGADYNE™ ULTRA VAC™ Smoke Evacuation Pencil	20
<a href="#">211010EC</a>	MEGADYNE™ ULTRA VAC™ Smoke Evacuation Pencil with 22mm connector	20
<a href="#">212009J</a>	MEGADYNE™ Attacha Vac Smoke Evacuation Shroud	25
<a href="#">212009EC</a>	MEGADYNE™ Attacha Vac Smoke Evacuation Shroud with 22mm connector	25
<a href="#">2190J</a>	Laparoscopic Tubing	25
<a href="#">2195J</a>	Speculum Tubing	25
<a href="#">2140J</a>	MEGADYNE™ Mega Vac Smoke Evacuator Connector	1
<a href="#">2145J</a>	22mm male to 10mm male Connector	1
<a href="#">2150J</a>	Universal Connector	1
<a href="#">2151J</a>	22mm ID Tube Connector	1
<a href="#">2155</a>	22mm female to 10mm male Connector	1

# Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## MEGADYNE™ Stainless Steel Blade Electrodes



Product Code	Description	Tip Length	EA/BX
<a href="#">0312</a>	Blade Electrode	2.5"	12
<a href="#">0312M</a>	Blade Electrode Modified	2.5"	12
<a href="#">0312A</a>	Blade Electrode	2.75"	12
<a href="#">0312AM</a>	Blade Electrode Modified	2.75"	12
<a href="#">0314A</a>	Blade Electrode	4"	12
<a href="#">0314AM</a>	Blade Electrode Modified	4"	12
<a href="#">0314</a>	Blade Electrode	6.5"	12
<a href="#">0314M</a>	Blade Electrode Modified	6.5"	12

## MEGADYNE™ Stainless Steel Needle Electrodes




Product Code	Description	Tip Length	EA/BX
<a href="#">0313</a>	Needle Electrode	2.75"	12
<a href="#">0313M</a>	Needle Electrode Modified	2.75"	12
<a href="#">0316AM</a>	Needle Electrode Modified	4"	12
<a href="#">0316</a>	Needle Electrode	6"	12
<a href="#">0316M</a>	Needle Electrode Modified	6"	12

## MEGADYNE™ Stainless Steel Ball Electrodes



Product Code	Description	Tip Length	EA/BX
<a href="#">0300</a>	Ball Electrode	5"	12
<a href="#">0315</a>	Ball Electrode	2"	12

## MEGADYNE™ Pencils With Stainless Steel Electrodes



Product Code	Description	EA/BX
<a href="#">0038</a>	Rocker switch pencil	50
<a href="#">0038H</a>	Rocker switch pencil with holster	50
<a href="#">0039</a>	Button switch pencil	50
<a href="#">0039H</a>	Button switch pencil with holster	50
<a href="#">FC39H</a>	Foot control pencil with holster	50
<a href="#">TC100</a>	Tip Cleaner	200

# Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## MEGADYNE™ Disposable Patient Return Electrodes and Accessories



Product Code	Description	Cord Length	EA/BX
<a href="#">0850C</a>	Single plate return electrode	9ft	50
<a href="#">0855C</a>	Dual plate return electrode	10ft	50
<a href="#">0855CL</a>	Dual plate return electrode	15ft	50
<a href="#">0855CN</a>	Dual plate return electrode	9ft	50
<a href="#">0865C</a>	Pediatric dual plate return electrode	10ft	50
<a href="#">0855</a>	Dual plate return electrode	No Cord	50
<a href="#">0870</a>	Return Electrode Cable, Single Style	10ft	1
<a href="#">0875</a>	Return Electrode Cable, Dual Style	10ft	1

# Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## ENDOPATH® Electrosurgery PROBE PLUS® II System



Product Code	Grip	Control	EA/BX
<a href="#">EPH02</a>	Pistol	Hand	6

<a href="#">EPH04</a>	Pencil	Hand	6
-----------------------	--------	------	---

## ENDOPATH® Electrosurgery PROBE PLUS® II Electrode Shafts



Product Code	Shaft Diameter	Shaft Length	Electrode	EA/BX
<a href="#">EPS01</a>	5mm	34cm	Hook	6
<a href="#">EPS02</a>	5mm	34cm	Spatula	6
<a href="#">EPS03</a>	5mm	34cm	Right angle	6
<a href="#">EPS04</a>	5mm	34cm	Curved dissector	6
<a href="#">EPS05</a>	5mm	29cm	Hook	6
<a href="#">EPS06</a>	5mm	29cm	Spatula	6
<a href="#">EPS07</a>	5mm	29cm	Right angle	6

## ENDOPATH® Electrosurgery PROBE PLUS® II Tubing



Product Code	Description	EA/BX
<a href="#">EPT01</a>	PROBE PLUS® II, dual-bottle irrigation tube set	6

<a href="#">EPT03</a>	PROBE PLUS® II, pressurized bag irrigation tube set with spike	6
-----------------------	--	---

# Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## MEGADYNE™ All-In-One Hand Control Unit



Product Code	Description	EA/BX
<a href="#">0055</a>	All-In-One Hand Control Unit	5

## MEGADYNE™ All-In-One Electrodes



Product Code	Tip Configuration	Length	EA/BX
<a href="#">0060</a>	L-Hook	32cm	5
<a href="#">0085</a>	J-Wire Hook	38cm	5
<a href="#">0062</a>	Curved Blade	32cm	5
<a href="#">0063</a>	L-Wire Hook	32cm	5
<a href="#">0080J</a>	L-Wire Hook	38cm	5
<a href="#">0064</a>	Standard Blade	32cm	5

## MEGADYNE™ All-In-One Electrosurgical Tip



Product Code	Tip Configuration	Tip Length	EA/BX
<a href="#">0066</a>	Standard Blade	2.5"	5

## MEGADYNE™ Suction Coagulators



Product Code	Description	French Size	EA/BX
<a href="#">004010J</a>	Foot control, 6" cannula, 10 ft cord	10 Fr	10
<a href="#">004025</a>	Foot control, 6" cannula, 10 ft cord	10 Fr	25
<a href="#">004125</a>	Hand control, 6" cannula, 10 ft cord	10 Fr	25
<a href="#">004225</a>	Hand control, 6" cannula, 10 ft cord	8 Fr	25
<a href="#">004325</a>	Hand control, 6" cannula, 10 ft cord	12 Fr	25

# Energy

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## MEGADYNE™ LLETZ Loop Electrodes Square



Product Code	Description	Color	EA/BX
<a href="#">0410</a>	LLETZ Loop 10 x 7 mm	Orange	12

## MEGADYNE™ LLETZ Loop Electrodes Round



Product Code	Description	Color	EA/BX
<a href="#">0420</a>	LLETZ Loop 10 x 5 mm	Purple	12
<a href="#">0430</a>	LLETZ Loop 15 x 6 mm	Gray	12
<a href="#">0440</a>	LLETZ Loop 20 x 10 mm	Red	12
<a href="#">0450</a>	LLETZ Loop 10 x 10 mm	Yellow	12
<a href="#">0460</a>	LLETZ Loop 15 x 12 mm	Green	12
<a href="#">0470</a>	LLETZ Loop 20 x 12 mm	White	12
<a href="#">0480</a>	LLETZ Loop 20 x 15 mm	Blue	12
<a href="#">0490</a>	LLETZ Loop 20 x 8 mm	Tan	12

## MEGADYNE™ Bipolar Forceps



Product Code	Description	Tip	Length	EA/BX
<a href="#">4010J</a>	Reusable Scoville-Greenwood Forcep	1.5mm	7"	1
<a href="#">4020J</a>	Reusable Jeweler Forcep	0.4mm	4"	1
<a href="#">4030J</a>	Reusable Cushing Forcep	1.5mm	6"	1
<a href="#">4040J</a>	Reusable Adson Forcep	1mm	4"	1
<a href="#">4050J</a>	Reusable Simkin Forcep	0.5mm	5"	1
<a href="#">4060J</a>	Reusable Gerald Forcep	1mm	7"	1
<a href="#">4000J</a>	Disposable Flying Lead Bipolar Cable		12'	10
<a href="#">4005J</a>	Disposable Fixed Lead Bipolar Cable		12'	10



# Hernia Repair & Fixation

Ethicon has been pioneering advances in hernia surgery to help restore patients back to their pre-hernia life. The Ethicon portfolio offers choice to help ensure you have the right solutions for the right procedure. Our innovative value-led solutions are meant to provide strong and safe hernia repair and are designed for good long-term patient outcomes.

**[Learn more about the portfolio on Ethicon.com](https://www.ethicon.com)**



# Hernia Repair & Fixation

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## PROLENE® Polypropylene Mesh



Product Code	Description	Size	EA/BX
<a href="#">PMH</a>	Square	6" x 6" (15cm x 15cm)	6
<a href="#">PMIL</a>	Rectangle	3" x 6" (8cm x 15cm)	6
<a href="#">PML</a>	Square	12" x 12" (30cm x 30cm)	3
<a href="#">PMLK</a>	Large pre-shaped keyhole	2.4" x 5.4" (6cm x 13.7cm)	6
<a href="#">PMSK</a>	Small pre-shaped keyhole	1.8" x 4.0" (4.6cm x 10cm)	6
<a href="#">PMXL</a>	Rectangle	2" x 12" (5cm x 30cm)	6
<a href="#">PMXS</a>	Rectangle	1" x 4" (2.54cm x 10cm)	6

## PROLENE® Soft Polypropylene Mesh



Product Code	Description	Size	EA/BX
<a href="#">SPMXS</a>	Rectangle	1" x 4" (2.54cm x 10cm)	6
<a href="#">SPMS</a>	Rectangle	2" x 4" (5cm x 10cm)	6
<a href="#">SPMIL</a>	Rectangle	3" x 6" (7.6cm x 15cm)	6
<a href="#">SPMH</a>	Square	6" x 6" (15cm x 15cm)	6
<a href="#">SPMLI</a>	Square	10" x 10" (25cm x 25cm)	3
<a href="#">SPMXXL</a>	Rectangle	14" x 12" (35.6cm x 30cm)	3
<a href="#">SPM3XL</a>	Square	19.6" x 19.6" (50 cm x 50 cm)	3

## ULTRAPRO® Partially Absorbable Lightweight Mesh




Product Code	Description	Size	EA/BX
<a href="#">UMS3</a>	Rectangle	2.4" x 4.3" (6cm x 11cm)	3
<a href="#">UMR3</a>	Rectangle	3" x 6" (7.6cm x 15cm)	3
<a href="#">UMM3</a>	Square	6" x 6" (15cm x 15cm)	3
<a href="#">UML1</a>	Square	12" x 12" (30cm x 30cm)	1

# Hernia Repair & Fixation

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## ULTRAPRO ADVANCED™



Product Code	Description	Size	EA/BX
<a href="#">UPA3612</a>	Rectangle	6cm x 12cm	3
<a href="#">UPA37615</a>	Rectangle	7.6cm x 15cm	3
<a href="#">UPA31015</a>	Rectangle	10cm x 15cm	3
<a href="#">UPA31515</a>	Square	15cm x 15cm	3
<a href="#">UPA1530</a>	Rectangle	15cm x 30cm	1
<a href="#">UPA3030</a>	Square	30cm x 30cm	1

## ULTRAPRO COMFORT PLUG™ Partially Absorbable Hernia Repair Device



Product Code	Plug Diameter	Onlay Mesh Size	EA/BX
<a href="#">UPLUG401</a>	40mm	7cm x 14cm	1
<a href="#">UPLUG403</a>	40mm	7cm x 14cm	3
<a href="#">UPLUG406</a>	40mm	7cm x 14cm	6
<a href="#">UPLUG551</a>	55mm	7cm x 14cm	1
<a href="#">UPLUG553</a>	55mm	7cm x 14cm	3
<a href="#">UPLUG556</a>	55mm	7cm x 14cm	6

## PROLENE® Polypropylene Hernia System



Product Code	Description	Underlay Patch Size	Onlay Patch Size	Connector Size	EA/BX
<a href="#">PHSM</a>	Medium	7cm	4.5cm x 10cm	1.9cm x 1.3cm	1
<a href="#">PSHM6</a>	Medium	7cm	4.5cm x 10cm	1.9cm x 1.3cm	3
<a href="#">PHSL</a>	Large	10cm	4.5cm x 10cm	1.9cm x 1.3cm	6
<a href="#">PHSL6</a>	Large	10cm	4.5cm x 10cm	1.9cm x 1.3cm	1
<a href="#">PHSE</a>	Extended	10cm	5.5cm x 12.8cm	1.9cm x 1.3cm	3
<a href="#">PHSE6</a>	Extended	10cm	5.5cm x 12.8cm	1.9cm x 1.3cm	6

# Hernia Repair & Fixation

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## ULTRAPRO® Hernia System



Product Code	Description	Underlay Patch Size	Onlay Patch Size	EA/BX
<a href="#">UHSL</a>	Large	2.4" x 4.7" (6cm x 12cm)	4.0" (-10cm)	3
<a href="#">UHSL6</a>	Large	2.4" x 4.7" (6cm x 12cm)	4.0" (-10cm)	6
<a href="#">UHSM</a>	Medium	2.4" x 4.7" (6cm x 12cm)	3.0" (7.5cm)	3
<a href="#">UHSM6</a>	Medium	2.4" x 4.7" (6cm x 12cm)	3.0" (7.5cm)	6
<a href="#">UHSOV</a>	Oval	2.4" x 4.7" (6cm x 12cm)	4.0" x 4.7" (-10cm x 12cm)	3

## PROCEED® Surgical Mesh



Product Code	Description	Size	EA/BX
<a href="#">PCDB1</a>	Rectangle	2" x 4" (5cm x 10cm)	1
<a href="#">PCDD1</a>	Rectangle	4" x 8" (10cm x 20cm)	1
<a href="#">PCDG1</a>	Oval	6" x 8" (15cm x 20cm)	1
<a href="#">PCDH1</a>	Oval	8" x 10" (20cm x 25cm)	1
<a href="#">PCDJ1</a>	Rectangle	8" x 12" (20cm x 30cm)	1
<a href="#">PCDL1</a>	Square	12" x 12" (30.5cm x 30.5cm)	1
<a href="#">PCDM1</a>	Square	6" x 6" (15cm x 15cm)	1
<a href="#">PCDN1</a>	Oval	4" x 6" (10cm x 15cm)	1
<a href="#">PCDR1</a>	Rectangle	3" x 6" C 7 .. 5cm x 15cm)	1
<a href="#">PCDT1</a>	Oval	10" x 13" (26cm x 34cm)	1
<a href="#">PCDW1</a>	Rectangle	10" x 14" (25cm x 35cm)	1

## PROCEED® Ventral Patch




Product Code	Description	Size	EA/BX
<a href="#">PVPS</a>	Circle	1.7" x 1.7" (4.3cm x 4.3cm)	2
<a href="#">PVPM</a>	Circle	2.5" x 2.5" (6.4cm x 6.4cm)	2

# Hernia Repair & Fixation

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ETH ICON PHYSIOMESH™ Open Flexible Composite Mesh Device



Product Code	Description	Size	EA/BX
<a href="#">OPHY1010</a>	≤3cm	10cm "ROUND"	1
<a href="#">OPHY1215</a>	3cm-5cm	12cmx15cm	1
<a href="#">OPHY1520</a>	3cm-≤6cm	15cmx20cm	1
<a href="#">OPHY1525</a>	<3cm-≤6cm	15cmx25cm	1
<a href="#">OPHY2030</a>	6cm-10cm	20cmx30cm	1
<a href="#">OPHY2536</a>	≥10cm-≤19cm	25cmx36cm	1

FlexHD® Structural Acellular Hydrated Dermis (allograft)



	Tissue Code	Thickness	Dimension
Thick (0.8mm-1.7mm)	<a href="#">471608</a>	0.4mm - 0.8mm	6cm x 8cm
	<a href="#">471812</a>	0.8mm -1.7mm	8cm x 12cm
	<a href="#">471016</a>	0.8mm -1.7mm	10cm x 16cm
	<a href="#">471216</a>	0.8mm -1.7mm	12cm x 16cm
	<a href="#">471620</a>	0.8mm -1.7mm	16cm x 20cm
	<a href="#">471202</a>	0.8mm -1.7mm	20cm x 20cm
	<a href="#">471224</a>	0.8mm -1.7mm	12cm x 24cm
	<a href="#">471225</a>	0.8mm -1.7mm	20cm x 25cm
Ultra Thick (1.8mm-4mm)	<a href="#">471230</a>	0.8mm -1.7mm	20cm x 30cm
	<a href="#">472812</a>	1.8mm-4mm	8cm x 12cm
	<a href="#">472016</a>	1.8mm-4mm	10cm x 16cm
	<a href="#">472216</a>	1.8mm-4mm	12cm x 16cm
	<a href="#">472620</a>	1.8mm-4mm	16cm x 20cm
	<a href="#">472224</a>	1.8mm-4mm	12cm x 24cm
	<a href="#">472225</a>	1.8mm-4mm	20cm x 25cm
	<a href="#">472230</a>	1.8mm-4mm	20cm x 30cm

# Hernia Repair & Fixation

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

FlexHD® Diamond XL Acellular Hydrated Dermis (allograft)



Product Code	Thickness	Dimensions	EA/BX
<a href="#">4D1331</a>	0.8mm-1.7mm	22cm x 30cm	1
<a href="#">4D1335</a>	0.8mm-1.7mm	24cm x 35cm	1
<a href="#">4D2331</a>	1.8mm-4mm	22cm x 30cm	1
<a href="#">4D2335</a>	1.8mm-4mm	24cm x 35cm	1

XCM BIOLOGIC® Tissue Matrix sterile (hydrated), Thick Average 1.5mm thick= 0.3mm



Product Code	Dimensions	Product Code	EA/BX
<a href="#">XM106.0612S</a>	6cm x 12cm	XM106.1620S	16cm x 20cm
<a href="#">XM106.0616S</a>	6cm x 16cm	XM106.1625S	16cm x 25cm
<a href="#">XM106.0812S</a>	8cm x 12cm	XM106.2020S	20cm x 20cm
<a href="#">XM106.0816S</a>	8cm x 16cm	XM106.2025S	20cm x 25cm
<a href="#">XM106.0820S</a>	8cm x 20cm	XM106.2030S	20cm x 30cm
<a href="#">XM106.1016S</a>	10cm x 16cm	XM106.2040S	20cm x 40cm
<a href="#">XM106.1020S</a>	10cm x 20cm	XM106.2540S	25cm x 40cm
<a href="#">XM106.1220S</a>	12cm x 20cm	XM106.2545S	25cm x 45cm

XCM BIOLOGIC® Tissue Matrix, Tailored Coverage™, Thick



Product Code	Dimensions
<a href="#">XM106.2535TCS</a>	25cm x 35cm

# Hernia Repair & Fixation

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ETHICON SECURESTRAP®  
Absorbable Strap Fixation  
Device



Product Code	Description	EA/BX
<a href="#">STRAP12</a>	5mm single use device with 12 absorbable straps	6
<a href="#">STRAP25</a>	5mm single use device with 25 absorbable straps	6

ETH ICON SECURESTRAP® Open  
Absorbable Strap Fixation Device



Product Code	Description	EA/BX
<a href="#">OPENSTRAP20</a>	Single use device with 20 absorbable straps	6

# Ligation

The Ethicon portfolio of endoscopic and open ligation devices are designed to provide optimal patient surgical experiences. Our deep understanding of the interaction between surgical tools and living tissue has enabled us to design solutions for confident occlusion.

**Learn more about the portfolio on [Ethicon.com](http://Ethicon.com)**





# Ligation

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

LIGAMAX™ 5 Endoscopic Clip Applier



Product Code	Shaft Diameter	Clip Size	Number of clips	EA/BX
<a href="#">EL5ML</a>	5mm	Medium/Large	15	3

10mm LIGACLIP® Endoscopic Rotating Multiple-Clip Applier



Product Code	Shaft Diameter	Clip Size	Number of clips	EA/BX
<a href="#">ER320</a>	10mm	Medium/Large	20	3

12mm LIGACLIP® Endoscopic Rotating Multiple-Clip Applier



Product Code	Shaft Diameter	Clip Size	Number of clips	EA/BX
<a href="#">ER420</a>	12mm	Large	20	3

LIGACLIP® Endoscopic Clip Applier



Product Code	Clip Size	EA/BX
<a href="#">EL214</a>	Medium	1
<a href="#">EL314</a>	Medium/Large	1
<a href="#">EL414</a>	Large	1

# Ligation

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## LIGACLIP® Multi-Clip Appliers (MCA)



Product Code	Shaft Length	Clip Size	Number of clips	EA/BX
<a href="#">MCS20</a>	9 3/8"	Small	20	6
<a href="#">MSM2</a>	9 3/8"	Medium	20	6
<a href="#">MCM20</a>	11"	Medium	20	6
<a href="#">MCM30</a>	11"	Medium	30	6
<a href="#">MCL20</a>	13"	Large	20	6

## LIGACLIP® Multi-Patient Single-Clip Applier - Blue



Product Code	Length	Handle Color	Clip Size	Specification	EA/BX
<a href="#">LX105</a>	5 3/4"	Blue	Small	-	1
<a href="#">LX107</a>	7 1/2"	Blue	Small	-	1
<a href="#">LX110</a>	10 1/2"	Blue	Small	-	1
<a href="#">LX130</a>	10 1/2"	Blue	Small	Right angle	1

## LIGACLIP® Multi-Patient Single-Clip Applier - Silver



Product Code	Length	Handle Color	Clip Size	Specification	EA/BX
<a href="#">LX205</a>	5 3/4"	Silver	Medium	-	1
<a href="#">LX207</a>	7 1/2"	Silver	Medium	-	1
<a href="#">LX210</a>	10 1/2"	Silver	Medium	-	1
<a href="#">LX220</a>	10 1/2"	Silver	Medium	Angled jaw	1
<a href="#">LX230</a>	10 1/2"	Silver	Medium	Right jaw	1

## LIGACLIP® Multi-Patient Single-Clip Applier - Green



Product Code	Length	Handle Color	Clip Size	Specification	EA/BX
<a href="#">LC307</a>	7 1/2"	Green	Medium/Large	-	1
<a href="#">LC310</a>	10 1/2"	Green	Medium/Large	-	1
<a href="#">LC320</a>	10 1/2"	Green	Medium/Large	-	1

# Ligation

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

LIGACLIP® Multi-Patient Single-Clip Applier - Yellow



Product Code	Length	Handle Color	Clip Size	Specification	EA/BX
<a href="#">LC407</a>	7 1/2"	Yellow	Large	-	1
<a href="#">LC410</a>	10 1/2"	Yellow	Large	-	1
<a href="#">LC420</a>	10 1/2"	Yellow	Large	Angled jaw	1
<a href="#">LC430</a>	10 1/2"	Yellow	Large	Right angle	1

LIGACLIP® Ligating Clip Cartridge Base



Product Code	Description	EA/BX
<a href="#">LC800</a>	LIGACLIP® Ligating Clip Cartridge Base	1

LIGACLIP® EXTRA Ligating Clips  
- For use with open and endoscopic single clip appliers



Product Code	Length	Handle Color	Clip Size	Specification	EA/BX
<a href="#">LT100</a>	Open single	Blue	Small	6	36
<a href="#">LT102</a>	Open single	Blue	Small	20	15
<a href="#">LT200</a>	Open/endo	Silver	Medium	6	36
<a href="#">LT202</a>	Open/endo	Silver	Medium	20	15
<a href="#">LT300</a>	Open/endo	Green	Medium/Large	6	18
<a href="#">LT400</a>	Open/endo	Yellow	Large	6	18

# Surgical Stapling

Reliable tissue repair is a critical factor for all surgery. That's why Ethicon has worked with surgeons worldwide to develop and refine our surgical staplers and cutters. Our surgical stapling portfolio addresses the issues of tissue thickness, staple line security, and tissue tension to provide you solutions that are designed to enhance control and patient safety.

**Learn more about the portfolio on [Ethicon.com](https://www.ethicon.com)**



# Surgical Stapling

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ECEHLON FLEX™ Powered Plus  
60mm Stapler



Product Code	Shaft Length	Maximum Firings	EA/BX
<a href="#">PCEE60A</a>	280mm	12	3
<a href="#">PSEE60A</a>	340mm	12	3
<a href="#">PLEE60A</a>	440mm	12	3

ECHELON FLEX™ Powered  
ENDOPATH® 45mm Stapler



Product Code	Shaft Length	Maximum Firings	EA/BX
<a href="#">PCE45A</a>	280mm	12	3
<a href="#">PSE45A</a>	340mm	12	3
<a href="#">PLE45A</a>	440mm	12	3

ECHELON FLEX™ Powered Vascular  
Stapler with Advanced Placement Tip



Product Code	Shaft Length	Maximum Firings	EA/BX
<a href="#">PVE35A</a>	320mm	12	3

ECHELON FLEX™ ENDOPATH®  
60mm Stapler



Product Code	Shaft Length	Articulation	Maximum Firings	EA/BX
<a href="#">SC60A</a>	280mm	Yes	12	3
<a href="#">EC60A</a>	340mm	Yes	12	3
<a href="#">LONG60A</a>	440mm	Yes	12	3

# Surgical Stapling

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ECHELON FLEX™ ENDOPATH®  
45mm Stapler



Product Code	Shaft Length	Articulation	Maximum Firings	EA/ BX
<a href="#">SC45A</a>	280mm	Yes	12	3
<a href="#">EC45A</a>	340mm	Yes	12	3
<a href="#">EC45AL</a>	440mm	Yes	12	3

ECHELON 60mm Reloads  
with Gripping Surface Technology



Product Code	Color Cartridge	Number of Staples	Tissue Thickness	Open Staple Height	Closed Staple Height	EA/ BX
<a href="#">GST60W</a>	White	6	Mesentery/ thin	2.6mm	1.0mm	12
<a href="#">GST60B</a>	Blue	6	Regular	3.6mm	1.5mm	12
<a href="#">GST60D</a>	Gold	6	Regular/ thick	3.8mm	1.8mm	12
<a href="#">GST60G</a>	Green	6	Thick	4.1mm	2.0mm	12
<a href="#">GST60T</a>	Black	6	Very thick	4.2mm	2.3mm	12

ECHELON 60mm Reloads  
- For use with all ECHELON  
60mm Endoscopic Staplers



Product Code	Color Cartridge	Number of Staples	Tissue Thickness	Open Staple Height	Closed Staple Height	EA/ BX
<a href="#">ECR60M</a>	Gray	6	Mesentery/ thin	2.0mm	.75mm	12
<a href="#">ECR60B</a>	Blue	6	Regular	3.5mm	1.5mm	12
<a href="#">ECR60D</a>	Gold	6	Regular/ thick	3.8mm	1.8mm	12
<a href="#">ECR60G</a>	Green	6	Thick	4.1mm	2.0mm	12

# Surgical Stapling

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ECHELON 45mm Reloads  
- For use with all ECHELON 4  
5mm Endoscopic Staplers



Product Code	Color Cartridge	Number of Staples	Tissue Thickness	Open Staple Height	Closed Staple Height	EA/ BX
<a href="#">ECR45M</a>	Gray	6	Mesentery/ thin	2.0mm	.75mm	12
<a href="#">ECR45W</a>	White	6	Vascular/ thin	2.5mm	1.0mm	12
<a href="#">ECR45B</a>	Blue	6	Regular	3.5mm	1.5mm	12
<a href="#">ECR45D</a>	Gold	6	Regular/ thick	3.8mm	1.8mm	12
<a href="#">ECR45G</a>	Green	6	Thick	4.1mm	2.0mm	12

ENDOPATH ECHELON™ Vascular  
White Reload for Advanced  
Placement Tip - 35mm



Product Code	Color Cartridge	Number of Staples	Tissue Thickness	Open Staple Height	Closed Staple Height	EA/ BX
<a href="#">VASECR35</a>	White	4	Vascular/thin	2.5mm	1.0mm	12

ENDOPATH® ETS 45mm Articulating  
Linear Cutter



Product Code	Shaft Length	Articulation	Reload Included	Maximum Firings	EA/ BX
<a href="#">ATS45</a>	340mm	Yes	No	8	3
<a href="#">ATS45NK</a>	340mm	Yes	No	8	3

# Surgical Stapling

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ENDOPATH® ETS 45mm Linear Cutters - For use with all ETS 45mm Endoscopic Cutters



Product Code	Color Cartridge	Number of Staples	Tissue Thickness	Open Staple Height	Closed Staple Height	EA/ BX
<a href="#">6R45M</a>	Gray	6	Mesentery/ thin	2.0mm	.85mm	12
<a href="#">TR45W</a>	White	6	Vascular/ thin	2.5mm	1.0mm	12
<a href="#">6R45B</a>	Blue	6	Regular	3.5mm	1.5mm	12
<a href="#">TR45G</a>	Green	4	Thick	4.1mm	2.0mm	12

PROXIMATE® Curved Intraluminal Stapler



Product Code	Shaft Length	Adjustable Closed Staple Height	Head Diameter	Lumen	Open leg Length	EA/ BX
<a href="#">CDH21A</a>	18cm	1.0-2.5mm	21mm	12.4mm	5.5mm	3
<a href="#">CDH25A</a>	18cm	1.0-2.5mm	25mm	16.4mm	5.5mm	3
<a href="#">CDH29A</a>	18cm	1.0-2.5mm	29mm	20.4mm	5.5mm	3
<a href="#">CDH33A</a>	18cm	1.0-2.5mm	33mm	24.4mm	5.5mm	3

PROXIMATE® Curved Intraluminal Stapler



Product Code	Shaft Length	Adjustable Closed Staple Height	Head Diameter	Lumen	Open leg Length	EA/ BX
<a href="#">ECS21A</a>	28cm	1.0-2.5mm	21mm	12.4mm	5.5mm	3
<a href="#">ECS25A</a>	28cm	1.0-2.5mm	25mm	16.4mm	5.5mm	3
<a href="#">ECS29A</a>	28cm	1.0-2.5mm	29mm	20.4mm	5.5mm	3
<a href="#">ECS33A</a>	28cm	1.0-2.5mm	33mm	24.4mm	5.5mm	3



# Surgical Stapling

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

PROXIMATE® ILS Circular Sizer Set



Product Code	Description	EA/BX
<a href="#">CSS</a>	Bowel sizers for the ENDOPATH® ILS Curved Staplers, the set includes one each of a 25mm, 29mm, and 33mm sizer	1

CONTOUR® Curved Cutter Stapler



Product Code	Cartridge Color	Number of Staple Rows	Tissue Thickness	Open Staple Height	Closed Staple Height	Reload Included	Maximum Firings	EA/BX
<a href="#">CS40B</a>	Blue	Regular	3.5mm	3.5mm	1.44mm	Yes	6	3
<a href="#">CS40G</a>	Green	Thick	4.7mm	4.7mm	2.0mm	Yes	6	3

CONTOUR® Curved Cutter Stapler Cartridge



Reload Code	Cartridge Color	Number of Staple Rows	Tissue Thickness	Open Staple Height	Closed Staple Height	EA/BX
<a href="#">CR40B</a>	Blue	4	Regular	3.5mm	1.44mm	3
<a href="#">CR40G</a>	Green	4	Thick	4.7mm	2.0mm	3

ETHICON ENDO-SURGERY™ Linear Cutter



Product Code	Shaft Length	Articulation	Reload Included	Maximum Firings	EA/BX
<a href="#">NTLC55</a>	55mm	6	1.5mm, 1.8mm, 2.0mm	12	3
<a href="#">NTLC75</a>	75mm	6	1.5mm, 1.8mm, 2.0mm	12	3

ETHICON ENDO-SURGERY™ Linear Cutter Cartridge



Product Code	Number of Stapled Rows	Selectable Closed Staple Height	EA/BX
<a href="#">SR55</a>	6	1.5mm, 1.8mm, 2.0mm	12
<a href="#">SR75</a>	6	1.5mm, 1.8mm, 2.0mm	12

# Surgical Stapling

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

PROXIMATE® Linear Cutter



Product Code	Cartridge Color	Number of Staple Rows	Tissue Thickness	Open Staple Height	Closed Staple Height	Reload Included	Maximum Firings	EA/BX
<a href="#">TVC55</a>	White	4	Vascular/thin	250mm	1.0mm	Yes	8	3
<a href="#">TLC55</a>	Blue	4	Regular	385mm	1.5mm	Yes	8	3
<a href="#">TCT55</a>	Green	4	Thick	450mm	2.0mm	Yes	8	3
<a href="#">TLC75</a>	Blue	4	Regular	385mm	1.5mm	Yes	8	3
<a href="#">TCD75</a>	Gold	4	Regular/thick	420mm	1.8mm	Yes	8	3
<a href="#">TCT75</a>	Green	4	Thick	450mm	2.0mm	Yes	8	3
<a href="#">TLC10</a>	Blue	4	Regular	385mm	1.5mm	Yes	8	3
<a href="#">TCT10</a>	Green	4	Thick	450mm	2.0mm	Yes	8	3

PROXIMATE® Linear Cutter Cartridge



Reload Code	Cartridge Color	Number of Staple Rows	Tissue Thickness	Open Staple Height	Closed Staple Height	EA/BX
<a href="#">TVR55</a>	White	4	Vascular/thin	2.50mm	1.0mm	3
<a href="#">TCR55</a>	Blue	4	Regular	3.85mm	1.5mm	3
<a href="#">TRT55</a>	Green	4	Thick	4.50mm	2.0mm	3
<a href="#">TCR75</a>	Blue	4	Regular	3.85mm	1.5mm	3
<a href="#">TRD75</a>	Gold	4	Regular/thick	4.20mm	1.8mm	3
<a href="#">TRT75</a>	Green	4	Thick	4.50mm	2.0mm	3
<a href="#">TCR10</a>	Blue	4	Regular/thick	3.85mm	1.5mm	3
<a href="#">TRT10</a>	Green	4	Thick	4.50mm	2.0mm	3

# Surgical Stapling

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## PROXIMATE® Reloadable Staplers (TL/TLH/TLV)



Product Code	Tissue Thickness	Maximum Closed Staple Height	Reload Included	Maximum Firings	EA/BX
<a href="#">TL30</a>	1.0mm-2.5mm	1.0mm	Yes	4	3
<a href="#">TL60</a>	1.0mm-2.5mm	1.0mm	Yes	4	3
<a href="#">TL90</a>	1.0mm-2.5mm	1.0mm	Yes	4	3
<a href="#">TLH30</a>	1.0mm-2.5mm	1.0mm	Yes	4	3
<a href="#">TLH60</a>	1.0mm-2.5mm	1.0mm	Yes	4	3
<a href="#">TLH90</a>	1.0mm-2.5mm	1.0mm	Yes	4	3
<a href="#">TLV30</a>	1.0mm	1.0mm	Yes	4	3

## TL Linear Stapler with Dial



Product Code	EA/BX
<a href="#">TR30</a>	3
<a href="#">TL6R</a>	3
<a href="#">TR90</a>	3

## TLH Linear Stapler with Dial



Product Code	EA/BX
<a href="#">TVR30</a>	12
<a href="#">TRH30</a>	12
<a href="#">TRH60</a>	12
<a href="#">TRH90</a>	12

# Surgical Stapling

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

PROXIMATE® Reloadable Staplers (TX)



Product Code	Cartridge Color	Number of Staple Rows	Tissue Thickness	Open Staple Height	Closed Staple Height	Reload Included	Maximum Firings	EA/ BX
<a href="#">TX30V</a>	White	3	Vascular	2.5mm	1.0mm	Yes	8	3
<a href="#">TX30B</a>	Blue	2	Regular	3.5mm	1.5mm	Yes	8	3
<a href="#">TX30G</a>	Green	2	Thick	4.8mm	2.0mm	Yes	8	3
<a href="#">TX60B</a>	Blue	2	Regular	3.85mm	1.5mm	Yes	8	3
<a href="#">TX60G</a>	Gold	2	Thick	4.8mm	1.8mm	Yes	8	3

PROXIMATE® Reloadable Staplers Cartridge (TX)



Product Code	Cartridge Color	Number of Staple Rows	Tissue Thickness	Open Staple Height	Closed Staple Height	EA/ BX
<a href="#">XR30V</a>	White	3	Vascular	2.5mm	1.0mm	12
<a href="#">XR30B</a>	Blue	2	Regular	3.5mm	1.5mm	12
<a href="#">XR30G</a>	Green	2	Thick	4.8mm	2.0mm	12
<a href="#">XR60B</a>	Blue	2	Regular	3.5mm	1.5mm	12
<a href="#">XR60G</a>	Green	2	Thick	4.8mm	2.0mm	12

PROXIMATE® Reloadable Staplers Cartridge (TX)




Product Code	Number of Rows	Closed Staple Height	Knife Diameter	Number of Staples	Staple Leg Length	EA/ BX
<a href="#">PPHO3</a>	2	0.75mm-1.5mm	24.4mm	28	4mm	3

# Uterine & Pelvic Surgery

Ethicon has worked with gynecologists to create innovative new solutions for a greater choice of procedure. Plus, we've aligned our portfolio to the surgical treatment pathway so you can ensure you have the right solutions for the right procedure. Enhanced technique combined with innovative technologies helps to improve surgical performance.

**Learn more about the portfolio on [Ethicon.com](http://Ethicon.com)**



# Uterine & Pelvic Surgery

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## GYNECARE INTERCEED® Absorbable Adhesion Barrier



Product Code	Description	EA/BX
<a href="#">4350</a>	GYNECARE INTERCEED® Absorbable Adhesion Barrier, 3" x 4" (7.6cm x 10.2cm)	10
<a href="#">4350XL</a>	GYNECARE INTERCEED® Absorbable Adhesion Barrier XL, 5" x 6" (12.7cm x 15.2cm)	10

## GYNECARE TVT ABBREVO® Continence System



Product Code	Description	EA/BX
<a href="#">TVTOML</a>	GYNECARE TVT ABBREVO®	1

## GYNECARE TVT EXACT® Continence System



Product Code	Description	EA/BX
<a href="#">TVTRL</a>	GYNECARE TVT EXACT®	1
<a href="#">810061</a>	Rigid catheter guide (reusable)	1

## GYNECARE TVT™ Obturator System tension-free support for incontinence



Product Code	Description	EA/BX
<a href="#">810081</a>	GYNECARE TVT™ Obturator System	1
<a href="#">810081L</a>	GYNECARE TVT™ Laser Cut Obturator System	1

## GYNECARE TVT™ Retropubic System tension-free support for incontinence



Product Code	Description	EA/BX
<a href="#">810041A</a>	GYNECARE TVT™ with abdominal guides	1
<a href="#">810041B</a>	GYNECARE TVT™ (mechanical cut mesh)	1
<a href="#">810041BL</a>	GYNECARE TVT™ (laser cut mesh)	1
<a href="#">830041</a>	GYNECARE TVT™ Undyed	3
<a href="#">830041B</a>	GYNECARE TVT™ (mechanical cut mesh)	3
<a href="#">830041BL</a>	GYNECARE TVT™ (laser cut mesh)	3
<a href="#">810051</a>	Introducer handle (reusable)	1
<a href="#">810061</a>	Rigid catheter guide (reusable)	1

# Uterine & Pelvic Surgery

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## GYNECARE VERSAPOINT™ Bipolar Electrosurgery System



Product Code	Description	EA/BX
<a href="#">00481</a>	Foot switch	1
<a href="#">00480</a>	Connector cable (Handpiece)	1
<a href="#">04995</a>	Universal Power Cord (U.S. and Canada)	1
<a href="#">01936</a>	Hysteroscope 4.0mm 30 degree wide angle lens	1
<a href="#">01938</a>	Hysteroscope 4.0mm 12 degree high magnification lens	1
<a href="#">01735</a>	Resectoscope hysteroscopic bridge/visual obturator	1
<a href="#">01939</a>	Resectoscope rotating outer sheath 27Fr	1
<a href="#">01934</a>	Resectoscope obturator - standard	1
<a href="#">01931</a>	Resectoscope working element - passive	1
<a href="#">00466</a>	Ball tip electrode	5
<a href="#">00467</a>	Twizzle tip electrode	5
<a href="#">00468</a>	Spring tip electrode	5
<a href="#">01950</a>	Bipolar electrode, resectoscopic - vaporizing tip CO degree)	5
<a href="#">01985</a>	Bipolar electrode, resectoscopic - angled loop tip	5
<a href="#">00463</a>	VERSAPOINT II Hysteroscopic Connector Cable	1
<a href="#">00464</a>	VERSAPOINT II Hysteroscopic Footswitch	1
<a href="#">00465</a>	VERSAPOINT II Hysteroscopic Electrosurgical Generator	1
<a href="#">01942</a>	VERSAPOINT II Resectoscopic System- Resectoscope Electrode 4.0mm Angled Loop	5

# Uterine & Pelvic Surgery

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ARTISYN® Y-Shaped Mesh



Product Code	Size	EA/BX
<a href="#">ARTY</a>	27cm x 5cm	1
<a href="#">ARTY5L</a>	27cm x 5cm	5

GYNECAREGYNEMESH® PS  
Nonabsorbable PROLENE®  
Soft Mesh



Product Code	Size	EA/BX
<a href="#">GPSL</a>	10cm x 15cm	1
<a href="#">GPSXL3</a>	25cm x 25cm	3



# Wound Closure

Secure skin closure is an integral part of nearly every surgical procedure. Ethicon has pioneered advances in skin closure that enable you to minimize tissue trauma, optimize healing, and improve patient outcomes. We offer a portfolio of innovative wound closure solutions designed for the specific needs of skin, fascia, and organ tissues.

**Learn more about the portfolio on [Ethicon.com](https://www.ethicon.com)**












# Wound Closure

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

STRATAFIX™ Spiral PDO Knotless Tissue Control Device - **Bidirectional dyed**



Cutting Edge Reverse 					
NEEDLE	SUTURE USP	COLOR	LENGTH	EACHES	CODE
 <b>FS-2</b> 19mm 3/8 Circle	5-0	Violet	7cm X 7cm	1	<a href="#">SXPDP2B421</a>
	4-0	Violet	7cm X 7cm	1	<a href="#">SXPDP2B422</a>
	4-0	Violet	14cm X 14cm	1	<a href="#">SXPDP2B423</a>
 <b>FS</b> 26mm 3/8 Circle	3-0	Violet	14cm X 14cm	1	<a href="#">SXPDP2B417</a>
	2-0	Violet	14cm X 14cm	1	<a href="#">SXPDP2B418</a>
	3-0	Violet	14cm X 14cm	1	<a href="#">SXPDP2B419</a>
	2-0	Violet	14cm X 14cm	1	<a href="#">SXPDP2B420</a>
 <b>CP-2</b> 26mm 1/2 Circle	0	Violet	24cm X 24cm	1	<a href="#">SXPDP2B415</a>
	1	Violet	24cm X 24cm	1	<a href="#">SXPDP2B416</a>
 <b>OS-6</b> 36mm 1/2 Circle	1	Violet	24cm X 24cm	1	<a href="#">SXPDP2B200</a>
	1	Violet	30cm X 30cm	1	<a href="#">SXPDP2B201</a>
 <b>OS-8</b> 40mm 1/2 Circle	1	Violet	36cm X 36cm	1	<a href="#">SXPDP2B202</a>








Tapercut 					
NEEDLE	SUTURE USP	COLOR	LENGTH	EACHES	CODE
 <b>V-4</b> 17mm 3/8 Circle	4-0	Violet	24cm X 24cm	1	<a href="#">SXPDP2B426</a>
 <b>V-26</b> 26mm 3/8 Circle	4-0	Violet	14cm X 14cm	1	<a href="#">SXPDP2B425</a>
	4-0	Violet	24cm X 24cm	1	<a href="#">SXPDP2B424</a>

# Wound Closure

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

STRATAFIX™ Spiral PDO Knotless  
Tissue Control Device -  
**Bidirectional dyed**







Taperpoint 					
NEEDLE	SUTURE USP	COLOR	LENGTH	EACHES	CODE
 <b>CT-3</b> 22mm 1/2 Circle	0	Violet	14cm X 14cm	1	<a href="#">SXPDP2B404</a>
 <b>SH</b> 22mm 1/2 Circle	2-0	Violet	10cm X 10cm	1	<a href="#">SXPDP2B413</a>
	2-0	Violet	14cm X 14cm	1	<a href="#">SXPDP2B414</a>
 <b>MH</b> 36mm 1/2 Circle	2-0	Violet	7cm X 7cm	1	<a href="#">SXPDP2B406</a>
	2-0	Violet	14cm X 14cm	1	<a href="#">SXPDP2B407</a>
	0	Violet	14cm X 14cm	1	<a href="#">SXPDP2B409</a>
	2-0	Violet	24cm X 24cm	1	<a href="#">SXPDP2B408</a>
	0	Violet	24cm X 24cm	1	<a href="#">SXPDP2B410</a>
	0	Violet	30cm X 30cm	1	<a href="#">SXPDP2B411</a>
	2-0	Violet	36cm X 36cm	1	<a href="#">SXPDP2B412</a>
 <b>CT-1</b> 36mm 1/2 Circle	0	Violet	14cm X 14cm	1	<a href="#">SXPDP2B401</a>
	0	Violet	24cm X 24cm	1	<a href="#">SXPDP2B402</a>
	0	Violet	30cm X 30cm	1	<a href="#">SXPDP2B403</a>
 <b>MO-4</b> 36mm 1/2 Circle	0	Violet	36cm X 36cm	1	<a href="#">SXPDP2B400</a>
 <b>CTX</b> 48mm 1/2 Circle	0	Violet	36cm X 36cm	1	<a href="#">SXPDP2B405</a>



# Wound Closure

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

STRATAFIX™ Spiral PGA-PCL  
Knotless Tissue Control Device  
**Bidirectional undyed**



Cutting Edge Reverse ▽					
NEEDLE	SUTURE USP	COLOR	LENGTH	EACHES	CODE
 <b>FS-2</b> 19mm 3/8 Circle	4-0	Undyed	7cm X 7cm	1	<a href="#">SXMD2B405</a>
	3-0	Undyed	7cm X 7cm	1	<a href="#">SXMD2B404</a>
	4-0	Undyed	14cm X 14cm	1	<a href="#">SXMD2B407</a>
	3-0	Undyed	14cm X 14cm	1	<a href="#">SXMD2B406</a>
	4-0	Undyed	30cm X 30cm	1	<a href="#">SXMD2B409</a>
	3-0	Undyed	30cm X 30cm	1	<a href="#">SXMD2B408</a>
 <b>FS-1</b> 24mm 3/8 Circle	4-0	Undyed	30cm X 30cm	1	<a href="#">SXMD2B150</a>
	3-0	Undyed	30cm X 30cm	1	<a href="#">SXMD2B410</a>
 <b>FS</b> 26mm 3/8 Circle	3-0	Undyed	14cm X 14cm	1	<a href="#">SXMD2B411</a>
	4-0	Undyed	30cm X 30cm	1	<a href="#">SXMD2B413</a>
	3-0	Undyed	30cm X 30cm	1	<a href="#">SXMD2B412</a>
 <b>CP-2</b> 26mm 1/2 Circle	2-0	Undyed	14cm X 14cm	1	<a href="#">SXMD2B414</a>




Taperpoint ⊙					
NEEDLE	SUTURE USP	COLOR	LENGTH	EACHES	CODE
 <b>RB-1</b> 17mm 1/2 Circle	4-0	Undyed	16cm X 16cm	1	<a href="#">SXMD2B403</a>
	3-0	Undyed	16cm X 16cm	1	<a href="#">SXMD2B402</a>
 <b>MH</b> 36mm 1/2 Circle	2-0	Undyed	14cm X 14cm	1	<a href="#">SXMD2B400</a>
	2-0	Undyed	36cm X 36cm	1	<a href="#">SXMD2B401</a>

# Wound Closure

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.





STRATAFIX™ Spiral Polypropylene  
Knotless Tissue Control Device  
**Bidirectional undyed**





Taperpoint 					
NEEDLE	SUTURE USP	COLOR	LENGTH	EACHES	CODE
 <b>CT-3</b> 22mm 1/2 Circle	0	Undyed	14cm X 14cm	1	<a href="#">SXPL2B401</a>
 <b>CT-1</b> 36mm 1/2 Circle	1	Undyed	24cm X 24cm	1	<a href="#">SXPL2B400</a>

STRATAFIX™ Spiral PGA-PCL  
Knotless Tissue Control Device  
**Unidirectional undyed**



Cutting Edge Reverse 					
NEEDLE	SUTURE USP	COLOR	LENGTH	EACHES	CODE
 <b>PS-2</b> 19mm 3/8 Circle	4-0	Undyed	30cm	1	<a href="#">SXMD1B105</a>
	3-0	Undyed	30cm	1	<a href="#">SXMD1B103</a>
	4-0	Undyed	60cm	1	<a href="#">SXMD1B106</a>
	3-0	Undyed	60cm	1	<a href="#">SXMD1B104</a>
 <b>PS-1</b> 24mm 3/8 Circle	3-0	Undyed	30cm	1	<a href="#">SXMD1B100</a>
	3-0	Undyed	60cm	1	<a href="#">SXMD1B101</a>
 <b>PSL</b> 30mm 3/8 Circle	2-0	Undyed	60cm	1	<a href="#">SXMD1B102</a>




Taperpoint 					
NEEDLE	SUTURE USP	COLOR	LENGTH	EACHES	CODE
 <b>RB-1</b> 17mm 1/2 Circle	4-0	Undyed	20cm	1	<a href="#">SXMD1B402</a>
	3-0	Undyed	20cm	1	<a href="#">SXMD1B403</a>
	2-0	Undyed	20cm	1	<a href="#">SXMD1B404</a>

# Wound Closure

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.



STRATAFIX™ Spiral PGA-PCL  
Knotless Tissue Control Device  
**Unidirectional undyed**




Taperpoint 					
NEEDLE	SUTURE USP	COLOR	LENGTH	EACHES	CODE
 <b>SH</b> 26mm 1/2 Circle	3-0	Undyed	20cm	1	<a href="#">SXMD1B405</a>
	2-0	Undyed	20cm	1	<a href="#">SXMD1B406</a>
 <b>CT-1</b> 36mm 1/2 Circle	3-0	Undyed	45cm	1	<a href="#">SXMD1B407</a>
	2-0	Undyed	45cm	1	<a href="#">SXMD1B408</a>
	0	Undyed	45cm	1	<a href="#">SXMD1B409</a>

STRATAFIX™ Spiral PDO Knotless  
Tissue Control Device  
**Unidirectional dyed**





Cutting Edge Reverse 					
NEEDLE	SUTURE USP	COLOR	LENGTH	EACHES	CODE
 <b>PS-2</b> 19mm 3/8 Circle	4-0	Violet	45cm	1	<a href="#">SXPD1B101</a>
	3-0	Violet	45cm	1	<a href="#">SXPD1B100</a>

Taperpoint 					
NEEDLE	SUTURE USP	COLOR	LENGTH	EACHES	CODE
 <b>CT-1</b> 36mm 1/2 Circle	2-0	Violet	20cm	1	<a href="#">SXPD1B400</a>
	2-0	Violet	30cm	1	<a href="#">SXPD1B401</a>

STRATAFIX™ Spiral  
Polypropylene Knotless  
Tissue Control Device  
**Unidirectional undyed**



Taperpoint 					
NEEDLE	SUTURE USP	COLOR	LENGTH	EACHES	CODE
 <b>SH-1</b> 22mm 1/2 Circle	2-0	Undyed	30cm	1	<a href="#">SXPL1B400</a>

# Wound Closure

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.



Product Code	Description	Adhesive Volume	EA/BX
<a href="#">DNX12</a>	Unique formulation with an easy-to-use applicator for precise control for wide and fine line applications	0.7mL	12

<a href="#">DNX6</a>	Unique formulation with an easy-to-use applicator for precise control for wide and fine line applications	0.7mL	6
----------------------	---	-------	---




Product Code	Description	Adhesive Volume	EA/BX
<a href="#">CLR222US</a>	2-octylcyanoacrylate topical skin adhesive applicator. Flexible, self-adhesive polyester mesh	22cm	2

<a href="#">CLR602US</a>	2-octylcyanoacrylate topical skin adhesive applicator. Flexible, self-adhesive polyester mesh tape dispenser	60cm	2
--------------------------	--	------	---



Product Code	Description	Adhesive Volume	EA/BX
<a href="#">DHVM12</a>	2-octylcyanoacrylate high viscosity formulation for small incisions and lacerations	0.36mL	12




Product Code	Description	EA/BX
<a href="#">EC11</a>	EN-3 Needle - O, 42" Suture	12

<a href="#">EC12</a>	EN-3 Needle - 2-O, 42" Suture	12
----------------------	-------------------------------	----

<a href="#">ESS15</a>	ENDO-HOLDER™ Cannula (for use with 5mm diameter trocar)	1
-----------------------	---	---


## Wound Closure

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.



Product Code	Description	EA/BX
<a href="#">EXIOG</a>	ETHIBOND® Extra Suture, ST-3 needle, O, 42"	12
<a href="#">JK10G</a>	Coated VICRYL® Suture, ST-3 needle, O, 42"	12
<a href="#">ZK10G</a>	PDS® II Suture, ST-3 needle, O, 42"	12

ENDOLOOP® Ligature made from PDS (polydioxanone) suture



Product Code	Description	EA/BX
<a href="#">EZIOG</a>	ENDOLOOP® Ligature made from PDS® II O 18" 12 (polydioxanone) Suture	12



Product Code	Description	EA/BX
<a href="#">EJIOG</a>	ENDOLOOP® Ligature made from VICRYL® O 18" 12 (polyglactin 910) Suture	12



Product Code	Description	EA/BX
<a href="#">E705R</a>	ENDOPATH® Needle Holder, reusable, 5mm, 30cm length	1
<a href="#">E705RL</a>	ENDOPATH® Needle Holder, reusable, 5mm, 45cm length	1
<a href="#">SRNHI</a>	Self Righting Needle Holder, reusable, 5mm, 30cm length	1
<a href="#">SRNHIL</a>	Self Righting Needle Holder, reusable, 5mm, 45cm length	1



Product Code	Description	EA/BX
<a href="#">KA200</a>	Sterile, one-piece clips molded from the polymer poly (p-dioxanone).	1
<a href="#">KA200L</a>	Sterile, one-piece clips molded from the polymer poly (p-dioxanone).	1
<a href="#">XC200</a>	Sterile, one-piece clips molded from the polymer poly (p-dioxanone).	6



# Wound Closure

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## Suture Assistant



Product Code	Needle Type	EA/BX
<a href="#">SW100</a>	N/A	4
<a href="#">SW110</a>	6" ETHIBOND EXCEL® (size 0) and EEN (canoe) needle	10
<a href="#">SW112</a>	6" ETHIBOND EXCEL® (size 2-0) and EEN (canoe) needle	10
<a href="#">SW120</a>	6" ETHIBOND EXCEL® (size 0) and SH needle	10
<a href="#">SW122</a>	6" ETHIBOND EXCEL® (size 2-0) and SH Needle	10
<a href="#">SW222</a>	6" VICRYL® (size 2-0) and SH Needle	10

## PROXIMATE® Rotating Head Skin Staplers



Product Code	Number of Staples	Staple Leg Length	Staple Wire Diameter	Staple Crown	Size	EA/BX
<a href="#">PRR35</a>	35	3.9	0.53	5.7	Regular	6
<a href="#">PRW35</a>	35	3.9	0.58	6.9	White	6

## PROXIMATE® Plus MD Skin Stapler



Product Code	Number of Staples	Staple Leg Length	Staple Wire Diameter	Staple Crown	Size	EA/BX
<a href="#">PMR35</a>	35	3.9	0.53	5.7	Regular	6
<a href="#">PMW35</a>	35	3.9	0.58	6.9	Wide	6

## PROXIMATE® PX Fixed Head Skin Staplers



Product Code	Number of Staples	Staple Leg Length	Staple Wire Diameter	Staple Crown	Size	EA/BX
<a href="#">PXR35</a>	35	3.9	0.53	5.7	Regular	6
<a href="#">PXW35</a>	35	3.9	0.58	6.9	Wide	6

## Staple Extractor




Product Code	Number of Staples	EA/BX
<a href="#">PSX</a>	Squeeze handle staple extractor	12

# Vascular Access & Infection Control

Infection Control has become a high priority for hospitals and throughout the continuum of care. Clinicians focus on fewer breaches of the skin and proper site maintenance which improves outcomes, lowers costs and increases patient satisfaction. Ethicon provides solutions for your vascular access and infection control needs.

**Learn more about the portfolio on [Ethicon.com](http://Ethicon.com)**



## Indications for Use

BIOPATCH® containing Chlorhexidine Gluconate is intended for use as a hydrophilic wound dressing that is used to absorb exudate and to cover a wound caused by the use of vascular and non-vascular percutaneous medical devices such as: IV catheters, central venous lines, arterial catheters, dialysis catheters, peripherally inserted coronary catheters, mid-line catheters, drains, chest tubes, externally placed orthopedic pins, and epidural catheters. It is also intended to reduce local infections, catheter-related blood stream infections (CRBSI), and skin colonization of microorganisms commonly related to CRBSI, in patients with central venous or arterial catheters.

# Vascular Access & Infection Control

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

BIOPATCH Protective Disk with CHG



Product Code	Size	French Size Range	Common Uses
<a href="#">4150</a>	1" disc (2.5cm) w/4mm center hole	6-12Fr	Central Lines, PICCs

BIOPATCH Protective Disk with CHG



Product Code	Size	French Size Range	Common Uses
<a href="#">4151</a>	3/4" disc (1.9cm) w/1.5mm center hole	<6Fr	Peripheral IVs, Huber Needles (ports), Arterial Lines, Extended Dwell PIVs, Midlines/ PICCS, Pins

BIOPATCH Protective Disk with CHG



Product Code	Size	French Size Range	Common Uses
<a href="#">4152</a>	1" disc (1.9cm) w/2.5mm center hole	13-20Fr	Dialysis, Drians, Sheaths, Cordis, VAD drive Lines

# Other Devices

Whatever your needs during the course of surgery, Ethicon is committed to providing you a solution. From laparoscopic hand instruments to gastric banding, drains, and reservoirs, you can count us to provide you with outstanding value, excellent service, and a broad range of product offerings.

**[Learn more about the portfolio on Ethicon.com](https://www.ethicon.com)**



## Other Devices

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

### ENDOPATH® Bipolar Forceps



Product Code	Shaft Diameter	Jaw	EA/BX
<a href="#">EBFO1</a>	5mm	Macro	5
<a href="#">EBFO2</a>	5mm	Micro	5

### ENDOPATH® Active Cords - for use with bipolar forceps

Product Code	Description	EA/BX
<a href="#">EBCO2</a>	Reusable coaxial pin connector bipolar cord, gray, 1 prong	1
<a href="#">EBCO5</a>	Active return cord for bipolar instruments fixed pin connector	1

### ENDOPATH® Scissors



Product Code	Description	EA/BX
<a href="#">5DCS</a>	Scissors, curved, with monopolar cautery, 5mm	6

### ENDOPATH® Graspers



Product Code	Description	EA/BX
<a href="#">10AG</a>	Anvil grasper	6
<a href="#">5DSG</a>	Grasper, 5mm	6

# Other Devices

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## ENDOPATH® Bipolar Forceps



Product Code	Description	EA/BX
<a href="#">SDCD</a>	Dissector, curved, with monopolar cautery, 5mm	6
<a href="#">BCD10</a>	Dissector, blunt tip, 10mm, 3 per pouch, 12 pouches per sales unit	36
<a href="#">BTDO5</a>	Dissector, blunt tip, 5mm, 3 per pouch, 12 pouches per sales unit	36

## ENDOPATH® Babcocks



Product Code	Description	EA/BX
<a href="#">SBB</a>	Babcock, 5mm	6
<a href="#">10BB</a>	Babcock, 10mm	6

## Purse String Clamp - Multi-patient use



Product Code	Description	EA/BX
<a href="#">EH40</a>	Reusable, purse string clamp	1

## ENDOPOUCH® Specimen Retrieval Bag



Product Code	EA/BX
<a href="#">POUCH</a>	6

## ENDOPATH® Electrosurgery PROBE PLUS® II System



Product Code	Grip	Control	Common Uses
<a href="#">EPHO2</a>	Pistol	Hand	6
<a href="#">EPHO4</a>	Pencil	Hand	6

# Other Devices

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

ENDOPATH® Electrosurgery  
PROBE PLUS® II Electrode Shafts



Product Code	Shaft Diameter	Shaft Length	Electrode	EA/BX
<a href="#">EPS01</a>	5mm	34cm	Hook	6
<a href="#">EPS02</a>	5mm	34cm	Spatula	6
<a href="#">EPS03</a>	5mm	34cm	Right angle	6
<a href="#">EPS04</a>	5mm	34cm	Curved dissector	6
<a href="#">EPS05</a>	5mm	29cm	Hook	6
<a href="#">EPS06</a>	5mm	29cm	Spatula	6
<a href="#">EPS07</a>	5mm	29cm	Right angle	6

ENDOPATH® Electrosurgery  
PROBE PLUS® II Tubing



Product Code	Description	EA/BX
<a href="#">EPT01</a>	PROBE PLUS® II, dual-bottle irrigation tube set	6
<a href="#">EPT03</a>	PROBE PLUS® II, pressurized bag irrigation tube set with spike	6

ETHICON ENDO-SURGERY®  
Gastric Band Dissector



Product Code	Description	EA/BX
<a href="#">RGBRSI</a>	Gastric Band dissector device	1

ETHICON ENDO-SURGERY®  
Gastric Calibration Tube



Product Code	Description	EA/BX
<a href="#">RGCT180</a>	Gastric calibration tube, 180 asymmetrical	1
<a href="#">RGCT360</a>	Gastric calibration tube, 360 asymmetrical	1

# Other Devices

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

## BLAKE® Cardio Connectors

Product Code	Description	EA/BX
<a href="#">BCC1</a>	BLAKE® Cardio Connector 1:1 - compatible with 19 and 24Fr round hubless drains	1
<a href="#">BCC2</a>	BLAKE® Cardio Connector 2:1 - compatible with 19 and 24Fr round hubless drains	1
<a href="#">BCC3</a>	BLAKE® Cardio Connector 3:1 - compatible with 19 and 24Fr round hubless drains	1

## BLAKE® Silicone Drains - Flat



Product Code	Material	EA/BX
<a href="#">2210</a>	7mm flat, 3/4 fluted	1
<a href="#">2211</a>	7mm flat, full fluted	1
<a href="#">2212</a>	7mm flat, full fluted, 3/16" trocar	1
<a href="#">2213</a>	10mm flat, 3/4 fluted	1
<a href="#">2214</a>	10mm flat, full fluted	1
<a href="#">2215</a>	10mm flat, full fluted, 3/16" trocar	1
<a href="#">2216</a>	7mm flat, 3/4 fluted, 3/16" trocar	1
<a href="#">2217</a>	10mm flat, 3/4 fluted, 3/16" trocar	1

## BLAKE® Silicone Drains - Hubless



Product Code	Material	EA/BX
<a href="#">2226</a>	Hubless, 10Fr round	1
<a href="#">2227</a>	Hubless, 10Fr round, 1/8" trocar	1
<a href="#">2228</a>	Hubless, 15Fr round	1
<a href="#">2229</a>	Hubless, 15Fr round, 3/16" trocar	1
<a href="#">2230</a>	Hubless, 19Fr round	1
<a href="#">2231</a>	Hubless, 19Fr round, 1/4" trocar	1
<a href="#">2232</a>	Hubless, 19Fr round, 1/4" bendable trocar	1
<a href="#">2233</a>	Hubless, 15Fr round, 3/16" bendable trocar	1
<a href="#">2234</a>	Hubless, 24Fr round	1



## Other Devices

Interested in ordering? Click on a [product code](#) and add it to your shopping cart on J&J Customer Connect. Requires customer login.

VICRYL® (polyglactin 910) Knitted Mesh  
VICRYL® (polyglactin 910) Woven Mesh



Product Code	Description	Size	EA/BX
<a href="#">VKML</a>	Knitted	12" x 12" (30cm x 30cm)	3
<a href="#">VKMM</a>	Knitted	6" x 6" (15cm x 15cm)	3
<a href="#">VWML</a>	Woven	12" x 12" (30cm x 30cm)	3
<a href="#">VWMM</a>	Woven	6" x 6" (15cm x 15cm)	3

VICRYL® (polyglactin 910) Mesh Bag

Product Code	Size	EA/BX
<a href="#">VM106K</a>	10" x 12" (26cm x 30cm)	1
<a href="#">VM210K</a>	17" x 12" (44cm x 30cm)	1
<a href="#">VM220K</a>	20" x 8" (50cm x 20cm)	1

J-VAC® Reservoirs - Sterile

Product Code	Description	EA/BX
<a href="#">2161</a>	150ml closed wound drainage reservoir	1
<a href="#">2162</a>	450ml closed wound drainage reservoir	1
<a href="#">2163</a>	300ml closed wound drainage reservoir	1

J-VAC® Bulb Suction Reservoirs  
- Sterile Bulbs

Product Code	Material	EA/BX
<a href="#">2160</a>	100cc bulb suction reservoir	1

J-VAC® Silicone Drains  
- Drain Adapters

Product Code	Material	EA/BX
<a href="#">2209</a>	J-VAC Drain Adapters 1/8"	1
<a href="#">2199</a>	J-VAC Drain Adapters 1/4"	1
<a href="#">2298</a>	J-VAC Drain Adapters 3/16"	1

# Products by Surgical Specialty

## Bariatric Surgery

Ethicon is committed to shaping the future of surgery in the fight against obesity and metabolic diseases and helping to elevate the standard of care through research and evidence, best-in-class education and training, innovating products and expanding patient access to care. Together, we're advancing surgical care, enabling patients to live longer, more fulfilling lives. For more information, [click here](#).


Ethicon offers a comprehensive portfolio of market-leading, innovative devices to meet the changing demands of today's most effective surgical treatments for severe obesity, including:



ECHELON FLEX™ GST System



HARMONIC ACE®+ 7 SHEARS with Advanced Hemostasis



STRATAFIX™ Knotless Tissue Control Devices

## Colorectal Surgery

Our commitment to colorectal surgery focuses on reducing risk factors that drive variability in outcomes and costs. By combining innovation with education and research, we enable more predictability in outcomes. For more information, [click here](#).

Ethicon is addressing the most pressing challenges of today's colorectal procedures through innovations such as:



ETHICON ENDO-SURGERY® Endoscopic Curved Intraluminal Circular Stapler



ECHELON FLEX™ GST System



HARMONIC ACE®+ 7 SHEARS with Advanced Hemostasis

# Products by Surgical Specialty

## Gynecologic Surgery

Our goal is to collaborate with the gynecologic community to help promote positive clinical outcomes. We offer solutions to support a number of gynecological surgical procedures – both open and laparoscopic – as well as a variety of techniques. For more information, [click here](#).


Ethicon provides a broad portfolio of devices and products to help meet your clinical and economic goals, including:



HARMONIC ACE®+ 7 SHEARS  
with Advanced Hemostasis



SURGICEL® Family  
of Absorbable Hemostats




STRATAFIX™ Knotless  
Tissue Control Devices

## Hernia Repair

Ethicon’s commitment to advancing hernia repair is reflected in our ongoing investments in evidence, innovative products, and best-in-class professional education. For more information, [click here](#).


Ethicon provides essential solutions inspired by surgeon and patient needs, including:



ETHICON PHYSIOMESH™  
Open Flexible Composite  
Mesh Device



ETHICON SECURESTRAP®  
Absorbable Strap Fixation  
Device



ULTRAPRO COMFORT  
PLUG™ Partially Absorbable  
Hernia Repair Device

# Products by Surgical Specialty

## Hepato-Biliary Surgery

We are committed to developing flexible solutions for every aspect of cancer procedures – from supporting education to funding research to delivering innovations that reduce the incidence and cost of managing bleeding and infection-related complications. For more information, [click here](#).

Ethicon offers a broad portfolio that delivers surgical precision and optimal hemostasis at every critical step of complex HPB procedures. Our portfolio includes such key innovations as:



HARMONIC™ HD 1000i Shears



SURGICEL® Family of Absorbable Hemostats



STRATAFIX™ Knotless Tissue Control Devices

## Thoracic Surgery

We are dedicated to procedural precision and have a renewed focus on bringing greater precision to every aspect of thoracic surgery. Through training, research, and innovative products, we provide tools to enable better outcomes in thoracic surgery. For more information, [click here](#).

Our portfolio of innovative products is designed to address the most pressing challenges of today’s thoracic procedures through innovation, including:



ECHELON FLEX™ Powered Vascular Stapler



ECHELON FLEX™ GST System



HARMONIC™ HD 1000i Shears

# Essential Product Information

## **SURGICEL Essential Product Information**

### **INDICATIONS**

SURGICEL® Absorbable Hemostat (oxidized regenerated cellulose) is used adjunctively in surgical procedures to assist in the control of capillary, venous, and small arterial hemorrhage when ligation or other conventional methods of control are impractical or ineffective. SURGICEL® ORIGINAL, SURGICEL® FIBRILLAR™ and SURGICEL® NU-KNIT® Hemostats can be cut to size for use in endoscopic procedures.

### **PRECAUTIONS**

Use only as much SURGICEL® Absorbable Hemostat as is necessary for hemostasis, holding it firmly in place until bleeding stops. Remove any excess before surgical closure in order to facilitate absorption and minimize the possibility of foreign body reaction.

In urological procedures, minimal amounts of SURGICEL® Absorbable Hemostat should be used and care must be exercised to prevent plugging of the urethra, ureter, or a catheter by dislodged portions of the product.

Since absorption of SURGICEL® Absorbable Hemostat could be prevented in chemically cauterized areas, its use should not be preceded by application of silver nitrate or any other escharotic chemicals.

If SURGICEL® Absorbable Hemostat is used temporarily to line the cavity of large open wounds, it should be placed so as not to overlap the skin edges. It should also be removed from open wounds by forceps or by irrigation with sterile water or saline solution after bleeding has stopped.

Precautions should be taken in otorhinolaryngologic surgery to assure that none of the material is aspirated by the patient. (Examples: controlling hemorrhage after tonsillectomy and controlling epistaxis.)

Care should be taken not to apply SURGICEL® Absorbable Hemostat too tightly when it is used as a wrap during vascular surgery (see Adverse Reactions).

### **ADVERSE REACTIONS**

"Encapsulation" of fluid and foreign body reactions have been reported.

There have been reports of stenotic effect when SURGICEL® Absorbable Hemostat has been applied as a wrap during vascular surgery.

Paralysis and nerve damage have been reported when SURGICEL® Absorbable Hemostat was used around, in, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasm.

Blindness has been reported in connection with surgical repair of a lacerated left frontal lobe when SURGICEL® Absorbable Hemostat was placed in the anterior cranial fossa.

Possible prolongation of drainage in cholecystectomies and difficulty passing urine per urethra after prostatectomy have been reported.

For more information, please consult your doctor or for product quality and technical questions, call 1-800-795-0012.

# Essential Product Information

## INSTAT Essential Product Information

### INDICATION

INSTAT MCH is recommended for use in surgical procedures (other than Urological and Ophthalmological surgery) as an adjunct to hemostasis when control of bleeding by ligature or conventional procedures is ineffective or impractical.

### PRECAUTIONS

As with other hemostatic agents, it is not recommended that INSTAT MCH be left in an infected or contaminated space.

Safety of this product has not been established in pregnant women; therefore, it should only be used when benefit to risk clearly warrants its use.

Another microfibrillar collagen agent has been reported to cause interference with the healing of skin edges when used in the closure of skin incisions, and to reduce the strength of methyl methacrylate adhesives used to attach prosthetic devices to bone surfaces.

INSTAT MCH is not intended to be used to treat systemic coagulation disorders.

The long term effects of leaving INSTAT MCH in situ are unknown.

### ADVERSE EVENTS

Adverse reactions reported with another microfibrillar collagen hemostatic agent that were possibly related to its use were adhesion formation, allergic reaction, foreign body reaction, and subgaleal seroma (report of a single case). The use of microfibrillar collagen in dental extraction sockets has been reported to increase the incidence of alveolalgia.

For more information, please consult your doctor or for product quality and technical questions, call 1-800-795-0012.

# Essential Product Information

## **SURGIFLO® Hemostatic Matrix Essential Product Information (Made from Absorbable Gelatin Sponge, USP)**

### **DESCRIPTION**

SURGIFLO® Hemostatic Matrix Kit is intended for hemostatic use by applying to a bleeding surface.

### **ACTIONS**

When used in appropriate amounts SURGIFLO® Hemostatic Matrix is absorbed completely within 4 to 6 weeks.

### **INTENDED USE/INDICATIONS**

SURGIFLO® Hemostatic Matrix, mixed with sterile saline or thrombin solution, is indicated in surgical procedures (other than ophthalmic) as an adjunct to hemostasis when control of bleeding by ligature or other conventional methods is ineffective or impractical.

### **CONTRAINDICATIONS**

Do not use SURGIFLO® Hemostatic Matrix in intravascular compartments because of the risk of embolization.

Do not use SURGIFLO® Hemostatic Matrix in patients with known allergies to porcine gelatin.

Do not use SURGIFLO® Hemostatic Matrix in closure of skin incisions because it may interfere with the healing of skin edges. This interference is due to mechanical interposition of gelatin and is not secondary to intrinsic interference with wound healing.

### **WARNINGS**

SURGIFLO® Hemostatic Matrix should not be used in the presence of infection and should be used with caution in contaminated areas of the body.

SURGIFLO® Hemostatic Matrix should not be used in instances of pumping arterial hemorrhage. SURGIFLO® Hemostatic Matrix will not act as a tampon or plug in a bleeding site.

SURGIFLO® Hemostatic Matrix should be removed from the site of application when used in, around, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasm because it may swell resulting in nerve damage.

Excess SURGIFLO® Hemostatic Matrix should be removed once hemostasis has been achieved.

The safety and effectiveness of SURGIFLO® Hemostatic Matrix for use in ophthalmic procedures has not been established.

SURGIFLO® Hemostatic Matrix should not be used for controlling post-partum intrauterine bleeding or menorrhagia.

The safety and effectiveness of SURGIFLO® Hemostatic Matrix has not been established in children and pregnant women.

The blue flexible applicator tip should not be trimmed to avoid exposing internal guidewire.

The white straight applicator tip should be trimmed away from the surgical area. Cut a square angle to avoid creating a sharp tip.

### **PRECAUTIONS**

Safe and effective use of SURGIFOAM® Sponge has been reported in a published neurologic retrospective study involving 1700 cases in Europe. Safe and effective use in neurosurgery has not been proven through randomized, controlled clinical studies in the United States.

SURGIFLO® Hemostatic Matrix is supplied as a sterile product and cannot be resterilized.

SURGIFLO® Hemostatic Matrix should not be used for packing unless excess product that is not needed to maintain hemostasis is removed.

SURGIFLO® Hemostatic Matrix may swell up to 20% upon contact with additional fluid.

SURGIFLO® Hemostatic Matrix should not be used in conjunction with autologous blood salvage circuits.

SURGIFLO® Hemostatic Matrix should not be used in conjunction with methylmethacrylate adhesives.

In urological procedures, SURGIFLO® Hemostatic Matrix should not be left in the renal pelvis or ureters to eliminate the potential foci for calculus formation.

# Essential Product Information

## **SURGIFLO® Hemostatic Matrix Essential Product Information (Made from Absorbable Gelatin Sponge, USP)**

### **ADVERSE EVENTS**

A total of 142 patients received SURGIFOAM® Sponge during a clinical trial comparing SURGIFOAM® Sponge to another absorbable gelatin sponge. In general, the following adverse events have been reported with the use of absorbable porcine gelatin-based hemostatic agents:

Gelatin-based hemostatic agents may serve as a nidus for infection and abscess formation and have been reported to potentiate bacterial growth.

Giant cell granulomas have been observed at implant sites when used in the brain.

Compression of the brain and spinal cord resulting from the accumulation of sterile fluid have been observed.

Multiple neurologic events were reported when absorbable gelatin-based hemostatic agents were used in laminectomy operations, including cauda equina syndrome, spinal stenosis, meningitis, arachnoiditis, headaches, paresthesias, pain, bladder and bowel dysfunction, and impotence.

The use of absorbable gelatin-based hemostatic agents during the repair of dural defects associated with laminectomy and craniotomy operations, has been associated with fever, infection, leg paresthesias, neck and back pain, bladder and bowel incontinence, cauda equina syndrome, neurogenic bladder, impotence, and paresis.

The use of absorbable gelatin-based hemostatic agents has been associated with paralysis, due to device migration into foramina in the bone around the spinal cord, and blindness, due to device migration in the orbit of the eye, during lobectomy, laminectomy, and repair of a frontal skull fracture and lacerated lobe.

Foreign body reactions, "encapsulation" of fluid, and hematoma have been observed at implant sites.

Excessive fibrosis and prolonged fixation of a tendon have been reported when absorbable gelatin-based sponges were used in severed tendon repair.

Toxic shock syndrome was reported in association with the use of absorbable gelatin-based hemostats in nasal surgery.

Fever, failure of absorption, and hearing loss have been observed when absorbable hemostatic agents were used during tympanoplasty.



# Essential Product Information

## **SURGIFLO® Hemostatic Matrix Kit Essential Product Information (Made from Absorbable Gelatin Sponge, USP) with Thrombin**

### **DESCRIPTION**

SURGIFLO® with Thrombin (SURGIFLO® Hemostatic Matrix Kit) is intended for hemostatic use by applying to a bleeding surface.

### **ACTIONS**

When used in appropriate amounts SURGIFLO® is absorbed completely within 4 to 6 weeks.

### **INTENDED USE/INDICATIONS**

SURGIFLO®, mixed with thrombin solution, is indicated in surgical procedures (other than ophthalmic) as an adjunct to hemostasis when control of bleeding by ligature or other conventional methods is ineffective or impractical.

### **CONTRAINDICATIONS**

Do not use SURGIFLO® in intravascular compartments because of the risk of embolization.

Do not use SURGIFLO® in patients with known allergies to porcine gelatin.

Do not use SURGIFLO® in closure of skin incisions because it may interfere with the healing of skin edges. This interference is due to mechanical interposition of gelatin and is not secondary to intrinsic interference with wound healing.

### **WARNINGS**

SURGIFLO® should not be used in the presence of infection and should be used with caution in contaminated areas of the body.

SURGIFLO® should not be used in instances of pumping arterial hemorrhage. SURGIFLO® will not act as a tampon or plug in a bleeding site.

SURGIFLO® should be removed from the site of application when used in, around, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasm because it may swell resulting in nerve damage.

Excess SURGIFLO® should be removed once hemostasis has been achieved.

The safety and effectiveness of SURGIFLO® for use in ophthalmic procedures has not been established.

SURGIFLO® should not be used for controlling post-partum intrauterine bleeding or menorrhagia.

The safety and effectiveness of SURGIFLO® has not been established in children and pregnant women.

The blue flexible applicator tip should not be trimmed to avoid exposing internal guidewire.

The white straight applicator tip should be trimmed away from the surgical area. Cut a square angle to avoid creating a sharp tip.

### **PRECAUTIONS**

Safe and effective use of SURGIFOAM® Sponge has been reported in a published neurologic retrospective study involving 1700 cases in Europe. Safe and effective use in neurosurgery has not been proven through randomized, controlled clinical studies in the United States.

SURGIFLO® is supplied as a sterile product and cannot be resterilized.

SURGIFLO® should not be used for packing unless excess product that is not needed to maintain hemostasis is removed.

SURGIFLO® may swell up to 20% upon contact with additional fluid.

SURGIFLO® should not be used in conjunction with autologous blood salvage circuits.

SURGIFLO® should not be used in conjunction with methylmethacrylate adhesives.

In urological procedures, SURGIFLO® should not be left in the renal pelvis or ureters to eliminate the potential foci for calculus formation.

# Essential Product Information

## **SURGIFLO® Hemostatic Matrix Kit Essential Product Information (Made from Absorbable Gelatin Sponge, USP) with Thrombin**

### **ADVERSE EVENTS**

A total of 142 patients received SURGIFOAM® Sponge during a clinical trial comparing SURGIFOAM® Sponge to another absorbable gelatin sponge. In general, the following adverse events have been reported with the use of absorbable porcine gelatin-based hemostatic agents:

Gelatin-based hemostatic agents may serve as a nidus for infection and abscess formation and have been reported to potentiate bacterial growth.

Giant cell granulomas have been observed at implant sites when used in the brain.

Compression of the brain and spinal cord resulting from the accumulation of sterile fluid have been observed.

Multiple neurologic events were reported when absorbable gelatin-based hemostatic agents were used in laminectomy operations, including cauda equina syndrome, spinal stenosis, meningitis, arachnoiditis, headaches, paresthesias, pain, bladder and bowel dysfunction, and impotence.

The use of absorbable gelatin-based hemostatic agents during the repair of dural defects associated with laminectomy and craniotomy operations, has been associated with fever, infection, leg paresthesias, neck and back pain, bladder and bowel incontinence, cauda equina syndrome, neurogenic bladder, impotence, and paresis.

The use of absorbable gelatin-based hemostatic agents has been associated with paralysis, due to device migration into foramina in the bone around the spinal cord, and blindness, due to device migration in the orbit of the eye, during lobectomy, laminectomy, and repair of a frontal skull fracture and lacerated lobe.

Foreign body reactions, "encapsulation" of fluid, and hematoma have been observed at implant sites.

Excessive fibrosis and prolonged fixation of a tendon have been reported when absorbable gelatin-based sponges were used in severed tendon repair.

Toxic shock syndrome was reported in association with the use of absorbable gelatin-based hemostats in nasal surgery.

Fever, failure of absorption, and hearing loss have been observed when absorbable hemostatic agents were used during tympanoplasty.

# Essential Product Information

## **EVITHROM®** Thrombin, Topical (Human) for Topical Use Only Lyophilized Powder for Solution

EVITHROM® is a topical thrombin indicated as an aid to hemostasis whenever oozing blood and minor bleeding from capillaries and small venules is accessible and control of bleeding by standard surgical techniques (such as suture, ligature or cautery) is ineffective or impractical. EVITHROM® may be used in conjunction with an Absorbable Gelatin Sponge, USP.

### **IMPORTANT SAFETY INFORMATION**

- For topical use only.
- Do not inject.
- Apply EVITHROM® on the surface of bleeding tissue only.
- The amount of EVITHROM® required depends upon the area of tissue to be treated and the method of application. In clinical studies, volumes up to 10 ml were used in conjunction with Absorbable Gelatin Sponge.
- Do not use for the treatment of severe or brisk arterial bleeding.
- Do not use in individuals known to have anaphylactic or severe systemic reaction to human blood products. Hypersensitivity reactions, including anaphylaxis, may occur.
- There is a potential risk of thrombosis if absorbed systemically.
- May carry a risk of transmitting infectious agents such as viruses and theoretically, the Creutzfeldt-Jakob disease (CJD) agent, despite manufacturing steps designed to reduce the risk of viral transmission.
- The most common adverse reactions during clinical trial (reported in at least 2% of subjects treated with EVITHROM®) were prolonged activated partial thromboplastin time, increased INR, decreased lymphocyte count, prolonged prothrombin time and increased neutrophil count.
- None of the patients treated with EVITHROM developed antibodies to human thrombin or to human Factor V/Va. The clinical significance of these findings is unknown.

For complete indications, contraindications, warnings, precautions, and adverse reactions, please reference full package insert.

# Essential Product Information

## **EVARREST® Fibrin Sealant Patch** **IMPORTANT SAFETY INFORMATION**

### INDICATIONS AND USAGE

EVARREST® is a fibrin sealant patch indicated for use with manual compression as an adjunct to hemostasis in adult patients undergoing surgery, when control of bleeding by conventional surgical techniques (such as suture, ligation, and cautery) is ineffective or impractical.

### LIMITATIONS FOR USE

- Cannot be used in place of sutures or other forms of mechanical ligation in the treatment of major arterial or venous bleeding.
- Not for use in children under one month of age.
- Laparoscopic and other minimally invasive surgeries where manual compression would be difficult to achieve.

### IMPORTANT SAFETY INFORMATION

For topical use only. Apply immediate manual compression over the entire surface of the patch and maintain contact pressure for 3 minutes to control the bleeding.

Do not apply intravascularly. This can result in life threatening thromboembolic events.

Do not use to treat bleeding from large defects in arteries or veins where the injured vascular wall requires conventional surgical repair and maintenance of vessel patency or where there would be persistent exposure of EVARREST® to blood flow and/or pressure during absorption of the product. Thrombosis can occur if absorbed systemically.

Do not use in individuals known to have anaphylactic or severe systemic reaction to human blood products. EVARREST® can cause hypersensitivity reactions including anaphylaxis.

Avoid application to contaminated areas of the body or in the presence of active infection. Infection can occur.

EVARREST® contains oxidized regenerated cellulose which adheres to bleeding surfaces. Inadvertent adhesions can occur.

Avoid use in, around, or in proximity to, foramina in bone or areas of bony confines where swelling may cause compression.

Use the least number of patches required to cover the entire bleeding area. Portions of excess patch material can become dislodged and migrate to other areas of the body.

Do not use more than eight 2x4 inch (5.1 x 10.2 cm) or more than four 4x4 inch (10.2 x 10.2 cm) patches.

Use in patients who have been previously exposed to EVARREST® has not been studied.

May carry a risk of transmitting infectious agents, e.g., viruses, the variant Creutzfeldt-Jakob disease (vCJD) agent and, theoretically, the Creutzfeldt-Jakob disease (CJD) agent.

The adverse reactions reported during clinical trials occurred in less than 1% of all cases and included deep venous thrombosis, pulmonary embolism, blood fibrinogen increase, anastomotic hemorrhage, post procedural and intra-abdominal hemorrhage, abdominal distension, anemia, gastrointestinal hemorrhage, thoracic cavity drainage, pleural effusion, abdominal abscess, ascites, localized intra-abdominal fluid collection, cardiac failure, operative hemorrhage, and ischemic bowel.

Pediatrics: Use in children under the age of one month may be unsafe or ineffective due to small size and limited ability to apply the patch as recommended.

Please see package insert for EVARREST® Full Prescribing Information.

To report SUSPECTED ADVERSE REACTIONS, contact ETHICON Customer Support Center at 1-877-384-4266 or FDA at 1-800-FDA-1088 or [www.fda.gov/medwatch](http://www.fda.gov/medwatch).

# Essential Product Information

## EVICEL® Fibrin Sealant (Human) IMPORTANT SAFETY INFORMATION

### INDICATIONS

EVICEL® Fibrin Sealant (Human) is indicated as an adjunct to hemostasis for use in patients undergoing surgery, when control of bleeding by standard surgical techniques (such as suture, ligation, or cautery) is ineffective or impractical.

### CONTRAINDICATIONS

Do not inject directly into the circulatory system. Intravascular application of EVICEL® may result in life-threatening thromboembolic events.

Do not use in individuals known to have anaphylactic or severe systemic reaction to human blood products.

Do not use for the treatment of severe or brisk arterial bleeding.

Do not use EVICEL® for spraying in endoscopic or laparoscopic procedures where the minimum recommended distance from the applicator tip to the target site cannot be ensured.

### WARNINGS AND PRECAUTIONS

Life-threatening gas embolism has occurred with the use of spray devices employing pressure regulator to administer fibrin sealants. These events appear to be related to the use of the spray device at pressures higher than recommended and/or at distances closer than recommended to the surface of the tissue. Follow labeled application instructions regarding pressure range and distance when using a spray device and monitor patients for the possibility of gas embolism.

Monitor changes in blood pressure, pulse oxygen saturation, and end-tidal CO<sub>2</sub> when spraying EVICEL® because of the possibility of gas embolism.

To reduce the risk of potentially life-threatening gas embolism, spray EVICEL® using only pressurized CO<sub>2</sub> gas at the recommended pressures and distances.

Use EVICEL® spray application only if it is possible to accurately judge the spray distance, especially during endoscopic or laparoscopic procedures. Apply as a thin layer.

Prior to applying EVICEL®, dry surface areas of the wound by standard techniques (e.g. intermittent application of compresses, swabs, use of suction devices). Prepare and administer EVICEL® according to the instructions and only with devices recommended for this product.

May carry a risk of transmitting infectious agents, e.g. viruses, the variant Creutzfeldt-Jakob disease (vCJD) agent and theoretically, the Creutzfeldt-Jakob disease (CJD) agent.

Most common adverse events reported in clinical trials (≥5%) are bradycardia, nausea, hypokalemia, insomnia, hypotension, pyrexia, graft infection, vascular graft occlusion, peripheral edema, and constipation.

### ADVERSE REACTIONS

The type and frequency of adverse events reported are consistent with events typically seen following surgery. Postsurgical adhesions may occur in the presence of GYNECARE INTERCEED®.

For complete indications, contraindications, warnings, precautions, and adverse reactions, please reference full package insert.

# Essential Product Information

## GYNECARE INTERCEED® ABSORBABLE ADHESION BARRIER

### INDICATIONS

GYNECARE INTERCEED® Absorbable Adhesion Barrier is indicated as an adjuvant in open Claparotomy) gynecologic pelvic surgery for reducing the incidence of postoperative pelvic adhesions after meticulous hemostasis is achieved consistent with microsurgical principles.

### CONTRAINDICATIONS

The use of GYNECARE INTERCEED® is contraindicated in the presence of frank infection. GYNECARE INTERCEED® is not indicated as a hemostatic agent. Appropriate means of achieving hemostasis must be employed.

### WARNINGS

The safety and effectiveness of GYNECARE INTERCEED® Adhesion Barrier in laparoscopic surgery or any procedures other than open (laparotomy) gynecologic microsurgical procedures have not been established.

Postoperative adhesions may be induced by GYNECARE INTERCEED® application if adjacent tissues (eg, ovary and tube) and structures are coapted or conjoined by the device, or if GYNECARE INTERCEED® is folded, wadded or layered. Postoperative adhesions may occur in the presence of GYNECARE INTERCEED® if meticulous hemostasis is not achieved prior to application. As with all foreign substances, GYNECARE INTERCEED® should not be placed in a contaminated surgical site.

### PRECAUTIONS

Use only a single layer of GYNECARE INTERCEED®, since multiple layers of packing or folding will not enhance the adhesion barrier characteristics and may interfere with the absorption rate of GYNECARE INTERCEED®. Care should be exercised in applying GYNECARE INTERCEED® to a pelvic organ not to constrict or restrict it. If the product comes in contact with blood prior to completing the procedure, it should be discarded, as fibrin deposition cannot be removed by irrigation and may promote adhesions formation. Ectopic pregnancies have been associated with fertility surgery of the female reproductive tract. No data exist to establish the effect, if any, of GYNECARE INTERCEED® on the occurrence of ectopic pregnancies. No adequate studies have been conducted in women who have become pregnant within the first month after exposure to GYNECARE INTERCEED®. No teratogenic studies have been performed. Therefore, avoidance of conception should be considered during the first complete menstrual cycle after use of GYNECARE INTERCEED®. The safety and effectiveness of using GYNECARE INTERCEED® in combination with other adhesion prevention treatments have not been clinically established. GYNECARE INTERCEED® is supplied sterile. As the material is not compatible with autoclaving or ethylene oxide sterilization, GYNECARE INTERCEED® must not be resterilized. Foreign body reactions may occur in some patients. Interactions may occur between GYNECARE INTERCEED® and some drugs used at the surgical site. Pathologists examining sites of GYNECARE INTERCEED® placement should be made aware of its usage and of the normal cellular response to GYNECARE INTERCEED® 'to facilitate proper evaluation of specimens'.

### ADVERSE REACTIONS

The type and frequency of adverse events reported are consistent with events typically seen following surgery. Postsurgical adhesions may occur in the presence of GYNECARE INTERCEED®.