



# CIDEX® Tray Systems

## INSTRUMENT HANDLING

CIDEX® Tray Systems have been designed by ASP to provide healthcare professionals with the ideal complementary accessories for manual cleaning/high-level disinfection. Manual processes have never been so easy.

## EXPERIENCE THE BENEFITS

CIDEX® Tray Systems offer a large choice of premium durable trays.

- **Easy and practical handling of medical instruments**
- **Ideal complementary accessories for manual disinfection**
- **5 shapes and sizes to meet every user preference**



# **CIDEX**® Tray Systems

## FEATURES & BENEFITS

- **FLEXIBLE UTILIZATION** – can be used for either manual cleaning or high-level disinfection stages.
- **EASY PERIODIC STERILIZATION** – heat resistant up to 134°C.
- **LONG LIFE TIME** – superior durability of glass filled polypropylene.
- **REDUCTION OF VAPOR AND FUMES ESCAPE** with their unique rounded lip design.
- **EASY CLEANING** with their wide radius corner.
- **EASY SOLUTION MONITORING** – reusable labels help monitor the ASP High-level disinfection solutions reuse.

## TECHNICAL INFORMATION

MEDICAL DEVICE CLASSIFICATION    Class I according to MDD 93/42/EEC

## ORDERING INFORMATION

REORDER NO.	DESCRIPTION AND DIMENSIONS	CONTENT
82010	344 mm x 177 mm x 132 mm 3.6 liters	1 tray, lid & liner
82016	503 mm x 186 mm x 132 mm 8.1 liters	1 tray, lid & liner
82027	740 mm x 220 mm x 132 mm 14.1 liters	1 tray, lid & liner
82032	449 mm diameter by 222 mm 17.5 liters	1 tray & lid
82076	594 mm x 381 mm x 141 mm 19.8 liters	1 tray, lid & liner



**FOR MORE INFORMATION, CONTACT YOUR LOCAL ADVANCED STERILIZATION PRODUCTS SALES REPRESENTATIVE**

[www.asppj.com/emea](http://www.asppj.com/emea)

### UK:

Advanced Sterilization Products  
A Division of J&J Medical Ltd.  
Pinewood Campus  
Nine Mile Ride  
Wokingham  
Berkshire RG40, 3EW, England  
T: +44 1344 871 081  
F: +44 1344 871 171

### Ireland:

Advanced Sterilization Products  
A Division of J&J Medical Ireland  
Airtown Road, Tallaght,  
Dublin 24, Ireland  
T: +353 1 466 5200  
F: +353 1 466 5340

### Middle East:

Advanced Sterilization Products  
A Division of J&J Medical Middle East  
Dubai Healthcare City, J&J building,  
3rd floor  
PO Box 505080, Dubai,  
United Arab Emirates  
T: +971 4 429 7235  
F: +971 4 429 7250

### South Africa:

Advanced Sterilization Products  
A Division of J&J Medical Pty.  
Ltd. SA  
PO Box 273, Midrand  
Halfway House  
1685  
South Africa  
T: +27 11 265 1120  
F: +27 11 265 1189

**ADVANCED STERILIZATION PRODUCTS**

Division of Cilag GmbH International  
a **Johnson & Johnson** company