

COMFORT.
CONVENIENCE.
SUPERIOR SHIELDING.



Protective eyewear and face shields offer a lightweight, economical eye protection solution against the splash and splatter of blood-borne pathogens and other potentially infectious materials as well as aerosol sprays.



FOR USE IN

- Prophylaxis
- Fillings
- Extractions
- Oral surgery
- Crowns and bridges
- Endodontic

ASSEMBLE 'N GO™ EYE SHIELD



- Separately available replacement lenses for maximized cost-effective replenishment
- Extremely lightweight and comfortable—weighs only 3/10 of an ounce!
- “Wrap-around” design provides both front and side eye protection
- Fits comfortably with a surgical mask and over prescription eyewear
- Frames are autoclavable and reusable
- Latex-free
- Patented design
- CE certified
- Meets requirements for both OSHA and EEC mandates

ASSEMBLE 'N GO™ FACE SHIELD



- Separately available replacement lenses for maximized cost-effective replenishment
- Lighter and more comfortable than other face shields—weighs less than 1.5 ounces!
- Designed to maximize airflow and prevent fogging
- Provides over-the-top, side and front face protection
- Disposable, optical grade, anti-static replacement shields
- Visors are autoclavable and reusable
- Includes a clear, snap-in nosepiece for support
- Latex-free
- Patented design
- CE certified
- Meets requirements for both OSHA and EEC mandates

ITEM NO.	DESCRIPTION
9210-100-6	Replacement lenses for use with 9452-10 frames (100 lenses)
9452-10	Combo Pack: Contains 5 assorted frames/ 25 eye shields
9452P-10	Combo Pack: Contains 5 pink frames/ 25 eye shields

ITEM NO.	DESCRIPTION
2210-40-6	Replacement shields for use with 2452-6 visors (40 shields)
2210-100-5	Replacement shields for use with 2452-6 visors (100 shields)
2452-6	Combo Pack: Contains 3 assorted visors/ 25 face shields

For a no-charge sample, call **1.800.521.1314** or visit **www.TIDIproducts.com**

See reverse side for Additional Eye Protection Products offered by TIDI Products.





TIDISHIELD GRAB 'N GO™ TOWER EYE SHIELDS DISPENSER

Part No. 9210A-100

25 assembled eye shields per dispenser
 Fits in Dispenser Holder Part No. 7003
 6.75" W x 17" H x 5.75" D



TIDISHIELD ASSEMBLE 'N GO™ EYE SHIELDS DISPENSER

Part No. 9209



100 frames and 100 lenses per dispenser
 Fits in Dispenser Holder Part No. 7001
 9.0" W x 8.25" H x 5.75" D

EYEWEAR ACCESSORIES

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
9810-50	Side Shield Clear	One size	50 Side Shields (25 pairs)	Disposable Side Shields, use with prescription glasses

FLIP 'N GO™ EYE SHIELD

- Perfect for patient or provider
- FLAT for compact storage, FLIP to wear
- Extremely lightweight and comfortable
- Wrap-around design for front and side eye protection
- Optical-grade, anti-static, lens




CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
9300FE-50	5 assorted colored frames/clear lens	One size fits most	50 (2 boxes of 25)	One piece disposable eyewear; stores flat and assembles with one simple step


PRODUCT ASSEMBLY INSTRUCTIONS

TidiShield™ Assemble 'n Go™ Shields


1 Slide frame (flat side down) through inside holes on lens.*



2 Slide frame through outside holes on lens.



3 Slide lens forward onto notches on frame.




Made in USA


*Images of Instructions are a trademark of TIDI Products, LLC

TidiShield™ Assemble 'n Go™ Face Shields

1 Insert nosepiece tab into visor slot, so nosepiece slants away from front of visor. (Nosepiece may be removed to wear over eyeglasses.)



2 Starting on one side of the visor, first mount the faceshield onto the visor pin (1), then insert the visor's support knuckle (2) through the faceshield slot.



Selection of mounting holes provides 3 options for adjusting the faceshield angle. Top hole will position the faceshield farther from chin, and bottom hole will position it closer to chin.

3 Repeat step 2 on the opposite side of the visor, pulling the faceshield tight and securing it to the opposite support knuckle and pin.

