



Chemotherapy Personal Protection Equipment (PPE)



P2 Personal Protection Gowns — Ultimate Protection

P2® **HiRisk™** Chemotherapy Gowns provide maximum protection while preparing or administering chemotherapy drugs.

- Closed back design for ultimate protection
- Meets OSHA, ASHP and ONS recommendations for wearing gowns when handling chemotherapy drugs
- Polyethylene coated front and sleeves
- Non-linting materials
- Long sleeves with knit cuffs
- Ties at waist with Velcro closure at neck



| CODE | DESCRIPTION | SIZE | COLOR | CASE QTY |
|------|-------------------------------|-----------|-------|----------------------------------|
| 8586 | P2® HiRisk™ Chemotherapy Gown | Universal | Blue | 10 gowns/poly bag; 30 gowns/case |
| 8587 | P2® HiRisk™ Chemotherapy Gown | X-Large | Blue | 10 gowns/poly bag; 30 gowns/case |
| 8588 | P2® HiRisk™ Chemotherapy Gown | XXL | Blue | 10 gowns/poly bag; 30 gowns/case |

P2 Personal Protection Gowns — Everyday Protection

P2® **HiRisk™** Impervious Gown and the **SafetyPlus™** Poly Coated Gown have an open back design for increased airflow.

- Open back design
- Non-linting materials
- Closed neck loop for easy donning and tear off
- Stirrup cuffs for wrist protection



#8576

- 100% embossed polyethylene
- High level of protection against hazardous materials



#8578


- White polyethylene coated polypropylene
- Light breathable materials for increased comfort

| CODE | DESCRIPTION | SIZE | COLOR | CASE QTY |
|------|----------------------------------|-----------|------------|-----------------------------|
| 8576 | P2® HiRisk™ Impervious Gown | Universal | Light Blue | 15 gowns/box; 75 gowns/case |
| 8578 | P2® SafetyPlus™ Poly Coated Gown | X-Large | White | 8 gowns/box; 40 gowns/case |

P2® Chemotherapy Prep Mat

P2® **SafetyPlus™** Chemotherapy Prep Mat is designed to contain the accidental spills that occur during chemotherapy drug preparation.

- Multiple layers of absorbency
 - > Top - specially designed wicking layer
 - > Middle - three ply absorbent layer
 - > Bottom - leak proof poly backing
- Poly bagged within carton to keep mat clean from lint, fiber and dust
- Economical option
- Latex free



| CODE | DESCRIPTION | SIZE | COLOR | CASE QTY |
|--------|---------------------------------------|-----------|------------|----------|
| 983930 | P2® SafetyPlus™ Chemotherapy Prep Mat | 13" x 18" | Light Blue | 250 |

P2® Gloves - Nitrile

P2® **HiRisk™ Nitrile** Gloves deliver superior protection for mixing and administering chemotherapy drugs and handling chemotherapy spills.

- Chemotherapy tested using ASTM D6978-05 standard
- 8 mils thick at fingertips
- Ambidextrous for quick and easy donning
- Textured for added grip
- Beaded cuff
- 12 inches long for added protection
- Convenient dispenser box



| CODE | DESCRIPTION | SIZE | COLOR | CASE QTY |
|----------|---------------------------|---------|------------|--|
| 940001-1 | P2® HiRisk™ Nitrile Glove | Small | Light Blue | 10 boxes of 50 gloves; 500 gloves/case |
| 940002-1 | P2® HiRisk™ Nitrile Glove | Medium | Light Blue | 10 boxes of 50 gloves; 500 gloves/case |
| 940003-1 | P2® HiRisk™ Nitrile Glove | Large | Light Blue | 10 boxes of 50 gloves; 500 gloves/case |
| 940004-1 | P2® HiRisk™ Nitrile Glove | X-Large | Light Blue | 10 boxes of 50 gloves; 500 gloves/case |

P2® Gloves - Latex

P2® **HiRisk™ Latex** Gloves provide the same maximum protection as the nitrile version when handling and administering chemotherapy drugs.

- Chemotherapy tested using ASTM D6978-05 standard
- 18 mils thick at fingertips
- Fitted for right and left hands, reduces hand fatigue
- Textured for added grip
- Beaded cuff
- 12 inches long for added protection
- Convenient dispenser box



Extra Long Protection: Measures 12 inches from tip to cuff for full wrist protection

| CODE | DESCRIPTION | SIZE | COLOR | CASE QTY |
|----------|-------------------------|---------|-------|--|
| 932477-1 | P2® HiRisk™ Latex Glove | Small | Blue | 10 boxes of 50 pair; 1,000 gloves/case |
| 932480-1 | P2® HiRisk™ Latex Glove | Medium | Blue | 10 boxes of 50 pair; 1,000 gloves/case |
| 932483-1 | P2® HiRisk™ Latex Glove | Large | Blue | 10 boxes of 50 pair; 1,000 gloves/case |
| 932486-1 | P2® HiRisk™ Latex Glove | X-Large | Blue | 10 boxes of 50 pair; 1,000 gloves/case |



P2® Cabinets | Stations

TIDI's offering of P2® Cabinets and Stations mount on walls, doors, and glass. They provide convenient point-of-use access to personal protective equipment (PPE) that help comply with OSHA guidelines for hospital preparedness. Other standards that are met are JCAHO, DOH, and Fire Safety requirements. These products eliminate the need for bulky carts and other obstructive storage units.

P2® Storage Units

- Keep personal protection and patient care equipment in a clean, organized, and accessible space
- Constructed with ABS plastic which is durable, mark-resistant, and easy to clean
- Storage units are all compatible with P2®'s Glove dispensers. Dispensers provide added space and easy accessibility to gloves

For Glass, Door, or Wall Mounting



Product #8468

P2® ICU Isolation Station

1 station/case | 20.5"H x 6.75"D x 17.75"W

- | | | |
|---|--|--|
| 1 | | glass-mount bracket with adhesive strips |
| 4 | | adjustable shelves |
| 2 | | pockets in door |
| 2 | | storage bins |
| 1 | | ambu bag shelf |
| 1 | | glove dispenser |



For door mounting, add item #8566

For Door Mounting



Product #8569

P2® Isolation Station

1 station/case
20.5"H x 6.75"D x 17.75"W

- | | | |
|---|--|--------------------|
| 4 | | adjustable shelves |
| 2 | | pockets in door |
| 2 | | storage bins |
| 1 | | ambu bag shelf |
| 1 | | glove dispenser |



For Wall or Glass Mounting



Product #8565

P2® Station

1 station/case
20.5"H x 6.75"D x 17.75"W

- 4 | adjustable shelves
- 2 | pockets in door
- 1 | storage bin
- 1 | ambu bag shelf
- 1 | station/case

For glass mounting,
add item #8461



For Wall Mounting



Product #8560

P2® Cabinet

1 cabinet/case
12.75"H x 5.5"D x 13.5"W

- 3 | adjustable shelves
- 2 | pockets in door
- 1 | station/case



Brackets



Product #8563

Vertical Extension Bracket for P2® Isolation Station

2 vertical extension brackets/case | 17.75"

Extension hardware accommodates taller hospital
doors so PPE is more easily reached by staff.

For use with door mount isolation station #8569



Product #8564

Horizontal Extension Bracket for P2® Isolation Station

1 horizontal extension bracket/case | 36"

Extension hardware accommodates most hospital
doors with built-in door damper slots.

For use with door mount
isolation station #8569




Product #8566

Door Hang Bracket For P2® Isolation Station

1 door hang bracket

Bracket converts #8468 and #8565 stations
to allow for door mounting.



Product #8461

P2® Glass-Mount Bracket

2/case | 3"H x 12"W

- 2 | glass-mount brackets
with adhesive strips

For use with wall mount station #8565

Glove Dispensers



Product #8555-SA

P2® Side Loaded Glove Dispenser

For use with any station or
cabinet. Can be stacked for
multiple glove dispensing.

10 glove brackets/case
7"H x 4.5"D x 12"W



Product #8562

P2® Glove Dispenser

For use with any station or
cabinet or as stand-alone
dispenser with double glove
bracket #8554-B.

10 glove dispensers/case
7"H x 4.5"D x 12"W



Product #8554-B

P2® Double Glove Bracket

Stand-alone bracket for
double glove dispensing
with dispenser #8562.

5 glove brackets/case
3.5"H x 4.5"D x 12"W

Pillowcases

Create a comfortable patient environment while providing reliable hygienic barrier protection.

- Spunbond Nonwoven
- Fabricel® Tissue/Poly-Back
- Tissue/Poly-Back



| CODE | DESCRIPTION | SIZE | COLOR | CASE QTY |
|--------|---|-----------|-------|----------|
| 980915 | Spunbond Nonwoven, Strong & Soft | 21" x 30" | White | 100 |
| 919353 | Fabricel® Tissue, Soft & Tear-Resistant, Leak-Proof | 21" x 30" | Blue | 100 |
| 919355 | Fabricel® Tissue, Soft & Tear-Resistant, Leak-Proof | 21" x 30" | White | 100 |
| 919356 | Fabricel® Tissue, Soft & Tear-Resistant, Leak-Proof | 21" x 30" | Mauve | 100 |
| 919365 | Tissue/Poly, Comfortable & Reliable | 21" x 30" | White | 100 |
| 919363 | Tissue/Poly, Comfortable & Reliable | 21" x 30" | Blue | 100 |

Biohazard Waste Bags

Assure proper disposal of infectious waste with these durable bags.

- Generous size options
- Red for easy identification



| CODE | DESCRIPTION | SIZE | COLOR | CASE QTY |
|------|-------------------------------------|-----------|-------|----------|
| 8653 | High Density, 12 Microns, 10 Gallon | 24" x 24" | Red | 250 |
| 8654 | High Density, 12 Microns, 16 Gallon | 34" x 37" | Red | 250 |
| 8656 | High Density, 17 Microns, 33 Gallon | 33" x 40" | Red | 250 |
| 8657 | High Density, 17 Microns, 40 Gallon | 40" x 48" | Red | 250 |



TIDI Products, LLC
570 Enterprise Drive
Neenah, WI 54956

Customer Service
800.521.1314

FAX 800.837.7770

www.tidiproductions.com