



# ACUTE CARE CONSUMABLES

Supporting Caregivers and  
Preventing Infections™



Call your TIDI Account  
Representative for  
more information -  
**1.800.521.1314** or visit  
**[www.tidiproductions.com](http://www.tidiproductions.com)**



## MAINTAIN STERILITY – DEFEND THE STERILE FIELD



## C-ARMOR®

C-Armor® is an expandable and collapsible sterile pouch that adheres to the patient drape above the sterile field line. It permits unlimited fluoroscopic imaging in the horizontal plane, maintains the sterile field and lies flat to the table when not in use, permitting unencumbered surgical access.

### C-ARMOR®

FOR USE WITH	CODE	ITEM	QUANTITY
C-Arms	5523	C-Armor®	Box of 12
	5536	C-Armor® Kit (includes Top C-Arm Cover)	Box of 12



## STERILE-Z® BACK TABLE COVER

Sterile-Z® Back Table Cover is designed to maintain sterility of surgical instruments and implants prior to use without risk of contaminating the sterile field upon removal. Sterile-Z separates from the center with its patent-pending, sterility-protecting “Z-Fold.”

### STERILE-Z® BACK TABLE COVER

FOR USE IN	CODE	ITEM	QUANTITY
Back Tables that have been prepared and will not be immediately used	5575	Sterile-Z® Back Table Cover	Box of 20
	5575-XL	Sterile-Z® Back Table Cover X-Large	Box of 20
	5582	Sterile-Z® Mayo Stand Cover	Box of 20



## STERILE-Z® PATIENT DRAPE

Sterile-Z® Patient Drape is designed to protect the intraoperative patient during 3-D radiographic imaging and maintain sterility during removal. Sterile-Z separates quickly following imaging with its patent-pending, sterility-protecting “Z-Fold.”

### STERILE-Z® PATIENT DRAPE

FOR USE IN	CODE	ITEM	QUANTITY
Intraoperative 3-D radiographic imaging	5574	Sterile-Z® Patient Drape	Box of 10

# MAINTAIN STERILITY – DEFEND THE STERILE FIELD




## NAVI-CRANI™

Navi-Crani™ is a sterile patient drape designed to reduce the risk of contamination during navigated craniotomy procedures. Navi-Crani™ is installed without cutting to fit over the navigation-system arm, reducing the risk for contamination.

### NAVI-CRANI™

FOR USE WITH	CODE	ITEM	QUANTITY
Navigation systems, including Brainlab, Stryker®, Medtronic	5568	Navi-Crani™ Patient Drape	Box of 8



## PURE-FIT™ O-ARM® COVER

Pure-Fit™ O-Arm® cover is easily deployed to maintain sterility throughout the use of the O-Arm®. It reduces the risk of sterile field breach with an intuitive fold and closure method for removing an O-Arm® from the surgical field.

### PURE-FIT™

FOR USE WITH	CODE	ITEM	QUANTITY
O-Arm®	5565	O-Arm® Tube Cover	Box of 20



## BACK TABLE CONTAINMENT COVER

### Reduce Contamination Risk

TIDIShield® Back Table Containment Cover fits standard back tables used for transport. It is ideal for any surgical procedure requiring extra instrumentation or instrumentation that will not fit on the case cart post-case.

### TIDISHIELD® BACK TABLE CONTAINMENT COVER

FOR USE WITH	CODE	ITEM	QUANTITY
Surgical procedures requiring extra instrumentation	5540	Back Table Containment Cover	Box of 20
	5540-XL	Back Table Containment Cover X-Large	Box of 20





## Protective Eyewear Dispensers Created Specifically for Convenience and Accessibility

### READILY AVAILABLE AT POINT-OF-USE FOR INCREASED COMPLIANCE

- Eye Shield Dispensers are easily portable so eye protection is always available
- Eyewear is ready to use—no need to search for lenses and frames
- Position the dispenser at point-of-use where gloves and masks are dispensed using a plastic dispenser holder (no-holes mounting using 3M Command™ Strips)



#### TIDISHIELD® GRAB 'N GO™ EYE SHIELDS DISPENSER (9210A-100)

25 assembled eye shields per dispenser, 4 dispensers per case  
Fits in Dispenser Holder 7003  
6.75" W x 17" H x 5.75" D

#### HORTON HEARS A WHO DISPENSER† (9210A-100HHW)



#### TIDISHIELD® EYE SHIELD DISPENSER HOLDER FOR TIDISHIELD® GRAB 'N GO™ (7003)

2 dispenser holders per case  
Can be mounted with 3M™ Command™ Strips  
6" W x 12" H x 7.5" D  
Metal holder for carts also available (7004)



#### TIDISHIELD® FLIP 'N GO™ EYE SHIELDS DISPENSER (9300FE-50)

50 individually bagged eye shields (2 boxes of 25)  
Fits in Dispenser Holder 7001  
9.25" W x 7.625" H x 4.625" D



#### TIDISHIELD® EYE SHIELD DISPENSER HOLDER FOR TIDISHIELD® FLIP 'N GO™ (7001)

2 dispenser holders per case  
Can be mounted with 3M™ Command™ Strips  
9.75" W x 8.375" H x 6.3125" D



### TIDISHIELD® FLIP 'N GO™ EYE SHIELD

- No assembly required
- FLAT for compact storage, FLIP to wear
- Extremely lightweight and comfortable
- Wrap-around design for front and side eye protection
- Optical-grade, anti-static lens



## Protection is Essential

Healthcare professionals work under hazardous and ever-changing conditions. What healthcare workers need is superior face protection from the splash and splatter of blood-borne pathogens and other potentially infectious materials. TIDI® Products offers a product perfect for use under any circumstances – TIDIShield® Grab 'n Go™ Disposable Face Shields.

### TIDISHIELD® GRAB 'N GO™ DISPOSABLE FACE SHIELDS

DESCRIPTION	CODE	DIMENSIONS	QUANTITY
Disposable, preassembled, optical-grade fog-free face shield with vented foam headband, adjustable elastic strap and over-the-top splash guard	<b>FS40EC</b>	Enhanced Coverage 12" wide x 9.5" long	40 (1 box of 40) Bulk
	<b>FS96EC</b>	Enhanced Coverage 12" wide x 9.5" long	95 (4 Dispensers of 24)
	<b>FS40FC</b>	Full Coverage 12" wide x 7.5" long	40 (1 box of 40) Bulk
	<b>FS96FC</b>	Full Coverage 12" wide x 7.5" long	95 (4 Dispensers of 24)

### MICROSCOPE DRAPES

(MOST POPULAR SIZES SHOWN - CALL TIDI® PRODUCTS FOR OTHER SIZES)

FOR USE WITH	CODE	DIMENSIONS	QUANTITY	LENS SIZE
Drapes for Zeiss® MD Microscope • Zeiss® OPMI® MD / CS / 11, 111 • ORL, PRO Magis, OPMI® Neuro • VISU 140, 150, 160, 200, 210 • Movena, Sensera, MKM • PROergo, Pentero, Moller-Wedel	<b>206550</b>	20.0"W x 65.0"L	20/case	50mm lens
	<b>206565</b>	20.0"W x 65.0"L	20/case	65mm lens
	<b>426568</b>	42.0"W x 65.0"L	10/case	68mm lens
	<b>421250</b>	42.0"W x 120.0"L	10/case	50mm lens
Drapes for Zeiss® OPMI® Microscope • Zeiss® OPMI® 1, 6, 9	<b>446565</b>	44.0"W x 65.0"L	10/case	65mm lens
	<b>441265</b>	44.0"W x 120.0"L	10/case	65mm lens
Drapes for Leica® Wild® Microscope • Leica® / Wild®	<b>541565</b>	54.0"W x 150.0"L	10/case	65mm lens
	<b>541568</b>	54.0"W x 150.0"L	10/case	68mm lens

### ULTRASOUND PROBE COVERS

POLYURETHANE, LATEX-FREE, STERILE WITH LATEX-FREE RUBBER BANDS AND GEL

CODE	DESCRIPTION	END CONFIG	FOLD CONFIG	QUANTITY
<b>5488</b>	US PROBE KIT 5.5" X 24" (14cm X 61cm)	Box End	Flat Fold	20/box
<b>5491</b>	US PROBE KIT 5.2" X 24" (12.7-5cm X 61cm) Tapered	Flat End	Flat Fold	20/box
<b>5492</b>	US PROBE KIT 5.5" X 36" (14cm X 91cm)	Box End	Telescope Fold	20/box
<b>5493</b>	US PROBE KIT 6" X 48" (15.2cm X 121.9cm)	Flat End	Flat Fold	20/box
<b>5529</b>	US PROBE KIT 3" X 24" (7.6cm X 61cm)	Box End	Telescope Fold	20/box
<b>5528</b>	US PROBE KIT 3" X 96" (7.6cm X 243.8cm)	Box End	Telescope Fold	20/box
<b>5527</b>	US PROBE KIT 4" X 24" (10.2cm X 61cm)	Box End	Telescope Fold	20/box
<b>5535</b>	US PROBE KIT 4" X 48" (10.2cm X 121.9cm)	Box End	Telescope Fold	20/box
<b>5530</b>	US PROBE KIT 4" X 96" (10.2cm X 243.8cm)	Box End	Telescope Fold	20/box
<b>5524</b>	US PROBE KIT 5" X 48" (12.7cm X 121.9cm)	Box End	Telescope Fold	20/box
<b>5525</b>	US PROBE KIT 6" X 24" (15.2cm X 61cm)	Box End	Telescope Fold	20/box
<b>5534</b>	US PROBE KIT 6" X 48" (15.2cm X 121.9cm)	Box End	Telescope Fold	20/box
<b>5502</b>	US PROBE KIT 6" X 96" (15.2cm X 243.8cm)	Box End	Telescope Fold	20/box



## OEM SPECIFIC C-ARM COVERS

CODE	PART DESCRIPTION	QUANTITY	APPLICATION
5429	9"/12" I.I. C-Arm Kit [Fits OEC 9600/9800, Siemens ARCADIS Avantec Siemens SIREMOBIL Compact L, Siemens ARCADIS Varic, Philips BV series, Ziehm]	20/case	Consists of (5427) C-Arm Drape with 5 mounting clips, (CFI-657T) I.I. Cover, (CFI-657T) X-Ray Tube Cover and a (5421) Footswitch Cover. Sterile.
5532	9"/12" I.I. C-Arm Kit [Fits OEC 9900]	20/case	Consists of C-Arm drape with clips, I.I. cover, X-Ray Tube Cover, and Footswitch Cover. Sterile.
5464	C-Arm Mini Kit [Fits Siemens SIREMOBILE ISO-C & Siemens ARCADIS Orbic]	20/case	Consists of (5461) C-Arm/I.I. Cover and (5462) X-Ray Tube Cover. Sterile.
5427	C-Bow Drape with Clips	30/case	Consists of C-Arm Drape with 5 mounting clips. Fits OEC 9600/9800 clips. Sterile.
5411	9" I.I. / X-Ray Tube Cover	60/case	30" x 32" Contoured Bandbag. Sterile, individually pouched.
5412	12" I.I. / X-Ray Tube Cover	30/case	36" x 36" Contoured Bandbag. Sterile, individually pouched.
5555	9" x 12" I.I. C-Arm Kit Ziehm Vision	12/case	Consists of 14cm Jett Wing C-Arm Drape; 12" x 18" Footswitch Cover; 30" x 28" Large Generator Cover; and 36" x 28" II Cover with Tape.

## ONE PIECE C-ARM COVERS

CODE	PART DESCRIPTION	QUANTITY	APPLICATION
5414	9" I.I. Universal Halfbag Cover	20/case	36" x 80" designed to fit the I.I. and midway down the C-Arm. Sterile.
5434	Mobile C-Arm Drape	20/case	Clear poly, 42" x 74" with poly straps. Sterile, individually pouched.
5454	Mobile C-Arm Drape	20/case	Clear poly, 54" x 78" with poly straps. Sterile, individually pouched.
5500	Nylon Storage Cover [111" x 131"]	1/case	Nylon. Protect expensive equipment with the C-Arm storage drape. Non-Sterile.

## MINI C-ARM DRAPES

CODE	PART DESCRIPTION	QUANTITY	APPLICATION
5402	Mini C-Arm Kit [Fits OEC Mini 6800/6600]	20/case	Kit consists of a (5403) MiniView C-Arm Drape, (5405) Keyboard Cover and a (5421) Footswitch Cover. Sterile, individually pouched.
5403	Mini C-Arm Drape [Fits OEC Mini 6800/6600]	25/case	Single piece drape to cover I.I., X-ray Tube and C-Arm. Sterile, individually pouched.
5423	Mini C-Arm Drape [Fluoroscanner/Hologic]	20/case	Single piece drape to cover I.I., X-ray Tube and C-Arm. Sterile, individually pouched.







## UROLOGY DISPOSABLES

UROLOGY DISPOSABLES ARE ALSO AVAILABLE IN NON-STERILE

CODE	PART DESCRIPTION	QUANTITY	APPLICATION
5000	Square Drain Bag	20/case	GE/OEC UroView 2000, Liebel Flarsheim, Mallinckrodt Hydrovision ED. Sterile.
5010	Urology Drain Collection Bag	10/case	Collects waste fluid when floor drain or suction is unavailable. Non-Sterile.
5313	Uro-Cather Bag	20/case	Fits all hoop style frames.
5416	Urology Disposables Pack	20/case	Fits OEC UroView 2600/2800. Includes 5 x 10 hand control cover, (5421) 24" x 20" footswitch cover, and (5422) 20 x 36 footswitch cover. Sterile.
5419	Hoop Drain Bag	20/case	GE/OEC UroView 2600/2800.
5442	Urology Drain Bag	20/case	Fits Siemens Uroskop System and Dornier.
5420	Hand Control Cover	25/case	5" x 10" rubber band closer. Bulk packaged. Sterile.



## MISCELLANEOUS DISPOSABLES

FOR USE WITH	CODE	DIMENSIONS	QUANTITY
 <p>Band Bags Clear poly, with a latex-free rubber band. Sterile, individually pouched.</p>	<b>CFI-600</b>	36"W x 28"D Blue	100/case
	<b>CFI-601</b>	36"W x 18"D Blue, CC	100/case
	<b>CFI-603</b>	5"W x 12"D Blue	50/case
	<b>CFI-664</b>	20"W x 15"D	100/case
	<b>CFI-651</b>	20"W x 20"D	100/case
	<b>CFI-661</b>	30"W x 30"D	100/case
	<b>CFI-667</b>	30"W x 36"D	100/case
	<b>CFI-654</b>	36"W x 20"D	100/case
	<b>CFI-655</b>	36"W x 30"D	100/case
	<b>CFI-657</b>	36"W x 36"D	100/case
	<b>CFI-879</b>	36"W x 40"D	100/case
	<b>CFI-761</b>	40"W x 30"D	100/case
	<b>CFI-827</b>	40"W x 40"D	100/case
	<b>CFI-668</b>	42"W x 36"D	100/case
 <p>Contour Band Bags Contoured shape minimizes excess material. Latex-free elastic band provides snug closure around equipment. Sterile.</p>	<b>CFI-6018</b>	18" Depth	100/case
	<b>CFI-6022</b>	22" Depth	100/case
	<b>CFI-6030</b>	30" Depth	50/case
 <p>Dome Bags Shower cap style bag, latex-free elastic opening. Sterile.</p>	<b>CFI-732</b>	12" Depth, 24" Expanded Width	100/case
	<b>CFI-762</b>	Bag, Clear Dome, 20"	100/case
	<b>CFI-786</b>	Bag, Clear Dome, 14"	100/case
	<b>CFI-749</b>	17" Depth, 34" Expanded Width	100/case
	<b>CFI-911</b>	18" Depth, 36" Expanded Width	100/case
	<b>CFI-821</b>	22" Depth, 44" Expanded Width	50/case
	<b>CFI-733</b>	24" Depth, 48" Expanded Width	50/case
	<b>CFI-756</b>	27" Depth, 54" Expanded Width	50/case
	<b>CFI-893</b>	26" Depth, 52" Expanded Width	50/case
	<b>5424</b>	Bag, Clear Dome, 27"	20/case
	<b>5425</b>	Bag, Clear Dome, 15"	20/case
 <p>Footswitch Covers Clear poly, non-sterile, various sizes, latex-free elastic cord.</p>	<b>5406</b>	17" x 15" [Fits with cinch lock]	100/case
	<b>5408</b>	Kit Consumables Pack, 9"	20/case
	<b>5422</b>	20" x 36" [Fits with cinch lock]	100/case
	<b>5421</b>	24" x 20" [Fits with cinch lock]	100/case
 <p>Camera Drapes</p>	<b>5431</b>	Camera Drape, 9" x 96", Clear polyethylene with perforated tip, Sterile.	25/case
	<b>5432</b>	Camera Drape, 5" x 96", Clear polyethylene with perforated tip, Sterile.	25/case
	<b>5433</b>	Blue Tip Camera Drape, 7" x 96", Clear polyethylene with elasticized fenestration to provide a secure fit around the hand piece, Sterile.	25/case
	<b>CFI-1500</b>	Blue Camera Drape, 5" x 96", Sterile.	20/case
	<b>CFI-1501</b>	Blue Ultrasound Probe Drape, 5" x 72", Sterile.	20/case
	<b>CFI-1504</b>	Clear Camera Drape, 7" x 96", Sterile.	20/case
 <p>Disposable Covers</p>	<b>CFI-702</b>	Set-Up Cover [44" x 70"] Sterile	25/case
	<b>CFI-708</b>	Mayo Stand Cover [24" x 58"] Sterile	25/case
	<b>5405</b>	Mini C Keyboard Cover [14" x 30 w/ tape] Sterile	25/case
	<b>5413</b>	IR Hand Control Cover [7" x 13" w tape closure] Sterile	25/case
	<b>5513</b>	Slush Drape, Polyurethane [66" x 44"] Sterile	24/case
	<b>5514</b>	Slush Drape, Polyurethane [44" x 44"] Sterile	24/case
	<b>5438</b>	X-Ray Cassette Drape [24.50" x 25"]	50/case
	<b>5505</b>	Spine Frame Pad Cover, Blue CSR [36" x 18"]	50/case
	<b>CFI-704</b>	Slush Drape, Polyethylene [54" x 54"] Sterile	25/case
	<b>CFI-788T</b>	Clear Drape [40" x 40" w/Tape]	50/case



## Sustainable Solutions

TIDI® Products is focused on reducing resource consumption, minimizing environmental impact and acting in an environmentally responsible manner, in order to provide effective solutions without compromising the well-being of future generations.

## Online Information

Explore our broad range of product areas and learn more about how our solutions can enhance the care experience. Visit **TIDIProducts.com**.

## Product Samples

Experience the softness of a gown, the absorbency of a towel, the coverage of a drape, or the performance of any TIDI® product with a no-charge sample pack.

Call **1.800.521.1314**, visit **TIDIProducts.com**, or e-mail us at **Excellence@TIDIProducts.com**.

## Value-Added Packaging

Newly designed packaging provides clear, concise, instant communication. Quickly and simply verify product identity and description for ease in ordering and selecting the right TIDI® products for your needs.



Certified in ISO 9001 and EN/ISO 13485

COPYRIGHT © 2016 TIDI PRODUCTS, LLC. ALL RIGHTS RESERVED.  
TIDISHIELD, C-ARMOR AND STERILE-Z REGISTERED TRADEMARKS ARE PROPERTY OF TIDI PRODUCTS, LLC. GRAB 'N GO, FLIP 'N GO, NAVI-CRANI, PURE-FIT, AND SUPPORTING CAREGIVERS AND PREVENTING INFECTIONS ARE TRADEMARKS OF TIDI PRODUCTS, LLC.  
FOR PATENT INFORMATION, SEE [GO.TIDIPRODUCTS.COM/PATENTS](http://GO.TIDIPRODUCTS.COM/PATENTS).

MADE IN USA LOGO IS REGISTERED TRADEMARK OF MADE IN THE USA BRAND, LLC. 3M COMMAND STRIPS IS REGISTERED TRADEMARK OF 3M. STRYKER REGISTERED TRADEMARK IS PROPERTY OF STRYKER CORPORATION. O-ARM IS A REGISTERED TRADEMARK OF MEDTRONIC, INC. ZEISS, OPMI, LEICA AND WILD ARE REGISTERED TRADEMARKS OF ECOLAB INC.



### TIDI Products, LLC

570 Enterprise Drive  
Neenah, WI 54956 USA

### Customer Service

800.521.1314  
+1 920.751.4036  
800.837.7770 FAX

[www.TIDIProducts.com](http://www.TIDIProducts.com)

