

# BEAUTIFIL Flow Plus<sup>®</sup>

**FINALLY, AN INJECTABLE HYBRID RESTORATIVE FOR ALL INDICATIONS**

**NOW IN TIPS –  
MORE SHADES!**

SHOFU

**F00  
Zero Flow**

Stackable

**F03  
Low Flow**

Self-leveling

## REVOLUTIONARY

**BEAUTIFIL Flow Plus** is the one flowable composite restorative approved for all classes of restoration – even for cusps and ridges! Unlike other flowables or Glass Ionomer (dentin replacements) which need a composite for the occlusal finish; **BEAUTIFIL Flow Plus** is a base, a liner, and a restorative all in one!

## A BASE, LINER & RESTORATIVE ALL IN ONE

- Indicated for class I, II, III, IV & V
- Stronger than many conventional composites
- Stackable & sculptable; stays put
- High radiopacity, greater than enamel
- Self-polishing and self-leveling





## RADIOPAQUE

BEAUTIFIL Flow Plus displays high radiopacity, greater than enamel which makes for accurate, quick and easy identification on an X-ray.

## CONTROLLED HANDLING

BEAUTIFIL Flow Plus combines the convenience of a traditional low flow and a unique Zero Flow that is stackable and sculptable; appropriate for all indications. It just stays put!

## AESTHETIC

BEAUTIFIL Flow Plus offers a variety of color stable chameleon-like shades that combine the light handling characteristics of both enamel and dentin into one material for natural looking restorations.



2.2g syringe / 0.21g x 20 tips

### BEAUTIFIL Flow Plus F00 (Zero Flow)

Shade	Item No.
A0.5	2000 / 2200
A1	2001 / 2201
A2	2002 / 2202
A3	2003 / 2203
A3.5	2004 / 2204
A4	2005 / 2205
A0.5O	2006 / 2211
A1O	2007 / 2212
A2O	2008 / 2213
B1	2062 / 2206
B2	2063 / 2207
C2	2064 / 2208
D2	2065 / 2209
INC	2010 / 2215
BW	2011 / 2216

### BEAUTIFIL Flow Plus F03 (Low Flow)

Shade	Item No.
A0.5	2013 / 2218
A1	2014 / 2219
A2	2015 / 2220
A3	2016 / 2221
A3.5	2017 / 2222
A4	2018 / 2223
A2O	2020 / 2229
A3O	2021 / 2230
B1	2066 / 2224
B2	2067 / 2225
C2	2068 / 2226
D2	2069 / 2227
INC	2023 / 2232
BW	2024 / 2233
MI	2026 / 2235
CV	2027 / 2236



### BEAUTIFIL Flow Plus/BeautiBond Kits

Contents: (6) 2.2g Beautifil Flow Plus syringes (1 - A1, 1 - A4, 2 - A2, 2 - A3), 50 BeautiBond unit dose, 50 micorbrushes

Description	Item No.
Zero Flow Kit	2000BF00
Low Flow Kit	2000BF03

### BEAUTIFIL Flow Plus Trial Kits\*

Description	Item No.
Pedo Kit	2000P
Standard Kit	2000S




# BEAUTIFIL Flow Plus

FINALLY, AN INJECTABLE HYBRID RESTORATIVE FOR ALL INDICATIONS


## STRONG

BEAUTIFIL Flow Plus exhibits higher compressive strength and lower wear than leading flowable restoratives. Flow Plus is even stronger than many conventional composites.



## FLUORIDE RELEASE

BEAUTIFIL Flow Plus has the proprietary GIOMER chemistry that delivers the unique ability to release and recharge fluoride as reported in the JADA 8-year study from the University of Florida Dental School. No failures or secondary caries – even after 8 years! At 13 years these same restorations had a retention rate of 66% and a secondary caries rate of only 3.27%.\*



\* 8-Year: J Am Dent Assoc. 2007 May;138(5):621-7  
13-Year: 2013 IADR Poster Presentation, Gordan et al.

### Additional BEAUTIFIL Products

BEAUTIFIL® II Description	Item No.
6 Color Set	1775
Tips 6 Color Set	1776
Cosmetic Set	1777
Tips 3 Color Set	1778

BEAUTIBOND® Description	Item No.
BeautiBond Unit Dose	1782
BeautiBond Bottle	1781



Shofu Dental Corporation  
1225 Stone Drive  
San Marcos, CA  
92078-4059, USA

\*Trial offer available while supplies last.