

page 1/6 BEAUTIFIL Flow Plus

Printing date: June 23, 2016

SECTION 1. Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

Trade Name:

BEAUTIFIL Flow Plus

1.2 Relevant identified uses of the substance or mixture and uses advised against

Relevant identified uses: Dental material Uses advised against: No further data

1.3 Details of the supplier of the safety date sheet

Company/Undertaking identification

Manufacturer's Name: SHOFU DENTAL GmbH Address: Am Brüll 17, 40878 Ratingen, Germany

Phone: +49 (0) 2102-8664-0 Fax: +49 (0) 2102-8664-64 E-Mail: info@shofu.de

Section in Charge: Quality Management & Regulatory Affairs

1.4 Emergency Telephone Number

+49-2102-8664-53 (SHOFU DENTAL GmbH) 24 hours / 7 days

SECTION 2. Hazards identification

2.1 Classification of the substance or mixture

CLASSIFICATION (EC 1272/2008)

Skin Sens. 1 H317 May cause an allergic skin reaction.

H318 Causes serious eye damage.

2.2 Label elements

LABEL IN ACCORDANCE WITH (EC) NO.1272/2008



HAZARD-DETERMINING COMPONENTS OF LABELLING Bis-GMA

SIGNAL WORD

Warning

HAZARD STATEMENTS

H317 May cause an allergic skin reaction.

H318 Causes serious eye damage.

PRECAUTIONARY STATEMENTS

P280 Wear protective gloves/protective clothing/eye protection/face

protection.

P261 Avoid breathing dust/fume/gas/mist/vapour s/spray.

P332+P313

If skin irritation occurs: Get medical advice/attention.

(Contd. on page2)



page 2/6
BEAUTIFIL Flow Plus

Printing date: June 23, 2016

(Contd. of page1)

P305+P351+P338

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

2.3 Other hazards

Results of PBT and vPvB assessment PBT: Not applicable. vPvB: Not applicable.

SECTION 3. Composition/information on ingredients

3.1 Chemical characterization: Mixtures

3.2 Description: Mixture of substances listed below with nonhazardous additions.

3.3 Dangerous components:

Bis-GMA Cas: 1565-94-2 15-25 %

EINECS: 216-367-7 Skin Sens. 1 H317, H318

3.4 Additional information: For the wording of the listed risk phrases refer to section 2.

SECTION 4. First-aid measures

4.1 Description of first aid measures

Eye contact: Rinse cautiously with water for several minutes. Remove contact lenses, if

present. and easy to do. If eye irritation persists, get medical advice/attention.

Skin contact: Wash immediately with soap and plenty of water. If on skin, skin irritation, get

medical advice/attention.

Ingestion: Rinse mouth. Get medical advice/attention.

Inhalation: Remove victim to fresh air and keep at rest in a position comfortable for breathing.

If symptom concerning breath goes out, call a POISON CENTER or doctor.

4.2 Most important symptoms and effects, both acute and delayed

No further relevant information available.

4.3 Indication of any immediate medical attention and special treatment needed

No further relevant information available.

SECTION 5. Fire-fighting measures

5.1 Extinguishing Media:

CO₂, Dry chemical, Foam

5.2 Special hazards arising from the substance or mixture:

In case of fire, do not breathe fumes and gases. (Irritation gases and fumes may emit.)

5.3 Advice for firefighters:

Wear fire protective cloth and self-contained breathing apparatus, if necessary.

SECTION 6. Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures: Avoid contact with eyes and skin.

(Contd. on page3)



page 3/6
BEAUTIFIL Flow Plus

Printing date: June 23, 2016

(Contd. of page2)

6.2 Environmental Precautions:

Send to approved treatment/disposal company or dispose according to local, state and federal regulations.

6.3 Methods and material for containment and cleaning Up:

Wipe up and discard in a suitable container.

6.4 Reference to other section:

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

SECTION 7. Handling and storage

7.1 Precautions for safe handling:

Avoid inhalation of grindings and prolonged skin contact with uncured resin.

7.2 Conditions for safe storage, including any incompatibilities:

Store in a cool and dark area. (Keep refrigerated when not in use)

7.3 Specific end use(s):

No further relevant information available

SECTION 8. Exposure controls/personal protection

8.1 Control parameters:

Ingredients with limit values that require monitoring at the workplace:

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

8.2 Exposure controls:

Respiratory Protection:

Not required.

Skin Protection: Hand Protection

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/the chemical mixture. Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation.

· Material of gloves

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

· For the permanent contact of a maximum of 15 minutes gloves made of the following materials are suitable:

(Contd. on page4)



page 4/6
BEAUTIFIL Flow Plus

Printing date: June 23, 2016

(Contd. of page3)

Butyl rubber, BR Nitrile rubber, NBR

Eye Protection: Safety goggles

SECTION 9. Physical and chemical properties

9.1 Information on basic physical and chemical properties

Appearance/Odor/Colour: Colored putty-like paste with sweet odor

Odour threshold Not determined. Ha Not determined. Melting point/freezing point Not determined. Not determined. **Boiling Point:** Flash point: Not determined. Not determined. Evaporation rate Flammability (solid, gas) Not applicable. Upper/lower flammability or explosive limits Not determined. Vapour pressure Not determined. Vapour density Not determined. Relative Density: 1.8 (water=1) Solubility: water solubility Insoluble.

Partition coefficient: n-octanol/water
Auto-ignition temperature
Decomposition temperature
Viscosity
Not determined.
Not determined.
Not determined.
Not determined.
Not applicable.
Oxidising properties
Not applicable.

9.2 Other information

No further relevant information available.

SECTION 10. Stability and reactivity

10.1 Reactivity:

Polymerization will occur with heat, light.

10.2 Chemical stability:

Stable under normal temperatures and pressures.

10.3 Possibility of hazardous reactions:

No dangerous reactions known.

10.4 Condition to Avoid:

Protect from exposure to the light.

Keep away from heat and direct sunlight.

10.5 Incompatible materials:

No further relevant information available.

10.6 Hazardous Decomposition Products:

None under normal conditions of storage and use.

(Contd. on page5)

Version Number 4 Revision date: June 2, 2016



page 5/6 BEAUTIFIL Flow Plus

Printing date: June 23, 2016

(Contd. of page4)

SECTION 11. Toxicological information

11.1 Information on toxicological effects:

Acute toxicity: Based on available data, the classification criteria are not met Skin corrosion/irritation: Based on available data, the classification criteria are not met.

Eye damage/irritation: H318 Causes serious eye damage.

Sensitization to the respiratory tract:

Based on available data, the classification criteria are not met.

Skin sensitization: Skin Sens. 1; H317 May cause an allergic skin reaction.

Germ cell mutagenicity/Genotoxicity:

Based on available data, the classification criteria are not met.

Carcinogenicity: Based on available data, the classification criteria are not met. Reproductive toxicity: Based on available data, the classification criteria are not met.

Effects on or via lactation: Lack of data.

Specific target organ toxicity (single exposure):

Based on available data, the classification criteria are not met.

Specific target organ toxicity (repeated exposure):

Based on available data, the classification criteria are not met.

SECTION 12. Ecological information

12.1 Toxicity:

No further relevant information available.

12.2 Persistence and degradability:

No further relevant information available.

12.3 Bioaccumulative potential:

No further relevant information available.

12.4 Mobility in soil:

No further relevant information available.

12.5 Results of PBT and vPvB assessment:

Not applicable.

12.6 Other adverse effects:

No further relevant information available.

SECTION 13. Disposal considerations

13.1 Waste treatment methods:

Dispose of contents/container to in accordance with local/regional/national/international regulations.

SECTION 14. Transport information

14.1 UN number: Void
14.2 UN proper shipping name: Void
14.3 Transport hazard class(es): Void
14.4 Packing group: Void

14.5 Environmental hazards: No further relevant information available.

14.6 Special precautions for user: Not applicable.

(Contd. on page6)

Revision date: June 2, 2016



page 6/6 BEAUTIFIL Flow Plus

Printing date: June 23, 2016

(Contd. of page5)

14.7 Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code: Not applicable.

SECTION 15. Regulatory information

- 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture:
 - EU REGULATIONS: See Section 2
 - Other regulations, limitations and prohibitive regulations:

The product is a medical device according to the EC-directive 93/42/EEC.

15.2 Chemical safety assessment:

A Chemical Safety Assessment has not been carried out.

SECTION 16. Other information

This product is intended for use by dental professionals. (instrument/material)

Relevant phrases:

H317 May cause an allergic skin reaction.

H318 Causes serious eye damage.

Abbreviations and acronyms:

EINESC: European Inventory of Existing Commercial Chemical Substances CAS: Chemical Abstracts Service (division of the American Chemical Society)

PBT: Persistent, Bioaccumulative and Toxic vPvB: very Persistent and very Bioaccumulative

Version Number 4 Revision date: June 2, 2016