

Evolve.

**GIOMER
TECHNOLOGY**
taking dentistry to the next level

Beautiful[®] II

A Nano-Hybrid Composite with
Fluoride Release & Recharge



Visit www.shofu.com or call 800.827.4638

BEAUTIFIL® II

Beautiful II creates indistinguishable restorations with optical characteristics of natural teeth. It has excellent versatility with natural shade reproduction that can be achieved with a chameleon effect, using a single shade that blends well with surrounding teeth making the restoration undetectable. Beautiful II is an aesthetic, universal nano-hybrid composite that releases and recharges fluoride and is approved for classes' I-V, perfect for direct restorations that require optimum aesthetics and biocompatibility. With superior handling and easy sculptability, Beautiful II will not slump or stick. After 13 years of placement, 61% of 41 restorations observed had either excellent or satisfactory ratings.¹

- Sustained fluoride release and recharge
- Natural fluorescence
- High level of radiopacity
- Superior handling – easy sculpting without slumping or sticking
- Exceptional marginal seal
- Anti-plaque affect minimizes plaque adhesion, inhibits bacterial colonization and plaque accumulation
- Fast and easy polishing creates outstanding surface luster that remains over time

BEAUTIFIL II Case Study



Dr. Frank Milnar, St. Paul, Minnesota

The Evolution of Giomer Technology.

As with all great technology, resin based composites continue to evolve and with this bioactive benefits have become increasingly important.

The Beautiful resin composite lines contain Shofu's proprietary Giomer technology, a bioactive surface pre-reacted glass (S-PRG) filler that actively releases six beneficial ions including fluoride.

Giomers unique ability to not only release but recharge fluoride helps to provide continual protection to the tooth. Giomer benefits:


- Decrease in acid production of cariogenic bacteria
- Formation of an acid-resistant layer
- An anti-plaque effect

8-Year: J Am Dent Assoc. 2007 May;138(5):621-7

13-Year: 2013 IADR Poster Presentation, Gordan et al.

References

1. University of Florida study. 3104 Clinical Evaluation of a Giomer Restorative System: Thirteen-Year Recall. V.V. Gordan, DDS, MS; P.K. Blaser; R.E. Watson, DDS, MAE; L. Sensi; D. McEdward, J.L. Riley III; I.A. Mjör, BDS, MSD, MS, Dr.odont. IADR General Session, March 2013



BEAUTIFIL® II Kits

Beautiful II 6 Color Set (Item No. 1775)

(6) 4.5g syringes in A2, A3, A3.5, B2, A3O, INC; Dura-White Stone, Super-Snap Singles, OneGloss Trial Set, paper pad

Beautiful II Tips 6 Color Set (Item No. 1776)

6 tip colors (0.25g x 120 tips) in A2, A3, A3.5, B2, A3O, INC; Dura-White Stone, Super-Snap Singles, OneGloss Trial Set, paper pad

Beautiful II Cosmetic Set (Item No. 1777)

(3) 4.5g syringes in A2, A3, INC; 2g Beautiful Opaquer-LO, 2g Beautiful Flow F02-A2, Dura-White Stone, Super-Snap Singles, OneGloss Trial Set, paper pad

Beautiful II Tips 3 Color Set (Item No. 1778)

3 tip colors (0.25g x 30 tips) in A2, A3, B2, OneGloss Trial Set

Beautiful Opaquers weight: 2grams

Beautiful Opaquer UO (Item No. 1363)

Beautiful Opaquer LO (Item No. 1364)