

Quala Dental Products

impression materials



Alginate Impression Material

This high algin impression material easily mixes to a smooth consistency producing shiny models with minimum flash. Providing more patient comfort, the thixotropic material does not flow from tray, also allowing for easier clean-up. Available in Spearmint flavor.

FEATURES

Easy wetting

Thixotropic

High algin content

High tear strength

No medicinal taste

BENEFITS

Mixes easier; more consistent; easier to clean up

Material will not flow from the tray, virtually eliminating discomfort

Smooth, shiny models with minimum flash

Removes from thin and undercut areas without tearing

Spearmint flavor; increased patient acceptance

	WORK TIME	INTRAORAL SET TIME
Regular Set	2:15	2:30
Fast Set	1:30	1:45

CODE	DESCRIPTION	UNIT
689400	Regular Set	1 lb canister
689410	Fast Set	1 lb canister

Distributed by:



#QSS032R2