

Quala Dental Products

Impression Materials

Alginate Substitute

For use with preliminary orthodontic models, provisionals, impressions for temporaries and many other applications. Quala Alginate Substitute material produces accurate impressions every time and the impression can be placed in disinfectant without damage.



FEATURES

Pre-mixed material

Excellent dimensional stability

Pleasant mint flavor

BENEFITS

No messy mixing of powder and liquid

Multiple pours can be made without loss of accuracy

Increases patient acceptance, helps to reduce gag reflexes during impressions

	WORK TIME	INTRAORAL SET TIME	SHORE "A"
Fast Set	1:00 minute	1:00 minute	50

CODE DESCRIPTION

Alginate Substitute

Q589510	50 mL Fast Set, Mint flavor <i>Contains: 8 x 50 mL cartridges</i>
Q589512	380 mL Fast Set, Mint flavor <i>Contains: 2 x 380 mL cartridges</i>

CODE	DESCRIPTION	QTY
Accessories		
Q-8020	Dynamic Mixing Tips, Yellow	50/bag
Q-8021	Bayonet Ring, Yellow	2/bag



Distributed by: