

Quala Dental Products

impression materials



High Performance Bite Registration

Mousse-like vinyl polysiloxane derivative material in a high performance cartridge dispensing system. Hard setting material with excellent tear strength. Maintains a 2-year shelf life if stored properly at or below 72°F and 50% humidity or less.

FEATURES

Nonflowing, thixotropic

Available in three different set times

Hard set

BENEFITS

No drip, sag or slump, stays in place on occlusal tooth surfaces

Reduces patient discomfort and speeds use of chair time

Prevents distortion in plastic, disposable bite trays

	WORK TIME	INTRAORAL SET TIME
Regular Set	2:00	1:00
Fast Set	0:30	1:00
Super Fast Set	0:30	0:30

Contains: 2 x 50 ml cartridges and 6 x 6.5mm mixing tips

CODE	DESCRIPTION
Q-9001HP	Regular Set
Q-9002HP	Fast Set
Q-9000HP	Super Fast Set
Q-9003HP	Vanilla Crème Fast Set
Q-9004HP	Vanilla Crème Super Fast Set
Q-9005HP	Rigid Bite - Fast Set, Unflavored (90+ Shore-A Hardness)



Quala
DENTAL PRODUCTS
www.quala.com

#QSS035R1

Distributed by: