

Quala Dental Products

Impression Materials

Clear Bite Registration High Performance

Clear vinyl polysiloxane derivative material in a high performance cartridge dispensing system. Clear Bite Registration can be used as a matrix material for temporaries. Maintains a 2-year shelf life.



FEATURES

Clear “see-through” material

Nonflowing, thixotropic

BENEFITS

Allows light to pass through in curing of temporaries, makes any bubbles or voids visible

No drip, sag or slump, stays in place on occlusal tooth surfaces

	WORK TIME	INTRAORAL SET TIME	SHORE “A”
Fast Set	1:00 minute	1:00 minute	60

Contains: 2 x 50 mL cartridges and 6 pink mixing tips

CODE	DESCRIPTION
Q-9006HP	Fast Set Kit

Distributed by: