

Quala Dental Products

infection control

Enzyme Ultrasonic Tablets

For general purpose ultrasonic cleaning, evacuation system cleaning, and instrument presoak. Effective, fast removal of blood, food particles, compounds, and other soils which are normally resistant to detergents alone. Two tablets make one gallon of solution. Saves storage space, and easy to use. Clean mint scent.



FEATURES

Compact, convenient package of 64 tablets dilutes to 32 gallons

Drop directly into ultrasonic cleaner

Clean mint scent

Contains fast-acting protease enzymes

Nonionic, low foaming, biodegradable

Also effective as Evacuation Cleaner

BENEFITS

Saves counter and storage space; No lifting of heavy liquids

No messy mixing of liquid concentrates; easy to use

Pleasant to work with; minimizes any odor from tank

Effective, fast removal of blood, food particles, compounds, and other soils which are normally resistant to detergents alone

Easy on instruments and the environment

One product for multiple needs

CODE

50036380

QTY

64 tablets/box

Distributed by:



#QSS047R1