



ETCH GEL

A 37% phosphoric acid blue gel that provides superb viscosity for accurate placement.

FEATURES

37% weight phosphoric acid

Thixotropic gel

Silica thickener

Blue gel

Syringe delivery system

BENEFITS


Optimally etches in 15 to 20 seconds

Will not slump or flow on a vertical surface, allowing greater control over the application

Easy to rinse off, with no remaining residue on tooth surface. Will not interfere with adhesive bonding

Excellent visibility in the oral environment

Provides accurate placement



CODE	DESCRIPTION
Q400042ML	4 x 2.0 mL Syringes with 20 tips
Q400122ML	Economy Kit, 12 x 2.0 mL Syringe with 50 tips
Q40030	Jumbo Kit, 2 x 30 mL Syringes, 50 tips, three 3 mL empty Syringes, 1 connector
Q400TP	23 Gauge (Blue) Tips, 25/pkg
Q400SR	3 mL empty delivery Syringes

Distributed by:



Quala
DENTAL PRODUCTS
www.quala.com

#QSS011R2