



EXTRAORAL DENTAL FILM

For panoramic and cephalometric applications. A high contrast film for use with Blue and Green emitting screens.



CODE	DESCRIPTION	QTY
18500-QU	5" x 12" Green Sensitive Medium Speed High Contrast	100/box
18520-QU	6" x 12" (15cm x 30cm) Green Sensitive Medium Speed High Contrast	100/box
18501-QU	5" x 12" Blue Sensitive, Full Speed	100/box
18521-QU	6" x 12" (15cm x 30cm) Blue Sensitive, Full Speed	100/box
18525-QU	8" x 10" Green Sensitive Medium Speed High Contrast	100/box
18530-QU	5" x 12" Duplicating	100/box
18531-QU	6" x 12" (15cm x 30cm) Duplicating	100/box

Distributed by:

