

Qnala Dental Products

infection control



Sterilization Pouches

Dual Indicator sterilization pouches feature process indicators that are printed internally and externally on each pouch and change color during sterilization to verify pouches have been processed. Pink arrows will turn brown when processed by steam. Blue arrows will turn gold/brown when EtO processed. These process indicators only measure temperature. Pouches are made from medical grade paper and feature blue tint film, solid seal and thumb notch.

FEATURES

Wide solid seal

Dual steam and gas process indicators

End seals

See-through four-ply film

Tack seal over perforation

Edge overlap

Clean, smooth opening

Steam absorbent paper

Available in an assortment of sizes

Internal and external indicators

BENEFITS

Protects against separation, puncture and instrument contamination

Arrows change color to indicate when processed

Prevent accidental opening and dust accumulation in storage

Lets you identify instrument set-up; helps reveal punctures

Prevents separation of sealing end, and prevents pouch from being compromised

Allows easy pouch opening to access contents

Assured with fiber orientation of fine, virgin medical grade paper

Allows steam through, yet remains strong even when wet

For all uses, from burs to small instrument cassettes

Assurance of process

Qnala Dental Products

infection control

CODE	SIZE	QTY
AD2440	2¼" x 4"	200/box
AD2790	2¾" x 9"	200/box
AD3554	3½" x 5¼"	200/box
AD3590	3½" x 9"	200/box
AD5410	5¼" x 10"	200/box

Bulk

LAD2790	2¾" x 9"	500/box
LAD3590	3½" x 9"	500/box
LAD5410	5¼" x 10"	500/box

Distributed by:



#QSS046R1