

# Quala Dental Products

## Infection Control

### Triple Enzymatic Evacuation Cleaner

Non-foaming formula cleans evacuation systems and deodorizes as it cleans. Triple enzyme formula (Protease; Amylase; Lipase) dissolves organic fluids and debris to restore maximum performance to system and hoses. Neutral pH formula leaves a clean, fresh scent.



#### FEATURES

Economical concentrate

Environmentally friendly, neutral pH

Easy dispensing

Safe to use every day

#### BENEFITS

One ounce mixes with warm water to make 32 oz. of solution

Safe for use with amalgam separators

Solution squeezes easily into pre-measured reservoir

Daily use improves cleaning efficiency

CODE	DESCRIPTION
Q182605	1 Quart (32 oz.) Concentrate

Distributed by: