



Surgical Products

BD Bard-Parker™ Protected Blade System

Protect yourself.
Protect your staff...
in the OR.

Protection Without Compromise

The BD Bard-Parker™ Protected Blade System has been designed with the weight and feel that clinicians prefer, with convenient and easy one-handed activation of the protective shield. The blade remains covered during both assembly and disposal to help prevent accidental contact with the blade.

Get a lock on safety

The system offers a unique locking shield that provides easy activation:

- Clear shield locks before and after use for protection
- Top activation button enables left and right-handed use
- Audible button click and tactile sensation confirm lock is in place
- Blade assembly and disassembly are quick and easy, no extra steps

No compromise to technique

The system provides all the quality you've come to expect from BD Bard-Parker™:

- BD Bard-Parker™ Stainless-Steel Blade
- Stainless-steel handle with conventional weight and feel
- Durable handle design for extended life

BD Bard-Parker™ Protected Blade System

Protected Blade System Metal Handles

PROD #	FOR USE WITH BLADE SIZE	HANDLE SIZE	PACKAGING
374030	10, 11, 15, 15C	3	5/box
374031	10, 11, 15, 15C	3L	5/box
374040	20-25	4	5/box
374041	20-25	4L	5/box

Protected Surgical Blades

PROD #	BLADE SIZE	PACKAGING
373910	10	50/box, 3/case
373911	11	50/box, 3/case
373915	15	50/box, 3/case
373920	20	50/box, 3/case
373921	21	50/box, 3/case
373922	22	50/box, 3/case
373923	23	50/box, 3/case
374305	15C	50/box, 3/case

BD Bard-Parker™
Stainless Steel Blade

Non-Slip Grip

Clear Protective Shield with integrated Blade

Locks when blade is covered –
Protection before and after use

Locks when blade is uncovered –
Security and control during use

One-handed activation – Ease of use

Transparent – Blade visible

Blade and shield are removed together –
Protection for blade disposal

Top Button

Tactile sensation – Confirms lock

Audible click – Confirms lock

Lock Indicator Line

Visual reference – Confirms lock

Centimeter Scale

Handle Size

Stainless Steel Handle

Conventional weight and feel –
Tactile sensitivity

Durable extended life



BD Bard-Parker™ Protected Disposable Scalpel


Featuring Locking Retractable Shield

Protect Yourself. Protect Your Staff. Designed to minimize the risk of surgical blade injuries, the BD Bard-Parker™ Protected Disposable Scalpel gives clinicians, throughout the hospital...

Protection Without Compromise

- **In Procedure:**
Permitting clinicians to focus first on the surgical site, rather than on the instrument.
- **In Handling:**
Providing instrument characteristics of weight, handling, and feel that clinicians prefer.
- **In Technique:**
Allowing clinicians to retract the protective shield, easily and safely with one hand.
- **In Reliability:**
Assuring quality and performance you have come to expect from the BD Bard-Parker™ brand.

PROD #	BLADE SIZE	PACKAGING
372610	#10	10/box, 10/case
372611	#11	10/box, 10/case
372612	#12	10/box, 10/case
372613	#12B	10/box, 10/case
372615	#15	10/box, 10/case
372620	#20	10/box, 10/case
372621	#21	10/box, 10/case
372622	#22	10/box, 10/case
372623	#23	10/box, 10/case



Surgical Products

BD Bard-Parker™ Safetylock™ Blades – Sterile BD Rib-Back™ Carbon Steel Blades, Individually Packaged

PROD #	BLADE SIZE	PACKAGING
371150	10	50/shelf pack, 150/case
371151	11	50/shelf pack, 150/case
371152	12	50/shelf pack, 150/case
371153	15	50/shelf pack, 150/case
371154	20	50/shelf pack, 150/case
371155	21	50/shelf pack, 150/case
371156	22	50/shelf pack, 150/case
371157	23	50/shelf pack, 150/case

BD Bard-Parker™ Sterile Rib-Back™ Carbon Steel Blades, Individually Packaged

PROD #	BLADE SIZE	PACKAGING
371110	10	50/shelf pack, 150/case
371111	11	50/shelf pack, 150/case
371112	12	50/shelf pack, 150/case
371115	15	50/shelf pack, 150/case
371120	20	50/shelf pack, 150/case
371121	21	50/shelf pack, 150/case
371122	22	50/shelf pack, 150/case
371123	23	50/shelf pack, 150/case

BD Bard-Parker Surgical Scalpels Sterile, Disposable Scalpels with Stainless Steel Blade

PROD #	BLADE SIZE	PACKAGING
371610	10	10/box, 100/case
371611	11	10/box, 100/case
371612	12	10/box, 100/case
371615	15	10/box, 100/case
371620	20	10/box, 100/case
371621	21	10/box, 100/case
371623	23	10/box, 100/case
371613	12B	10/box, 100/case

BD Bard-Parker™ Sterile Stainless Steel Blades, Individually Packaged

PROD #	BLADE SIZE	PACKAGING
371210	10	50/shelf pack, 150/case
371700	10A	50/shelf pack, 150/case
371211	11	50/shelf pack, 150/case
371212	12	50/shelf pack, 150/case
371712	12B	50/shelf pack, 150/case
371215	15	50/shelf pack, 150/case
371716	15C	50/shelf pack, 150/case
371220	20	50/shelf pack, 150/case
371221	21	50/shelf pack, 150/case
371222	22	50/shelf pack, 150/case
371223	23	50/shelf pack, 150/case

BD Bard-Parker™ Rib-Back™ Carbon Steel Blades Nonsterile, Individually Packaged, in 25 Strips of 6

PROD #	BLADE SIZE	PACKAGING
371310	10	150/case
371311	11	150/case
371315	15	150/case
371320	20	150/case
371321	21	150/case
371322	22	150/case
371325	25	150/case
371340	60 Autopsy (fits #8 Handle only)	150/case
371312	12	150/case
371323	23	150/case
371324	24	150/case
371325	31	150/case

BD Bard-Parker™ Surgical Blade Handles

PROD #	FOR USE WITH BLADE SIZE	PACKAGING
371030	3/10 through 15	5/case
371031	3L/10 through 15	5/case
371040	4/20 through 23	5/case
371041	4L/20 through 23	5/case
371050	5/10 through 15	5/case
371060	6/20 through 23	5/case
371070	7/10 through 15	5/case
371080	8/Autopsy Handle (for Blade 60)	5/case
371090	9/10 through 15	5/case



Antiseptic Surgical Scrubs and Surgical Accessories

BD™ E-Z Scrub Impregnated Preoperative Scrub Brushes and Accessories

PROD #	DESCRIPTION	PACKAGING
377479	BD E-Z Scrub with 4% CHG (brush/sponge)	10 dispensers of 30 ea.
374089	BD E-Z Scrub with BD Ultradex™ Solution; 3% Chloroxylenol (PCMX) with emollients (brush/sponge)	10 dispensers of 30 ea.
371901	Bracket	1 ea.

BD Chest Drainage

Heimlich Chest Drain Valve-Sterile

PROD #	DESCRIPTION	PACKAGING
373460	Drainage of the Pleural Cavity	10/case

