

BD Vacutainer® Plus Plastic K₂EDTA Tubes



BD Reference #
367842

BD Reference #
367899

BD Reference #
368589

A comprehensive offering of Pink Closures for easy identification

These tubes are designed with distinguishable Pink Closures to easily differentiate Hematology from Blood Bank Specimens.

Each has a specially designed “crossmatch” label for required American Association of Blood Banks (AABB) patient identification.

All BD Vacutainer® Plus Plastic K₂EDTA Tubes may be used for routine immunohematology testing and blood donor screening.*

BD Vacutainer® Plus Plastic K₂EDTA Tubes

BD Reference Number	Draw Volume (mL)	Tube Size (mm)	Closure Type	Packaging
367842	2	13 x 75	BD Hemogard™	100 Box/ 1,000 Case
367899	6	13 x 100	BD Hemogard™	100 Box/ 1,000 Case
368589	10	16 x 100	Rubber Stopper	100 Box/ 1,000 Case



Helping all people
live healthy lives

To learn more about these products or any other BD Vacutainer® specimen collection products, educational materials, or services offered by **BD Diagnostics – Preanalytical Systems**, please contact your local BD Sales Consultant today.

You can also contact us via:

BD Global Technical Services at 1.800.631.0174 or e-mail
vacutainer_techservices@bd.com

BD Customer Service at 1.888.237.2762

Or visit us anytime online at www.bd.com/vacutainer

**The performance characteristics of these tubes have not been established for immunohematology testing in general; therefore, users must validate the use of these tubes for their specific assay-instrument/reagent system combinations and specimen storage conditions.*