



BD veterinary products

Devices you can trust for every patient's journey through life

Providing safe and innovative products is at the core of BD, and we are committed to being your partner in the life-long relationships you build with each patient.



Improving the veterinary experience on every level

Whether it's a routine exam or a surgical procedure, one of the keys to successful veterinary medicine is confidence. The importance of confidence is magnified since pets have no voice but can sense trepidation or insecurity. In turn, pet owners pick up on their pets' anxiety. This may result in discomfort in the exam room or surgical suite. BD designs and manufactures medical devices that enable accurate dosing and are easy to use. We understand that you and your veterinary staff feel secure in your medical knowledge. Now, you can feel secure in the execution of that care.

For more than 100 years, the BD name has been synonymous with safe, reliable products. We've built a reputation the world over for consistent, dependable quality. Our expertise extends from products used for wellness exams and surgery to specialized areas of use such as oncology and diabetes. All products are designed with two things in mind—patient care and ease of use for the caregiver. That's why medical professionals trust BD as a partner for every patient's journey through life.

Contact your authorized BD veterinary distributor
or our logistics partner.



Table of contents



Wellness products

4

Hypodermic products

Blood collection products

Disposal solutions

Surgical products¹

7

Vascular access products

Anesthesia products

Antimicrobial solutions

Delivery devices for diabetes

9

Diabetes care products

Oncology products²

10

Biopsy products

The BD PhaSeal™ system

1. Additional BD veterinary surgical products include: hypodermic syringes and needles, blood collection products, sharps disposal solutions

2. Additional BD veterinary oncology products include: hypodermic syringes and needles, blood collection products, spinal needles, sharps disposal solutions

Wellness products

Designed to help pets live healthier lives

Routine injections. Vaccinations. Blood collection. BD has a syringe, needle and disposal system for every injection application. Our philosophy is that routine procedures should be just that, routine. So it's not surprising that with this as our priority, BD hypodermic products are designed to provide consistent quality and ease of use for the entire veterinary team.

Hypodermic products

- BD™ syringes are a precision drug delivery device that are an essential part of your daily patient care toolkit, which you can use across a range of critical applications and departments.
- BD syringes provide smooth plunger movement ensured by silicone lubrication of the stopper.
- BD syringes help clinicians provide accurate medication doses due to the clear barrel with contrasted scale markings.



BD syringes with needles



Reference no.	Product description	Packaging
305664	1 mL, slip-tip syringe with 25 G x 5/8 in. needle, regular bevel	100 / box, 600 / case
309623	1 mL TB syringe with 27 G x 1/2 in. needle	100 / box, 800 / case
309625	1 mL TB syringe with 26 G x 3/8 in. needle, intradermal bevel	100 / box, 800 / case
309626	1 mL TB syringe with 25 G x 3/8 in. needle, regular bevel	100 / box, 800 / case
309569	3 mL, BD Luer-Lok™ tip syringe with 22 G x 3/4 in. needle, regular bevel	100 / box, 800 / case
309570	3 mL, BD Luer-Lok tip syringe with 25 G x 3/8 in. needle, regular bevel	100 / box, 800 / case
309572	3 mL, BD Luer-Lok tip syringe with 22 G x 1 in. needle, regular bevel	100 / box, 800 / case
309578	3 mL, BD Luer-Lok tip syringe with 20 G x 1 in. needle, regular bevel	100 / box, 800 / case
301628	3 mL, slip-tip syringe with 22 G x 3/4 in. needle, regular bevel	100 / box, 800 / case

BD PrecisionGlide™ needles



Reference no.	Product description	Packaging
305122	25 G x 5/8 in.	100 / box, 1000 / case
305142	22 G x 3/8 in.	100 / box, 1000 / case
305175	20 G x 1 in.	100 / box, 1000 / case
305195	18 G x 1 in.	100 / box, 1000 / case

BD veterinary syringes with needles VS600

Reference no.	Product description	Packaging
305660	3 mL, BD Luer-Lok syringe with 25 G x 5/8 in. needle	600 syringes / case
305661	3 mL, BD Luer-Lok syringe with 22 G x 1 in. needle	600 syringes / case
305662	3 mL, BD Luer-Lok syringe with 22 G x 3/8 in. needle	600 syringes / case
305663	3 mL, slip-tip syringe with 22 G x 3/4 in. needle	600 syringes / case
305664	1 mL, slip-tip syringe with 25 G x 5/8 in. needle	600 syringes / case

BD syringes

Reference no.	Product description	Packaging
309659	1 mL TB syringe, slip tip	200 / box, 1600 / case
309628	1 mL syringe, BD Luer-Lok tip	100 / box, 800 / case
309656	3 mL, slip-tip	200 / shelf carton
309657	3 mL, BD Luer-Lok tip	200 / shelf carton
309646	5 mL, BD Luer-Lok tip	125 / shelf carton
309647	5 mL, slip-tip	125 / shelf carton
302995	10 mL, BD Luer-Lok tip	100 / box, 400 / case
303134	10 mL, slip-tip	100 / box, 400 / case
302830	20 mL syringe, BD Luer-Lok tip	48 / box, 192 / case
302831	20 mL, slip-tip	48 / box, 192 / case

BD oral syringes

Reference no.	Product description	Packaging
305217	1 mL, clear	100 / bag, 500 / case
305220	3 mL, clear	100 / bag, 500 / case
305207	1 mL, amber	100 / box, 500 / case
305210	3 mL, amber	100 / box, 500 / case



Blood collection products

Quality results begin with a quality specimen. BD has engineered needles and tubes to make blood and specimen collection easier for medical professionals. The BD Vacutainer® Safety-Lok™ blood collection set is simple, easy to use and engineered for safety. The safety mechanism can be activated immediately after the blood draw and helps protect against needlestick injury. BD Microtainer® and BD Vacutainer tubes bring confidence to caregivers with clinically proven efficacy and performance.

BD Vacutainer® Safety-Lok™ blood collection sets



- The BD Vacutainer Safety-Lok shield provides safety, yet requires no change in technique when collecting blood from patients for routine testing.
- The slim, easy-to-handle design provides control during needle insertion.
- Large wing size allows positioning during use.
- Minimized coiling of tubing facilitates effective handling from insertion to removal of needle.
- Blood flashback is clearly viewed through a translucent shield.

Reference no.	Product description	Packaging
367287	with Luer, 21 G x ¾ in., 7 in. tube, green	50 / box, 200 / case
367281	with Luer, 21 G x ¾ in., 12 in. tube, green	50 / box, 200 / case
367296	without Luer, 21 G x ¾ in., 12 in. tube, green	50 / box, 200 / case
367292	with Luer, 23 G x ¾ in., 7 in. tube, light blue	50 / box, 200 / case
367283	with Luer, 23 G x ¾ in., 12 in. tube, light blue	50 / box, 200 / case
367297	without Luer, 23 G x ¾ in., 12 in. tube, light blue	50 / box, 200 / case
367294	with Luer, 25 G x ¾ in., 7 in. tube, royal blue	50 / box, 200 / case
367285	with Luer, 25 G x ¾ in., 12 in. tube, royal blue	50 / box, 200 / case
367298	without Luer, 25 G x ¾ in., 12 in. tube, royal blue	50 / box, 200 / case
368652	with pre-attached holder, 21 G x ¾ in., 12 in. tube, green	25 / box, 200 / case
368653	with pre-attached holder, 23 G x ¾ in., 12 in. tube, light blue	25 / box, 200 / case

BD Microtainer® tubes with BD Microgard™ closure



- BD Microtainer blood collection tubes for capillary blood collection are non-evacuated and have a textured inner surface to improve blood flow and mixability.
- Printed fill lines ensure correct blood volume during blood collection for testing.
- The wide-diameter BD Microtainer tubes with BD Microgard™ closure feature an integrated collection scoop for improved mixing ability.

Reference no.	Product description	Packaging
365974	Plastic, 250–500 µL, lavender, K ₂ EDTA	50 / box, 200 / case
365965	Plastic, 200–400 µL, green, lithium heparin	50 / box, 200 / case
365985	Plastic, 400 – 600 µL, mint green, lithium heparin and gel for plasma separation	50 / box, 200 / case
365967	Plastic, 400 – 600 µL, gold, clot activator and gel for serum separation	50 / box, 200 / case
365978	Plastic, 400 – 600 µL, gold / amber, clot activator and gel for serum separation	50 / box, 200 / case
365963	Plastic, 250 – 500 µL, red, no additive	50 / box, 200 / case

BD Vacutainer® tubes with Hemogard™ closure

- The efficacy, performance and ease of use of BD Vacutainer plastic tubes is clinically proven.
- The BD Vacutainer Hemogard™ closure is a safety-engineered closure, featuring an integrated shield that minimizes human contact with blood on the stopper or tube rim and effectively guards against blood splashing and splatter.
- A variety of different additives, tube sizes and draw volumes are available to meet the workflow needs of all laboratories for sample evaluation.

Reference no.	Product description	Packaging
367841	Plastic, 13 x 75, 2.0 mL lavender, spray-coated K ₂ EDTA	100 / box, 1,000 / case
367856	Plastic, 13 x 75, 3.0 mL lavender, spray-coated K ₂ EDTA	100 / box, 1,000 / case
363083	Plastic, 13 x 75, 2.7 mL light blue, buffered sodium citrate 0.109M	100 / box, 1,000 / case
367986	Plastic, 13 x 100, 5.0 mL gold, SST™ tube clot activator and gel	100 / box, 1,000 / case

BD Vacutainer tubes with conventional stopper

Reference no.	Product description	Packaging
367981	Plastic, 13 x 75, 3.5 mL red / gray, SST tube clot activator and gel	100 / box, 1,000 / case
366667	Plastic, 13 x 75, 3.0 mL green, lithium heparin	100 / box, 1,000 / case
366668	Plastic, 13 x 75, 3.0 mL red, silicone-coated	100 / box, 1,000 / case
367987	Plastic, 16 x 100, 7.5 mL red / gray, SST tube clot activator and gel	100 / box, 1,000 / case

Disposal solutions

Our BD™ disposal solutions meet Occupational Safety and Health Administration (OSHA) Bloodborne Pathogens Standard CFR 1910.1030. BD is also committed to reducing environmental impact. Our BD Recykleen™ products are made from recycled materials while still maintaining the high quality you desire.



BD Recykleen™ patient / exam room

Reference no.	Product description
305053	3 gal next-generation counterbalanced door—20% recycled plastic, pearl
305054	3 gal next-generation counterbalanced door—20% recycled plastic, red
305076	3 gal plug cap, yellow
305493	3 gal plug cap, yellow

Chemotherapy

Reference no.	Product description
305076	3 gal plug cap, yellow
305493	3 gal plug cap, yellow

Patient / exam room

Reference no.	Product description
305517	5.4 qt counterbalanced door, red
305426	5.4 qt counterbalanced door, red
305443	5.4 qt side entry, red
305444	5.4 qt side entry, pearl
305425	5.4 qt side entry, pearl
305422	2 gal counterbalanced door, pearl
305435	2 gal counterbalanced door, red
305423	3 gal counterbalanced door, pearl
305436	3 gal counterbalanced door, red

Surgical products

Confidence with every procedure

Hypodermic needles and syringes. Catheter systems. Blood collection products. Anesthesia and spinal needles. BD excels in the devices you need to make surgery not only successful, but safe for both patient and veterinary staff. From routine spay and neuter procedures to more complex amputation or orthopedic surgeries, BD has a pet-focused portfolio. We offer delivery devices for every stage of the surgical process from anesthesia to post-surgery pain medication.

Vascular access products

Our catheters are made with BD Vialon™ biomaterial that softens up to 70% in the vein, enabling longer dwell times* and reducing the chance of mechanical phlebitis by up to 50%.

Catheters made with BD Vialon biomaterial are:

- Unique with properties that allow it to become softer and more pliable during use*
- Designed to resist kinking
- More flexible than FEP polymer catheter material

BD Insyte™ IV catheters with BD Vialon™ biomaterial



Reference no.	Product description	Packaging
381212	24 G x 0.75 in. yellow	50 / box, 200 / case
381223	22 G x 1.00 in. blue	50 / box, 200 / case
381233	20 G x 1.00 in. pink	50 / box, 200 / case
381234	20 G x 1.16 in. pink	50 / box, 200 / case
381237	20 G x 1.88 in. pink	50 / box, 200 / case
381244	18 G x 1.16 in. green	50 / box, 200 / case
381247	18 G x 1.88 in. green	50 / box, 200 / case
381254	16 G x 1.16 in. gray	50 / box, 200 / case
381257	16 G x 1.77 in. gray	50 / box, 200 / case
381267	14 G x 1.75 in. orange	50 / box, 200 / case

BD Insyte-W™ IV catheters with BD Vialon biomaterial and wings



Reference no.	Product description	Packaging
381312	24 G x 0.75 in. yellow	50 / box, 200 / case
381323	22 G x 1.00 in. blue	50 / box, 200 / case
381333	20 G x 1.00 in. pink	50 / box, 200 / case
381334	20 G x 1.16 in. pink	50 / box, 200 / case
381337	20 G x 1.88 in. pink	50 / box, 200 / case
381344	18 G x 1.16 in. green	50 / box, 200 / case
381347	18 G x 1.88 in. green	50 / box, 200 / case
381354	16 G x 1.16 in. gray	50 / box, 200 / case
381357	16 G x 1.77 in. gray	50 / box, 200 / case

BD Insyte-N™ IV catheters with BD Vialon biomaterial

Reference no.	Product description	Packaging
381211	24 G x 0.56 in. yellow	50 / box, 200 / case
381311	24 G x 0.56 in. yellow with wings	50 / box, 200 / case

*Compared to an FEP catheter

Anesthesia products

For the safety of both pets and the veterinary team, anesthesia is used for many procedures to help keep pets still so they don't hurt themselves. The BD Quincke spinal needles are available in gauges from 18 to 27 and in lengths from 1 in. to 7 in.

BD Quincke spinal needles, sterile, single use

Reference no.	Product description	Packaging
405081	27 G x 3 ½ in., gray	25 / box, 100 / case
405164	26 G x 3 ½ in., tan	25 / box, 100 / case
405170	25 G x 3 in., blue	25 / box, 100 / case
405180	25 G x 3 ½ in., blue	25 / box, 100 / case
400106	23 G x 3 ½ in., turquoise	25 / box, 100 / case
405161	22 G x 1 ½ in., black	25 / box, 100 / case
405074	22 G x 2 ½ in., black	25 / box, 100 / case
405171	22 G x 3 in., black	25 / box, 100 / case
405181	22 G x 3 ½ in., black	25 / box, 100 / case
405162	20 G x 1 ½ in., yellow	25 / box, 100 / case
405071	20 G x 2 ½ in., yellow	25 / box, 100 / case
405172	20 G x 3 in., yellow	25 / box, 100 / case
405182	20 G x 3 ½ in., yellow	25 / box, 100 / case
405174	18 G x 3 in., pink	25 / box, 100 / case
405184	18 G x 3 ½ in., pink	25 / box, 100 / case



Antimicrobial solutions

BD E-Z Scrub™ surgical scrub brushes are designed to provide maximum clinical efficacy.



BD E-Z Scrub™ surgical scrub brushes

Reference no.	Product description	Packaging
371073	BD E-Z Scrub with 4% CHG (brush / sponge)	30 / box, 300 / case
371163	BD E-Z Scrub with 3% chloroxylenol (PCMX) (brush / sponge)	30 / box, 300 / case
372053	BD E-Z Scrub with povidone iodine, 1% available iodine (brush / sponge)	30 / box, 300 / case
372453	BD E-Z Scrub with povidone iodine, 0.5% available iodine (brush / sponge)	30 / box, 300 / case
371603	BD E-Z Scrub dry, no detergent (sterile, brush / sponge)	30 / box, 300 / case
371901	BD E-Z Scrub brush stainless steel dispenser bracket	1 bracket / case

Additional BD veterinary products include:

- Hypodermic syringes and needles
- Blood collection products
- Disposal solutions

Delivery devices for diabetes

Discover a pet-sized approach to diabetes

As pet obesity rates rise, so does the prevalence of diabetes in cats and dogs. One in 100 dogs and one in 200–250 cats has diabetes.^{3,4} With this many pets and pet owners affected, it should be no surprise that BD has made it a priority to design specific devices to make this disease manageable for veterinarians and pet owners. For example, our U-40 syringe is carefully designed to make insulin administration comfortable for pets, minimizing squirming and potential injuries during the injection process.

Diabetes care products

- Cost effective for clinics so you can pass the savings along to clients.
- Only veterinary-specific insulin syringe in the U.S. with red scale marks to indicate use for U-40 insulin, which may help avoid misdosing with other insulin syringes.
- Half-unit scale marks designed to reduce the risk of over- or underdosing.
- Manufactured in Holdrege, Nebraska, U.S.
- Pet-specific packaging helps drive pet owners back into the clinic for refill purchases.



BD pen needles

Reference no.	Product description	Packaging
328203	BD Ultra-Fine™ original pen needle 29 G x 12.7 mm	100 / pack

BD U-40 insulin syringe

Reference no.	Product description	Packaging
323000	BD veterinary insulin syringe ¾ mL, 29 G x 12.7 mm	100 / pack
323001	BD veterinary insulin syringe ½ mL, 29 G x 12.7 mm	100 / pack

Diabetes swabs

Reference no.	Product description	Packaging
326895	BD alcohol swabs	100 / pack

Diabetes accessories

Reference no.	Product description	Packaging
328235	BD Safe-Clip™ needle clipping device	1
328233	BD Magni-Guide™	1

Sharps containers

Reference no.	Product description	Packaging
323488	BD sharps disposal by mail, 1.4 qt	1

BD U-100 insulin syringe

Reference no.	Product description	Packaging
324702	BD insulin syringe ¾ mL, 29 G x 12.7 mm	200 / pack
328431	BD Ultra-Fine insulin syringe ¾ mL, 30 G x 12.7 mm	100 / pack
328438	BD Ultra-Fine short-needle insulin syringe ¾ mL, 31 G x 8 mm	100 / pack
328440	BD Ultra-Fine short-needle ½ unit scale insulin syringe ¾ mL, 31 G x 8 mm	100 / pack
329461	BD Lo-Dose™ insulin syringe ½ mL, 28 G x 12.7 mm, in blister pack	
324703	BD insulin syringe ½ mL, 29 G x 12.7 mm	200 / pack
329410	BD insulin syringe 1.0 mL, 28 G x 12.7 mm	100 / pack
329420	BD Micro-Fine™ insulin syringe 1.0 mL, 28 G x 12.7 mm, in blister pack	100 / pack

BD Microtainer contact-activated lancet

Reference no.	Product description	Packaging
366592	Needle: 30 G x 1.5 mm, purple	200 / box, 2,000 / case
366593	Needle: 21 G x 1.8 mm, pink	200 / box, 2,000 / case
366594	Blade: 1.5 mm width x 2.0 mm, blue	200 / box, 2,000 / case

3. Nelson RE. Canine diabetes mellitus. In: Ettinger SJ, Feldman EC, eds. Textbook of Veterinary Internal Medicine 7th ed. St. Louis, MO: Saunders; 2010: 1782–1816.

4. Armstrong PJ. Diabetes mellitus in cats: Risk factors and treatment. Apr 01, 2009. Dvm360.com. Accessed 29 Dec 2016 at <http://veterinarycalendar.dvm360.com/diabetes-mellitus-cats-risk-factors-and-treatment-proceedings>

Oncology products

Making the tough stuff a little less tough for you

When a pet is undergoing something as serious as a tissue or bone marrow biopsy, there's enough stress. Both pet owners and pets are nervous and they look to you to instill confidence in them. They need to know they are getting the best care. The devices you use in these situations need to feel comfortable in your hand and be easy to administer—so there's no added anxiety in the room.

Biopsy products*

Bone marrow biopsy products

- Jamshidi™ and Illinois biopsy products are ergonomically designed to help you and your staff accurately capture and retrieve samples from cats and dogs.
- Sterile trays include prepping, draping, local anesthetic and the procedural components all in one.



Reference no.	Product description	Packaging
TJC3513	13 G x 3.5 in. T-handle needle	10 / case
TJC4011	11 G x 4 in. T-handle needle	10 / case
TJC4008	8 G x 4 in. T-handle needle	10 / case

Illinois sternal / iliac bone marrow aspiration needles

Reference no.	Product description	Packaging
TIN3018	T-handle 18 G x 3 in. (adjustable length 1–3 1/16 in. [26–79 mm])	10 / case

Standard trays with T-handle Jamshidi™ needle

Reference no.	Product description	Packaging
BCT3411	Tray with T-handle Jamshidi needle 11 G x 4 in.	10 / case
BCAT4511	Tray with T-handle Jamshidi needle 11 G x 4 in. and T-handle Illinois handle 15 G x 3 in.	10 / case

Soft tissue biopsy products

- BD offers a range of needles and trays that can be used for veterinary procedures such as kidney, liver, lung, prostate, breast and thyroid biopsies.
- Our proven brand, Temno™, is designed to offer procedural flexibility, accurate penetration and precise control.



Temno biopsy devices

Reference no.	Product description	Packaging
T206	20 G x 6 cm	5 / case
T209	20 G x 9 cm	5 / case
T2011	20 G x 11 cm	5 / case
CT2011	20 G x 11 cm, with introducer	5 / case
T2015	20 G x 15 cm	5 / case
CT2015	20 G x 11 cm, with introducer	5 / case
T2020	20 G x 20 cm	5 / case
CT2020	20 G x 20 cm, with introducer	5 / case
T186	18 G x 6 cm	5 / case
T189	18 G x 9 cm	5 / case
T1811	18 G x 11 cm	5 / case
CT1811	18 G x 11 cm, with introducer	5 / case



*Biopsy products are for single use only

The BD PhaSeal™ system

Our proprietary BD PhaSeal™ system is a leakproof and airtight closed-system drug transfer device. It mechanically prohibits the transfer of environmental contaminants into the system⁵ and also mechanically prohibits the escape of hazardous drugs, vapors or aerosols.

BD PhaSeal protector drug vial access device

Reference no.	Product description	Packaging
515100	protector (P14); equalizing capacity of 20 mL air; for vials with 31 mm necks	50 / box, 200 / transport box
515102	protector (P21); equalizing capacity of 20 mL air; for vials with 20 mm necks	50 / box, 200 / transport box
515104	protector (P28); equalizing capacity of 60 mL air; for vials with 28 mm necks	25 / box, 100 / transport box
515105	protector (P50); equalizing capacity of 20 mL air; for vials with 50 mm necks	25 / box, 100 / transport box

BD PhaSeal injector syringe safety device

Reference no.	Product description	Packaging
515003	Injector Luer lock (N35)	50 / box, 200 / transport box

BD PhaSeal IV line access

Reference no.	Product description	Packaging
515200	connector Luer lock (C35)	50 / box, 200 / transport box
515202	connector Luer lock (C45)	50 / box, 200 / transport box
515304	Y-site connector (C80); length 21 cm, priming volume 1.6 mL	30 / box, 120 / transport box

Additional BD veterinary products include:

- Hypodermic syringes and needles
- Blood collection products
- Spinal needles for anesthesia
- Disposal solutions



5. For up to 168 hours, within an ISO Class V environment following aseptic technique.



bd.com/animalhealth

BD, 1 Becton Drive, Franklin Lakes, NJ, 07417, USA

bd.com

