

A full line of restorative solutions to fit your unique needs.

| Versatility & Esthetics | | Efficiency & Simplicity | |
|--|--|---|---|
| 3M™ Filtek™ Supreme Ultra Universal Restorative <ul style="list-style-type: none"> • Exceptional esthetics • Excellent polish retention • Exceptional handling • Outstanding strength for anterior and posterior • Wear resistance similar to enamel* • Natural-looking fluorescence • True nanotechnology • 36 shades, 4 opacities | 3M™ Filtek™ Supreme Ultra Flowable Restorative <ul style="list-style-type: none"> • Excellent esthetics • Excellent adaptation • Excellent polish retention • Excellent wear resistance • True nanotechnology • 12 shades | 3M™ Filtek™ One Bulk Fill Restorative <ul style="list-style-type: none"> • One-step placement up to 5 mm • Increased opacity for improved esthetics • Excellent adaptation • Excellent handling • Excellent wear resistance • Stress relief • True nanotechnology • 5 shades | 3M™ Filtek™ Bulk Fill Flowable Restorative <ul style="list-style-type: none"> • One-step placement up to 4 mm • Low shrinkage stress • Easy adaptation • Excellent wear resistance • 4 shades |

*In vitro data on file.

We have you covered for every class of direct restoration.

| | Versatility & Esthetics | | Efficiency & Simplicity | |
|---|--|---|---|---|
| | 3M™ Filtek™ Supreme Ultra Universal Restorative (Incremental up to 2 mm) | 3M™ Filtek™ Supreme Ultra Flowable Restorative (Incremental up to 2 mm) | 3M™ Filtek™ One Bulk Fill Restorative (Bulk up to 5 mm) | 3M™ Filtek™ Bulk Fill Flowable Restorative (Incremental up to 4 mm) |
| Class I  | ✓ | ✓ (Liner/Base) | ✓ (up to 4 mm) | ✓ (Liner/Base) |
| Class II  | ✓ | ✓ (Liner/Base) | ✓ (up to 5 mm) | ✓ (Liner/Base) |
| Class III  | ✓ | ✓ | | ✓ |
| Class IV  | ✓ | | | |
| Class V  | ✓ | ✓ | | ✓ |
| Minimally Invasive | ✓ | ✓ | ✓ | ✓ |
| Stress Bearing | ✓ | | ✓ | |

✓ = recommended use

Refer to Instructions for Use for full list of indications.



www.3M.com/Filtek

3M Oral Care
2510 Conway Avenue
St. Paul, MN 55144-1000 USA
Phone 1-800-634-2249
Web 3M.com/dental

3M Canada
Post Office Box 5757
London, Ontario N6A 4T1
Canada
Phone 1-888-363-3685

3M and Filtek are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. © 3M 2017. All rights reserved.
70-2013-0807-2