

Miltex[®]
Redefining Excellence



GENERAL SURGICAL INSTRUMENTS
...the most popular patterns!



...for GENERAL INSTRUMENTS

SCISSORS

Miltex recognizes and meets a variety of clinical needs in the careful construction of our Surgical Scissors. Whether Dissecting, Cutting Tissue, Suture or Dressing material or performing any number of important functions, our finely crafted scissors perform at a level equal to the surgical task.

FORCEPS and RING-HANDLE FORCEPS

Miltex Forceps are for the most delicate procedural requirements as well as routine yet important surgical needs. It may be instrumentation for exacting Plastic Surgery with just the right handle or a simple "pick-up". Miltex is able to fulfill either demand with just the right instrument.

HEMOSTATS

Miltex uses not only the very best material but the finest craftsmanship to build into our instruments the surgical "feel" that is essential to achieve ideal hemostasis. Our hemostatic forceps are distinguished from all others by an unrelenting and complete commitment to creating an instrument which surpasses the most critical evaluation.

NEEDLE HOLDERS

A close examination of Miltex Needle Holders, either Standard or Tungsten Carbide, will reveal meticulously crafted surgical tools that will perform precisely as intended. Smoothly beveled box locks and jaw edges protect delicate suture material from snags or cutting. The jaw surface is milled to ensure a secure grasp and to prevent needle slippage or rotation. They are, very simply, the best.



SCISSORS

STANDARD (STD)



TUNGSTEN CARBIDE (TC)



SUPERCUT (SC)



STD	TC	SC	Description
5-122	5-122TC	5-SC-122	Mayo Scissors 5-1/2 Curved
5-120	5-120TC	5-SC-120	Mayo Scissors 5-1/2 Straight
5-126	5-126TC	5-SC-126	Mayo Scissors 6-3/4 Curved
5-124	5-124TC	5-SC-124	Mayo Scissors 6-3/4 Straight

5-182	5-182TC	5-SC-182	Metzenbaum Scissors 7 Curved
5-181	5-181TC		Metzenbaum Scissors 7 Straight
5-284		5-SC-284	Metzenbaum Scissors 5 Curved Delicate
5-283		5-SC-283	Metzenbaum Scissors 5 Straight Delicate

5-26	5-26TC		OR Scissors 5-1/2 Straight B/B
5-16	5-16TC	5-SC-16	OR Scissors 5-1/2 Straight S/B
5-6	5-6TC		OR Scissors 5-1/2 Straight S/S
5-36	5-36TC		OR Scissors 5-1/2 Curved S/S
5-46	5-46TC	5-SC-46	OR Scissors 5-1/2 Curved B/B
5-56	5-56TC		OR Scissors 5-1/2 Curved B/B

5-306	5-306TC	5-SC-306	Iris Scissors 4-1/2 Curved
5-304	5-304TC	5-SC-304	Iris Scissors 4-1/2 Straight

9-100			Spencer Stitch Scissors 3-1/2
9-108			Stitch Scissors 4-1/2 Angle, Delicate
9-102			Littauer Stitch Scissors 4-1/2
9-104			Littauer Stitch Scissors 5-1/2

5-514			Lister Bandage Scissors 5-1/2
5-516			Lister Bandage Scissors 7-1/4

5-800			Bandage & Utility Scissors 8-1/2
5-998			Bandage & Utility Scissors 6-1/2
5-1000			Bandage & Utility Scissors 7-1/2

5-122



5-182



5-26



5-306



9-102



5-514



5-800



Miltex® offers a large selection of our SuperCut Scissors as well as Standard and Tungsten Carbide.

SuperCut & Tungsten Carbide Brochure

FORCEPS STANDARD (STD)



STD Description

6-118	Adson Dressing Forceps 4-3/4 (serrated)
6-120	Adson Tissue Forceps 4-3/4 (1x2)
6-121	Adson Tissue Forceps 4-3/4 (1x2 and cross serrated)
6-123	Adson Tissue & Suture Forceps 4-3/4 (1x2 with tying platform)
6-124	Brown-Adson Tissue Forceps 4-3/4 (7x7 side grasping)
6-112	Hudson Dressing Forceps 4-3/4 (serrated)
6-114	Hudson Tissue Forceps 4-3/4 (1x2)
6-4	Dressing Forceps 4-1/2 (serrated)
6-6	Dressing Forceps 5 (serrated)
6-8	Dressing Forceps 5-1/2 (serrated)
6-40	Tissue Forceps 4-1/2 (1x2)
6-42	Tissue Forceps 5 (1x2)
6-44	Tissue Forceps 5-1/2 (1x2)
6-104	Semkin Dressing Forceps 5 (serrated)
18-1107	Swiss Cilia Forceps (suture)

6-118



6-123



6-124



6-114



6-4



6-40



6-104



18-1107



RING HANDLE FORCEPS STANDARD (STD)



STD Description

16-6	Allis Tissue Forceps 6, 4x5 teeth
16-8	Allis Tissue Forceps 6, 5x6 teeth
7-504	Backhaus Towel Clamp 3-1/2
7-506	Backhaus Towel Clamp 5-1/4
7-600	Foerster Sponge Forcep 9-1/2, Straight (serrated)
7-610	Foerster Sponge Forcep 9-1/2, Curved (serrated)
7-620	Bozeman Sponge Forceps 10-1/2, Straight
7-622	Bozeman Sponge Forceps 10-1/2, Double Curved
7-624	Bozeman Sponge Forceps 10-1/2, Double Curved, One Large Finger Ring

16-8




7-504



7-600



7-620



In addition to a broad selection of Stainless Steel Instruments, Miltex® also manufactures Surgical Disposables of exceptional quality.


Disposables Brochure

HEMOSTATS

STANDARD (STD)



STD	Description
7-4	Halsted Mosquito Forceps 5 Curved
7-2	Halsted Mosquito Forceps 5 Straight
7-26	Hartman Mosquito Forceps 3-1/2 Curved
7-24	Hartman Mosquito Forceps 3-1/2 Straight
7-44	Crile Forceps 5-1/2 Curved
7-42	Crile Forceps 5-1/2 Straight
7-38	Kelly Forceps 5-1/2 Curved
7-36	Kelly Forceps 5-1/2 Straight




NEEDLE HOLDERS

STANDARD (STD) TUNGSTEN CARBIDE (TC)



STD	TC	Description
8-40	8-40TC	Baumgartner Needle Holder 5-1/2
8-8	8-8TC	Halsey Needle Holder 5 (smooth)
	8-8A-TC	Halsey Needle Holder 5 (serrated)
8-4	8-4TC	Derf Needle Holder 4-3/4
8-50	8-50TC	Crile-Wood Needle Holder 6
8-6	8-6TC	Webster Needle Holder 5 (smooth)
8-15	8-15TC	Olsen-Hegar Needle Holder 5-1/2 (with suture scissors)
8-16	8-16TC	Olsen-Hegar Needle Holder 6-1/2 (with suture scissors)
8-42	8-42TC	Mayo-Hegar Needle Holder 5
8-44	8-44TC	Mayo-Hegar Needle Holder 6
8-46	8-46TC	Mayo-Hegar Needle Holder 7



The Miltex® Surgical Catalog illustrates one of the largest selections of General Surgical Instruments of any manufacturer.

Surgical Catalog



www.miltex.com

toll-free phone 866-854-8300

toll-free fax 866-854-8400

phone 717-840-9335

fax 717-840-9347

**Please call Customer
Service today for additional
information, an in-office
demonstration, to place an
order, or to receive other
Miltex catalogs.**

