

Integra® Miltex®

Sterilization Containers

Limit uncertainty by protecting
your instrument investment



INTEGRA®
LIMIT UNCERTAINTY

Summary of Integra Sterilization Container Product Features

Gaskets within retention plates and container perimeter ensure a secure seal of the container, lid and filter.

Snap on filter plate with a simple locking mechanism that helps provide easy change of disposable filters.

Round disposable filters with centering hole to assist in proper alignment.

Lid and Bottom manufactured from anodized aluminum alloy. Easy to maintain, corrosion resistant, lightweight yet durable, and conducts heat for material drying period.

4 ID Tag Windows (2 front & 2 back)

ID tags available in 6 colors and can be engraved with up to 18 characters per line.

Low profile latch easily and securely closes lid over container bottom.

Tamper Evident Lock – helps to identify a sterility breach.

Fold-down handles extend and stop at 90° to aid lifting.

Corner design allows for same size container stacking in storage.

2 Indicator Card Windows (1 front & 1 back). Indicator card has visible bar for clear indication that steam or gas sterilization has been achieved.

Note: Baskets are constructed with solid wire mesh to allow ultrasonic cleaning; handles can be positioned to assist with a secure grasp for removal under sterile conditions.

One-Half Size (11¼" x 11¾" x 4⅛") Containers

- Vacuum Sterilizer Containers - Solid/Perforated Bottom
- Gravity Sterilizer Containers - Perforated Bottom

Height	Solid Bottom Full Container	Solid Bottom Only	Perforated Bottom Full Container	Perforated Bottom Only
4⅛"	730400	730400B	730410	730410B
5¼"	730401	730401B	730411	730411B
6⅞"	730402	730402B	730412	730412B
7¾"	730403	730403B	730413	730413B
10¼"	730404	730404B	730414	730414B

Note: Complete container assembly only available in silver.
See below for other color options with lids and bottoms sold separately.



One-Half Size Colored Lids (bottoms sold separately)

Includes Filter Retention Plate and Latches

Color	Lid Only
Silver	730420
Red	730422
Gold	730423
Green	730424
Black	730425



One-Half Size Wire or Perforated Baskets

Height	Wire Basket	Perforated Basket
2"	740400	750400
3"	740401	750401
4"	740402	750402
Lid	-	750420



Wire Basket



Perforated Basket

Note: The wire baskets come with 'screw-on' feet and the perforated baskets come with 'snap-on' feet.

Three-Quarter Size (18¼" x 11¼" x 4¼") Containers

- Vacuum Sterilizer Containers - Solid/Perforated Bottom
- Gravity Sterilizer Containers - Perforated Bottom

Height	Solid Bottom Full Container	Solid Bottom Only	Perforated Bottom Full Container	Perforated Bottom Only
4¼"	730500	730500B	730510	730510B
5¼"	730501	730501B	730511	730511B
6"	730502	730502B	730512	730512B



Three-Quarter Size Colored Lids
(bottoms sold separately)

Includes Filter Retention Plate and Latches

Color	Lid Only
Silver	730520
Red	730522
Gold	730523
Green	730524
Black	730525




Three-Quarter Size Wire or Perforated Baskets

Height	Wire Basket	Perforated Basket
2"	740500	750500
3"	740501	750501
4"	740502	750502
Lid	-	750520



Wire Basket



Perforated Basket

Full-Size (23¼" x 11¼" x 4⅛") Containers

- Vacuum Sterilizer Containers - Solid/Perforated Bottom
- Gravity Sterilizer Containers - Perforated Bottom

Height	Solid Bottom Full Container	Solid Bottom Only	Perforated Bottom Full Container	Perforated Bottom Only
4⅛"	730600	730600B	730610	730610B
5¼"	730601	730601B	730611	730611B
5⅞"	730602	730602B	730612	730612B
7¾"	730603	730603B	730613	730613B
10¼"	730604	730604B	730614	730614B

Note: Complete container assembly only available in silver.
See below for other color options with lids and bottoms sold separately.



Full-Size Colored Lids (bottoms sold separately)

Includes Filter Retention Plate and Latches

Color	Lid Only
Silver	730620
Red	730622
Gold	730623
Green	730624
Black	730625




Full-Size Wire or Perforated Baskets

Height	Wire Basket	Perforated Basket
2"	740600	750600
3"	740601	750601
4"	740602	750602
Lid	-	750620



Wire Basket



Perforated Basket

Wide-Body (23½" x 15¾" x 4¾") Containers

- Vacuum Sterilizer Containers - Solid/Perforated Bottom
- Gravity Sterilizer Containers - Perforated Bottom

Height	Solid Bottom Full Container	Solid Bottom Only	Perforated Bottom Full Container	Perforated Bottom Only
4¾"	730701	730701B	730711	730711B
7"	730702	730702B	730712	730712B



Wide-Body Lid (bottoms sold separately)


Includes Filter Retention Plate and Latches

Color	Lid Only
Silver	730220




Wide-Body Wire Baskets

Height	Wire Basket
2"	740700
4"	740701



Instrument Stringers and Holders

Catalog Number	Description	Catalog Number	Description
3-70	Stringer 2½" x 8" Side Wide	3-86	Stringer 2½" x 6" Open
3-71	Stringer 2½" x 10" Side Wide	3-87	Stringer 2½" x 8" Open
3-72	Stringer 2½" x 12" Side Wide	3-88	Stringer 2½" x 10" Open
3-74	Stringer 2½" x 16" Side Wide	3-89	Stringer 2½" x 12" Open
3-75	Stringer 2½" x 6" W/Lock	3-90	Center Lock Stringer 2½" x 6"
3-76	Stringer 2½" x 8" W/Lock	3-91	Center Lock Stringer 2½" x 8"
3-77	Stringer 2½" x 10" W/Lock	3-92	Center Lock Stringer 2½" x 10"
3-78	Stringer 2½" x 12" W/Lock	3-93	Center Lock Stringer 2½" x 12"
3-79	Stringer 2½" x 16" W/Lock	3-94	Center Lock Stringer 2½" x 16"
3-80	Stringer 2½" x 18" W/Lock	3-95	Center Lock Stringer 2½" x 18"
3-81	Stringer 5" x 8" W/Lock	3-96	Closed U Stringer 2½" x 6"
3-82	Stringer 5" x 10" W/Lock	3-97	Closed U Stringer 2½" x 8"
3-83	Stringer 5" x 12" W/Lock	3-98	Closed U Stringer 2½" x 10"
3-84	Stringer 5" x 16" W/Lock	3-99	Closed U Stringer 2½" x 12"
3-85	Stringer 5" x 18" W/Lock		



Identification Tags

Each tag allows for 2 lines of custom text, 18 characters per line.

Color	Catalog Number
Silver	775100
Blue	775101
Red	775102
Gold	775103
Green	775104
Black	775105



Snap-on Basket Tags (with key ring)


Each small tag allows for 3 lines of custom text, 32 characters per line.
Each large tag allows for 5 lines of custom text, 32 characters per line.

Color	Small	Size	Large	Size
Silver	775120	2⅜" x 1½"	775130	2⅜" x 1¼"
Gold	775121	2⅜" x 1½"	775131	2⅜" x 1¼"
Red	775122	2⅜" x 1½"	775132	2⅜" x 1¼"
Blue	775123	2⅜" x 1½"	775133	2⅜" x 1¼"
Green	775124	2⅜" x 1½"	775134	2⅜" x 1¼"
Black	775125	2⅜" x 1½"	775135	2⅜" x 1¼"




Tamper Evident Locks and Indicator Cards

Catalog Number	Description
775400	Tamper Evident Locks, single use, blue, box of 1000
775300	Indicator Cards, box of 250




Tamper Evident Lock




Indicator Card

Retention Plate and Filters

Catalog Number	Description
775212INT	Retention Plate for Paper Filters, Round
775200	Disposable Filters for Regular Containers PKG/1000



Retention Plate




Disposable Filters

Instrument Rack

- 9 slots
- 5mm - 10mm diameter silicone bars
- Compatible with the 8" Full-Size container
- The rack fits inside the 2" or 3" wire or perforated basket, allowing ample storage underneath for other instruments and auxiliary items

Catalog Number	Description
790010	Instrument Rack



Silicone Holders

Catalog Number	Description	Size
775630	Silicone Holder #1	8mm x 26mm x 40mm
775631	Silicone Holder #2	8mm x 26mm x 40mm
775632	Silicone Holder #3	8mm x 26mm x 40mm
775633	Silicone Holder #4	10mm
775634	Silicone Holder #5	5mm
775650	Silicone Holder Clamps, 4 slot, plastic	-



#1



#2



#3




#4



#5



Silicone Mats

Catalog Number	Description
776400	One-Half Size Silicone Mat
776500	Three-Quarter Size Silicone Mat
776600	Full-Size Silicone Mat



Partition Sheets and Holding Pins

Catalog Number	Description	Size
775610	Partition Sheet	2" x 7/8" (50mm x 20mm)
775611	Partition Sheet	5" x 7/8" (130mm x 20mm)
775612	Partition Sheet	9 7/8" x 7/8" (225mm x 20mm)
775613	Partition Sheet	18 1/8" x 7/8" (460mm x 20mm)
775600	Holding Pin for Wire Baskets	
775601	Holding Pin for Partition Sheets	



Holding Pin for Wire Baskets



Holding Pin for Partition Sheets

Holding Clamps

Catalog Number	Size Range
775620	1/8" to 1/2" (7mm-12mm)
775621	1/2" to 5/8" (12mm-16mm)
775622	5/8" to 1" (16mm-26mm)
775623	1" to 1 1/2" (26mm-36mm)
775624	1 1/2" to 1 3/4" (36mm-45mm)



Screw-on Feet

Catalog Number	Description
780FT100	Plastic Foot for Baskets, Pkg/1

Tolerances in product specifications, due to the hand craftsmanship of the instruments, are subject to minor variations in dimensions and finishes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

integralife.com/contact

