


Limit uncertainty with quality instrumentation
from the brand you trust

McGivney Hemorrhoidal Ligator

- For use with non-natural rubber latex (NRL) o-rings
- Lightweight model with offset handle helps to provide a controlled grip and may provide better viewing access
- The shaft with ligating drum can be rotated 360 degrees for correct positioning and readily disassembled for thorough cleaning
- Stainless Steel loading cone facilitates placement of the o-ring bands onto the ligating drum



For Use With Non-NRL O-Rings


Reference	Descriptionw
28-154A	McGivney Hemorrhoidal Ligator Set for non-NRL applications - includes McGivney Hemorrhoidal Ligator 8" shaft (28-154B), McGivney Grasping Forceps 7-1/2" (28-152).
28-154B	McGivney Hemorrhoidal Ligator improved model for non-NRL applications - 8" (20.3 cm) working length. Supplied with Loading Cone (28-156) and 100 Latex-free O-Rings (28-155).
28-152	McGivney Hemorrhoidal Grasping Forceps, 7-1/2" (19.1 cm), angled shanks. Specially designed for use with McGivney Ligators.
28-155	Non-NRL O-Rings compatible with 28-154B Ligator only. 100 rings per bag. Provided Non-Sterile.
28-156	Loading Cone compatible with 28-154B Ligator only.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Always refer to the appropriate instructions for use for complete clinical instructions.
- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America
USA 800-654-2873 ▪ 888-980-7742 fax
International +1 609-936-5400 ▪ +1 609-750-4259 fax
integralife.com/contact



Integra York PA, Inc.
589 Davies Drive
York, PA 17402 ▪ USA