


Integra™ Miltex®

Gynecology Biopsy Punch Forceps

Limit uncertainty with
precision cutting performance
and optimal biopsy integrity



INTEGRA™
LIMIT UNCERTAINTY



589 Davies Drive

York, PA 17402

toll-free phone 866-854-8300

toll-free fax 866-854-8400

phone 717-840-9335

fax 717-840-9347

email customerservice.miltex@integralife.com


integralife.com/integra-miltex

Integra™ Miltex®

Premium Cervical Biopsy Forceps

- German stainless steel and craftsmanship
- Designed to facilitate a clean cut and help reduce crush artifact during biopsies
- Locking mechanism protects the jaws during processing and storage
- Cutting stop limits depth of tissue sample

Reference	Description	Bite Size/Shape	Jaws	Length
301445WL	Townsend Biopsy Forceps with Lock	4 mm x 2 mm, mini bite	single tooth on lower jaw	7-3/4" (19.7 cm)
301443WL	Burke (Modified Baby Tischler) Biopsy Forceps with Lock	5 mm x 3 mm, oval	single tooth on both jaws	7-3/4" (19.7 cm)
301450WL	Baggish Biopsy Forceps with Lock	6.5 mm x 4 mm, oval	two teeth on lower jaw	7-3/4" (19.7 cm)
301442WL	Tischler Biopsy Forceps with Lock	7 mm x 3 mm, oblong	single tooth on both jaws	9-1/4" (23.5 cm)



* Bites approximate to size and shape. Tolerances in product specifications, due to the hand craftsmanship of the instruments, are subject to minor variations in dimensions and finishes.

Integra™ Miltex®

Premium Cervical Biopsy Forceps

- German stainless steel and craftsmanship
- Designed to facilitate a clean cut and help reduce crush artifact during biopsies
- Locking mechanism protects the jaws during processing and storage
- Cutting stop limits depth of tissue sample

Reference	Description	Bite Size/Shape	Jaws	Length
301462WL	Baby Tischler Biopsy Forceps with Lock	4 mm x 2.3 mm, mini bite	single tooth on both jaws	7-3/4" (19.7 cm)
301461WL	Coppleson Biopsy Forceps with Lock	7.5 mm x 3 mm, tapered rectangle	single tooth on both jaws	7-3/4" (19.7 cm)
301460WL	Appex Biopsy Forceps with Lock	7 mm x 3 mm, bullet shape	angled up tip	7-3/4" (19.7 cm)
301463WL	Kevorkian-Younger Biopsy Forceps with Lock	8 mm x 3 mm, rectangular	teeth on lower jaw	7-3/4" (19.7 cm)
301464WL	Tischler-Kevorkian Biopsy Forceps with Lock	9.5 mm x 3mm, rectangular	single tooth on lower jaw	7-3/4" (19.7 cm)



301462WL

301461WL

301460WL

301463WL

301464WL

* Bites approximate to size and shape. Tolerances in product specifications, due to the hand craftsmanship of the instruments, are subject to minor variations in dimensions and finishes.

Integra™ Miltex®

Premium Uterine Biopsy Forceps

- German stainless steel and craftsmanship
- No cutting stop feature to maximize depth of tissue sample
- Curved bites and angled shanks facilitate easier access for hard-to-reach biopsy sites

Reference	Description	Bite Size/Shape	Jaws	Shanks	Length
30-1400	Gaylor Uterine Biopsy Forceps	5 mm, round	curved, interfitted cups	straight	8-1/2" (21.6 cm)
30-1405	Thoms-Gaylor Uterine Biopsy Forceps	5 mm, round	curved, interfitted cups	angled	8-1/4" (21.0 cm)
30-1410	Thoms-Gaylor Uterine Biopsy Forceps	5 mm, round	curved, interfitted cups	angled	9-1/4" (23.5 cm)
30-1440	Gellhorn Uterine Biopsy Forceps	8 mm x 5 mm, oval	curved	angled	9-1/2" (24.1 cm)



30-1400



30-1405



30-1410



30-1440
Top View




30-1440
Bottom View

* Bites approximate to size and shape. Tolerances in product specifications, due to the hand craftsmanship of the instruments, are subject to minor variations in dimensions and finishes.

- German stainless steel and craftsmanship
- Shank design allows the physician to maximize pressure applied during sampling
- Angle and shape of Schubert bite pattern facilitates easier access for hard-to-reach biopsy sites
- Cutting stop limits the depth of tissue sample

Reference	Description	Bite Size/Shape	Jaws	Shanks	Length
30-1425	Schubert Uterine Biopsy Forceps	9 mm x 5 mm, triangular	curved	angled	8" (20.3 cm)
30-1420	Schubert Uterine Biopsy Forceps	9 mm x 5 mm, triangular	-	straight	9-1/2" (24.1 cm)
30-1430	Schubert Uterine Biopsy Forceps	9 mm x 5 mm, triangular	-	angled	9-1/2" (24.1cm)
30-1435	Schubert Uterine Biopsy Forceps	9 mm x 5 mm, triangular	curved	S-shape	10-3/4" (27.3 cm)



30-1425

30-1420


30-1430

30-1435

* Bites approximate to size and shape. Tolerances in product specifications, due to the hand craftsmanship of the instruments, are subject to minor variations in dimensions and finishes.

- German stainless steel and craftsmanship
- Longer, thinner shafts improve site visualization and are useful in the treatment of bariatric patients
- Cutting stop limits depth of tissue sample
- Angle and shape of Wittner bite pattern facilitates easier access for hard to reach biopsy sites

Reference	Description	Bite Size/Shape	Jaws	Length
30-1480	Miltex Uterine Biopsy Forceps	6 mm, round	non interfitted cups	9" (22.9 cm)
30-1485	Eppendorfer Uterine Biopsy Forceps	8.5 mm x 4.5 mm, oblong	no teeth	8-1/2" (21.6 cm)
30-1487	Wittner Uterine Biopsy Forceps	8 mm x 3.5 mm, tapered oblong	straight tip, teeth on lower jaw	8-1/2" (21.6 cm)
30-1486	Wittner Uterine Biopsy Forceps	8 mm x 3.5 mm, tapered oblong	curved tip, teeth on lower jaw	8-1/2" (21.6 cm)
30-1482	Kevorkian-Younger Uterine Biopsy Forceps	8 mm x 3.5 mm, rectangular	teeth on lower jaw	9-1/2" (24.1 cm)



* Bites approximate to size and shape.
Tolerances in product specifications, due to the hand craftsmanship of the instruments, are subject to minor variations in dimensions and finishes.

- Shafts are interchangeable with pistol grip handle
- Shafts rotate 360 degrees allowing easier access to target tissue
- Longer, thinner shafts improve site visualization and are useful in the treatment of bariatric patients
- Two-piece design facilitates cleaning

Reference	Description	Bite Size/Shape	Jaws	Length
30-1492	Eppendorfer-Krause Rotating Biopsy Forceps with Pistol Grip Handle	6.5 mm x 4.5 mm, oval	-	9" (22.9 cm)
30-1493	Eppendorfer-Krause Biopsy Forceps (Shaft Only)	6.5 mm x 4.5 mm, oval	-	9" (22.9 cm)
30-1494	Tischler Rotating Biopsy Forceps with Pistol Grip Handle	7 mm x 3 mm, oblong	single tooth on both jaws	9" (22.9 cm)
30-1495	Tischler Biopsy Forceps (Shaft Only)	7 mm x 3 mm, oblong	single tooth on both jaws	9" (22.9 cm)
30-1490	Kevorkian-Younger Rotating Biopsy Forceps with Pistol Grip Handle	8 mm x 3.5 mm, rectangular	teeth on lower jaw	9" (22.9 cm)
30-1491	Kevorkian-Younger Rotating Biopsy Forceps (Shaft Only)	8 mm x 3.5 mm, rectangular	teeth on lower jaw	9" (22.9 cm)
28-400	Pistol Grip Spring Handle for Rotating Biopsy Forceps	-	-	-



* Bites approximate to size and shape. Tolerances in product specifications, due to the hand craftsmanship of the instruments, are subject to minor variations in dimensions and finishes.


Integra™ Miltex®

Standard Grade Cervical and Uterine Biopsy Forceps

- Crafted from surgical grade German stainless steel
- Designed to facilitate operating room procedural needs
- Locking mechanism protects the jaws during processing and storage*

Reference	Description	Bite Size/Shape	Jaws	Shanks	Length
MH301443WL*	Burke (Modified Baby Tischler) Biopsy Forceps with Lock	4.5 mm x 3 mm, oval	single tooth on both jaws	-	7-3/4" (19.7 cm)
MH301442WL*	Tischler Biopsy Forceps with Lock	7 mm x 3 mm, oblong	single tooth on both jaws	-	8" (20.3 cm)
MH30-1405	Thoms-Gaylor Uterine Biopsy Forceps	5 mm, round	curved, interfitted cups	angled	8-1/2" (21.6 cm)
MH30-1485	Eppendorfer Uterine Biopsy Forceps	8.5 mm x 4 mm, oblong tapered	no teeth	-	8-1/2" (21.6 cm)
MH30-1482	Kevorkian-Younger Uterine Biopsy Forceps	8 mm x 3.5 mm, rectangular	teeth on lower jaw	-	9-1/2" (24.1 cm)

*Available on MH301443WL and MH301442WL only



* Bites approximate to size and shape. Tolerances in product specifications, due to the hand craftsmanship of the instruments, are subject to minor variations in dimensions and finishes.

Standard Colposcopy Set

Reference	Description	Bite Size/Shape	Jaws/Tips	Length
30-15	Graves Vaginal Speculum, Medium	-	-	1-3/8" (3.5 cm) x 4" (10.2 cm)
30-55	Pederson Vaginal Speculum, Narrow Blades, Medium	-	-	7/8" (2.2 cm) x 4" (10.2 cm)
30-951	Emmett Uterine Tenaculum Hook, Style 1	-	sharp, half-curved tip	9" (22.9 cm)
30-1350	Kogan Endocervical Speculum with Gauge & Set Screw	-	fenestrated, 7mm x 25 mm	9-1/2" (24.1 cm)
30-1382	Kevorkian-Younger Endocervical Curette with Basket	-	loop, 3 mm x 12 mm	12" (30.5 cm)
301442WL	Tischler Biopsy Forceps with Lock	7mm x 3mm, oblong	single tooth on both jaws	9-1/4" (23.5 cm)



* Sets may be customized to meet procedural needs

* Tolerances in product specifications, due to the hand craftsmanship of the instruments, are subject to minor variations in dimensions and finishes.

Delicate Colposcopy Set

Reference	Description	Bite Size/Shape	Jaws/Tips	Length
30-10	Graves Vaginal Speculum, Small	-	-	3/4" (1.9 cm) x 3" (7.6 cm)
30-50	Pederson Vaginal Speculum, Narrow Blades, Small	-	-	1/2" (1.3 cm) x 3" (7.6 cm)
30-951	Emmett Uterine Tenaculum Hook, Style 1	-	sharp, half-curved tip	9" (22.9 cm)
30-1351	Miltex-Kogan Endocervical Speculum with Ratchet & Gauge	-	fenestrated, 3 mm x 25 mm	11" (27.9 cm)
30-1383	Kevorkian-Younger Endocervical Curette without Basket	-	loop, 3 mm x 12 mm	12" (30.5 cm)
301443WL	Burke (Modified Baby Tischler) Biopsy Forceps with Lock	5mm x 3mm, oval	single tooth on both jaws	7-3/4" (19.7 cm)



* Sets may be customized to meet procedural needs

* Tolerances in product specifications, due to the hand craftsmanship of the instruments, are subject to minor variations in dimensions and finishes.

Integra™ Miltex® Instrument Care & Cleaning Products

Instrument Prep Enzyme Foam

Reference	Description
3-760	24 oz. (0.71 liter) spray bottles (case of 12)*



3-760

EZ-Zyme® All Purpose Enzyme Cleaner

- Easy to Use
- Neutral pH
- Multi-use as an instrument presoak, ultrasonic cleaner, laundry presoak, and evacuation system cleaner
- Space-saving or bulk options

Reference	Description
3-750	32 packets (each packet, 3/4 oz., 22 ml) 1 box makes 32 gallons of solution
3-755	1-gallon bottle (3.79 liter), makes 170 gallons of solution



3-755



3-750

Surgical Instrument Cleaner

- Neutral pH
- Biodegradable
- Rinses free, leaves no residue
- Phosphate-free

Reference	Description
3-720	8 fl. oz. (0.24 liter) bottles (case of 12), makes 32 gallons of solution*
3-725	1-gallon bottle (3.79 liter), makes 512 gallons of solution



3-725



3-720

*Contact Distributor for availability of individual units

Surgical Instrument Cleaner & Lubricant

- Cleaner and brighter instruments
- Removes stains, blood and mineral deposits
- Biodegradable
- Concentrate, lower cost per use

Reference	Description
3-710	1 Gallon, 128 fl. oz. (3.79 liter) (see label for concentration rates)

Note: when using this product, no additional lubrication is required



3-710

Instrument Lubricant

- Formulated to extend the working life of surgical instruments
- Helps prevent corrosion
- Formulated with gentle, non-aggressive ingredients
- Frees up box locks - lubricates moving parts of surgical instruments

Reference	Description
3-705	1 Gallon, 128 fl. oz. (3.79 liter) (see label for concentration rates)



3-705

Spray Lube (Surgical Instrument Lubricant)

Reference	Description
3-700	8 fl. oz. (0.24 liter) bottles (case of 12)*



3-700

*Contact Distributor for availability of individual units

Surgical Instrument Stain Remover

- Effective on instruments, autoclaves, ultrasonic cleaners or any stainless steel, tile or plated surface
- Will not harm instrument surfaces
- Easy, fast, convenient to use
- Economical

Reference	Description
3-740	3 oz. (85 grams), (case of 12)*



3-740

*Contact Distributor for availability of individual units

Nylon & Stainless Steel Brushes

- Nylon bristles for effective cleaning of instrument box locks
- Stainless steel bristles for cleaning serrations
- Autoclavable



3-1000



3-1001



21-379



23-1600



23-1602

Reference	Description	Type	Length
3-1000	Instrument Cleaning Brushes	nylon bristles (set of 3)	7-1/4" (18.4 cm)
3-1001	Instrument Cleaning Brushes	stainless steel bristles (set of 3)	7-1/4" (18.4 cm)
21-379	Stainless Wire Brush	-	6-1/4" (15.9 cm)
23-1600	Cleaning Brushes for Trachea and Laryngectomy Tubes	small brush for tubes sizes 2 To 5	8" (20.3 cm)
23-1602	Cleaning Brushes for Trachea and Laryngectomy Tubes	large brush for tubes sizes 6 To 12	6" (15.2 cm)

*Contact Distributor for availability of individual units

Instrument Care Kit System

A great way to try various Integra™ Miltex® instrument care products!

Reference	Description
3-800	Spray Lube (3-700), Surgical Instrument Cleaner (3-720), Surgical Instrument Stain Remover (3-740), Instrument Care DVD (3-810), Nylon & Stainless Steel Brushes (3-1000, 3-1001)



3-800

Instrument Care Training

A 25-minute Instructional DVD assists practitioners and staff to properly care for instruments and to keep them in top working condition. The DVD offers step-by-step guidance on handling, cleaning, rinsing, inspecting, lubricating and sterilizing surgical instruments.

The step-by-step poster (SURGHIPPOST) provides easy instructions to clean and prepare instruments for sterilization.

Try a variety of our popular cleaning products with a H.I.P. sample box.

Reference	Description
3-810	Instrument Care DVD
SURGHIPPOST	H.I.P. Poster
SURGHIPBOX	H.I.P. Sample Box

SURGHIPBOX



SURGHIPPOST

3-810



For more information or to place an order, please contact:
Integra ▪ 589 Davies Drive, York, PA 17402
866-854-8300 USA ▪ 717-840-2763 outside USA ▪ 717-840-9335 fax
integralife.com/integra-miltex

Integra and the Integra logo are trademarks of, and Miltex is a registered trademark of,
Integra LifeSciences Corporation or its subsidiaries.
©2011 Integra LifeSciences Corporation. All rights reserved

0112/6/A/1 SURGGYNBPFROC PM11262

