

PacEndo™ Vital-Ice™ Pulp Vitality Refrigerant Spray


Availability: In Stock



Short Description

- A quick way to check cold sensitivity.
- Significantly colder than ethyl chloride or ice.
- For accuracy use with micro applicators or cotton pellets.
- Directional nozzle for easy placement on micro applicator or cotton pellet.

Description



Vitality testing with PacEndo™ Vital-Ice™ is a critical part of the endodontic diagnosis, in addition to a complete dental exam with exposed radiographs. Stimulating the nerve endings of the pulp with a cold spray is a effective way to test for vitality.

The results of this test will indicate the vitality (existing blood flow) of the tooth before starting the procedure. Note: pain relief drugs may cause a false reading that the tooth is vital but does not respond to cold sensation. Adjacent metal restorations can also give a false reaction to the cold test.

- A quick way to check cold sensitivity.
- Significantly colder than ethyl chloride or ice.
- For accuracy use with micro applicators or cotton pellets.
- Directional nozzle for easy placement on micro applicator or cotton pellet.



Research
Text...

