

TPC

Mirage

Chair Mounted Operatory System



#MP2015



Model #MP2000 with touchpad



Model #MP2000 with touchpad and LCD

MirageOperatorySystem

Model 4000 Dental Chair

Standard Features

Two operator-programmable position, one preset exit position
Electronic foot switch preposition / auto return switch
Fuse / surge protection
Toe guard protection
Synchronized seat and toe board tilt
All metal construction
Thin back with lumbar support
Double articulating headrest
Safety shut-off plate
Pivoting armrest
Seamless asepsis upholstery
Input voltage (110 V or 220 V)- specify
5 year warranty

Mirage Chair Mounted Delivery System

Standard Features

Asepsis automatic control for three handpieces with syringe
Individual handpiece pressure adjustments
Non-retracting water coolant with on/off
Individual water coolant flow controls
Air coolant flow control
Handpiece tubing flush system
Ultra soft handpiece tubing
4 hole straight gray asepsis tubing with Midwest connector
Stainless steel tray with non-skid pad
Disc foot control wet/dry toggle with 7" of smooth tubing
Equipment utility center with asepsis housing
Vacuum canister with solids collector screen
Straight light post for unit mount dental light
Pre-wiring for fiber optic light source
Air-actuated counterbalanced arm lock

Mirage Post Mount Utility Center

Standard Features

Gravity drain swing away porcelain cuspidor
Adjustable timed bowl rinse with push-button cup filler
Self contained water system with city/ bottle selector
Water outlet with flow control
Air outlet with flow control
Vacuum canister with solids collector screen

Assistant's Instrumentation

Standard Features

Telescoping arm with movable assistant's instrument holder
One autoclavable HVE with asepsis tubing
One autoclavable saliva ejector with tubing

Dental Operatory Light

L700-1 Advance Operatory Light, 12,000 lux / 24,000 lux, 4,500° Kelvin

L825 Lucent Operatory Light (sensor), 29,000 lux / 17,000 lux, 4,600° Kelvin

L800 Lucent Operatory Light (no sensor), 29,000 lux / 17,000 lux, 4,600° Kelvin

Factory Options

Additional 3 way syringe
Additional HVE with tubing
Power optic system, ISO B (5-hole)
Power optic system, ISO C (6-pin)
Built-in piezo scaler with transformer
Super thin x-ray viewer- unit mounted
Standard LCD monitor bracket
Deluxe LCD monitor bracket
Mirage LCD monitor bracket
15" standard multi media LCD monitor
Plastic utility shelf- (grey, white or beige)-post mounted
Touch pad-unit head mounted

Item Number

#TPC3WAY+T
#P2335+TB
#H7000
#H7660
#800N
#XV100/150
#WIC800
#WIC825
#WIC840
#M015
#US100
#TP-2000

Packing Dimensions

Chair: 54" x 34" x 31"
Cushion: 51" x 9" x 23"
Unit: 51" x 21" x 17"
PMU: 23" x 19" x 21"
Light: 37" x 17" x 15"

Upholstery Selections:

Mauve DU-03
Desert Rose DU-06
Purple Grey DU-07
Dusty Jade DU-11
Wedgewood DU-14
Grey Teal DU-15
Taupe DU-16
Greystone DU-19
Black DU-100
Cobalt DU-116
Atlantis DU-201
Autumn DU-202
Tumbleweed DU-210
Cinnamon DU-213
Stormy DU-220

