

PROVISION™ SAFETY EYEWEAR

The Insight to Safeguard Your Eyesight

Palmero's Safety Eyewear is a sight to behold. From broad range of protective products consisting of full face shields and goggles to more procedure specific systems like laser safety and bonding glasses, our protective eyewear offers the ultimate in comfort, convenience, and affordability. Give your eyesight the protection it deserves with Palmero Healthcare Safety Eyewear.

The Three C's of Safety Eyewear Compliance

Certification + Clarity + Comfort = Compliance

CERTIFIED Meets industry specifications such as ANSI.Z87.1, OSHA, CDC Guidelines to provide wraparound protection








CLARITY High level of visibility provided by anti-fog coating, scratch resistant technologies, and superior lens designs

COMFORT Broad range of options to provide an optimal fit for all-day security and comfort

COMPLIANCE Willingness and discipline to use safety eyewear for everyone, in every procedure, every day



Eyewear Collections at a Glance

Traditional Safety Eyewear					Disposables	
Wraparound	Goggles	Bifocals	Bonding	Laser	Face Shields	Eyewear
						
*7 styles *Ergonomic *Lightweight *Full & Minis **Adjustability	*4 styles *Wear over prescription glasses *Side shields *Fog-free **Adjustability	*4 styles *Up to +3.0 diopter **Adjustability	*7 styles *Red and orange polycarbonate lens options **Adjustability	*3 styles *Yellow and green lens *All include case & lanyard	*2 styles, NEW Snap-On or Foam Elastic Band style *Full & half sizes *Anti-fog, Anti-static	*1 frame, 4 lens color options *Snap-On feature *Anti-fog, Anti-static
**Adjustability = telescopic/ratcheting feature on selected styles						

To view Palmero full eyewear line, visit www.palmerohealth.com and go to CATALOG on the navigation menu for online access or download version.

PROVISION™ SAFETY EYEWEAR

Find the right fit of safety eyewear to protect and enhance your practice, staff, and patients!

Wraparound Collection

ProVision™ Econo Wraps™ & Mini Econo Wraps™

- Sleek, feather-weight polycarbonate wraps, universal fit
- Available in 2 colors and sizes
- Meets ANSI Z87.1-2003 standard

Ref#	Frame/Lens	Type	Facial Fit
3601C	Clear/Clear	Clinician/Patient	Small/Narrow/Medium
3607C	Clear/Clear Mini	Clinician/Patient	Small/Narrow/Medium
3601G	Grey/Grey	Patient	Small/Narrow/Medium
3607G	Grey/Grey Mini	Patient	Small/Narrow/Medium

ProVision™ Econo Wraps™ 10-Pack Dispenser

- 10-Pack Bulk Dispenser with the same features as the full size Econo Wraps™

Ref#	Frame/Lens	Type	Facial Fit
3601CB	Clear/Clear	Clinician/Patient	Small/Narrow/Medium
3601GB	Grey/Grey	Patient - NEW	Small/Narrow/Medium

ProVision™ Contour Wraps™

- Lightweight polycarbonate frame with telescoping arms ensure a great contoured fit
- Meets ANSI Z87.1-2003 standard

Ref#	Frame/Lens	Type	Facial Fit
3552	Black/Grey	Patient	Small/Narrow/Medium
3553	Black/Clear	Clinician/Patient	Small/Narrow/Medium

ProVision™ Bad Dogs™

- Designed to provide an uninterrupted peripheral view and ultimate protection
- Meets ANSI Z87.1-2003 standard

Ref#	Frame/Lens	Type	Facial Fit
3620	Clear/Clear	Clinician/Patient	Small/Narrow/Medium

ProVision™ Chic™

- Designed to fit smaller, narrow facial structure with colorful rubber tipped temples and soft rubber nose bridge
- Meets ANSI Z87.1-2003 standard

Ref#	Frame/Tips/Lens	Type	Facial Fit
3760CP	Clear/Pink/Clear	Clinician/Patient	Small/Narrow/Medium



PROVISION™ SAFETY EYEWEAR

Wraparound Collection

ProVision™ See-Breeze

- Ultra-lightweight frames with sculpted dual lenses and telescoping arms for a comfortable fit
- Great for children or adults with smaller faces
- Meets ANSI Z87.1-2003 standard

Ref#	Frame/Lens	Type	Facial Fit
3560BL	Black/Clear	Clinician/Patient	Universal
3560B	Blue/Clear	Clinician/Patient	Universal
3561BL	Black/Grey	Patient	Universal
3560PL	Platinum/Clear	Clinician/Patient	Universal
3560R/3561R	Clear/Grey	Replacement Lenses	----

ProVision™ Tech Specs™ Safety

- Lightweight design that provides excellent protection and comfort
- Contoured frame style hugs the face with rubberized (non-latex) attachments on the nosepiece and ear rests
- Meets ANSI Z87.1-2003 standard

Ref#	Frame/Lens	Type	Facial Fit
3710C	Grey/Clear	Clinician/Patient	Small/Narrow/Medium
3711C	Blue/Clear	Clinician/Patient	Small/Narrow/Medium
3710G	Grey/Grey	Patient	Small/Narrow/Medium

Goggles Collection

ProVision™ End-Fog

- Extra anti-fog coating sets End-Fog apart from other economy eyewear
- Secure snug fit and optically correct, crystal-clear lenses for those special procedures
- Meets ANSI Z87.1-2003 standard

Ref#	Frame/Lens	Type	Facial Fit
3556C	Clear/Clear	Clinician/Patient	Universal
3556P	Purple/Clear	Clinician/Patient	Universal

ProVision™ Eyesavers™

- Ultra-light, economical, durable no-frills eye protection
- Hinged temples allow for convenient storage when not in use
- Meets ANSI Z87.1-2003 standard

Ref#	Frame/Lens	Type	Facial Fit
3S	Grey/Grey	Patient	Universal
17S	Clear/Clear	Clinician/Patient	Universal
1968	Clear/Clear, 10-Pack	Clinician/Patient	Universal
18S	Green/Green	Patient	Universal



PROVISION™ SAFETY EYEWEAR

Bifocals Collection

ProVision™ Bifocal

- Fog-free, dual curvature mono-lens with integrated side shields
- Fully adjustable, lens pivots with telescopic arms
- Meets ANSI Z87.1-2003 standard

Ref#	Diopter	Frame/Lens	Type	Facial Fit
3701B	1.5	Black/Clear	Clinician	Universal
3701C	2.0	Black/Clear	Clinician	Universal
3701D	2.5	Black/Clear	Clinician	Universal
3701E	3.0	Black/Clear	Clinician	Universal

ProVision™ Tech Specs™ Bifocal

- Lightweight design that provides excellent protection and comfort
- Contoured frame style hugs the face with rubberized (non-latex) attachments on the nosepiece and ear rests
- Meets ANSI Z87.1-2003 standard

Ref#	Diopter	Frame/Lens	Type	Facial Fit
3720B	1.5	Black/Clear	Clinician	Small/Narrow
3720C	2.0	Black/Clear	Clinician	Small/Narrow
3720D	2.5	Black/Clear	Clinician	Small/Narrow
3720E	3.0	Black/Clear	Clinician	Small/Narrow

ProVision™ Cool Wraps™ Bifocal

- Sleek, lightweight wraparound design
- Over 95% UVA and UVB protection
- Meets ANSI Z87.1-2003 standard

Ref#	Diopter	Frame/Lens	Type	Facial Fit
3730A	1.0	Clear/Clear	Clinician	Small/Narrow
3730B	1.5	Clear/Clear	Clinician	Small/Narrow
3730C	2.0	Clear/Clear	Clinician	Small/Narrow

Eyewear Accessories

Dynamic Disposables Eyewear Cords

- Disposable, flexible plastic eyewear cord
- Eliminate lost time searching for misplaced eyewear
- Allows eyewear to hang conveniently around neck when not in use
- Available in plastic or string style cord, 24" length cords

Ref#	Description	Quantity
3913	Black plastic cords	10-Pack
3904	Multicolor classic string style cords	5-Pack



Eyewear not included

PROVISION™ & DYNAMIC DISPOSABLES SAFETY EYEWEAR

Bonding Collection

ProVision™ Bonding Wraps™

- Durable lightweight eyewear with sleek styling and a close fit
- Meets ANSI Z87.1-2003 standard, filters over 95% of blue, UVA & UVB light

Ref#	Frame/Lens	Type	Facial Fit
3600	Red bonding wraps	Clinician/Patient	Small/Narrow/Medium

Dynamic Disposables Bonding Frames and Lenses

- Disposable lens with reusable ergonomically frame eyewear that can be worn over prescription glasses, filters over 95% of blue, UVA & UVB light

Ref#	Description	Type	Facial Fit
3901	Black frames 10-Pack	Clinician/Patient	Universal
3906/3906B	Orange Bonding 10/100-Pack	Clinician/Patient	Universal

ProVision™ Econo Bonding

- One piece design with detachable side shields
- Meets ANSI Z87.1-2003 standard, filters over 95% of blue, UVA & UVB light

Ref#	Frame/Lens	Type	Facial Fit
1100	Red Econo Bonding	Clinician/Patient	Universal

ProVision™ Specs™ Bonding

- Stylish, impact-and-stress resistant material
- Meets ANSI Z87.1-2003 standard, filters over 95% of blue, UVA & UVB light

Ref#	Frame/Lens	Type	Facial Fit
3712	Grey/Red Bonding	Clinician/Patient	Small/Medium

Laser Collection

ProVision™ Laser Safety Eyewear

- Diode Alexandrite design with green polycarbonate filter lens
- Optical density >5+ from 800-980 nm (3588DA)/800-830 nm (3589YG & 3590YG) and 50% (3588DA) / 68% (3589YG & 3590YG) visible light transmission
- Rubber nose pad and rubberized temple bars provide comfortable grip while the wrap around design ensures your eyes are protected
- All models include case and lanyard
- Meets ANSI Z87-1 & Z136-1 standard for laser protection

Ref#	Frame/Lens	Type	Facial Fit
3588DA	Silver & Black/Green	Clinician	Universal
3589YG	Black & Yellow/Yellow-Green	Clinician	Universal
3590YG	Black/Yellow-Green	Clinician	Universal



PROVISION™ SAFETY EYEWEAR

Find the right fit to enhance your practice, staff, and patients!

Eyewear Maintenance

ProVision™ Anti-Fog Solution

- Compatible with all eyewear, face shields, mirrors, loupes, and scopes
- Resists scratches, static, and stains
- Foam formula yields 5 times the applications compared to a traditional pump or squeeze bottle
- Available in 8 oz. and 1.7 oz. bottles

Ref#	Description	Product Qty	Case Qty
3533	ProVision Anti-Fog Solution 8 oz.	1 bottle	16 bottles
3536	ProVision Anti-Fog Solution 1.7 oz.	1 bottle	24 bottles



ProVision™ Crystal Clear Optical Cleaning Wipes

- Convenient pre-moistened optical wipes
- Contains Ultra Crystal Anti-Fog Clear Technology providing anti-fog, anti-static lens cleaning, and polishing solution
- Cleans dirt, grease, oils, and blood-borne pathogens from eyewear

Ref#	Description	Product Qty
3534	ProVision Crystal Clear Optical Cleaning Wipes 160 count, 6" x 6 3/4"	1 canister



ProVision™ Lens Cleaning Station

- Convenient eyewear and lens cleaning station
- Designed for countertop placement or wall mounting with hanging hardware
- Station includes 8 oz. pump spray bottle (same formulated liquid to our Optical Cleaning Wipes) and 600 optical equipment dry wipes

Ref#	Description	Product Qty
3535	ProVision Lens Cleaning Station 7 3/4" x 8" x 5"	1 station



For more information, visit palmerohealth.com,
call 800.344.6424 or email customerservice@palmerohealth.com.