

OCTOBER 9, 2018
F0178



MORE INVENTIVE LESS INVASIVE

2018 U.S.A. CATALOG

INDEX.....	3
PIEZO ULTRASONIC SURGERY	
PIEZOTOME CUBE and Implant Center II LED.....	4
PIEZOTOME II and PIEZOTOME Tip Kits.....	5
PIEZOTOME II and PIEZOTOME Tips.....	6
PIEZOTOME II Handpieces, Cords and Accessories.....	7
PIEZOTOME Handpieces, Cords and Accessories.....	8
PIEZO ULTRASONIC SCALERS & TIPS	
NEWTRON, NEWTRON B.LED, NEWTRON B.LED XS & Accessories.....	9
SP NEWTRON & SP NEWTRON LED & Accessories.....	10
NEWTRON Tips on One Page.....	11
NEWTRON Tips.....	12
NEWTRON Tip Kits.....	13
Accessories for P5, P5 LED, P5XS, P5XS LED & P MAX NEWTRON XS.....	14
AIR-N-GO EASY and AIR POLISHING ACCESSORIES	
AIR-N-GO EASY Kits.....	15
AIR-N-GO Powders.....	16
AIR-N-GO Accessories.....	17
AIR MAX, PROPHY MAX Accessories.....	18
CURING LIGHTS	
Mini LED Black, SuperCharged, Ortho2 & Mini LED.....	19
OEM Mini LED CURING LIGHT & Accessories for all MINI LED LIGHTS.....	20
GINGIVAL RETRACTION PASTES	
EXPASYL and HEMOSTASYL.....	21
RISKONTROL & BIOSPLINT	
RISKONTROL Tips and Adapters.....	22
BIOSPLINT Ribbons and Accessories.....	23
ACTEON IMAGING	
X-Mind Trium.....	24
X-Mind Unity and SOPIX Inside.....	25
PSPIX ² and SPIX ² Accessories.....	26
SOPIX ² and SOPIX ² Accessories.....	27
SoproCARE and SoproLIFE Cameras.....	28
SOPRO 717First and SOPRO 617 Cameras.....	29
SOPRO Docking Stations and Accessories.....	30
SOPRO OEM Docking Stations (linking cables and DOCK U-USB).....	31
ACTEON One-Year Warranty.....	32

PLEASE NOTE

ACTEON North America is pleased to provide you with this publication outlining the company's products and associated pricing; however, possession, control, publication or other method of acquisition of this information does not constitute the "holder" as an authorized distributor. Additionally, any reproduction, duplication, facsimile or copy of this information without the expressed written and/or verbal approval from a corporate officer of ACTEON North America is strictly prohibited. Due to product updates, changes and/or other modifications, all pricing outlined in this publication may be subject to change without prior notification. In order to verify pricing, please contact ACTEON North America during normal business hours at 800-289-6367.



PIEZOTOME CUBE

REFERENCE	PRODUCT	DESCRIPTION
F50108	Piezotome Cube	Piezotome Cube device with power cord, one Piezotome CUBE handpiece, torque wrench, foot pedal, irrigation bracket, and 5 sterile disposable irrigation lines
F12816	CUBE LED Handpiece	Autoclavable handpiece and cord with LED light ring for use with Piezotome CUBE device
F50109	Foot pedal	Foot pedal for use with Piezotome CUBE device
F50150	Dynamometric wrench	Torque wrench for safe application of tips on the Piezotome CUBE LED handpiece
F50129	Bracket	Pole for holding irrigation bag on the Piezotome CUBE
F50133	Handpiece support	Autoclavable handpiece holder for Piezotome CUBE handpiece
F50116	Clips x10	Autoclavable clips that hold the irrigation line against Piezotome CUBE handpiece cord



IMPLANTCENTER™ 2

REFERENCE	PRODUCT	DESCRIPTION
F27168- CA	Implant Center 2 LED	Delivered with i-Surge LED motor, i-Surge LED motor cord, power cord, a STARTER KIT with LED progressive Foot Pedal, 2 brackets, 2 handpiece holders and 2 sterilizable Surgical Contra Angle irrigation lines, 30 perforators, an Extraction II pack and Bein Air LED (20:1) Surgical Contra Angle (20:1) - (Extraction II Pack includes Piezotome LED Handpiece and Cord, Extraction2 Kit, Torque Wrench & Sterilization Box)
F27168	Implant Center 2 LED	Delivered with i-Surge LED motor, i-Surge LED motor cord, power cord, a STARTER KIT with LED progressive Foot Pedal, 2 brackets, 2 handpiece holders and 2 sterilizable irrigation lines, 30 perforators, and Extraction II pack *Contra Angles not supplied
F87521	Specialist Upgrade Tip Kit	Bone Surgery 2 Tip Kit, Sinus Lift 2 Tip Kit, IntraLift 2 Tip Kit and Clinical Book
F27210	i-Surge LED Motor	I-SURGE LED MOTOR for use with the Implant Center 2 LED only *Contra Angles not supplied
F27211	i-Surge LED Motor Cord	Cord for use with the i-Surge LED Motor on Implant Center 2 LED only *Contra Angles not supplied
F27170	Progressive Foot Pedal	Progressive Foot Pedal for use with Implant Center 2
E27197	Metal Carrying Case	Carrying Case for Implant Center 2, Piezotome 2 or Piezotome Solo LED

SINUS LIFT 2 KIT



AUTOCLAVABLE STAND



PIEZOTOME 2 TIP KITS



REFERENCE	PRODUCT	DESCRIPTION
F87321	NEWTRON - Perio Kit	Tips - H1, H2L, H2R, H3, H4L, H4R, Autoclavable Stand, Universal Wrench (These tips will only work with the NEWTRON Handpiece)
F87325	NEWTRON - Micro Retro Kit	Tips - P14D, P15LD, P15RD, Autoclavable Stand, Universal Wrench (These tips will only work with the NEWTRON Handpiece)
F87326	NEWTRON Retro Kit	Tips - S12-70D, AS6D, ASLD, ASRD, Autoclavable Stand, Universal Wrench (These tips will only work with the NEWTRON Handpiece)
F87509	Bone Surgery 2 Kit	Tips - BS1-2, BS2L-2, BS2R-2L, BS4-2, BS5-2, BS6-2, Autoclavable Stand, Universal Wrench (These tips will only work with the PIEZOTOME LED HP)
F87519	Sinus Lift 2 Kit	Tips - SL1-2, SL2-2, SL3-2, SL4-2, SL5-2, Autoclavable Stand, Universal Wrench (These tips will only work with the PIEZOTOME LED HP)
F87528	Essential 2 Kit	Tips - BS1S, BS4-2, LC2-2, SL1-2, SL2-2, SL3-2, Autoclavable Stand, Universal Wrench (These tips will only work with the PIEZOTOME LED HP)
F87536	Intralift 2 Kit	Tips - TKW1-2, TKW2-2, TKW3-2, TKW4-2, TKW5-2, Autoclavable Stand, Universal Wrench (These tips will only work with the PIEZOTOME LED HP)
F87546	Extraction 2 Kit	Tips - LC1-2, LC1 90-2, LC2-2, LC2L-2, LC2R-2, Ninja-2, Autoclavable Stand, Universal Wrench (These tips will only work with the PIEZOTOME LED HP)
F87554	Crown Extension 2 Kit	Tips - BS6-2, CE1-2, CE2-2, CE3-2, Autoclavable Stand, Universal Wrench (These tips will only work with the PIEZOTOME LED HP)
F87567	Crest Splitting 2 Tip Kit	Tips - CS1-2, CS2-2, CS3-2, CS4-2, CS5-2, CS6-2, Autoclavable Stand and Universal Wrench (These tips will only work with the PIEZOTOME LED HP)
N88449	Autoclavable Stand For Tips	Silver Metal Autoclavable Stand - Holds 6 NEWTRON or Piezotome tips (Tips not included)

SINUS LIFT KIT



AUTOCLAVABLE STAND



PIEZOTOME TIP KITS



REFERENCE	PRODUCT	DESCRIPTION
F87309	Bone Surgery Kit	Tips - BS1, BS2L, BS2R, BS4, BS5, BS6, Autoclavable Stand, Universal Wrench (These tips will only work with the Piezotome handpiece)
F87319	Sinus Lift Kit	Tips - SL1, SL2, SL3, SL4, SL5, Autoclavable Stand, Universal Wrench (These tips will only work with the Piezotome handpiece)
F87336	SL Intralift Kit	Tips - TKW1, TKW2, TKW3, TKW4, TKW5, Autoclavable Stand, Universal Wrench (These tips will only work with the Piezotome handpiece)
F87346	Extraction Kit	Tips - LC1, LC1 90, LC2, LC2L, LC2R, Ninja, Autoclavable Stand, Universal 1,691 Wrench (These tips will only work with the Piezotome handpiece)
F87354	Crown Extension Kit	Tips - BS6, CE1, CE2, CE3, Autoclavable Stand, Universal Wrench (These tips will only work with the Piezotome handpiece)
F87321	NEWTRON - Perio Kit	Tips - H1, H2L, H2R, H3, H4L, H4R, Autoclavable Stand, Universal Wrench (These tips will only work with the NEWTRON Handpiece)
F87325	NEWTRON - Micro Retro Kit	Tips - P14D, P15LD, P15RD, Autoclavable Stand, Universal Wrench (These tips will only work with the NEWTRON Handpiece)
F87326	NEWTRON Retro Kit	Tips - S12-70D, AS6D, ASLD, ASRD, Autoclavable Stand, Universal Wrench (These tips will only work with the NEWTRON Handpiece)
N88449	Autoclavable Stand For Tips	Silver Metal Autoclavable Stand - Holds 6 NEWTRON or Piezotome tips (Tips are not included)



PIEZOTOME 2 TIPS



PIEZOTOME TIPS

REFERENCE	DESCRIPTION
BONE SURGERY II TIPS - These tips will only work with the Piezotome II Handpiece	
F87525	BS1-2 Slim Bone Surgery Tip (9mm)
F87527	BS1-2 Long Bone Surgery Tip (15mm)
F87557	BS1-RD Rounded Saw Tip Tip (40mm)
F87502	BS2L-2 Bone Surgery Tip
F87503	BS2R-2 Bone Surgery Tip
F87504	BS4-2 Bone Surgery Tip
F87505	BS5-2 Bone Surgery Tip
F87506	BS6-2 Bone Surgery Tip
EXTRACTION II TIPS - These tips will only work with the Piezotome II Handpiece	
F87507	LC1 - 2 Extraction Tip
F87541	LC1 90 - 2 Extraction Tip
F87542	LC2 - 2 Extraction Tip
F87543	LC2L - 2 Extraction Tip
F87544	LC2R - 2 Extraction Tip
F87545	Ninja - 2 Extraction Tip
INTRALIFT II TIPS - These tips will only work with the Piezotome II Handpiece	
F87531	TKW1 - 2 Intralift Tip
F87532	TKW2 - 2 Intralift Tip
F87533	TKW3 - 2 Intralift Tip
F87534	TKW4 - 2 Intralift Tip
F87535	TKW5 - 2 Intralift Tip
F87538	TKW6 - 2 Intralift Tip
SINUS LIFT II TIPS - These tips will only work with the Piezotome II Handpiece	
F87511	SL1 - 2 Sinus Lift Tip
F87512	SL2 - 2 Sinus Lift Tip
F87513	SL3 - 2 Sinus Lift Tip
F87514	SL4 - 2 Sinus Lift Tip
	SL5 - 2 Sinus Lift Tip
CROWN EXTENSION II TIPS - These tips will only work with the Piezotome II Handpiece	
F87506	BS6-2 Bone Surgery Tip
F87551	CE1 - 2 Crown Extension Tip
F87552	CE2 - 2 Crown Extension Tip
F87553	CE3 - 2 Crown Extension Tip
CREST SPLITTING II TIPS - These tips will only work with the Piezotome II Handpiece	
F87561	CS1 - 2: Crest Splitting Tip
F87562	CS2 - 2: Crest Splitting Tip
F87563	CS3 - 2: Crest Splitting Tip
F87564	CS4 - 2: Crest Splitting Tip
F87565	CS5 - 2: Crest Splitting Tip
F87566	CS6 - 2: Crest Splitting Tip

REFERENCE	DESCRIPTION
BONE SURGERY TIPS - These tips will only work with the Piezotome Handpiece	
F87301	BS1 Bone Surgery
F87302	BS2L Bone Surgery
F87303	BS2R Bone Surgery
F87304	BS4 Bone Surgery
F87305	BS5 Bone Surgery
F87306	BS6 Bone Surgery
EXTRACTION TIPS - These tips will only work with the Piezotome Handpiece	
F87307	LC1 Extraction Tip
F87341	LC1 90 Extraction Tip
F87342	LC2 Extraction Tip
F87343	LC2L Extraction Tip
F87344	LC2R Extraction Tip
F87345	Ninja Extraction Tip
INTRALIFT TIPS - These tips will only work with the Piezotome Handpiece	
F87331	TKW1 Intralift Tip
F87332	TKW2 Intralift Tip
F87333	TKW3 Intralift Tip
F87334	TKW4 Intralift Tip
F87335	TKW5 Intralift Tip
SINUS LIFT TIPS - These tips will only work with the Piezotome Handpiece	
F87311	SL1 Sinus Lift
F87312	SL2 Sinus Lift
F87313	SL3 Sinus Lift
F87314	SL4 Sinus Lift
F87315	SL5 Sinus Lift
CROWN EXTENSION TIPS - These tips will only work with the Piezotome Handpiece	
F87306	BS6 Bone Surgery Tip
F87351	CE1 Crown Extension Tip
F87352	CE2 Crown Extension Tip
F87353	CE3 Crown Extension Tip

PIEZOTOME 2 & IMPLANT CENTER 2 Handpieces, Cords and Accessories



PIEZOTOME LED
Handpiece and Cord




I-Surge
LED Motor



NEWTRON
LED Handpiece

REFERENCE	PRODUCT	DESCRIPTION
F12815	Piezotome LED HP & Cord	Piezotome 2, Piezotome Solo & Implant Center 2 LED Handpiece and Cord: Non-disconnectable assembly (For use with PIEZOTOME II Tips)
F12609	NEWTRON LED Handpiece	Autoclavable NEWTRON LED Handpiece (For use with NEWTRON LED Handpiece Cord - Not supplied)
F12605	Light Ring	Light ring for the NEWTRON LED Handpiece
E12674	LED Light Guide	Light guide for the NEWTRON LED Handpiece
E12681	Grey Handpiece Nose	Grey Handpiece nose for the NEWTRON LED Handpiece 30 18
F12830	NEWTRON LED HP Cord for Piezotome II	Piezotome II NEWTRON LED Handpiece Cord (For use with NEWTRON LED Handpiece - Not supplied)
F12821	Light Ring	Light ring for the Piezotome LED Handpiece
E12812	LED Light Guide	Light guide for the Piezotome LED Handpiece
F27210	I-Surge LED Motor I	I-Surge LED Micro-Motor (Motor Only) for use with Implant Center II LED (Contra-Angles are not supplied)
F27211	I-Surge LED Cord	I-Surge LED Micro-Cord (Cord Only) for use with Implant Center II LED (Contra-Angles are not supplied)
E27113	Bracket	Irrigation Bracket for use with Implant Center and I-Surge
F02463	Motor Holder	Motor / Handpiece holder
F57350	Easytorque Wrench	Easytorque Wrench For Piezotome Tips only
E13125	Power Cord	110 Volt Power Cord for Implant Center 1& 2, Piezotome 1, 2 & Solo LED
E57353	Irrigation Clip	Irrigation Clip for the Piezotome (Single)
F57404	Progressive Foot Pedal	Progressive Foot Pedal for use with Piezotome 2

REFERENCE	PRODUCT	DESCRIPTION
F57374	Irrigation Kit (Autoclavable)	Kit includes 1 autoclavable irrigation line and 10 disposable perforators. (For use with all current SATELEC Surgical Devices)
F57370	Irrigation Kit (Disposable)	Kit includes 10 sterile disposable lines and cartridge. (For use with all current SATELEC Surgical Devices)
F57372	Perforators (Disposable)	Disposable perforators for use with the F57374 Irrigation lineSterilization Box,
F87523	Piezotome 2 Pack	Piezotome LED and cord and Easytorque Dynamometric Wrench. (This handpiece will only work with Piezotome II tips)
F87524	NEWTRON 2 Pack	Sterilization Box, NEWTRON LED Handpiece, NEWTRON LED Cord, Easytorque Wrench
F87520	NEWTRON Perio 2 Pack	Sterilization Box, NEWTRON LED Handpiece, NEWTRON LED Cord, Autoclavable Stand, Tips H1, H2L, H2R, H3, H4L, H4R and Easytorque Wrench
F87540	Extraction II Pack	Sterilization Box, Piezotome 2 Handpiece, Piezotome 2 Cord, Extraction 2 Kit and Easytorque Wrench (These tips will only work with the Piezotome LED Handpiece)
F27156	Sterilization Box	Piezotome II Sterilization box: can hold one handpiece with cord, 2 Tip Kits and 1 Easytorque Wrench (items are not included)



PIEZOTOME II Pack



NEWTRON Perio Pack



PIEZOTOME I & IMPLANT CENTER I Handpieces, Cords and Accessories

PIEZOTOME Handpiece




NEWTRON Handpiece




I-Surge Motor

REFERENCE	PRODUCT	DESCRIPTION
F57321	Piezotome HP + Cord	Non-disconnectable assembly (For use with Piezotome Tips)
F12281	NEWTRON Handpiece	Autoclavable Handpiece: GII 2-Tone Grey (Cord not supplied)
F57305	NEWTRON Cord	Sterilizable and disconnectable cord (NEWTRON Handpiece not supplied)
F27123	I-Surge Motor	I-Surge Micro-Motor (Motor Only) for use with Ic1 and Ic2 Non-LED (Contra-Angles are not supplied)
F27124	I-Surge Cord	I-Surge Micro-Cord (Cord Only) Contra-Angles are not supplied
E27113	Bracket	Irrigation Bracket for use with Implant Center and I-Surge
E57325	Bracket	Irrigation Bracket for the Piezotome
E57348	Hook	Irrigation Bottle Bracket Hook
F02463	Motor Holder	Motor / Handpiece holder
F57350	Easytorque Wrench	Easytorque Wrench (For Piezotome Tips only)
E13125	Power Cord	110 Volt Power Cord for Implant Center, I-Surge and Piezotome
E57353	Irrigation Clip	Irrigation Clip for the Piezotome (Single)

REFERENCE	PRODUCT	DESCRIPTION
F27120	Implant Center Pedal	Pedal for use with the Implant Center
E57341	Piezotome Pedal	Foot Pedal for use with the Piezotome
F57374	Irrigation Kit (Autoclavable)	Kit includes 1 autoclavable irrigation line and 10 disposable perforators (For use with all current SATELEC Surgical Devices)
F57370	Irrigation Kit (Disposable)	Kit includes 10 sterile disposable lines and cartridge (For use with all current SATELEC Surgical Devices)
F57372	Perforators (Disposable)	10 Disposable perforators for use with the F57374 Irrigation line
F87323	Piezotome Pack	Sterilization Box, Piezotome Handpiece, Piezotome Cord and Easytorque Wrench
F87320	NEWTRON Perio Pack	Sterilization Box, NEWTRON Handpiece, NEWTRON Cord, Auto Stand, Tips H1, H2L, H2R, H3, H4L, H4R and Easytorque Wrench
F87324	NEWTRON Pack	Sterilization Box, NEWTRON Handpiece, NEWTRON Cord and Easytorque Wrench
F27156	Sterilization Box	Piezotome II Sterilization box: can hold one handpiece with cord, 2 Tip Kits and 1 Easytorque Wrench (items are not included)



PIEZOTOME II Pack



NEWTRON Perio Pack

NEWTRON



NEWTRON

B.LED



NEWTRON

B.LED XS



New Reference #	Product	Description
F60216	NEWTRON	NEWTRON unit delivered with NEWTRON Handpiece, US Power Supply
F60211	NEWTRON with Tip Kit	NEWTRON unit delivered with NEWTRON Handpiece, intro kit (#1, 10Z, H3 tips, 2 Blue Torque Wrenches, 1 Green Torque Wrench, 1 Autoclavable Plastic Box, US Power Supply
F61116	NEWTRON B.LED	NEWTRON B.LED unit delivered with NEWTRON Slim B.LED Handpiece with White LED Ring, US Power Supply
F61117	NEWTRON B.LED with Tip Kit	NEWTRON B.LED delivered with NEWTRON Slim B.LED Handpiece with White LED Ring, Intro Kit (#1, 10Z, H3 tips, 1 Blue Torque Wrenches, 1 Green Torque Wrench, Autoclavable plastic box , US Power Supply
F62118	NEWTRON B.LED XS	NEWTRON B.LED XS unit delivered with NEWTRON Slim B.LED Handpiece with White LED Ring, US Power Supply
F62117	NEWTRON B.LED XS with Tip Kit	NEWTRON B.LED XS delivered with NEWTRON Slim B.LED Handpiece with White LED Ring, Intro kit (#1, 10Z, H3 tips, 2 Blue Torque Wrenches, 1 Green Torque Wrench, Autoclavable Plastic Box , US Power Supply

Pure NEWTRON Accessories



Reference	Product	Description
F12901	NEWTRON Slim Handpiece	NEWTRON Slim Handpiece (For use with NEWTRON Booster sold prior to 1/1/2016)
F12905	NEWTRON Slim B LED HP w White LED Ring	NEWTRON Slim B.LED Handpiece with WHITE LED ring (For use with NEWTRON B.LED, XS B.LED devices only)
F12900	NEWTRON Slim B LED HP w Blue LED Ring	NEWTRON Slim B.LED Handpiece with blue LED ring; Compatible with NEWTRON B.LED and NEWTRON XS B.LED devices for use with FLAG fluorescent liquid assistance guide (For use with NEWTRON B.LED, XS B.LED devices only)
F12915	LED Ring - White	WHITE LED Ring for NEWTRON Slim B.LED Handpiece
F62201	LED Ring - Blue	Blue LED Ring for NEWTRON Slim B.LED handpiece (to illuminate FLAG)
E12904	Nose Cone	Dark Grey nose cone for NEWTRON Slim & NEWTRON Slim B.LED
E12903	Light Guide	LED Light Guide for NEWTRON Slim B.LED Handpiece
E15002	Medium "O" Ring	Medium "O" Ring for NEWTRON Slim & NEWTRON Slim B.LED Handpieces
F62004	Tank 300ml	Tank 300ml for NEWTRON B.LED XS
F62005	Tank 500ml	Tank 500ml for NEWTRON B.LED XS
F62209	FLAG 1.5ML Vials	FLAG 5, 1.5 ML vials. Fluorescent liquid assistance guide adheres to plaque and is illuminated by the blue light of the NEWTRON Slim B.LED handpiece
F62207	FLAG 100ML Bottle	FLAG 100 ML bottle. Fluorescent liquid assistance guide adheres to plaque and is illuminated by the blue light of the NEWTRON Slim B.LED handpiece
F62206	FLAG Package	FLAG Package 100 ML bottle with Blue LED ring. Fluorescent liquid assistance guide adheres to plaque and is illuminated by the blue light of the NEWTRON Slim B.LED.



OEM PIEZO ULTRASONIC UNITS

Reference	Product	Description
F55292	Piezo Ultrasonic Module	SP NEWTRON universal w/ connector, AHP, HP Tubing, Potentiometer. Note: Does not include Air-Electric switch or power supply. (Items are available separately below)
F55292S	Piezo Ultrasonic Module with Hygiene Tips	Reference F55292 with Autoclavable Box, (2) Blue Torque Wrenches, (1) Green Torque Wrench Tips: # 1, 10Z, H3. Note: Does not include Air-Electric switch or power supply.
F55297	Module	SP NEWTRON Module with connector (10 & 4 pin)
F55287	Module	SP NEWTRON Module with connector (6 & 2 pin)
F12281	NEWTRON Handpiece	Handpiece GII two tone grey autoclavable for SP NEWTRON
F00195	Handpiece Tubing	Handpiece Tubing- Light grey silicone
F00197	Handpiece Tubing	Handpiece Tubing- Light grey silicone with connector (2 pin)
F00184	Handpiece Tubing	Handpiece Tubing Surf 6 (61")
F00188	Handpiece Tubing	Handpiece Tubing Surf 6 (81")
F01153	Handpiece Tubing	Handpiece Tubing Surf 6 GREY3 (61")
F01154	Handpiece Tubing	Handpiece Tubing Surf 6 GREY3 (81")

OEM PIEZO ULTRASONIC UNIT with LED LIGHTING!

Reference	Product	Description
F55521	Piezo Ultrasonic LED Module	SP NEWTRON LED (Universal), w/ LED handpiece, HP tubing and Potentiometer. Does not include Air-Electric switch or power supply. (Items are available separately below)
F55521S	Piezo Ultrasonic LED Module with Hygiene Tips	Reference F55521 with Autoclavable Box, (2) Blue Torque Wrenches, (1) Green Torque Wrench Tips: # 1, 10Z, H3. Note: Does not include Air-Electric switch or power supply. (Items available separately below)
X55508	LED Module	SP NEWTRON LED Module with connector
F12609	NEWTRON LED	Autoclavable LED Handpiece
F12605	Light Ring	Light ring for the LED Handpiece
E12674	LED Light Guide	Light guide for the NEWTRON LED Handpiece
F00631	LED HP Tubing	Handpiece Tubing- Light grey silicone for use with SP NEWTRON LED
F00616	LED HP Tubing	Handpiece Tubing LED Surf6 (61")
F00629	LED HP Tubing	Handpiece Tubing LED Surf6 (81")
F00632	LED HP Tubing	Handpiece Tubing LED Grey3 (61")
F00633	LED HP Tubing	Handpiece Tubing LED Grey3 (81")

OEM PIEZO ULTRASONIC UNIT ACCESSORIES

Reference	Product	Description
DTHPU101-36C13M	Transformer 40 V DC, 100W	Power supply for SP NEWTRON + Mini LED + Built-in Sopro Dock
E00596-OEM	Air-Electric Switch	Air-Electric Switch- For automatic operation
F55213K	Potentiometer Complete	Potentiometer, wire harness, knob and sticker
F55213	Potentiometer w/harness	Potentiometer with wire harness
F00524	Potentiometer	Power control for SP NEWTRON & SP NEWTRON LED
M00101	Air Bleed Pack	Air Bleed Pack for SP NEWTRON
M00100	Water Adapter	Water Adapter for SP NEWTRON
M00103	Water Solenoid Kit	Water Solenoid Kit for installation w/o water block
F00526	Knob and Sticker	Knob and Color-Coded Sticker for the SP NEWTRON Potentiometer

Reference	Description
SCALING TIPS	
F00245	#1S Longer & thinner version of the #1 Tip
F00246	#1 Universal Scaling Tip
F00247	#2 (Heavy Calculus Tip)
F00248	#3 (Stain Removal Tip)
F00253	10P (Shallow Pockets Tip)
F00359	10X (Interproximal Scaling Tip)
F00254	10Z (Subgingival Scaling Tip)
PERIO FINE TIPS	
F02170	PFU (Perio Fine Universal)
F02171	PFL (Perio Fine Left)
F02172	PFR (Perio Fine Right)
CONSERVATION and RESTORATION	
F02250	PM1 (1/4 Round-Tip Diamond-Coated Tip)
F02251	PM2 (1/4 Round-Tip Diamond-Coated Tip)
F02252	PM3 (1/4 Round-Tipped Smooth Tip)
F02253	PM4 (1/4 Conical-Tipped Diamond Coated Tip)
CONSERVATION and RESTORATION	
F02254	PMS1 (1/4 Round-Tip Diamond Coated Tip)
F02255	PMS2 (1/4 Round-Tip Diamond Coated Tip)
F02256	PMS3 (1/4 Round-Tipped Smooth Tip)
CONSERVATION and RESTORATION	
F02021	PMV1 (Diamond Coated Tip)
F02022	PMV2 (Half Ball Diamond Tip)
F02023	PMV3 (Half Ball Diamond Tip)
F02024	PMV4 (Half Ball Polishing Tip)
F02025	PMV5 (Half Ball Polishing Tip)
F02026	PMV6 (Non-Diamond Polishing Tip)
CAVITY PREPARATION	
F02040	EX1 (Ball Diamond Tip)
F02041	EX2 (Half Ball Diamond Tip)
F02042	EX3 (Half Ball Diamond Tip)
F02044	EX-L (Left Half Ball Diamond Tip)
F02043	EX-R (Right Half Ball Diamond Tip)
PERIO TIPS & ROOT PLANING	
F00366	H1 (Universal Diamond Perio Tip)
F00367	H2L (Left Diamond Perio Tip)
F00368	H2R (Right Diamond Perio Tip)
F00369	H3 (Universal Curette Tip)
F00114	H4L (Left Curette Tip)
F00115	H4R (Right Curette Tip)
PERIO-PRECISION TIPS	
F00090	P2L (Perio Precision Left Tip)
F00091	P2R (Perio Precision Right Tip)
BIOFILM DISRUPTION TIPS (BDR)	
F01001	TK1-1S (Short BDR Tip)
F01004	TK1-1L (Long BDR Tip)
F02162	TK2-1L (Left BDR Tip)
F02161	TK2-1R (Right BDR Tip)
PERIO-SOFT TIPS (Carbon Composite Tips for Implants)	
F00702	PH1 (4 pack) (Universal Carbon Tip)
F00705	PH2L (4 pack) (Left Carbon Tip)
F00706	PH2R (4 pack) (Right Carbon Tip)
F00702	PH1 (4 pack) (Universal Carbon Tip)

Reference	Description
IMPLANT PROTECT TIPS (Titanium Tips for Implants)	
F02121	IP1 Universal Implant Cleaning Tip
F02122	IP2L Left - for use on medium implant threads
F02123	IP2R Right - for use on medium implant threads
F02124	IP3L Left - for use on narrow implant threads
F02125	IP3R Right - for use on narrow implant threads
ENDODONTIC K-Files (4 Files per package)	
F43710	K10/21mm(4 pk)
F43712	K10/25mm(4 pk)
F43715	K15/21mm(4 pk)
F43717	K15/25mm(4 pk)
F43725	K25/21mm(4 pk)
F43727	K25/25mm(4 pk)
F43730	K30/21mm(4 pk)
F43732	K30/25mm(4 pk)
IRRISAFE TIPS (Irrigating Tips)	
F43807	K20/21mm (4pk) Irrigating Files
F43808	K20/25mm (4pk) Irrigating Files
F43805	K25/21mm (4pk) Irrigating Files
F43806	K25/25mm (4pk) Irrigating Files
ENDODONTIC TIPS	
F88009	SO4 Gutta Percha Condensing Tip
F88013	ET20D Diamond Coated version of the ET 20 Tip
F88012	ET40 40mm Canal Probing Tip
F88014	ET40D Diamond Coated version of the ET40 Tip
ENDO SUCCESS RETREATMENT TIPS	
F88017	ET18D Heavy Course Diamond Tip
F88011	ET20 20mm Canal Probing Tip
F88018	ET25 Titanium Niobium Alloy Tip
F88021	ET25S Short Titanium Niobium alloy Tip
F88022	ET25L Long Titanium Niobium alloy Tip
F88020	ETBD Ball Diamond Tip
ENDO SUCCESS CAP TIPS	
F88181	CAP1 Micro-blade Mini Tip - Non Active End
F88182	CAP2 Micro-blade Mini Tip - Active End
F88183	CAP3 Micro-blade Mini Tip - Active Extremity
RETRO & MICRO RETRO TIPS (Diamond Coated)	
F00118	S12/70D Back Action Tip
F00106	P14D Universal Micro Retro Tip
F00107	P15LD Left Micro Retro Tip
F00108	P15RD Right Micro Retro Tip
APICAL SURGERY TIPS (Diamond Coated)	
F00065	AS3D Tip - 3mm working length
F00079	AS6D Tip - 6mm working length
F00067	AS9D Tip - 9mm working length
F00080	ASLD Tip - Left - 3mm working length
F00081	ASRD Tip - Right - 3mm working length
LOOSENING TIPS	
F00113	C20 (Inlay / Onlay Tip)
F00249	5AE (Crown Removal Tip)
F88019	ETPR (Loosening Crowns and Posts)



**NEWTRON Intro
Kit F00932**



**NEWTRON TIP KITS
Implant Protect
Kit F02120**



**Perfect Margin
Veneers Kit
F00739**

Reference	Product	Description
F00932	NEWTRON Intro Kit	Tips - #1, 10Z, H3 with 2 Blue Torque Wrenches, 1 Green Torque Wrench and Autoclavable Box
F00939	Perio Precision Kit	Tips - #TK1-1S, P2L, PR with 3 Green Torque Wrench and Autoclavable Box
F00900	Scaling Kit	Tips - #1, #2, 10X, 10Z with 4 Blue Wrenches and Autoclavable Box
F00902	Perio Kit	Tips - H3, H4L, H4R with 3 Green Wrenches and Autoclavable Box
F00904	Bio-Film Disruption Kit BDR	Tips - TK1- 1S /TK1- 1L / TK2- 1R / TK2- 1L, 4 Green Wrenches and Autoclavable Box
F00906	Perio Soft Kit	Tips - 4 pk PH1, 4 pk PH2L, 4 pk PH2R, 3 Black Wrenches and Autoclavable Box
F02120	Implant Protect Kit	Tips - IP1, IP2L, IP2R, IP3L, IP3R, 1 Universal Wrench & Autoclavable Stand
F00738	Perfect Margin Rounded Kit	Tips - PM1, PM2, PM3, PM4, 1 Universal Wrench & Autoclavable Stand
F00736	Perfect Margin Shoulder Kit	Tips - PMS1, PMS2, PMS3, PM4, 1 Universal Wrench & Autoclavable Stand
F02020	Perfect Margin Veneers Kit	Tips - PMV1, PMV2, PMV3, PMV4, PMV5, PMV6, 1 Universal Wrench & Autoclavable Stand
F00739	Excavus Kit	Tips - EX1, EX2, EX3, EXR, EXL, 1 Universal Wrench & Autoclavable Stand



**Endo Success
CAP**



**Endo Success
Apical Surgery
F00069**



**Endo Success
Retreatment Kit
F00737**

Reference	Product	Description
F88180	ENDO Success CAP Kit	Tips - CAP1, CAP2, CAP3, 1 Universal Wrench and Autoclavable Stand
F00069	ENDO Success Apical Surgery Kit	Tips - AS3D, AS6D, AS9D, ASLD, ASRD, 1 Universal Wrench and Autoclavable Stand
F00737	ENDO Success Retreatment Kit	Tips - ET18D, ETBD, ET20, ET25S, ET25, ETPR, Univ. Wrench & Auto. Stand
F00910	Retro Kit	Tips - S12/70D, AS6D, ASLD, ASRD 4 Yellow Wrenches and Autoclavable Box
F00912	Micro Retro Kit	Tips - P14D, P15LD, P15RD, 3 Yellow Wrenches and Autoclavable Box
F00914	Loosening Kit	Tips - #5AE with 1 Orange Wrench and Autoclavable Box

ACCESSORIES FOR THE P5, P5LED, P5XS, P5XS LED and P MAX NEWTRON XS UNITS



Reference	Product	Description
F12281	NEWTRON Handpiece	Autoclavable Handpiece: GII 2-Tone (Grey) NEWTRON
F58705	Sterijoint	Sterijoint for NEWTRON Handpieces
F12609	NEWTRON LED Handpiece	Autoclavable LED Handpiece for P Max XS, P5 LED, P5 XS LED, SP NEWTRON LED and Piezotome 2 / Implant Center 2
E12681	Grey Handpiece Nose	Grey Handpiece nose for the NEWTRON LED Handpiece
F12605	Light Ring	Light ring for the NEWTRON LED Handpiece
E12674	LED Light Guide	Light guide for the NEWTRON LED Handpiece
F12640	LED Sterijoint	Sterijoint for the NEWTRON LED Handpieces
F00412	LUX Blub	Bulb for the LUX Handpiece
F12304	"O" Ring with Tool	Handpiece 3 gasket "O" rings and tool
E15019	Large "O" ring	Large "O" Ring for Handpiece
E15052	Small "O" ring	Small "O" Ring for Handpiece
F58136	Sterikit P5XS	Includes metal base, autoclavable bracket, Sterijoint and 3 universal single-use irrigation lines
F58225	Sterikit P5XS LED	Includes metal base, autoclavable bracket, LED Sterijoint and 3 universal single-use irrigation lines
E15121	"O" Ring Kit	"O" Ring Kit for P Max XS and P5 NEWTRON XS Irrigation Tanks



Reference	Product	Description
F12351	Autoclavable Box (Round)	Clear Cover GII (Empty) Autoclavable Box - Holds four tips total - Requires Satelec Torque wrenches to hold the tips (Tips not included)
F81320	Wrench (round)	Green Torque Color-Coded Wrench
F81321	Wrench (round)	Yellow Torque Color-Coded Wrench
F81322	Wrench (round)	Blue Torque Color-Coded Wrench
F81323	Wrench (round)	Orange Torque Color-Coded Wrench
F12350	Autoclavable Box	Clear Cover GII (Empty) Autoclavable Box - Holds four tips total - Requires Satelec plastic wrenches to hold the tips (Tips not included)
F81300	Wrench	Green Plastic Color coded Wrench
F81301	Wrench	Yellow Plastic Color coded Wrench
F81302	Wrench	Blue Plastic Color coded Wrench
F81303	Wrench	Orange Plastic Color coded Wrench
F81305	Wrench	Black Plastic Color coded Wrench - For carbon tips only
N88449	Autoclavable Stand for NEWTRON Tips	Silver Metal Autoclavable Stand - Holds 6 NEWTRON or Piezotome tips (Tips are not included)
F00143	Autoclavable box	Metal Blue Autoclavable Box - Holds 6 tips (Tips are not included)
F00406	Wrench	Universal Wrench (metal)



**Air n Go Easy
Air Polishing**

Reference	Product	Description
F10130	AIR-N-GO Easy KAVO	AIR-N-GO Easy with KAVO adapter. ANG Easy handpiece, (1) Kavo Multiflex adapter, (2) Supra 120° nozzles, (2) Autoclavable Front Bodies, (1) Intro kit of 10 multiflavor sticks of Classic Powder and (2) Pearl Powder sticks, (1) Maintenance kit (silicone grease syringe, 3 brushes, spare o-rings, syringe, cannula and 2 probes)
F10131	AIR-N-GO Easy SIRONA	AIR-N-GO Easy with SIRONA adapter. ANG Easy handpiece, (1) Sirona adapter, (2) Supra 120° nozzles, (2) Autoclavable Front Bodies, (1) Intro kit of 10 multiflavor sticks of Classic Powder and (2) Pearl Powder sticks, (1) Maintenance kit (silicone grease syringe, 3 brushes, spare o-rings, syringe, cannula and 2 probes)
F10132	AIR-N-GO Easy W & H	AIR-N-GO Easy with W & H adapter. ANG Easy handpiece, (1) W&H adapter, (2) Supra 120° nozzles, (2) Autoclavable Front Bodies, (1) Intro kit of 10 multiflavor sticks of Classic Powder and (2) Pearl Powder sticks, (1) Maintenance kit (silicone grease syringe, 3 brushes, spare o-rings, syringe, cannula and 2 probes)
F10133	AIR-N-GO Easy BIEN AIR	AIR-N-GO Easy with BIEN AIR adapter. ANG Easy handpiece, (1) Bien Air adapter, (2) Supra 120° nozzles, (2) Autoclavable Front Bodies, (1) Intro kit of 10 multiflavor sticks of Classic Powder and (2) Pearl Powder sticks, (1) Maintenance kit (silicone grease syringe, 3 brushes, spare o-rings, syringe, cannula and 2 probes)
F10135	AIR-N-GO Easy MIDWEST	AIR-N-GO Easy with MIDWEST adapter. ANG Easy handpiece, (1) Midwest adapter, (2) Supra 120° nozzles, (2) Autoclavable Front Bodies, (1) Intro kit of 10 multiflavor sticks of Classic Powder and (2) Pearl Powder sticks, (1) Maintenance kit (silicone grease syringe, 3 brushes, spare o-rings, syringe, cannula and 2 probes)



air N GO
easy **Accessories**

Reference	Product	Description
F10120	AIR-N-GO Perio Easy Ki	AIR-N-GO Perio Easy nozzle, (1) 160g bottle Perio glycine powder
F10125	AIR-N-GO Supra	AIR-N-GO Supra nozzle
F10126	AIR-N-GO Perio Easy Nozzle	AIR-N-GO Perio Easy nozzle
F10127	AIR-N-GO Perio Nozzle	AIR-N-GO Perio nozzle
F10128	AIR-N-GO Perio Maintenance Nozzle	AIR-N-GO Perio Maintenance nozzle
F10121	AIR-N-GO Easy O-ring	AIR-N-GO Easy o-rings kit: includes silicone grease syringe, (3) brushes,
E10207	AIR-N-GO Front Body	Autoclavable Front Body for AIR-N-GO Easy Handpiece



Bottles of Powder

AIR-N-GO
Powder

Reference	Product	Description
F10250	AIR-N-GO: Propphy Powder- Peppermint (4-pk bottles)	AIR-N-GO: Sodium Bicarbonate Propphy Powder - Peppermint (4-pk - 250g Bottles)
F10251	AIR-N-GO: Propphy Powder- Lemon (4-pk bottles)	AIR-N-GO: Sodium Bicarbonate Propphy Powder - Lemon (4-pk - 250g Bottles)
F10252	AIR-N-GO: Propphy Powder- Raspberry (4-pk bottles)	AIR-N-GO: Sodium Bicarbonate Propphy Powder - Raspberry (4-pk - 250g Bottles)
F10253	AIR-N-GO: Propphy Powder- Cola (4-pk bottles)	AIR-N-GO: Sodium Bicarbonate Propphy Powder - Cola (4-pk - 250g Bottles)
F10254	AIR-N-GO: Propphy Powder- Neutral (4-pk bottles)	AIR-N-GO: Sodium Bicarbonate Propphy Powder - Neutral (4-pk - 250g Bottles)
F10255	AIR-N-GO: Propphy Powder- Variety Pack (4-pk bottles)	AIR-N-GO: Sodium Bicarbonate Propphy Powder - Variety Pack (4-pk - 250g Bottles) (1- Peppermint, 1- Lemon, 1- Raspberry & 1- Cola)
F10262	AIR-N-GO: Propphy Powder- Pearl (4-pk bottles)	AIR-N-GO: Calcium Carbonate Propphy Powder - Pearl (4-pk - 250g Bottles)
F10270	AIR-N-GO: Propphy Powder- PERIO (3-pk bottles)	AIR-N-GO: Glycine Propphy Powder - PERIO ONLY (3-pk - 160g Bottles)



Boxes of 50 x 20g Packets



AIR-N-GO
Powder

F10256	AIR-N-GO: Propphy Powder- Peppermint: Box of 50 (20g) Sticks	AIR-N-GO: Sodium Bicarbonate Propphy Powder - Peppermint: Box of 50 (20g) Sticks
F10257	AIR-N-GO: Propphy Powder- Lemon: Box of 50 (20g) Sticks	AIR-N-GO: Sodium Bicarbonate Propphy Powder - Lemon: Box of 50 (20g) Sticks
F10258	AIR-N-GO: Propphy Powder- Raspberry: Box of 50 (20g) Sticks	AIR-N-GO: Sodium Bicarbonate Propphy Powder - Raspberry Box of 50 (20g) Sticks
F10259	AIR-N-GO: Propphy Powder- Cola: Box of 50 (20g) Sticks	AIR-N-GO: Sodium Bicarbonate Propphy Powder - Cola Box of 50 (20g) Sticks
F10260	AIR-N-GO: Propphy Powder- Neutral: Box of 50 (20g) Sticks	AIR-N-GO: Sodium Bicarbonate Propphy Powder - Neutral: Box of 50 (20g) Sticks
F10263	AIR-N-GO: Propphy Powder- Pearl: Box of 50 (20g) Sticks	AIR-N-GO: Calcium Carbonate Propphy Powder - Pearl: Box of 50 (20g) Sticks



**Air n Go
Blue Tank**

airN GO (Accessories)



**Air n Go
Green Tank**

Reference	Product	Description
F10276	AIR-N-GO: Cap w/button	AIR-N-GO: Cap with button for powder tank
F10213	AIR-N-GO: Blue Tank- SUPRA	AIR-N-GO: Blue Tank- For SUPRA-Gingival Polishing
F10207	AIR-N-GO: Probe (.5mm) SUPRA	AIR-N-GO: Probe (.5mm) SUPRA
F10208	AIR-N-GO: Probe (1mm) SUPRA	AIR-N-GO: Probe (1mm) SUPRA
F10205	AIR-N-GO: KIT of (11) Rubber Joints	AIR-N-GO: KIT of (11) Rubber Joints
F10216	AIR-N-GO: Green Tank- SUB/PERIO	AIR-N-GO: Green Tank- SUB-GINGIVAL / PERIO
F10206	AIR-N-GO: PERIO Probe	AIR-N-GO: PERIO Probe



**Air n Go Adaptor
MIDWEST**

Reference	Product	Description
F10240	AIR-N-GO: Adaptor- KAVO Multiflex	AIR-N-GO: Adaptor- KAVO Multiflex
F10241	AIR-N-GO: Adaptor- SIRONA	AIR-N-GO: Adaptor- SIRONA
F10242	AIR-N-GO: Adaptor- W&H	AIR-N-GO: Adaptor- W&H
F10243	AIR-N-GO: Adaptor- BIEN AIR	AIR-N-GO: Adaptor- BIEN AIR
F10245	AIR-N-GO: Adaptor- MIDWEST	AIR-N-GO: Adaptor- MIDWEST

PROPHY PENS AND NOZZLES FOR AIR MAX AND PROPHY MAX UNITS



<i>Reference</i>	<i>Product</i>	<i>Description</i>
F10554	Prophy Pen	Prophy Pen (Bayonet) - Newest Style - Blue
F10557	Prophy Pen	Prophy Pen (Bayonet) - Newest Style - White
F10552	Titanium Nozzle	Titanium Prophy Pen Nozzle (0.6mm diameter)
F10555	Titanium Nozzle	Titanium Prophy Pen Nozzle (0.7mm diameter)

ACCESSORIES FOR AIR MAX AND PROPHY MAX UNITS



<i>Reference</i>	<i>Product</i>	<i>Description</i>
F59750	Irrigation Kit	Peristaltic Pump, 2 - 400ml containers (Autoclavable), Zip Tool, Sterijoint, 3 Pump Irrigation Lines, Hook & Irrigation Line
F59751	LUX Irrigation Kit	Peristaltic Pump, 2 - 400ml containers (Autoclavable), Zip Tool, LUX Sterijoint, 3 Pump Irrigation Lines, Hook & Irrigation Line
F59752A	Irrigation Lines	10 pack (double-lines for Max line units)
F59756A	Irrigation Lines	10 pack (lines for IV Needle with new Max line units)
F59761	Clips	Kit of 3 Clips for Irrigation Lines
E59764	Bracket	Bracket for IV Needle Line source
E59774-KIT	Irrigation Tanks 400ml	400ml Autoclavable Single Tank with cap
E59760	Zip Tool	Zip Tool for irrigation line installation
F12620	LUX Sterijoint	LUX Sterijoint GII for P Max LUX and Prophy Max LUX Handpieces
F58705	Sterijoint GII	Sterijoint GII for P Max and Prophy Max Handpieces
F10389	Water Filter Cartridge	Replacement water filter cartridge
F10386	Water Filter	Water Filter, complete
F32125	Water Line	Water line for Max line units
F32126	Air Line	Air line for Max line units
F00412	LUX Bulb	Bulb for P Max LUX Handpiece



Mini LED
BLACK

**Works with
Three "AA" Batteries**

Reference	Product	Description
F05400	Mini LED BLACK	1,250 mWcm2 Output - Mini LED Black Curing Light Handpiece, 7.5mm Light Guide, Shield and 3 x AA Batteries
F05409	Light Guide	7.5mm Light Guide for Mini LED Black



MINILED
— SUPERCHARGED

**3 Second
curing for
most composites**



MINILED
— ORTHO 2

**Designed
for
Orthodontics**



MINILED
— STANDARD

**High-power -
with no
heat build-up**

Reference	Product	Description
F05214	Mini LED SuperCharged	2,200 mWcm2 Output, (Modes: Full, Pulse & Ramp) Handpiece, Battery, Charger, Power Supply, Light Guide 7.5mm, Oval Light Shield (Cordless)
F05217	Handpiece	Mini LED SuperCharged Handpiece
F05220	Mini LED Ortho2	3,000 mWcm2 Output (Modes 4, 8, 16 & 32 seconds) Handpiece, Battery, Charger, Power Supply, Light Guide 5.5mm, Oval Light Shield (Cordless)
F02611	Mini LED	1,250 mWcm2 Output (Modes: Full, Pulse & Ramp) Handpiece, Battery, Charger, Power Supply, Light Guide 7.5mm, Oval Light Shield (Cordless)
F02544	Power Supply	Power supply with US plug for Mini LED
F02510	Battery Charger	Mini LED Battery Charger Base
F02520	Battery	Mini LED Battery Lithium Ion
F02530	Handpiece	Mini LED Handpiece for Battery operated unit

MINI LED LIGHT GUIDES

**(All of the light guides listed below work with all of the
Mini LED Lights except the Mini LED Auto Focus)**



Reference	Product	Description
F02648	Light Guide	Light Guide 7.5mm
F02652	Light Guide	Light Guide, 5.5mm
F02658	Light Guide	5.5mm Black Light Guide
F05407	Light Shield	Light Shield for Mini LED Black, Mini LED Supercharged, Mini LED Ortho2
F02555	Light Shield	Oval Light Shield for 7.5mm and 5.5mm

OEM LED CURING LIGHTS

(Easily integrated into your dental unit...)



Reference	Product	Description
F02720	Mini LED OEM	Mini LED Module, Light Grey Cord, Handpiece, Light Guide 7.5mm, Protective Light Shield (Fast, Pulse & Ramp) Note - Requires a 24 VAC Transformer - see below.
F02710	Handpiece	Mini LED Handpiece OEM
F05205	Mini LED Supercharged OEM	2,200 mWcm2 Output, (Modes: Full, Pulse & Ramp) Handpiece, Module, Light Grey Cord, Light Guide 7.5mm, Oval Light Shield) Note - Requires a 24 VAC Transformer - see below
F05218	Handpiece	Mini LED SuperCharged OEM Handpiece
F02721	Tubing	Mini LED Handpiece Tubing OEM - light grey
F02326	Tubing	Mini LED Handpiece Tubing - Serf 6 - 9 foot
F02730	Module	Mini LED Module OEM
DTHPU101-36C13M	Transformer 40 V DC, 100W	Power supply for SP NEWTRON + Mini LED + Built-in Sopro Dock

MINI LED LIGHT GUIDES

(All of the light guides listed below work with all of the Mini LED Lights except the Mini LED Auto Focus)



Reference	Product	Description
F02648	Light Guide	Light Guide 7.5mm
F02652	Light Guide	Light Guide, 5.5mm
F02658	Light Guide	5.5mm Black Light Guide
F02555	Light Shield	Oval Light Shield for 7.5mm and 5.5mm



294100



261030



260900



261040

EFFECTIVE AND ATRAUMATIC SULCULAR OPENING

Reference	Product	Description
261030		Expasyl Capsule - 20 Pack (Standard)
261001		Expasyl Capsule - 20 Pack (Strawberry)
294100		Expasyl Starter Kit - (Gun + 6 Capsules + 12 Applicator Tips)
261040		Expasyl Applicator Tips - Box of 40
261005		Expasyl Applicator Tips - Box of 100
260900		Expasyl Applicator Gun (Alone)



246200

Reference	Product	Description
246200		Hemostasyl Intro Kit - (2 Syringe + 40 Tips)
246050		Hemostasyl Syringe - 2 Pack
246100		Hemostasyl Applicator Tips - Box of 40

Ris<Control[®]
MINT

200291



Ris<Control[®]
Disposable tips for a air-water syringe



WATER SYRINGE TIPS & ADAPTORS

Reference	Product	Description
201700	Riskontrol MINT Tips	Box of 250 MINT tips (disposable)
201664	Riskontrol Perfect Tip	Box of 250 WHITE tips (disposable)
201665	Riskontrol Perfect Tip	Box of 250 BLUE tips (disposable)
201668	Riskontrol Perfect Tip	Box of 250 PINK tips (disposable)

Reference	Product	Description
200195	Syringe Adapter	Adapter ADEC II "NEW" Fits Riskontrol Perfect only
200158	Syringe Adapter	Adapter BELMONT
200152	Syringe Adapter	Adapter DCI
200164	Syringe Adapter	Adapter Dental EZ
200178	Syringe Adapter	Adapter Siemens 4000
200154	Syringe Adapter	Adapter Sirona
200171	Syringe Adapter	Adapter KAVO
200150	Syringe Adapter	Adapter Luzzani



201660

BIO SPLINTTM

SPLINTING RIBBONS

Reference	Product	Description
201661	BioSplint Ribbons- Refill	BioSplint Ribbons- Refill: 8 Splinting Ribbons



X MIND
trium

*Panoramic, Cephalometric and CBCT
capable X-ray Imaging device for
Demanding dental professionals.*

Reference	Product	Description
W0900225	XMT Pan Upgradeable R	XMT Pan Upgradeable R (AIS SW)
W0900226	XMT Pan Upgradeable L	XMT Pan Upgradeable L (AIS SW)
W0900003	XMT Pan R Ceph L	XMT Pan R Ceph L (Pan Upgradeable R W0900225 + Left Ceph UPG W0900019, AIS SW)
W0900009	XMT Pan L Ceph R	XMT Pan L Ceph R (Pan Upgradeable L W0900226 + Right Ceph UPG W0900018, AIS SW)
W0900232	XMT Pan L Ceph L	XMT Pan L Ceph L (Pan Upgradeable L W0900226 + Left Ceph UPG W0900019, AIS SW)
W0900230	XMT Pan R Ceph R	XMT Pan R Ceph R (Pan Upgradeable R W0900225 + Right Ceph UPG W0900018, AIS SW)
W0900004	XMT Pan 3D R	XMT Pan 3D R (AIS SW, 2 3DApp Standard, 5 3DApp Basic)
W0900010	XMT Pan 3D L	XMT Pan 3D L (AIS SW, 2 3DApp Standard, 5 3DApp Basic)
W0900006	XMT Pan 3D R Ceph L	XMT Pan 3D R Ceph L (Pan 3D R W0900004 + L Ceph UPG W0900019 - AIS SW, 2 3DApp Standard, 5 3DApp Basic)
W0900012	XMT PAN 3D L Ceph R	XMT Pan 3D L Ceph R (Pan 3D L W0900010 + R Ceph UPG W0900018 - AIS SW, 2 3DApp Standard, 5 3DApp Basic)
W0900233	XMT Pan 3D L Ceph L	XMT Pan 3D L Ceph L (Pan 3D L W0900010 + L Ceph UPG W0900019 - AIS SW, 2 3DApp Standard, 5 3DApp Basic)
W0900231	XMT PAN 3D R Ceph R	XMT Pan 3D R Ceph R (Pan 3D R W0900004 + R Ceph UPG W0900018 - AIS SW, 2 3DApp Standard, 5 3DApp Basic)
W0900227	XMT UPG 3D	XMT UPG 3D (AIS SW, 2 3DApp Standard, 5 3DApp Basic)
W0900018	XMT UPG Ceph R	XMT UPG Ceph R
W0900019	XMT UPG Ceph L	XMT UPG Ceph L
W0900236	XMT SELF STANDING	XMT SELF STANDING BASE (FOR WORKING TRIUM)

X-Mind Trium Accessories

Reference	Product	Description
W0900110	XMT Ear Protection	XMT Ear C2026 Protection (280pcs)
W0900080	XMT Chin Rest, Bite Block. TMJ	XMT Chin Rest, Bite Block, TMJ
W0900142	XMT Sliding Bite Blocks	XMT Sliding Bite Blocks (10pcs)
W0900113	XMT Bite Protection	XMT Bite Protection (500 pcs)
W1000012	AIS 3D App Basic 1 Electronic License	AIS 3D App Basic 1 Electronic License
W1000013	AIS 3D App Basic 5 Electronic Licenses	AIS 3D App Basic 5 Electronic Licenses
W1000007	AIS 3D APP Standard 1 Dongle License	AIS 3D APP Standard 1 Dongle License
W1000015	AIS 3D APP Standard 5 Dongle Licenses	AIS 3D APP Standard 5 Dongle Licenses
W1000016	AIS 3D APP Premium 1 Dongle License	AIS 3D APP Premium 1 Dongle License
W1000017	AIS 3D APP Premium 5 Dongle Licenses	AIS 3D APP Premium 5 Dongle Licenses
39800751	License DICOM Print	License DICOM Print
W0900153	XMT License AIS DICOM Package (No DICOM Print)	XMT License AIS DICOM Package (No DICOM Print)
W0900154	XMT First Aid Kit	XMT First Aid Kit (Incl most common parts)



X MIND
unity

Digital X-Ray Systems (For use with Sopix Inside Sensors)

Reference	Product	Description	Top Mount
HXMU-TS	XMU TOP 40 SHORT CONE	DC Tube Head, Scissor Arm, Timer, 2 Stud Wall Mount, Short Cone, Top Mount Extender Arm - 16"	
HXMU-TM	XMU TOP 80 SHORT CONE	DC Tube Head, Scissor Arm, Timer, 2 Stud Wall Mount, Short Cone, Top Mount Extender Arm - 31"	
HXMU-TL	XMU TOP 110 SHORT CONE	DC Tube Head, Scissor Arm, Timer, 2 Stud Wall Mount, Short Cone, Top Mount Extender Arm - 43"	

Reference	Product	Description	Bottom Mount
HXMU-BS	XMU BOTT 40 SHORT CONE	DC Tube Head, Scissor Arm, Timer, 2 Stud Wall Mount, Short Cone, Bottom Mount Extender Arm - 16"	
HXMU-BM	XMU BOTT 80 SHORT CONE	DC Tube Head, Scissor Arm, Timer, 2 Stud Wall Mount, Short Cone, Bottom Mount Extender Arm - 31"	
HXMU-BL	XMU BOTT 110 SHORT CONE	DC Tube Head, Scissor Arm, Timer, 2 Stud Wall Mount, Short Cone, Bottom Mount Extender Arm - 43"	



Reference	Product	Description
S 802 6500	SOPIX INSIDE - USB SENSOR SIZE 1	SOPIX Inside (Size 1) Digital Sensor - Direct USB2 (For use with X Mind Unity Only)
S 802 6501	SOPIX INSIDE - USB SENSOR SIZE 2	SOPIX Inside (Size 2) Digital Sensor - Direct USB2 (For use with X Mind Unity Only)



W0800021

X-Mind Unity Accessories



W0800034

Reference	Product	Description
W0800018	XMU SHORT CONE	SHORT CONE: 8" Ø60 (20cm)
W0800019	XMU LONG CONE	LONG CONE: 12" Ø60 (31cm)
W0800020	XMU SQ CONE	SQUARE SHORT CONE: 8" Ø60 (20cm)
W0800021	XMU REMOTE EXPOSURE SWITCH	REMOTE EXPOSURE SWITCH
W0800034	XMU EXTERNAL LIGHT 115V	XMU EXTERNAL LIGHT - 115V



Phosphor Plate Dental X-Ray System

PSPIX²

Reference	Product	Description
S 702 0011	PSPIX - Intra Oral Scanner	Includes 4 imaging plates (2 x size 1, 2 x size 2), 200 bag & cover for imaging plates (100 x size 1, 100 x size 2), 1 imaging plate storage box, 1 stylus for touchscreen, 1 power supply 100-240 VAC, 1 Ethernet cable, SOPRO Imaging software for MAC and PC with license

PSPIX Accessories



Reference	Product	Description
990215	IMAGING PLATE SIZE 0 (X2)	STANDARD IMAGING PLATE - SIZE 0 22 x 31 mm (Set of 2 pieces)
990216	IMAGING PLATE SIZE 1 (X2)	STANDARD IMAGING PLATE- SIZE 1 24x 40 mm (Set of 2 pieces)
990217	IMAGING PLATE SIZE 2 (X2)	STANDARD IMAGING PLATE - SIZE 2 31 x 41 mm (Set of 2 pieces)
990218	IMAGING PLATE SIZE 3 (X2)	STANDARD IMAGING PLATE - SIZE 3 27 x 54 mm (Set of 2 pieces)
700546	BAGS & COVERS S0 - 300	BAG & COVER FOR IMAGING PLATE SIZE 0 (set of 300 pieces)
700547	BAGS & COVERS S1 - 300	BAG & COVER FOR IMAGING PLATE SIZE 1 (set of 300 pieces)
700548	BAGS & COVERS S2 - 300	BAG & COVER FOR IMAGING PLATE SIZE 2 (set of 300 pieces)
700549	BAGS & COVERS S3 - 300	BAG & COVER FOR IMAGING PLATE SIZE 3 (set of 300 pieces)
700550	BAGS & COVERS START UP KIT- 100 ea	BAG & COVER FOR IMAGING PLATE START UP KIT (set of 100 pieces of size 1 and 100 pieces of size 2)
990219	IP STORAGE BOX PSPIX	IMAGING PLATE STORAGE BOX (Delivered without imaging plates)
005132	Stylus	STYLUS FOR PSPIX TOUCHSCREEN
608019	Auto. insertion guide	Autoclavable insertion guide (Top)
608018	Autoclavable receptacle	Autoclavable receptacle (Bottom)
700342L	Imaging Plate Bag: Size 0	Protective Bag for Imaging Plate (250 per pack) Size 0
700343L	Imaging Plate Bag: Size 1	Protective Bag for Imaging Plate (250 per pack) Size 1
700344L	Imaging Plate Bag: Size 2	Protective Bag for Imaging Plate (250 per pack) Size 2
700345L	Imaging Plate Bag: Size 3	Protective Bag for Imaging Plate (250 per pack) Size 3
700195	PSPIX Positioning Set - Bitewing	Positioning Set for PSPIX Imaging Plates (Includes 5 holders for bitewing)
700196	PSPIX Positioning Set - Anterior & Posterior	Positioning Set for PSPIX Imaging Plates (2 for Anterior & 2 for Posterior)
990243	PSPIX Cleaning Plate	Cleaning plate for use with PSPIX



SOPIX²

Digital X-Ray Sensors

Reference	Product	Description
S 802 0006	SOPIX ² - Size 1	SOPIX ² (Size 1) Digital Sensor - Direct USB2
S 802 0005	SOPIX ² - Size 2	SOPIX ² (Size 2) Digital Sensor - Direct USB2

SOPIX Inside and SOPIX² ACCESSORIES



926039



700042



403026

Reference	Product	Description
926039	Sensor Holder	SOPIX Inside or SOPIX ² Digital Sensor Holder
700406	Sopro Tabs	Positioning Tabs - Box of 500
700042	SOPIX Positioning Set	SOPIX Inside or SOPIX ² Positioning Set (5 sets with positioning Ring for: Bitewing, Ant, Post, Endo Ant, Endo Post)
403026	Sensor Sheaths: Size 1	Sheaths for Size 1 SOPIX Inside or SOPIX ² Digital Sensor (500 per pack)
403027	Sensor Sheaths: Size 2	Sheaths for Size 2 SOPIX Inside or SOPIX ² Digital Sensor (500 per pack)

*** ORANGE Font Indicates 30% Dealer Discount**

SOPRO CARE

INTRA-ORAL CAMERA

**for detection of occlusal caries, plaque, calculus
and to expose gingival inflammation**



Reference	Product	Description
S 950 0002	SoproCARE - Camera Handpiece	Intraoral camera used for detection of occlusal caries, plaque, calculus and to expose gingival inflammation. Includes SoproCARE Handpiece, Handpiece Holder & 4 Intraoral Tips
926063	Handpiece holder	Universal Handpiece Holder for SOPRO Cameras
990135	SoproCARE: Intra Oral Tip	Intra Oral Tip for SoproCARE & SoproLIFE: (4 per pack)
403017	Sleeves for SoproCARE & others	Sheaths for SoproCARE, SoproLIFE and SOPRO 717 First Cameras (500 sheaths in a box)
S 595 1145	Dock USB 2.0	DOCK USB 2.0 This allows the connection of a Sopro Camera directly to the USB 2.0 port of any computer (Single Cable Docking Station)


Daylight Mode



Perio Mode



Cario Mode



SOPRO LIFE

Light Induced Fluorescence Evaluator

INTRA-ORAL CAMERA with CARIES DETECTION




Reference	Product	Description
S 900 0002	SoproLIFE - Camera Handpiece	Light Induced Fluorescence Evaluator: Caries Detection and Treatment Intraoral Camera. Includes SoproLIFE Handpiece, Handpiece Holder & 4 Intra Oral Tips
926063	Handpiece holder	Universal Handpiece Holder for SOPRO Cameras
990135	SoproLIFE: Intraoral Tip	Intraoral Tip for SoproCARE & SoproLIFE: (4 per pack)
990136	SoproLIFE: Light Shield	Light Shield for SoproLIFE: (2 per pack)
403017	Sleeves for SoproLIFE & others	Sheaths for SoproCARE, SoproLIFE and SOPRO 717 First Cameras (500 sheaths in a box)
S 595 1145	Dock USB 2.0	DOCK USB 2.0 This allows the connection of a Sopro Camera directly to the USB 2.0 port of any computer (Single Cable Docking Station)


Daylight Mode



Diagnosis Mode



Treatment Mode





Reference	Product	Description
S 717 0101	Sopro 717 First Camera Handpiece	SOPRO 717 First Camera Handpiece. Includes Sopro 717 First camera handpiece, 10 protective sleeves and Handpiece Holder.
926063	Handpiece holder	Universal Handpiece Holder for SOPRO Cameras
403017	Sleeves for Sopro 717 & others	Sheaths for SoproCARE, SoproLIFE and SOPRO 717 First Cameras (500 sheaths in a box)
S 595 1145	Dock USB 2.0	DOCK USB 2.0 This allows the connection of a Sopro Camera directly to the USB 2.0 port of any computer (Single Cable Docking Station)



State of the seal of the amalgam



Micro fissure with lesion



Dental Cavity Preparation



Infiltrated occlusal groove



Reference	Product	Description
S 617 0002	Sopro 617 Camera Handpiece	SOPRO 617 Camera Handpiece Includes Sopro 617 Camera Handpiece, 10 protective sleeves and Handpiece Holder.
926063	Handpiece Holder	Universal Handpiece Holder for SOPRO Cameras
403023	Sleeves for Sopro 617	SHEATHS FOR SOPRO 617 CAMERA (500 sheaths in a box)
S 595 1145	Dock USB 2.0	DOCK USB 2.0 This allows the connection of a Sopro Camera directly to the USB 2.0 port of any computer (Single Cable Docking Station)

Fixed-focus and a large depth of field



From portrait ←

→ to a single tooth.

SOPRO CAMERA DOCKING STATIONS



Dock USB 2.0

Reference	Product	Description
S 595 1145	Dock USB 2.0	DOCK USB 2.0 This allows the connection of a Sopro Camera directly to the USB 2.0 port of any computer (Single Cable Docking Station)

SOPRO CAMERA ACCESSORIES

Reference	Product	Description
SOPRO HC5	Sleeves for Sopro 595	SHEATHS FOR SOPRO 595 CAMERA (500 sheaths in a box)
009002	Foot Pedal	Camera Foot Pedal for M-USB Dock
009009	Foot Pedal	USB Foot Pedal
002003	Composite cable	VIDEO CABLE With RCA/RCA connector. Length : 5 ft
002004	S-Video cable	S-VIDEO CABLE With YC/YC connector. Length : 5 ft
700012	S-Video cable	S-VIDEO CABLE With YC/YC connector. Length : 15 ft
002081	2M USB Cable	USB CABLE 6 1/2ft for Dock MUSB reference S 717 1001
002091	5M USB Cable	USB CABLE 16 1/2ft for Dock MUSB reference S 717 1001
002148	5M Mini USB Cable	5M Mini USB Cable for U USB (Dock to PC connection)
123540	8M USB Cable	USB CABLE 26 ft for Dock MUSB reference S 717 1001
922050	Linking cable for M-USB dock	LINKING CABLE 8ft FOR THE DOCK MUSB
922051	Linking cable for M-USB dock	LINKING CABLE 15ft FOR THE DOCK MUSB



SOPRO OEM

DOCKING STATIONS without memory

SOPRO Dock U-USB

Reference	Product	Description
S 717 1507	Dock U-USB - Light Grey	BUILT-IN DOCKING STATION WITHOUT MEMORY: Features USB 2.0 Outputs for Dental Unit - LIGHT GREY colored tubing - One built-in docking station with one 21 foot linking cable with LIGHT GREY color tubing & one USB cable with mini USB connection <i>(Does not work with SOPRO 717 SoproShade Cameras)</i>
S 717 1508	Dock U - USB - Grey	BUILT-IN DOCKING STATION WITHOUT MEMORY: Features USB 2.0 Outputs for Dental Unit - GREY colored tubing - One built-in docking station with one 21 foot linking cable with GREY color tubing & one USB cable with mini USB connection <i>(Does not work with SOPRO 717 SoproShade Cameras)</i>
997051	Dock U-USB - Dock Only	BUILT-IN DOCKING STATION WITHOUT MEMORY: Features USB 2.0 Outputs for Dental Unit

LINKING CABLES FOR THE U USB Docking Stations

Reference	Product	Description
002148	5 M Mini USB Cable	Mini USB CABLE 16 1/2ft for Dock U USB <i>(This cable only works with the U USB DOCK)</i>
922138	Surf6 linking cable U USB	SURF6 LINKING CABLE FOR DOCK U USB w/swivel connector (2m/6ft) w/clip <i>(This cable only works with the U USB DOCK)</i>
922139	Surf6 linking cable U USB	SURF6 LINKING CABLE FOR DOCK U USB w/swivel connector (7m/23ft) w/clip <i>(This cable only works with the U USB DOCK)</i>
922140	Grey3 linking cable U USB	GREY3 LINKING CABLE FOR DOCK U USB w/swivel connector (2m/6ft) w/clip <i>(This cable only works with the U USB DOCK)</i>
922141	Grey3 linking cable U USB	GREY3 LINKING CABLE FOR DOCK U USB w/swivel connector (7m/23ft) w/clip <i>(This cable only works with the U USB DOCK)</i>
922143	Light Grey Linking Cable U USB	LIGHT GREY LINKING CABLE 21ft FOR BUILT-IN U USB DOCK SOPRO camera linking cable to build into the arm of the dental unit. <i>(This cable only works with the U USB DOCK)</i>
922163	Grey Linking Cable U USB	GREY LINKING CABLE 21ft FOR BUILT-IN U USB DOCK SOPRO camera linking cable to build into the arm of the dental unit. <i>(This cable only works with the U USB DOCK)</i>



SOPRO OEM

**DOCKING STATIONS with memory
(VIDEO and USB 2.0 OUTPUT)**

SOPRO
Dock MU USB

Reference	Product	Description
997002	Dock MU-USB dock only	BUILT-IN DOCKING STATION WITH MEMORY: Features Analog and USB 2.0 Outputs for Dental Unit.
S 717 1034	Dock MU-USB - White	BUILT-IN DOCKING STATION WITH MEMORY: Features Analog and USB 2.0 Outputs for Dental Unit - WHITE colored tubing - One built-in docking station with Image Memory, one 21 foot linking cable with WHITE color tubing, 1 S-Video Cable, 1 Composite Cable, 1 USB 2.0 Cable and SOPRO Imaging Disc.
S 717 1035	Dock MU-USB - Grey	BUILT-IN DOCKING STATION WITH MEMORY: Features Analog and USB 2.0 Outputs for Dental Unit - GREY colored tubing - One built-in docking station with Image Memory, one 21 foot linking cable with GREY color tubing, one S-Video Cable, one Composite Cable, one USB 2.0 Cable and SOPRO Imaging Disc.
S 717 1036	Dock MU-USB - Light Grey	BUILT-IN DOCKING STATION WITH MEMORY: Features Analog and USB 2.0 Outputs for Dental Unit - LIGHT GREY colored tubing - One built-in docking station with Image Memory, one 21 foot linking cable with LIGHT GREY color tubing, one S-Video Cable, one Composite Cable, one USB 2.0 Cable and SOPRO Imaging Disc.
S 717 1047	Dock MU-USB Without Sheathing on Camera Linking Cable	BUILT-IN DOCKING STATION WITH MEMORY: Features Analog and USB 2.0 Outputs for Dental Unit. Includes handpiece holder.

SOPRO OEM

LINKING CABLES FOR THE MU USB and MU VIDEO Docking Stations

Reference	Product	Description
922112	Surf 6 linking cable	SURF6 LINKING CABLE 15 ft FOR BUILT-IN DOCK SOPRO camera linking cable to build into the arm of the dental unit
922117	Surf 6 linking cable	SURF6 LINKING CABLE FOR DOCK MU USB w/swivel connector (2m/6ft) w/clip
922118	Surf 6 linking cable for MU USB	SURF6 LINKING CABLE FOR DOCK MU USB w/swivel connector (7m/23ft) w/clip
922126	Grey3 linking cable MU USB	GREY3 LINKING CABLE FOR DOCK MU USB w/swivel connector (2m/6ft) w/clip
922127	Grey3 linking cable MU USB	GREY3 LINKING CABLE FOR DOCK MU USB w/swivel connector (7m/23ft) w/clip
922067	White linking cable	WHITE LINKING CABLE 21ft FOR BUILT-IN DOCK SOPRO camera linking cable to build into the arm of the dental unit
922070	Grey linking cable	GREY LINKING CABLE 21ft FOR BUILT-IN DOCK SOPRO camera linking cable to build into the arm of the dental unit
922073	Light grey linking cable	LIGHT GREY LINKING CABLE 21ft FOR BUILT-IN DOCK SOPRO camera linking cable to build into the arm of the dental unit



RELAX. WE'VE GOT YOU COVERED

ACTEON® now offers a One-Year, No-Questions-Asked Warranty on all its devices.* If anything goes wrong within the first year, we will fix it for free.**

* Excludes all consumable products including but not limited to: ultrasonic tips, wrenches, imaging plates, prophylaxis powder, air-water syringes, curing light guides.

** ACTEON reserves the right to change this policy without notice or modify accordingly due to malicious or deliberate damage.

124 Gaither Drive | Suite 140
Mount Laurel | NJ | 08054 | USA
Tel. 1-800-289-6367
Fax. 1-856-222-4726
www.acteongroup.com
info@acteonusa.com

MORE INVENTIVE LESS INVASIVE

