

Product: AIR N GO CLASSIC PEPPERMINT POWDER

SDS No. AC 10-20-1

Created on: 15/01/2011

Updated on: 17/10/2013

1. IDENTIFICATION OF THE SUBSTANCE / MIXTURE AND OF THE COMPANY / UNDERTAKING

1.1 PRODUCT identification

AIR N GO CLASSIC PEPPERMINT POWDER

1.2 Relevant use of mixture

Powder for dental polishing

1.3 Information concerning safety data sheet supplier

SATELEC

Address

Z.I. du Phare - 17 avenue Gustave Eiffel - 33700 MERIGNAC - FRANCE

e-mail

olivier.petit@acteongroup.com

Phone:

+ 33 (0)5 56.34.06.07

Fax:

+ 33 (0)5 56.34.92.92

EMERGENCY CALL NUMBER

<http://echa.europa.eu/web/guest/support/helpdesks/national-helpdesks/list-of-national-helpdesks>

2. HAZARDS IDENTIFICATION

2.1 Classification of substance or mixture

In conformity with the directives 67/548/EEC, 1999/45/EC and adaptations of them

The product is not classified as dangerous in accordance with Directives 67/548/EEC and 1999/45/EC and/or Regulation (EC) 1272/2008 (CLP) and subsequent amendments.

2.2 Label elements

Symbol(s)

None

Indication(s) of danger

Wear mask, gloves and suitable goggles

Risk phrase(s)

Keep out of reach of children

Avoid contact with eyes

Safety advices

In case of contact with eyes, rinse immediately with plenty of

Water for several minutes

Contains

Sodium bicarbonate, Peppermint essential oil

2.3 Other hazards

Handling of this substance can potentially generate dust formation can cause mechanical irritation of eyes, skin, nose and throat

3. COMPOSITION /INFORMATION ON INGREDIENTS

Chemical nature

Powder

Chemical name (Concentration %)				
CAS Nb	EC Nb	Index	Classification 67/548/EEC	Classification 1272/2008/CE (CLP)
Peppermint essential oil (< 1%)				
8006-90-4	308-770-2		Xi R36/38 N R51/53	Skin irritation 2 H315 Skin sensitization 1 H317 Aquatic Chronic 2H411

4. FIRST AID MEASURES

4.1 Description of the first aid measures

Inhalation

Remove to fresh air. In case of unconsciousness place patient stably in side position for transportation. Seek medical advice.

Contact with skin

Remove contaminated clothing. Wash immediately with plenty of water. If irritation persists, seek medical attention.

Contact with eyes

Irrigate copiously with clean, fresh water at least 15 minutes keeping the eyelids opened. Seek medical advice.

Ingestion

Rinse mouth with water. Do not drink, do not induce vomiting, seek medical advice.

4.2 Most important symptoms and effects, both acute and delayed

Not available

4.3 Indications of possible immediate medical attention and particular treatment needed

Not available

Product: AIR N GO CLASSIC PEPPERMINT POWDER

SDS No. AC 10-20-1

Created on: 15/01/2011

Updated on: 17/10/2013

5. FIREFIGHTING MEASURES

5.1 Extinguishing media

Suitable extinguishing media All are suitable

5.2 Special hazards resulting from substance or mixture Avoid breathing vapors

5.3 Advice to firefighters Use self-contained breathing apparatus

6. ACCIDENTAL RELEASE MEASURES

6.1 Individual precautions, protection equipment and emergency procedures

Basic precautions Wear mask, gloves and suitable goggles.

Means of protection Consult section 8.

6.2 Environmental precautions Do not dump into drains.

6.3 Methods and material of confinement and cleaning up

Vacuum or a broom collect spillage. Place everything in a labeled waste container.
Elimination by a licensed waste collection company.

6.4 Reference to other sections

Consult section 8 for personal protection

Consult section 13 for disposal conditions

7. HANDLING AND STORAGE

7.1 Precautions for safe handling

Keep out of reach of children. Avoid contact with skin and eyes.

Do not swallow. Avoid breathing dust.

Avoid release to the environment.

Wash your hands thoroughly after handling.

Wash soiled clothing before reuse.

Comply with conditions of use.

7.2 Conditions for safe storage, including any incompatibilities

Conditions of storage

Keep the product in its original pack, tightly closed and protected from exposure to light.

Keep in a cool place at a temperature between 5 - 25°C.

7.3 Special end use(s)

Product for professional dental use

8. EXPOSURE CONTROLS / PERSONAL PROTECTION

8.1 Control parameters None according our knowledge

8.2 Exposure controls

Respiratory protection Not necessary in normal conditions of use

Hand protection Wear appropriate impervious gloves.

Eye/face protection Wear protective goggles.

Skin and body protection Wear safety clothing.

Wash hands after use

Do not eat, drink or smoke during use.

9. PHYSICAL AND CHEMICAL PROPERTIES

9.1 Information on basic physical and chemical properties

Appearance Powder

Color white

Odour Peppermint

pH Not applicable

Solubility In water

Product: AIR N GO CLASSIC PEPPERMINT POWDER

SDS No. AC 10-20-1

Created on: 15/01/2011

Updated on: 17/10/2013

10. STABILITY AND REACTIVITY

- | | |
|---|--|
| 10.1 Reactivity | Not available |
| 10.2 Chemical stability | Chemically stable under normal conditions of use. |
| 10.3 Possibility of hazardous reactions | No hazardous reactions are foreseeable in normal condition of use and storage. |
| 10.4 Conditions to avoid | Heat, humidity. Do not mix with other products. |
| 10.5 Incompatible materials | None according to our knowledge |
| 10.6 Hazardous decomposition products | None under normal conditions of storage and use |

11. TOXICOLOGICAL INFORMATION

- | | |
|---|--|
| 11.1 Information on toxicological effects | May cause irritation effects similar to other dusts.
Possibility of allergy/sensitization to Peppermint Essential Oil |
|---|--|

12. ECOLOGICAL INFORMATION

- | | |
|--|-----------------------------------|
| 12.1 Toxicity | Not available |
| 12.2 Persistence and degradability | Not available |
| 12.3 Bioaccumulation potential | Not available |
| 12.4 Mobility in the soil | Not available |
| 12.5 Results of PBT and vPvB assessments | Not applicable |
| 12.6 Other adverse effects | Do not release to the environment |

13. DISPOSAL CONSIDERATIONS

- | | |
|------------------------------|---|
| 13.1 Waste treatment methods | Do not dump residue into drain. The chemicals must be eliminated in keeping with applicable regulations.
Packaging of chemical products must be eliminated in keeping with applicable regulations. |
|------------------------------|---|

14. TRANSPORT INFORMATION

- . The product isn't dangerous according the requirements of the regulations applicable to the transportation

15. REGULATORY INFORMATION

- | | |
|--|--|
| 15.1 Special regulations/legislation for substance or mixture in terms of safety, health and environment | Medical device in class IIa according to directive 93/42/EEC |
| 15.2 Chemical safety evaluation | Not available |

16. OTHER INFORMATION

Successive modifications:

Version 1 of 15/01/2011 : Created

Version 2 of 17/10/2013 : Changes according to CLP 1272/2008/EC

Abbreviations : /



SAFETY DATA SHEET
According to regulation (EC) No. 453/2010 annex I

Product: AIR N GO CLASSIC PEPPERMINT POWDER

SDS No. AC 10-20-1

Created on: 15/01/2011

Updated on: 17/10/2013

Recommended restrictions on use

CAUTION: Product intended exclusively for professional dental use

Other information

This sheet complements the technical sheets but does not replace them.

The information that is contained herein is based on the state of our knowledge related to the product concerned at the date of issue and is given in good faith. Moreover, the user's attention is drawn to the possible risks incurred by using the product for any other use than that for which it was intended.

With the present, the user is not dispensed to know and to apply the whole regulating text of activity. He should be responsible for all precautions relating to the use of his product.

The whole regulation prescription mentioned is simply to help the user to fulfill the obligations in the use of hazardous product. This enumeration must not be considered as exhaustive. The user is not exempted to secure than those cited concerning the detention and the manipulation of product for which he is alone responsible.


Olivier PETIT

Quality System Manager

End of the document.