



MATERIAL SAFETY DATA SHEET
According to 91/155/EEC and 2001/58/CE
Product : Classic Powder Lemon
MSDS No AC 10-20-2 ind B 15/01/2011

1. IDENTIFICATION/ PREPARATION OF THE COMPANY / COMPANY

TRADE NAME	Classic Powder Lemon
USE OF THE SUBSTANCE	Powder for Dental Polishing
PREPARATION	
COMPANY	SATELEC
	Z.I. du Phare - BP 30216 33708 MERIGNAC Cedex - FRANCE
Tel:	+ 33 (0)5 56.34.06.07
Fax:	+ 33 (0)5 56.34.92.92
EMERGENCY INFORMATION	+ 33 (0)5 56.34.06.07

2. COMPOSITION / INFORMATION ON INGREDIENTS

PREPARATION

<i>Components</i>	Sodium Hydrogenocarbonate	CAS : 144-55-8	Concentration : < 97 %
	Sodium Saccharinate	CAS : 82385-42-0	Concentration : < 3%
	Lemon Oil	CAS n° 8008-56-8-	Concentration : < 1%%
	Hydrophobic Silica	CAS n 68611-44-9	Concentration : < 0,6%
	Tricalcic Phosphate	CAS n° 7758-87-4	Concentration : < 0,6%

3. HAZARDS IDENTIFICATION

PRINCIPAL HAZARDS

<i>Potential health effects</i>	No specific hazards when the material is handled in accordance with instructions of use.
<i>Physical and chemical hazards</i>	No specific hazards when the material is handled in accordance with instructions of use.
<i>Environment effects</i>	No environment hazards when the material is handled in accordance with official regulations

4. FIRST AID MEASURES

<i>Eyes contact</i>	Wash immediately with running water.
<i>Skin contact</i>	Wash immediately with running water and soap. Remove contaminated clothing.

5. FIRE-FIGHTING MEASURES

Extinguishing media

- *appropriate* All fire-extinguishers are suitable.

Special risks Not considered to be a fire or an explosive hazard.

Individual intervention methods Use any means suitable for extinguishing surrounding fire.

6. ACCIDENTAL RELEASE MEASURES

Personal protection

- *elementary precautions* Wear gloves and appropriated glasses
- *means of protection* Refer to section 8

Environmental protection measure

Do not discard into sewers or streams

Cleaning procedure

Pick up with appropriate products and store for a late disposal according to official regulations.

Eliminate in an authorized center.

7. HANDLING AND STORAGE

HANDLING

Technical measures Keep out of the reach of children.

STORAGE

Technical measures Store in the original packaging

Storage conditions Keep container tightly closed between each use. Store in a dry place.

8. EXPOSURE CONTROLS / PERSONAL PROTECTION

EXPOSURE CONTROLS

Technical measures No particular technical measures if rules of hygiene and safety are respected.

Personal protection

Hands Safety gloves
Eyes Appropriate goggles
Skin Security clothing

Hygiene measures Do not eat, drink or smoke when using

9. PHYSICAL AND CHEMICAL PROPERTIES

APPEARANCE

Physical state Powder
Color White
Odor Lemon
pH Not applicable
Solubility Soluble in Water

10. STABILITY AND REACTIVITY

STABILITY

Stable under ordinary conditions of use and storage.

REACTIVITY

Conditions to avoid Heat, moisture.
Hazardous decomposition None to our knowledge



MATERIAL SAFETY DATA SHEET
According to 91/155/EEC and 2001/58/CE
Product :Classic Powder Lemon
MSDS No AC 10-20-2 ind B 15/01/2011

11. TOXICOLOGICAL INFORMATION

May cause irritation effects similar to other dusts.
Possibility of allergy to Lemon Oil.

12. ECOLOGICAL INFORMATION

ECOTOXICITY DATA

No information found
Do not discard into streams

13. DISPOSAL CONSIDERATIONS

There is no homogeneous regulation within the EC regarding the disposal of chemical products and waste. Dispose of in accordance with local, state and federal regulations. Do not burn empty packaging not decontaminated. Dispose of as hazardous product.

14. TRANSPORT INFORMATION

International regulations

Not subject to transport regulations

15. REGULATORY INFORMATION

Labelling according to EC Directives

Not restricted.

16. OTHERS INFORMATIONS

Product use

Product exclusively reserved for use professional dental use

This data sheet comes in addition to the instructions for use but it does not replace it. The above information is based on our present day knowledge and relates solely to the safety requirements of the product.

The purchaser / user should nevertheless comply with the regulation of the product and its therapeutic use. Should he not comply with it, may the liability of the manufacturer NOT be challenged in any way.

CAUTION: Use of the product for any other purpose than the one it was aimed for is at user's own risks.

Olivier Petit
Quality System Manager