


Filtek™ Z250

Universal Restorative



All-around versatility. 10 years of proven performance.

Summary of advantages

- **Long-term clinical history:** performance is proven by 10 years of clinical service
- **Excellent handling:** resists slumping, and virtually no sticking to instruments
- **Outstanding versatility:** useful for anterior and posterior restorations
- **Speed and convenience:** easy-to-use, direct-delivery capsule cures in just 20 seconds
- **Natural esthetics:** available in 15 shades matching tooth structure
- **Radiopaque**

Filtek™ Z250 Universal Restorative is a trusted restorative solution.


Since it was introduced, dentists around the world have trusted Filtek Z250 restorative to create strong, long-lasting, natural-looking restorations. Any question about how well it performs has been answered by millions of success stories.

The classic in hybrid technology meets a variety of needs with beautiful results.

Filtek Z250 restorative has grown to exceed **35 million restorations globally** per year. That's because dentists know they can consistently achieve excellent results in both anterior and posterior situations, due to the restorative's excellent performance.

- **Esthetics** include a wide range of shades and it's easy to finish and polish.

- **Excellent marginal adaptation** makes it easy to shape.
- **Proven physical properties** include reduced shrinkage for less post-operative sensitivity, low wear and high fracture toughness.




Results of a single notch fracture toughness test proves Filtek™ Z250 Universal Restorative is more fracture-resistant than leading composites tested.

Source: 3M ESPE internal data



Filtek™ Z250 Universal Restorative receives top ratings from The Dental Advisor.

A respected publication dedicated to improving patient care through education and research, recalled and reviewed more than 560 restorations at nine years. The study results proved Filtek Z250 restorative offers remarkable wear resistance and overall strength. In addition, the restorations remained smooth and color-stable to earn a “very good to excellent” rating for esthetics. The study confirms Filtek Z250 restorative is suited for both anterior and posterior use.



Filtek™ Z250 Universal Restorative rating at 9 years.

Source: The Dental Advisor, June 2008

Indications:

- Direct anterior and posterior restorations
- Core buildups
- Splinting
- Indirect restorations including, inlays, onlays, and veneers

Ordering Information



Item # Product Information

REFILLS

| Syringe | Capsule | Shade | Syringe | Capsule | Shade |
|----------|----------|-------|---------|---------|------------------|
| 6020A1 | 6021A1 | A1 | 6020B3 | 6021B3 | B3 |
| 6020A2 | 6021A2 | A2 | 6020C2 | 6021C2 | C2 |
| 6020A3 | 6021A3 | A3 | 6020C3 | 6021C3 | C3 |
| 6020A3.5 | 6021A3.5 | A3.5 | 6020C4 | 6021C4 | C4 |
| 6020A4 | 6021A4 | A4 | 6020D3 | 6021D3 | D3 |
| 6020B0.5 | 6021B0.5 | B0.5* | 6020I | 6021I | Incisal |
| 6020B1 | 6021B1 | B1 | 6020UD | 6021UD | Universal Dentin |
| 6020B2 | 6021B2 | B2 | | | |

Syringe refills contain 1 – 4g syringe. Capsule refills contain 20 – 0.2g capsules.

*3M ESPE shade and not a part of the Vitapan® Classical Shade Guide.

Customer Care Center: 1-800-634-2249 www.3MESPE.com

3M ESPE

Dental Products

3M Center
Building 275-2SE-03
St. Paul, MN 55144-1000
U.S.A.
1-800-634-2249

3M Canada

Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-888-363-3685

Please recycle. Printed in U.S.A.
© 3M 2008. All rights reserved.
70-2013-0158-0

3M, ESPE, Adper and Filtek are trademarks of 3M or 3M ESPE AG. Used under license in Canada. Charisma, XRV, Prodigy, TPH and Tetric Ceram are not trademarks of 3M.