

MICRO-TOUCH[®] Elite[®]

Essential Series

Vinyl, powder-free examination glove for use in low-risk clinical applications

This non-sterile, synthetic glove is produced with a stretchy, polyvinyl chloride formulation providing a more comfortable fit than traditional vinyl gloves. They are an alternative for latex-sensitive individuals.

KEY FEATURES AND BENEFITS:

- Unique formulation for a more comfortable and secure fit than traditional vinyl gloves
- Smooth finish offers enhanced tactile sensitivity
- Protects from latex Type I allergy

Recommended For

- Brief, low-risk examination procedures
- Protection from Type I latex allergy in HCW's or patients

FEATURED TECHNOLOGIES



MICRO-TOUCH[®] Everyday, Everywhere Protection

MICRO-TOUCH[®] premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

MICRO-TOUCH[®] Elite[®]

Essential Series

PRODUCT INFORMATION

Material	Vinyl (Polyvinyl Chloride)	Not Made From Natural Rubber Latex	Yes
Color	Cream	Cuff Length	Standard
Powder Content	Powder-Free	External Glove Surface	Smooth
Freedom From Holes	1.5 AQL	Palm Thickness (mm/mil)	0.11 / 4.3
Finger Thickness (mm/mil)	0.10 / 3.9	Allergy Prevention	Latex (Type I)
Available Sizes	S, M, L, XL	Tested For Use With Chemotherapy Drugs	No
Sterile	No	Antistatic	Not Tested

ORDERING INFORMATION

Size		S	M	L	XL
Ansell Product Code		3091	3092	3093	3094

PACKAGING AND STORAGE

Packaging	100 gloves per box 10 boxes per carton/case
Shelf Life	3 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition

STANDARDS AND CERTIFICATIONS

EN 420:2003 + A1:2009, FDA21 CFR 177-2600, ISO 10993-10:2002, ISO 9001:2008, ISO 13485:2003

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email:
insidesalesus@ansell.com

Latin America & Caribbean
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil Tel: +55-11-3356-3100
Brazil Email:
luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email:
apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email:
protection@ap.ansell.com