

# GAMMEX<sup>®</sup> Non-Latex PI Green

Surgical Gloves

## Polyisoprene, skin-friendly surgical glove ideal for a broad range of surgical procedures

A next-generation, skin-friendly polyisoprene (PI) glove free of chemical accelerators known to cause Type IV glove allergies and sensitivities.<sup>1</sup> This standard non-latex glove delivers all the performance characteristics of natural rubber latex without its ensuing risks. This glove is tested for use with chemotherapy drugs widely used in hospitals.

### KEY FEATURES AND BENEFITS:

- Softer and 'stretchier' PI glove featuring superior comfort compared to Ansell's latex gloves
- Anti-slip, smooth external finish for secure instrument handling and grip
- Tested for use with chemotherapy drugs and a wide range of chemicals<sup>2</sup>

1. Data on file: HRIPT (Human Repeat Insult Patch Testing) conducted by an independent laboratory in the U.S.

2. Contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations.

### Recommended For

- All surgeries, especially high risk procedures in orthopaedics, trauma, and cardiac
- Protection from chemical accelerators such as thiurams, carbamates and mercaptos known to cause chemical Type IV allergies
- Tested for use with chemotherapy drugs
- Having passed the Mouse Embryo Assay (MEA) test, may be used for In Vitro Fertilization (IVF) procedures

### FEATURED TECHNOLOGIES



### GAMMEX<sup>®</sup> Innovation that Protects

Driven by innovation, GAMMEX<sup>®</sup> surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

# GAMMEX<sup>®</sup> Non-Latex PI Green

Surgical Gloves

## PRODUCT INFORMATION

|  |   |
|--|---|
| Material                               | PI (Polyisoprene)                                     |
| Powder Content                         | Powder-Free   |
| Color                                  | Green   |
| Thickness                              | Standard  |
| Cuff Style                             | Beaded with SUREFIT™ Technology                       |
| External Glove Surface                 | Smooth with Micro-Textured Finish                     |
| Internal Glove Surface                 | Polymer Coated with DERMASHIELD™ Technology           |
| Allergy Prevention                     | Latex (Type I)  |
| Grip Level                             | Moderate  |
| Double Gloving Recommendation          | Outer glove or underglove                             |
| Tested For Use With Chemotherapy Drugs | Yes, in accordance with ASTM D6978 and US FDA cleared |

## ORDERING INFORMATION

| Size                | 5.5      | 6        | 6.5      | 7         | 7.5      | 8        | 8.5      | 9        |
|---------------------|----------|----------|----------|-----------|----------|----------|----------|----------|
| Ansell Product Code | 20685255 | 20685260 | 20685265 | 206852570 | 20685275 | 20685280 | 20685285 | 20685290 |

## PACKAGING AND STORAGE

|                      |   |
|----------------------|---|
| Packaging            | 50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case                                    |
| Shelf Life           | 3 years   |
| Storage Instructions | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition |

## STANDARDS AND CERTIFICATIONS

|                                   |  |
|-----------------------------------|--|
| Product Standards                 | ASTM D3577, ASTM D6978, ASTM D7160, ASTM F1670, ASTM F739, EN 420:2003 + A1:2009 |
| Quality / Environmental Standards | EN556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001                          |

Contact your Ansell representative for ordering or more information.

North America  
US Tel: 800 952 9916  
CA Tel: 1-844-494-7854  
Email: insidesalesus@ansell.com

Central & South America  
Tel: +52(442) 296 20 50  
Email: csllac@ansell.com  
Brazil  
Tel: +55-11-3356-3100  
Email: luvas.medicas@ansell.com

Europe, Middle East & Africa  
Tel: +32 (0) 2 528 74 00  
Email: info@ansell.eu

Asia Pacific  
Tel: +603 8310 6688  
Email: apac.medical@ansell.com

Australia & New Zealand  
Tel: +61 3 9270 7270  
Email: protection@ap.ansell.com